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[www.cousin-surgery.com](http://www.cousin-surgery.com)

Made in France

**Caution: Federal law (USA) restricts this device to sale,  
distribution and use by or on the order of a physician.**

This release is the last update of the instructions of use and replace the previous edition

## Sterile spinal dynamic posterior stabilization device

### **1-DESCRIPTION**

The sterile spinal dynamic posterior stabilization device BDyn is intended to restore the stabilization of the non-cervical segment of the spine by preserving the anatomical lordosis and the deadening of the intervertebral joint. It is composed of the BDyn spinal shock absorber fixed on the vertebrae thanks to the dedicated screws and can be replaced according to the decision of the surgeon by a dedicated replacement rigid rod.

### **2-MATERIALS**

The B-DYN Device components are manufactured from medical grade titanium TA6V (ISO 5832-3; ASTM F 136), silicone, long-term implantable polyurethane.

*BDyn implant:* Titanium alloy Ti6Al4V ELI (ISO 5832-3; ASTM F 136), Polydimethyl siloxane (PDMS) MED 4770 unrestricted Polycarbonate urethane (PCU) Bionate® II 80A

*BDyn ancillaries:* Stainless steel, Titanium, Polypropylene

No human or animal origin - Non resorbable.

### **3-INDICATIONS**

The sterile spinal dynamic posterior stabilization device BDyn is intended for posterior stabilization from thoracic vertebrae T10 to sacrum S1 with or without bone graft for the following indications:

- Degenerative intervertebral disc disease and/or articular facets confirmed by further examinations
- Spinal canal stenosis
- Degenerative spondylolisthesis grade 1
- Segmental instability

### **4-PERFORMANCE**

The BDyn spinal shock absorber of the sterile spinal dynamic posterior stabilization device BDyn consists of a metallic cylindrical hollow part containing elastomer components made of silicone elastomers (polydimethyl siloxane (PDMS)) and long-term implantable polycarbonate urethane (PCU) which are bending out under the effect of a metallic piston rod connected with the vertebra of the treated segment by the pedicular screws tested and/or approved by the company COUSIN BIOTECH. The elastomer components assure the absorption of the mechanical loads applied on the intervertebral joint in compression, pulling, flexion-traction and lateral flexion. Some configurations enable fixation across several spinal segments.

### **5-CONTRAINDICATIONS**

- Active infectious process or significant risk of infection (immunocompromise)
- Signs of local inflammation
- Fever or leukocytosis
- Morbid obesity
- Pregnancy
- Mental illness
- Grossly distorted anatomy caused by congenital abnormalities
- Any other medical or surgical condition which would preclude the potential benefit of spinal implant surgery, such as the presence of congenital abnormalities, elevation of the sedimentation rate unexplained by other diseases, elevation of the white blood count
- Suspected or documented metal allergy or intolerance
- Any case where the implant components selected for use would be too large or too small to achieve a successful result
- Any patient having inadequate tissue coverage over the operative site or inadequate bone stock or quality
- Any patient in which implant utilization would interfere with anatomical structures or expected physiological performance
- Any patient unwilling to follow postoperative instructions
- Any case not describe in the indications
- Traumas (i.e. fracture or dislocation)
- Abnormal curvatures (i.e. scoliosis and/or hyper lordosis)
- Tumors.
- Spondylolisthesis grade 2 and more
- Pseudarthrosis and/or failed previous fusion
- Severe bone resorption, osteomalacia, severe osteoporosis

### **6-UNDESIRABLE SIDE EFFECTS**

All of the possible adverse events associated with spinal surgery and without instrumentation are possible:

- Infection
- Pseudomeningocele, fistula, breach dura, persistent CSF leakage, meningitis
- Loss of neurological function, sensorial and/or motor, including complete or incomplete paralysis, dysesthesias, hyperesthesia, anesthesia, paresthesia, appearance of radiculopathy, and/or the development or continuation of pain, numbness, neuroma, spasms, sensory loss, tingling sensation, and/or visual deficits
- *Cauda equina* syndrome, neuropathy, transient or permanent neurological deficits, paraplegia, paraparesis, reflex deficits, irritation, arachnoïditis, and/or muscle loss
- Urinary retention or loss of bladder control or other types of urological system compromise

- Scar formation possibly causing by a neurological compromise or compression around nerves and/or pain
- Fracture, microfracture, resorption, damage or penetration of any spinal bone (including the sacrum, pedicles, and/or vertebral body) and/or bone graft or bone graft harvest site at, above and/or below the level of surgery
- Herniated *nucleus pulposus*, disc disruption or degeneration at, above, or below the level of surgery, canal adjacent stenosis
- Non-union or pseudarthrosis, delayed union. Mal union
- Cessation of any potential growth of the operated portion of the spine
- Loss of or increase in spinal mobility or function
- Inability to perform the activities of daily living
- Bone loss or decrease in bone density
- Graft donor site complications including pain, fracture, or wound healing problems
- Ileus, gastritis, bowel obstruction or loss of bowel control or other types of gastrointestinal system compromise
- Hemorrhage, hematoma, occlusion, seroma, edema, hypertension, embolism, stroke, excessive bleeding, phlebitis, wound necrosis, wound dehiscence, damage to blood vessels, or other types of cardiovascular system compromise
- Reproduction system compromise, sterility, sexual dysfunction
- Development of respiratory problems, e.g. pulmonary embolism, atelectasis, bronchitis, pneumonia, etc
- Change in mental status
- Death

All of the possible adverse events associated with spinal surgery with instrumentation are possible. A listing of potential adverse events linked to the medical device includes, not limited to:

- Early or late loosening of any or all of the components
- Disassembly, bending and/or breakage of any or all of the components (screw breakage)
- Foreign body (allergic) reaction to implants, debris, corrosion products (from crevice, fretting, and/or general corrosion), including metallosis, tumor formation and/or autoimmune disease
- Pressure on the skin from components parts with inadequate tissue coverage over the implant possibly causing skin penetration, irritation, fibrosis, neurosis, and/or pain
- Tissue or nerve damage caused by improper positioning and placement of implants or instruments
- Post-operative change in spinal curvature, loss of correction, height, and/or reduction

## **7-PRECAUTIONS FOR USE**

The BDyn devices are delivered sterile, sterilized by ethylene oxide.

Before any use, inspect the integrity of the packaging and device (including peelable pouches).

Do not use in the event of deterioration of the labels and/or the device and/or the packaging.

Do not use if the device is out of date.

The sterile spinal dynamic posterior stabilization device BDyn is only intended to adult human female or male more than 40 Kg.

COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned.

The installation of BDyn device can only be made with pedicular screws tested and/or approved by the company COUSIN BIOTECH.

The eventual substitution of the BDyn spinal shock absorber has to be carried out only with the dedicated replacement rigid rods.

You should never use implants made of stainless steel and titanium alloy in the same construct.

A successful result is not always achieved in every surgical case. This fact is especially true in spinal surgery where many extenuating circumstances may compromise the results. This device system is not intended to be the sole means of spinal support. No spinal implant can withstand body loads without the support of bone. In this event, bending, loosening, disassembly and/or breakage of the device(s) will eventually occur.

Preoperative and operating procedures, including precise knowledge of suitable surgical techniques, and proper selection of the good reference of the device adapted to the patient and its narrow setting up are important considerations in the successful use of the device by the surgeon.

The BDyn device must be implanted only by a qualified surgeon, having knowledge in the use of the product and who has the knowledge of the anatomy, spinal surgery, pedicle screws fixation technique and specific BDyn device surgical technique

*Postoperative precautions*

### **MRI safety**

The implants are composed of non-ferromagnetic materials and present a geometry non susceptible to generate induced currents. Moreover, as they are fixed to bone /or/ tissues, they are unlikely to be mobilised. A priori they can be considered compatible with an MRI scan. Their safety, in particular in terms of heating and migration of implant has been evaluated through bibliographic data by comparison with data available on devices with similar composition, shape and use. This evaluation concluded to a safety use for MRI scan between 1.5 and 3 Tesla.

As a precautionary measure, it is recommended to avoid MRI scans within the 48h of the implant placement, and to inform the person in charge of the scan of the recent implant placement, if such examination is essential.

It has to be noted that the devices which present a high contrast with the biological environment can generate « artifacts » that has to be taken into account for the perfect execution and interpretation of imaging exams.

For this purpose, it has to be recommended to the patient who has this implant to warn as far as possible the concerned health professionals (radiologists and radiology operators) about the presence of this implant before these exams.

### **IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE**

In accordance with the labelling of this product, BDyn implants are for single use; they must not in any case be reused or re-sterilised (potential risks include, but are not limited to, loss of sterility of the products, risk of infection, loss of product efficacy, recurrence, etc.)

## **8-SURGICAL TECHNIQUE**

The placement of the spinal dynamic posterior stabilization device BDyn begins by the preparation of the pedicles of the vertebrae and by the pedicle aim and the fixation of the pedicle screws.

When tapping the vertebrae to prepare the hole before the insertion of the pedicle screw, flexion efforts on the tap are to be avoided. It must stay in the lined-up in order to make as straightest a pre-hole as possible to ensure as secure a screw placement as possible, with a straight guidance way. Thus, when tapping the vertebrae in order to stay in the line-up and to create a straight pre-hole.

Viewed in the sagittal plane, the pedicle screws must be as parallel as possible in order to let the biggest available space to place the cylindrical body of BDyn spinal shock absorber and facilitate its insertion between the head of the polyaxial pedicle screws. It is preferable to be sure that the adequate implant, with an appropriate size is used with the instrumentation which suit. A 2 mm space between the top of the cylindrical body of the BDyn spinal shock absorber and the lower part of the screw head must be respected to insure the shock absorption stock necessary for the good working of the BDyn device. For that, the surgeon can use the trial device RCBANTD50U to check if there is enough space to place the BDyn spinal shock absorber. The mobile piston rod of the BDyn spinal shock absorber must be fixed in the head of the screw concentrically to the cylindrical body.

The surgery ends by the final tightening by using the dedicated torque limiting of the plugs on the head of the polyaxial pedicle screws by using systematically the suitable antirotator.

It is important to manipulate correctly the implants. The components of the implants to be arched do not have to be against-arched in the opposite direction. In any case they have to be chipped. These operations could effectuate the concentration of internal stresses which could become the place of an eventual failure of the implant.

## **9-PRESENTATION AND WARNING ABOUT ANCILLARIES**

- The BDyn ancillaries are Class I medical devices, intended for temporary use and are re-usable.

Unlike the BDyn implants, the ancillaries included in the instrument set are furnished non sterile.

Before use, the ancillaries must be:

- cleaned by the appropriate mode of cleaning
- sterilized in an autoclave by moist heat

After use, the ancillaries have to be cleaned and decontaminated according to a suitable protocol especially by taking into account the reduction of the risk of transmission of unconventional transmissible agents - UTA.

In case of a patient suspected of having Creutzfeld-Jacob Disease (CJD), perform suitable decontamination.

In the event of a confirmed case of CJD, incinerate the ancillaries.

## **10-DECONTAMINATION AND CLEANING OF BDYN ANCILLARIES**

Decontamination and cleaning are performed under the responsibility of the healthcare centers.

The following prescribed methods and materials must be used to reduce the risk of transmission of UTA, (French health ministry - DGS/R13/2011/449 dated 01/12/2011).

This step has to take place before **the first utilization** and **immediately after the utilization** to avoid adhesion of particles or dry secretion to the instruments

**The detachable instruments have to be dismantled**

**Products advised for the cleaning**

- an enzymatic detergent neutral
- Or products adapted for the cleaning (neutral or soft alkaline) of surgical instruments in compliance with the regulation in effect.

**- Warning:** Don't use corrosive factor or caustic cleaning product

**2 process are possible**

### **A- Chart Automated cleaning in a disinfector washer with manual Pre-cleaning:**

Step	Step Description	Step Instruction	Accessories	Duration
Cleaning Steps	1	Contamination Removal Rinse product room temperature running tap water removing any visible organic material with assistance of a soft bristle brush	- Tap water - Soft bristle brush ( <b>Do not utilize metal cleaning brushes</b> )	Until all visible soil is removed
	2	Prepare a detergent solution at ambient temperature (15-25°C / 59-77°F). Afterwards lumina, threads, joints and gaps of the instruments have to be flushed with a water jet pistol for a minimum time of 10 seconds for each position.	- Detergent - Water jet pistol	Until product is visually clean
	3	Clean the instrument in the ultrasonic tray including the solution. <b>Be careful every unclench connector or screws and bolts by the vibrations will be reassured.</b>	- Ultrasonic tray	Minimum duration 10 minutes, this duration can vary according product

Cleaning Steps	4	Soak	Soak the instrument in tap water during 30 seconds. After, ultrasonic treatment the lumina, threads, joints and gaps of the instrument have to be flushed again with a water jet pistol for a minimum time of 10 seconds for each position.	- Reserve osmosis water	Minimum 30 seconds of soak
	5	Automated Washer	Place entire disassembled device into the automated washer	- Automated Washer - Disassembled Device	Minimum total cycle time: 39 minutes

Instructions of cleaning in washer / disinfecter			
Cycles parameters	Time	Minimum temperature	Type of detergent/water
Pre-cleaning	2 minutes	<45°C / <113°F	Tap water
Cleaning	5 minutes	55°C / 131°F	Detergent
Draining	2 minutes	<45°C / <113°F	Tap water
Thermic draining	5 minutes	90°C / 194°F	Reserve osmosis water
Drying	25 minutes	> 70°C / > 158°F	Not applicable

Finishing Steps	6	Final Rinse	<b>Rinse</b> carefully with <b>distilled water</b>	- Distilled water	Minimum 1 minute
	7	Final Drying	<b>Dry</b> devise utilizing medical quality filtered air	- Medical quality filtered air	Until product is visually dry
	8	Tidying	The instruments have to be placed <b>in the right position</b> inside the trays as described <b>on the position map fixed</b> on the tray for the sterilization.		

**B- Chart Just manual cleaning :**

	Step	Step Description	Step Instruction	Accessories	Duration
Cleaning Steps	1	Contamination Removal	Rinse product room temperature running tap water removing any visible organic material with assistance of a soft bristle brush	- Tap water - Soft bristle brush ( <b>Do not utilize metal cleaning brushes</b> )	Until all visible soil is removed
	2	Drying	Dry the device utilizing a dry non-shedding wipe. Medical quality filtered air may be utilized if available	- Non-shedding wipe - Medical quality filtered compressed air	Until product is visually dry
Disinfection Steps	3	Disinfection Application	Prepare a low-foaming neutral enzymatic detergent solution, using tap water (15-25°C / 59-77°F). Soak the instrument in the open position (if possible) during minimum 1 minute. During the immersion, actuate mobile parts of the instrument minimum 3 times (if possible), so that the detergent can access all parts of the instrument.	- Detergent	Minimum duration 1 minute, this duration can vary according product
	4	Drying	After the soaking time, take out the instrument and wipe it with a disposable cloth. Then, place the instrument in a new bath of enzymatic detergent solution using warm tap water.	- Disposable cloth - Detergent	Until product is visually dry
	5	Manual Disinfection	Brush all surface of the instrument using a soft brush for 2 minutes. Adjust brushing time if needed. Actuate mobile parts of the instrument 3 times (when applicable) and use a syringe to access all difficult areas. Use a volume of 60 mL of the detergent solution.	- Soft bristle brush - Detergent - Syringe	Manual cleaning time duration is complete when the device's surface, joints, and crevices have been manually cleaned
	6	Rinse and finition	Soak the instrument in reverse osmosis water for 1 minute. Use a syringe and 60mL of reverse osmosis water to access all difficult areas. Repeat soaking 2 additional times using fresh water	- Reserve osmosis water - Syringe	Minimum 1 minute
	7	Final Rinse	<b>Rinse</b> carefully with <b>distilled water</b>	- Distilled water	Minimum 1 minute
	8	Final Drying	<b>Dry</b> devise utilizing medical quality filtered air	- Medical quality filtered air	Until product is visually dry
	9	Tidying	The instruments have to be placed <b>in the right position</b> inside the trays as described <b>on the position map fixed</b> on the tray for the sterilization.		

**Note:** In case of **suspicion of Creutzfeld-Jakob Disease (CJD)**, incinerate the ancillaries and the fixation pliers

In case of non respect of the instructions mentioned above, the healthcare center will have to apply a validated equivalent cleaning process for which it will be fully responsible.

## **11-STERILIZATION OF BDYN ANCILLARIES**

Sterilization is performed under the responsibility of the healthcare centers.

The following prescribed method and materials must be used to reduce the risk of transmission of UTA, (French health ministry - DGS/R13/2011/449 dated 01/12/2011).

In case of non respect of the instructions mentioned above, the healthcare center will have to apply a validated equivalent cleaning process for **which it will be fully responsible**.

1- The sterilization in autoclave is to make in a **specific container** (fenced basket, Tray, Plastic tray placed in a peelable pouches autoclavable).

**NB: It isn't recommended** to realize the sterilization with peelables pouches **autoclavables** without a container Ancillaries can pierce the films components the pouches.

2-A sterilization **BY AUTOCLAVE** is advised following this protocol:

Steam sterilization parameters	Values
Cycle Type	Pre-vacuum
Set Point Temperature	134°C / 273°F
Cycle Time	Minimum 3 minutes
Dry Time	Minimum 30 minutes
Cool Time (inside and outside the autoclave chambers)	Minimum 40 minutes

3 - Remark: After each cycle of cleaning/sterilization, insure of all the good working of the ancillaries: verify the integrality and the correct operation, of the locking systems (clipping, prehension) without excessive play.

4- If the ancillaries utilization is different, after sterilization stock them in a tray on a peelable pouches.

In case of non respect of the instructions mentioned above, the healthcare center will have to apply a validated equivalent cleaning process for which it will be fully responsible.

## **12-STORAGE PRECAUTIONS**

- The BDyn devices must be stored in a clean, dry and tempered place under atmospheric pressure, away from sunlight and rays in its original packaging. Ensure that the box is not dropped, hit or crushed.

- The ancillaries must be stored in a dedicated container or equivalent packaging to prevent any deterioration.

## **13-EXPLANTATION AND ELIMINATION OF DEVICES**

Explantation and handling should be done following recommendations of ISO 12891-1:2015 « Implants for surgery - Retrieval and analysis of surgical implants » Part 1: « Retrieval and Handling ».

Any explanted device must be send back, for analysis, following the current protocol. This protocol is available on demand to COUSIN BIOTECH. It is important to note that any implant that should not have been cleaned and disinfected before expedition must be contained in a sealed package.

The elimination of explanted medical device must be conducted in accordance with standards in the country for the disposal of infectious hazards waste.

The elimination of a non-implanted device is not the subject of specific recommendations.

The surgeon will decide the necessity of the resection of the BDyn, its replacement or the rachis fusion according to the known and recognized techniques.

## **14-INFORMATION REQUEST AND CLAIMS**

Following its quality policy, COUSIN BIOTECH is committed to make every effort to produce and supply a high quality medical device. However, if a health professional (client, user, prescriber ...) had a claim or cause of dissatisfaction with a product in terms of quality, safety or performance, he must inform COUSIN BIOTECH as soon as possible. In case of failure of an implant or if it contributed to cause serious adverse effects for the patient, the healthcare center must follow the legal procedures in his country, and inform COUSIN BIOTECH in the shortest time. For any correspondence, thank you to specify the reference, batch number, the coordinates of a reference and a comprehensive description of the incident or claim. Brochures, documentation and surgical technique are available on request from COUSIN BIOTECH and its distributors. If further information is needed or required, please contact your COUSIN BIOTECH representative or distributor.

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## Patient information leaflet BDyn

Item	Patient information leaflet
1	a. <b>The name of the device :</b> BDyn b. <b>The model of the device :</b> RCBDYMD55U, RCBDYSD55U
2	a. <b>The intended purpose of the device :</b> To restore the stabilization of the intervertebral thoracic, lumbar and sacral segments by preserving the anatomical lordosis and the deadening of the intervertebral joint. b. <b>The kind of patient on whom the device is intended to be used :</b> Adult human and female or male more than 40 kg.
3	<b>Special operating instructions for the use of the device :</b> The BDyn device must be implanted only by a qualified surgeon, having knowledge in the use of the product and who has the knowledge of the anatomy, spinal surgery, pedicle screws fixation technique and specific BDyn device surgical technique.
4	a. <b>The intended performance of the device :</b> The BDyn spinal shock absorber of the sterile spinal dynamic posterior stabilization device BDyn consists of a metallic cylindrical hollow part containing elastomer components made of silicone elastomers (polydimethyl siloxane (PDMS)) and long-term implantable polycarbonate urethane (PCU) which are bending out under the effect of a metallic piston rod connected with the vertebra of the treated segment by the pedicular screws tested and/or approved by the company COUSIN BIOTECH. The elastomer components assure the absorption of the mechanical loads applied on the intervertebral joint in compression, pulling, flexion-traction and lateral flexion. Some configurations enable fixation across several spinal segments.

This leaflet was revised in March 2023

Page 58/64

**Headquarters & Factory:** COUSIN BIOTECH 8 Rue Abbé Bonpains - 59117 Wervicq-Sud – France **Tel** +33 (0) 3 20 14 41 20 –  
**Fax** +33 (0) 3 20 14 40 13  
[www.cousin-surgery.com](http://www.cousin-surgery.com) - **e-mail** f.pelletier@cousin-surgery.com

**Item Patient information leaflet**

4	<p>b. <b>Undesirable side effects :</b> Like any implantable medical device, this implant is susceptible to generate possible undesirable side effects that may lead to reoperation: Infection, Pseudomeningocele, fistula, breach dura, persistent CSF leakage, meningitis, loss of neurological function, sensorial and/or motor, including complete or incomplete paralysis, dysesthesias, hyperesthesia, anesthesia, paresthesia, appearance of radiculopathy, and/or the development or continuation of pain, numbness, neuroma, spasms, sensory loss, tingling sensation, and/or visual deficits, cauda equina syndrome, neuropathy, transient or permanent neurological deficits, paraplegia, paraparesis, reflex deficits, irritation, arachnoiditis, and/or muscle loss, urinary retention or loss of bladder control or other types of urological system compromise, scar formation possibly causing by a neurological compromise or compression around nerves and/or pain, fracture, microfracture, resorption, damage or penetration of any spinal bone (including the sacrum, pedicles, and/or vertebral body) and/or bone graft or bone graft harvest site at, above and/or below the level of surgery, herniated nucleus pulposus, disc disruption or degeneration at, above, or below the level of surgery, canal adjacent stenosis, Non-union or pseudarthrosis, delayed union, mal union, Cessation of any potential growth of the operated portion of the spine, Loss of or increase in spinal mobility or function, Inability to perform the activities of daily living, Bone loss or decrease in bone density, Graft donor site complications including pain, fracture, or wound healing problems, Ileus, gastritis, bowel obstruction or loss of bowel control or other types of gastrointestinal system compromise, hemorrhage, hematoma, occlusion, seroma, edema, hypertension, embolism, stroke, excessive bleeding, phlebitis, wound necrosis, wound dehiscence, damage to blood vessels, or other types of cardiovascular system compromise, Reproduction system compromise, sterility, sexual dysfunction, Development of respiratory problems, e.g. pulmonary embolism, atelectasis, bronchitis, pneumonia, etc, Change in mental status, Death</p>
5	<p><b>When BDyn Should Not be Used :</b> Active infectious process or significant risk of infection (immunocompromise), signs of local inflammation, Fever or leukocytosis, Morbid obesity, Pregnancy, Mental illness, Grossly distorted anatomy caused by congenital abnormalities, Any other medical or surgical condition which would preclude the potential benefit of spinal implant surgery, such as the presence of congenital abnormalities, elevation of the sedimentation rate unexplained by other diseases, elevation of the white blood count Suspected or documented metal allergy or intolerance, Any case where the implant components selected for use would be too large or too small to achieve a successful result, Any patient having inadequate tissue coverage over the operative site or inadequate bone stock or quality, Any patient in which implant utilization would interfere with anatomical structures or expected physiological performance, Any patient unwilling to follow postoperative instructions, Any case not describe in the indications, Traumas (i.e. fracture or dislocation), Abnormal curvatures (i.e. scoliosis and/or hyper lordosis), Tumors, Spondylolisthesis grade 2 and more, Pseudarthrosis and/or failed previous fusion, Severe bone resorption, osteomalacia, severe osteoporosis</p>

This leaflet was revised in March 2023

Page 59/64

**Headquarters & Factory:** COUSIN BIOTECH 8 Rue Abbé Bonpains - 59117 Wervicq-Sud – France **Tel** +33 (0) 3 20 14 41 20 – **Fax** +33 (0) 3 20 14 40 13  
[www.cousin-surgery.com](http://www.cousin-surgery.com) - **e-mail** f.pelletier@cousin-surgery.com

**Item    Patient information leaflet**

6	<b>MRI Information :</b> The implants are composed of non-ferromagnetic materials and present a geometry non susceptible to generate induced currents. Moreover, as they are fixed to bone /or/ tissues, they are unlikely to be mobilised. A priori they can be considered compatible with an MRI scan. Their safety, in particular in terms of heating and migration of implant has been evaluated through bibliographic data by comparison with data available on devices with similar composition, shape and use. This evaluation concluded to a safety use for MRI scan between 1.5 and 3 Tesla. As a precautionary measure, it is recommended to avoid MRI scans within the 48h of the implant placement, and to inform the person in charge of the scan of the recent implant placement, if such examination is essential. It has to be noted that the devices which present a high contrast with the biological environment can generate « artifacts » that has to be taken into account for the perfect execution and interpretation of imaging exams. For this purpose, it has to be recommended to the patient who has this implant to warn as far as possible the concerned health professionals (radiologists and radiology operators) about the presence of this implant before these exams.
7	<b>Precautions :</b> The BDyn devices are delivered sterile, sterilized by ethylene oxide. Before any use, inspect the integrity of the packaging and device (including peelable pouches). Do not use in the event of deterioration of the labels and/or the device and/or the packaging. Do not use if the device is out of date. The sterile spinal dynamic posterior stabilization device BDyn is only intended to adult human female or male more than 40 Kg. COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned. The installation of BDyn device can only be made with pedicular screws tested and/or approved by the company COUSIN BIOTECH. The eventual substitution of the BDyn spinal shock absorber has to be carried out only with the dedicated replacement rigid rods. You should never use implants made of stainless steel and titanium alloy in the same construct. A successful result is not always achieved in every surgical case. This fact is especially true in spinal surgery where many extenuating circumstances may compromise the results. This device system is not intended to be the sole means of spinal support. No spinal implant can withstand body loads without the support of bone. In this event, bending, loosening, disassembly and/or breakage of the device(s) will eventually occur. Preoperative and operating procedures, including precise knowledge of suitable surgical techniques, and proper selection of the good reference of the device adapted to the patient and its narrow setting up are important considerations in the successful use of the device by the surgeon. The BDyn device must be implanted only by a qualified surgeon, having knowledge in the use of the product and who has the knowledge of the anatomy, spinal surgery, pedicle screws fixation technique and specific BDyn device surgical technique Postoperative precautions

This leaflet was revised in March 2023

Page 60/64

**Headquarters & Factory:** COUSIN BIOTECH 8 Rue Abbé Bonpains - 59117 Wervicq-Sud – France **Tel** +33 (0) 3 20 14 41 20 –  
**Fax** +33 (0) 3 20 14 40 13  
[www.cousin-surgery.com](http://www.cousin-surgery.com) - **e-mail** f.pelletier@cousin-surgery.com

Item Patient information leaflet	
<b>8</b>	<p><b>Materials included in the device :</b> The BDyn Device components are manufactured from medical grade titanium TA6V (ISO 5832-3; ASTM F 136), silicone, long-term implantable polyurethane.</p> <p><b>BDyn implant:</b> Titanium alloy Ti6Al4V ELI (ISO 5832-3; ASTM F 136), Polydimethyl siloxane (PDMS) MED 4770 unrestricted Polycarbonate urethane (PCU) Bionate® II 80A</p>
<b>9</b>	<p><b>Reporting a serious incident :</b></p> <p>If you encounter any serious incident with your Bdyn implant you should report it immediately to your doctor. You should also report the incident to COUSIN BIOTECH via the link <a href="mailto:serviceclients@cousin-surgery.com">serviceclients@cousin-surgery.com</a> or <a href="mailto:materiovigilance@cousin-biotech.com">materiovigilance@cousin-biotech.com</a></p> <p>and to the Therapeutic Goods Administration using the following website: <a href="http://www.tga.gov">http://www.tga.gov</a></p>
<b>10</b>	<p><b>Legal Manufacturer &amp; Manufacturing Facility :</b> COUSIN BIOTECH, 8 Rue Abbé Bonpains, 59117 Wervicq-Sud, France</p>

This leaflet was revised in March 2023

Page 61/64

**Headquarters & Factory:** COUSIN BIOTECH 8 Rue Abbé Bonpains - 59117 Wervicq-Sud – France **Tel** +33 (0) 3 20 14 41 20 –  
**Fax** +33 (0) 3 20 14 40 13  
[www.cousin-surgery.com](http://www.cousin-surgery.com) - **e-mail** f.pelletier@cousin-surgery.com

## NOTES:

## **NOTES :**

	en	CE mark and identification number of Notified Body. Product conforms to essential requirements of the Medical device directive 93/42/EEC
	fr	Marque CE et numéro d'identification de l'organisme notifié. Produit conforme aux exigences essentielles de la directive 93/42/CEE
	de	CE-Kennzeichnung und Identifikationsnummer der benannten Stelle. Produkt entspricht grundlegenden Anforderungen der Richtlinie für Medizinische Geräte 93/42/EWG
	it	Marchio CE e numero di identificazione dell'ente notificato. Prodotto conforme ai requisiti essenziali della direttiva sui dispositivi medici 93/42/EEC
	es	Marca CE y número de identificación del organismo notificado. El producto está conforme con los requerimientos de la directiva 93/42/EEC sobre aparatos médicos.
	pt	Marca CE e número de identificação de Organismo Notificado. O produto obedece aos principais requisitos da Diretiva 93/42/CEE
	nl	CE-markering en identificatienummer van de aangemelde instantie. Product voldoet aan de essentiële eisen van de richtlijn 93/42/EEG
	pl	Znak CE i numer identyfikacyjny uprawnionego organu. Produkt zgodny z zasadniczymi wymaganiami dyrektywy 93/42/CEE
	ru	Знак CE и идентификационный номер аккредитованного органа по сертификации. Продукция соответствует основным требованиям директивы по изделиям медицинского назначения 93/42/EEC
ελ	Σήμανση CE και αριθμός αναγνώρισης του κοινοποιημένου φορέα. Το προϊόν συμμορφώνεται με τις βασικές απαιτήσεις της Οδηγίας 93/42/EEC για τα αιστροτεχνολογικά προϊόντα	

en	- Symbols used on labelling	pt	- Símbolos usados na etiqueta		<a href="http://www.cousin-biotech.com/ifu">www.cousin-biotech.com/ifu</a>
fr	- Symboles utilisés sur l'étiquette	nl	- Op de etikettering gebruikte symbolen		
de	- Bei Etiketten verwandete Symbole	pl	- Symbole użyte na etykiety		
it	- Simboli utilizzati sull'etichetta	ru	- Символы, используемые на маркировке		
es	- Símbolos utilizados en el etiquetaje	ελ	- Σύμβολα που χρησιμοποιούνται στις ετικέτες		IFU can be downloaded through the QR code. A hard copy can be sent within 7 days on request by email to <a href="mailto:ifurequest@cousin-biotech.com">ifurequest@cousin-biotech.com</a> or by using the order form on our website.

	LOT	REF		
en	Batch number	Reference on the brochure	Caution (See instructions for use)	Manufacturer
fr	Numéro de lot	Référence du catalogue	Attention, voir notice d'instructions	Fabricant
de	Chargennummer	Verweis auf der Broschüre	Vorsicht (Siehe Gebrauchsanweisung)	Hersteller
it	Número di lotto	Referimento del catalogo	Attenzione (Vedere manuale istruzioni)	Produttore
es	Número de lote	Referencia al folleto	Precaution (Veanse las instrucciones de uso )	Fabricante
pt	Número do lote	Referência do catálogo	Atenção (Consultar nota de instruções)	Fabricante
nl	Batchnummer	Referentie in de brochure	Vorsichtigkeit (Zie gebruiksinstructies)	Fabrikant
pl	Numer partii	Symbol katalogowy	Uwaga, patrz instrukcja obsługi	Producent
ru	Номер партии	Ссылка на брошюру	Смотрите инструкцию по применению	Производитель
ελ	Αριθμός παρτίδας	Αριθμός καταλόγου	Δείτε τις οδηγίες χρήσης	Κατασκευαστής

	<p>en Do not re-use      fr Ne pas réutiliser      de Nicht wiederverwenden      it Non riutilizzare      es No volver a utilizar      pt Não reutilizar      nl Niet hergebruiken      pl Nie stosować ponownie      ru Не предназначен для повторного использования      ελ Να μην επαναχρησιμοποιείται</p>		<p>en Use before: year and month      fr Utiliser jusque : année et mois      de Verbrauchen bis : Jahr und Monat      it Utilizzare entro e non oltre: anno e mese      es Úsese antes de: año y mes      pt Utilizar até : ano e mês      nl Te gebruiken vóór: jaar en maand      pl Stosować do: rok i miesiąc      ru Срок годности: год и месяц      ελ Ημερομηνία λήξης: έτος και μήνας</p>
	<p>en Keep in a dry place      fr Conserver au sec      de Trocken lagern      it Conservare in un luogo asciutto      es Manténgase en un lugar seco      pt Manter em lugar seco      nl Op een droge plaats bewaren      pl Przechowywać w suchym miejscu      ru Хранить в сухом месте      ελ Φυλάξτε το προϊόν μακριά από την υγρασία</p>		<p>en Do not use if the packaging is damaged      fr Ne pas utiliser si l'emballage est endommagé      de Nicht verwenden, wenn Packung beschädigt ist      it Non utilizzare se la confezione è danneggiata      es No usar si el embalaje está dañado      pt Não usar se a embalagem estiver danificada      nl Niet gebruiken indien de verpakking beschadigd is      pl Nie używać, jeśli opakowanie jest uszkodzone      ru Не использовать в случае повреждения упаковки      ελ Να μην χρησιμοποιείται εάν η συσκευασία έχει υποστεί κάποια βλάβη</p>
	<p>en Do not re-sterilize      fr Ne pas restériliser      de Nicht erneut sterilisieren      it Non risterilizzare      es No volver a esterilizar      pt Não reesterilizar      nl Niet opnieuw steriliseren      pl Nie sterylizować ponownie      ru Не подлежит повторной стерилизации      ελ Να μην επαναστεριώνεται</p>		<p>en Keep away from sunlight      fr A stocker à l'abri de la lumière du soleil      de Vor Sonnenlicht geschützt lagern      it Tenere lontano dalla luce      es Manténgase oscuras      pt Manter ao abrigo da luz      nl Vrij van licht houden      pl Przechowywać w miejscu chronionym przed światłem      ru Хранить в темном месте вдали от солнечных лучей      ελ Προφύλαξτε το προϊόν από το φως</p>

<b>STERILE</b>	R : Polyaxial pedicle screws, Replacements rigid rods	<b>STERILE</b>	EO : BDyn spinal shock absorber
en	Sterilized by gamma radiation.	en	Sterilized by ethylene oxide
fr	Produit stérile. Méthode de stérilisation : irradiation	fr	Produit stérile. Méthode de stérilisation : oxyde d'éthylène
de	Sterilisierungsmethode : Gammastrahlung.	de	Mit Ethylenoxid sterilisiert
it	Il prodotto è sterile. Sterilizzazione : raggi Gamma.	it	Sterilizzato all'ossido di etilene
es	Esterilizado por rayos gamma.	es	Esterilizado con óxido de etileno
pt	Produto estéril. Método de esterilização: irradiação	pt	Produto estéril. Método de esterilização: óxido de etileno
nl	Sterilizováno gamma zářením	nl	Gesteriliseerd met ethylenoxide
pl	Sterylizowane przy użyciu promieni gamma.	pl	Sterylizowany tlenkiem etylenu
ru	Стерилизовано гамма-излучением.	ru	Стерилизовавшись с помощью окиси этилена
ελ	Στείρω προϊόν. Μέθοδος αποστείρωσης: ακτινοβολία	ελ	Στείρω προϊόν. Μέθοδος αποστείρωσης: αιθυλενοξείδιο