STERILE SELF-LOCKING POLYMER CERCLAGE SYSTEM

en	Instructions for use	Page	2
fr	Instructions d'utilisation	Page	9
de	Gebrauchsanweisung	Seite	16
es	Instrucciones de uso	Pagina	23
it	Istruzioni per l'uso	Pagina	30
fi	Käyttöohjeet	Sivu	37
sv	Bruksanvisning	Sida	44

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NOT278/230627

Version du 27/06/2023

Caution : Federal law (USA) restricts this device to sale, distribution and use by or on the order of a physician.

This release is the last update of the instructions of use and replace the previous edition

STERILE SELF-LOCKING POLYMER CERCLAGE SYSTEM

SINGLE USE STERILE PRODUCT

1-DESCRIPTION

The "Sterile self-Locking Polymer Cerclage system" consists of a braided cerclage polymer cable with metal connector plate and self-locking mechanism.

2- MATERIALS

For the prosthesis

- Cable: polyethylene terephthalate multifilament braided cable, non resorbable, implanted
- Connector: Titanium TA6V, non resorbable, implanted
- Semi-rigid end : polyethylene terephthalate multifilament braided cable, non resorbable, not implanted

For the ancillary

- Tensioner : Stainless steel

Origin neither human nor animal - Biocompatible. Non resorbable.

3-INDICATIONS

- Repair of long bone fractures due to trauma or reconstruction;
- Reattachment of the greater trochanter in total hip arthroplasty, surface replacement arthroplasty, or other procedures involving trochanteric osteotomy

4-PERFORMANCE

The role of the cerclage device (with or without bone plate) is to hold the fragments of the bone in position till the bone heals.

In details, it works according three sequences:

- Ease of cerclage device setting around the bone (soft and conformable materials)
- Ease of hand pre-setting during initial fracture reduction to stabilise multi-fragmented bone elements
- Reliable tightening of the bone during healing process

5-CONTRAINDICATIONS

- Skeletally immature patients;
- Active or suspected infection, either systemic or localized, in or around the implant site;
- Patient conditions, mental or neurological, that would tend to impact the patient's ability to follow physician's instructions during the post-operative healing phase;
- Compromised vascularity that would inhibit adequate blood supply to the fracture or operative site;
- Insufficient quality or quantity of bone, which would inhibit rigid device fixation;
- Any disease affecting the support and fixation of the prosthesis;
- Distant foci of infection, such as genitourinary, pulmonary, skin or other sites which may spread to the implant site. The foci of infection should be treated prior to, during, and after implantation

6-UNDESIRABLE SIDE EFFECTS

- Cable failure, early or late loosening, wear, and change in position of one or more components
- Malunion of fracture / fracture displacement due to inadequate reattachment and/or early weight bearing;
- Progressive bone resorption and osteolysis
- Infection
- Neurological injury
- Vascular injuries due to the passer or the cable
- Foreign body reaction
- Allergy to a component of the implant
- Migration of the cerclage material
- Pain

7-POST-OPERATIVE PROCEDURES

- The patient should be released from the hospital with complete instructions and warnings regarding further exercises and therapies coming from medical staff. The patient should be encouraged to report any unusual changes to his/her physician.
- The patient should be cautioned to govern his/her activity level. Postoperative therapy should be structured to prevent excessive loading.

8-PRECAUTIONS

- The patient must be advised of the short- and long-term limitations of the procedure and the need to protect the implant from full weight bearing until adequate fixation and healing have occurred. Excessive activity and trauma affecting the reattachment and/or fixation have been implicated in failure of this procedure by loosening, fracture, and/or wear of the implants. Loosening of the components can result in increased production of wear particles, as well as damage to the bone, making successful revision surgery more difficult. The patient should be cautioned to limit activities, protect the bone from unreasonable stresses, and to follow the instructions of the physician with respect to follow-up care and treatment. The patient should be warned of the surgical risks and made aware of possible adverse effects. The patient should be warned that the device cannot and will not replicate the flexibility, strength, reliability or durability of normal healthy bone, that the implant can break or become damaged, particularly as a result of strenuous activity or trauma, and that the device has a finite expected service life and may need to be replaced in the future.
- Transient bacteremia can occur in daily life. Dental manipulation, endoscopic examination, and other minor surgical procedures have been associated with transient bacteremia. To prevent infection at the implant site, it may be advisable to use antibiotic prophylaxis before and after such procedures.
- An implant should never be re-used. Even though an implant may appear undamaged, previous handling and in-service stresses may have created imperfections that would reduce the service life of the implant.
- Patients should be questioned regarding sensitivity to metals. If there is a question pertaining to the patient's tolerance for titanium, appropriate testing should be performed
- The following conditions, singularly or concurrently, tend to impose severe loading on the affected bone, thereby placing the patient at higher risk for failure: obesity, heavy labor, active sports participation, high general activity level, likelihood of falls, alcohol or drug addiction, and other disabilities
- Cerclage cable breakage is reported in the literature. Possible causes of breakage include, but are not limited to, over tensioning the cable, fretting, and delayed or non-union of the fracture, osteotomy, or fusion site. In many instances, adverse effects may be clinically related rather than implant related. The surgical and post-operative management of the patient must be carried out with due consideration for all existing conditions of increased risk. Mental attitudes or personality disorders resulting in patient's failure to adhere to his/her physician's orders might delay post-operative recovery and aggravate adverse effects.
- Over-tensioning can create the potential for the cable to cut into the bone. Under-tensioning can cause inadequate fixation and increased fatigue stress.
- Improper placement, positioning and fixation of these devices may result in unusual local inflammation. Proper handling of any implant is mandatory.
- Use care in utilizing the Cable Passer and other instrumentation to avoid damaging neurovascular structures. Exercise caution when positioning and tightening the Cerclage Cables to avoid impinging or tethering neurovascular structures and causing attendant significant sequelae
- Prior to surgical use, a visual inspection of each implant for possible imperfections should be performed. Damaged or altered implants may produce undesirable stresses and cause defects, which could lead to failure.
- This device is delivered sterile (gamma sterilisation). Before any use, inspect the integrity of packaging and device (of which pelables sachets). Do not use in the event of deterioration of the device and/or the packaging.
- Do not use if the device is out of date.
- This device must be implanted only by a qualified surgeon trained on the use of the product (knowledge of anatomy and orthopaedic surgery). The surgeon is to be thoroughly familiar with the surgical procedure, instruments and implant characteristics prior to performing surgery.

MRI SAFETY:

The implants are composed of non ferromagnetic materials and present a geometry non susceptible to generate induced currents. Moreover, as they are fixed to bone/or /tissues, they are unlikely to be mobilised. A priori they can be considered compatible with an MRI scan. Their safety, in particular, in terms of heating and migration of implant has been evaluated through bibliographic data by comparison with data available on devices with similar composition, shape and use. This evaluation concluded to a safety use for MRI scan between 1.5 and 3 Tesla.

As a precautionary measure, it is recommended to avoid MRI scans within the 48h of the implant placement, and to inform the person in charge of the scan of the recent implant placement, if such examination is essential.

It has to be noted that the devices which present a high contrast with the biological environment can generate « artifacts » that has to be taken into account for the perfect execution and interpretation of imaging exams.

For this purpose, it has to be recommended to the patient who has this implant to warn as far as possible the concerned health professionals (radiologists and radiology operators) about the presence of this implant before these exams.

WARNING ABOUT ANCILLARIES

The Self locking polymer cerclage system ancillaries are Class I medical devices, intended for temporary use and are re-usable. Unlike the Self locking polymer cerclage system prosthesis, ancillaries are furnished unsterile.

Before first use

Ancillaries must be:

- Cleaned in the appropriate way
- Sterilised in an autoclave

DECONTAMINATION AND CLEANING OF ANCILLARIES

Decontamination and cleaning are performed under the responsibility of the healthcare centers.

The following prescribed methods and materials must be used to reduce the risk of transmission of Non Conventional Transmissible Agent (NCTA), (French health ministry - DGS/R13/2011/449 dated 01/12/2011).

This step has to take place before **the first utilization** and **immediately after the utilization** to avoid adhesion of particles or dry secretion to the instruments.

Products advised for the cleaning

- an enzymatic detergent neutral
- Or products adapted for the cleaning (neutral or soft alkaline) of chirurgial instruments in compliance with the regulations in effect.

- **Warning:** Don't use corrosive factor or caustic cleaning product

2 process are possible

A- Chart Automated cleaning in a disinfectant washer with manual Pre-cleaning:

	Step	Step Description	Step Instruction	Accessories	Duration
Cleaning Steps	1	Contamination Removal	Rinse product room temperature running tap water removing any visible organic material with assistance of a soft bristle brush	<ul style="list-style-type: none"> - Tap water - Soft bristle brush <p>(Do not utilize metal cleaning brushes)</p>	Until all visible soils removed
	2	Pre-soak	Prepare a detergent solution at ambient temperature (15-25°C / 59-77°F). Afterwards lumina, threads, joints and gaps of the instruments have to be flushed with a water jet pistol for a minimum time of 10 seconds for each position.	<ul style="list-style-type: none"> - Detergent - Water jet pistol 	Until product is visually clean
	3	Ultrasonic cleaning	Clean the instrument in the ultrasonic tray including the solution.	<ul style="list-style-type: none"> - Ultrasonic tray 	Minimum duration 10 minutes, this duration can vary according product
	4	Soak	Soak the instrument in tap water during 30 seconds. After, ultrasonic treatment the lumina, threads, joints and gaps of the instrument have to be flushed again with a water jet pistol for a minimum time of 10 seconds for each position.	<ul style="list-style-type: none"> - Reverse osmosis water 	Minimum 30 seconds of soak
	5	Automated Washer	Place device into the automated washer	<ul style="list-style-type: none"> - Automated Washer 	Minimum total cycle time: 39 minutes

Instructions of cleaning in washer / disinfectant			
Cycles parameters	Time	Minimum temperature	Type of detergent/water
Pre-cleaning	2 minutes	<45°C / <113°F	Tap water
Cleaning	5 minutes	55°C / 131°F	Detergent
Draining	2 minutes	<45°C / <113°F	Tap water
Thermic draining	5 minutes	90°C / 194°F	Reverse osmosis water
Drying	25 minutes	> 70°C / > 158°F	Not applicable

Finishing Steps	6	Final Rinse	Rinse carefully with distilled water	- Distilled water	Minimum 1 minute
	7	Final Drying	Dry device utilizing medical quality filtered air	- Medical quality filtered air	Until product is visually dry
	8	Tidying	The instruments have to be placed in a box (non specific) for the sterilization.		

B- Chart just manual cleaning :

	Step	Step Description	Step Instruction	Accessories	Duration
Cleaning Steps	1	Contamination Removal	Rinse product room temperature running tap water removing any visible organic material with assistance of a soft bristle brush	- Tap water - Soft bristle brush (Do not utilize metal cleaning brushes)	Until all visible soils removed
	2	Drying	Dry the device utilizing a dry non-shedding wipe. Medical quality filtered air may be utilized if available	- Non-schedding wipe - Medical quality filtered compressed air	Until product is visually dry
Disinfection Steps	3	Disinfection Application	Prepare a low-foaming neutral enzymatic detergent solution, using tap water (15-25°C / 59-77°F). Soak the instrument in the open position (if possible) during minimum 1 minute. During the immersion, actuate mobile parts of the instrument minimum 3 times (if possible),so that the detergent can access all parts of the instrument.	- Detergent	Minimum duration 1 minute, this duration can vary according product
	4	Drying	After the soaking time, take out the instrument and wipe it with a disposable cloth. Then, place the instrument in a new bath of enzymatic detergent solution using warm tap water.	- Disposable cloth - Detergent	Until product is visually dry
	5	Manual Disinfection	Brush all surface of the instrument using a soft brush for 2 minutes. Adjust brushing time if needed. Actuate mobile parts of the instrument 3 times (when applicable) and use a syringe to access all difficult areas. Use a volume of 60 mL of the detergent solution.	- Soft bristle brush - Detergent - Syringe	Manual cleaning time duration is complete when the device's surface, joints, and crevices have been manually cleaned
	6	Rinse and finition	Soak the instrument in reverse osmosis water for 1 minute. Use a syringe and 60mL of reverse osmosis water to access all difficult areas. Repeat soaking 2 additional times using fresh water	- Reserve osmosis water - Syringe	Minimum 1 minute

Disinfection Steps	7	Final Rinse	Rinse carefully with distilled water	- Distilled water	Minimum 1 minute
	8	Final Drying	Dry devise utilizing medical quality filtered air	- Medical quality filtered air	Until product is visually dry
	9	Tidying	The instruments have to be placed in a box (non specific) for the sterilization		

Note: In case of **suspicion** of **Creutzfeld-Jakob Disease (CJD)**, incinerate the ancillaries and the fixation pliers
 In case of no respect of the instructions mentioned above, the healthcare center will have to apply a validated equivalent cleaning process for which it will be fully responsible.

STERILIZATION OF ANCILLAIRES

Sterilization is performed under the responsibility of the healthcare centers.
 The following prescribed methods and materials must be used to reduce the risk of transmission of Non Conventional Transmissible Agent (NCTA), (French health ministry - DGS/R13/2011/449 dated 01/12/2011).

1- The sterilization in autoclave is to make in a **non-specific container** (fenced basket, Tray, Plastic tray placed in a peelable pouches autoclavable).

NB: It **isn't recommended** to realize the sterilization with peelables pouches **autoclavable** without a container. Ancillaries can pierce the films components the pouches.

2-A sterilization **BY AUTOCLAVE** is advised following this protocol:

Steam sterilization parameters	Values	Values
Cycle Type	Pre-vacuum	Pre-vacuum
Set Point Temperature	134°C / 273°F	132°C / 269.6°F
Cycle Time	Minimum 3 minutes	Minimum 4 minutes
Dry Time	Minimum 20 minutes	Minimum 20 minutes
Cool Time (inside and outside the autoclave chambers)	Minimum 40 minutes	Minimum 40 minutes

3 -Remark: After each cycle of cleaning/sterilization, insure of all the good working of the ancillaries: verify the integrity and the correct operation without excessive play.

4- If the ancillaries utilization is different, after sterilization stock them in a tray on a peelable pouches.
 In case of no respect of the instructions mentioned above, the healthcare center will have to apply a validated equivalent sterilization process for which it will be fully responsible.

STORAGE PRECAUTIONS

Ancillaries must be stored in their container or equivalent packaging to prevent any deterioration.
 Keep in a dry place away from sunlight, at room temperature, in its original packaging.

9-IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE

As specified on the product's labeling, the device is for single use only. It can not be re-used and/or re-sterilised (potential risks would be and are not limited to loss of the product's sterility, to risk of infection, to loss of the product's efficiency, recurrence)

10-DESCRIPTION OF THE SURGICAL PROCEDURES:

STEP 1. CERCLAGE SETTING

It is important to note the direction of introduction of the Cable Passer around the bone as this affects the orientation of the cable and thus the orientation of the tensioning instrument relative to the user and the incision. Position the appropriate Cable Passer around the bone such that the distal end of the cannula emerges on the operating surgeon's side of the bone.

Prior to Cerclage fixation, Remove carefully soft tissues around bone (or bone fragments) to secure a direct hard contact between polymer cerclage device and bone surface.

Pass the Cable Passer around the bone from posterior to anterior. Insert the semi-rigid end into the tip of the Cable Passer until the wire is seen exiting from the shaft of the Cable Passer

Use care in passing the cable passer to avoid damage to soft tissues, bone surface and neurovascular structures (including muscle insertion line)

Remove carefully Cable Passer while the semi-rigid end remains in place, wrapped around the bone.

Wrap Polymer Cerclage (connector plate side) all the way through around the bone.

Remove entirely the semi-rigid end.

Pull through completely the tensioning end.

Avoid wrapping the cables over sharp metal or bone graft surfaces.

Pass the fixation loop into the connector. The auto-fix mechanism of the connector will keep the connection while moving or rotating the cable.

Maintain cerclage wires parallel. If twisted wires, unlock the connector, untwist wires and restart connecting to the connector.

STEP 2. HAND PRE-SETTING

Fracture Reduction and Hand pre-setting

For fracture reduction, polymer cerclage must be used complementary to bone clamps. Bone clamps will address proper bone alignment and reduction.

If needed, several cables could be used. The order of cable insertion depends upon the fracture pattern, particular use or situation. The cables on each end of the plate may be wrapped around the bone to stabilize the bone while the final reduction is carried out.

The cables can be loosely tensioned (pre-setting) by hand as the self-locking mechanism will hold the tension.

The cable can be softly hand re-tight as required.

Be aware of sharp age bone structure that could damage cerclage system.

STEP 3. FINAL TIGHTENING

Insert the tensioning-end into the Tensioner (first through the front nose then second into the rotational axis)

Position Tensioner in direct contact with buckle eyelet for stability while tensioning.

Note: for optimal tensioning, it is recommended to maintain tensioner well aligned with wire splice (reduce angle).

Rotate manually the Tensioner axis until desired compression load is achieved.

Self-locking Mechanism will hold tension to allow re-tightening as desired.

To release the Tensioning instrument, De-Rotate (counterclockwise) the axis. That will allow to verify compression load around the bone.

As the self-locking mechanism hold the Tension, you can re-tight the cerclage as needed.

NOTE: We recommend 3 Times Tightening process to ensure stable desired compression load

Do not cut the tensioning end until complete operation is achieved.

After removal of the tensioning instrument to cut off at 2 cm length.

NOTE: The cable may be placed over smooth metal implants, such as bone plates.

NOTE: Complementary to bone plate, clinical literature recommends to adjust minimum of two cerclages devices (at both extremities) before cerclage final tensioning. For long bone fracture, it is described up to four cerclages to improve stability and secure healing process.

NOTE: With bone plate with screws, it is recommended to check cerclage final tensioning (tightening) after screws fixation and retight if it's necessary. This sequence will allow full operational performance for cerclage.

NOTE: Considering its specific design and for optimal effectiveness, it's recommended to place the polymer cerclage at the exterior surface of the bone plate however the cerclage could be placed under the bone plate avoiding proximity to screw drilling and the final tightening of the cerclage must be done prior to the final tightening of the bone plate. The titanium connector mustn't be put under the bone plate.

NOTE: Damaged instruments may lead to improper cable tension or implant position and result in implant failure. Thorough familiarity with the surgical technique is essential to ascertain their proper working condition.

NOTE: It is recommended to tension the self-locking polymer cerclage system with the tensioner developed by COUSIN BIOTECH. The use of another tensioner must be approved by the company COUSIN BIOTECH before use.

EXPLANTATION STEP:

After bone fusion, the sterile self locking polymer cerclage system needn't to be explanted. The device can be implanted for life. However in case of explantation, the cerclage device can be cut with blade and after releasing it from possible adhesions and fibrosis around the cerclage. The extraction of the titanium connector must be controlled during the explantation procedure and ensure that it is extracted with the cerclage. The cerclage should be cut with preference outside the splice area and loop to facilitate the extraction of the cerclage with the titanium connector.

11-STORAGE PRECAUTION

Use extreme care in handling and storage of implant components. Cable and connector plate must be handled with care. Twisting, kinking, cutting, notching, or scratching the braided surface may reduce the strength, fatigue resistance and/or wear characteristics of the implant system. These, in turn may induce internal stresses that are not obvious to the eye and may lead to fracture of the component. To be stored in a dry place away from light and at room temperature in its original packaging.

12-EXPLANTATION AND ELIMINATION OF THE PRODUCT

Devices should be retrieved and handled according to the recommendations of ISO 12891-1:2015 "Implants for surgery - Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling".

Any explanted device must be sent back for analysis, following the current protocol. This protocol is available on request from COUSIN BIOTECH. It is important to note that any implant that must not be cleaned or disinfected prior to despatch must be contained in sealed packaging.

The removed medical device must be disposed of in accordance with the country's standards for the disposal of infectious waste.

There are no specific recommendations for the disposal of a non-implanted device.

13-INFORMATION REQUEST AND COMPLAINTS

Following its quality policy, COUSIN BIOTECH is committed to making every effort to produce and supply a high-quality medical device. However, if a health professional (client, user, prescriber, etc.) has a complaint or cause for dissatisfaction with a product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In case of failure of an implant or if it has contributed to causing a serious adverse reaction in the patient, the health centre must follow the legal procedures in that country and inform COUSIN BIOTECH immediately.

For any correspondence, please specify the reference, batch number, contact person's details and a comprehensive description of the incident or complaint.

Brochures, documentation and surgical technique are available on request from COUSIN BIOTECH and its distributors. If further information is needed or required, please contact your COUSIN BIOTECH representative or distributor.

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