BIOMESH®

CA.B.S'Air® Composite

STERILE INTRAPERITONEAL NON-RESORBABLE PARIETAL REINFORCEMENT IMPLANT

en	<u>Instructions for use</u>	Page	2
fr	<u>Notice d'instructions</u>	Page	5
de	<u>Gebrauchsanweisung</u>	Seite	8
it	<u>Istruzioni per l'uso</u>	Pagina	11
es	<u>Instrucciones de uso</u>	Página	14
pt	<u>Instruções de utilização</u>	Página	17
ελ	<u>Οδηγίες χρήσης</u>	Σελίδα	20
da	<u>Brugsvejledning</u>	Side	23
fi	<u>Käyttöohjeet</u>	Sivu	26
nl	<u>Gebruiksaanwijzing</u>	Pagina	29
ro	<u>Instrucțiuni de utilizare</u>	Pagină	32
tr	<u>Kullanım Talimatları</u>	Sayfa	35
sv	<u>Bruksanvisning</u>	Sida	38

**COUSIN BIOTECH**

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NOT188 / 230530

Version du 30/05/2023

Caution : Federal law (USA) restricts this device to sale,
distribution and use by or on the order of a physician.

This release is the last update of the instructions of use and replace the previous edition

STERILE INTRAPERITONEAL NON-RESORBABLE PARIETAL REINFORCEMENT IMPLANT

BIOMESH® CA.B.S.'AIR® COMPOSITE

EXPANDABLE PARIETAL PROSTHESIS WITH INFLATABLE BALLOON FOR VISCERAL SURGERY

SINGLE USE STERILE PRODUCT**DESCRIPTION**

The BIOMESH® CA.B.S.'Air® Composite prosthesis consists of:

- A parietal prosthesis made of an ePTFE and polypropylene plug, two or four polyethylene terephthalate with biocompatible colorant (D&C Green #6) and ePTFE fixation threads crimped with stainless steel needles (see on package).
- A balloon which allows the deployment of the parietal prosthesis (except for CABST1R05T).

IMPLANTED MATERIALS:

Polypropylene, Polyethylene terephthalate, ePTFE, Biocompatible colorant (D&C Green #6).

MATERIALS IN CONTACT WITH THE PATIENT DURING THE SURGERY:

Stainless steel, Poly vinyl chloride.

Origin neither human nor animal - Non-absorbable.

INDICATIONS:

Repair and parietal reinforcement for umbilical and ventral hernias.

PERFORMANCES

The BIOMESH® CA.B.S.'Air® Composite prosthesis has a non-adhesive side, which is designed to avoid or reduce visceral adherences during intra peritoneal use or extra peritoneal use in the presence of a peritoneal defect for the mini invasive surgical treatment of umbilical and ventral hernias.

The BIOMESH® CA.B.S.'Air® Composite is biocompatible and non resorbable. They have the advantage of having a very high resistance to suture.

If necessary, the surgeon has the possibility of adding two additional points of fixing, according to his appreciation.

CONTRAINDICATIONS

Do not use in the following cases:

- Allergy to any of the components
- Infected site
- Pregnancy
- Growing children
- Anticoagulant therapy

UNDESIRABLE SIDE EFFECTS

Like any implantable medical device, this implant could generate possible undesirable effects:

- | | | | |
|-------------------|----------------------|---------------------|----------------------------|
| - Discomfort/Pain | - Adhesion formation | - Hematoma | - Foreign Body Reaction |
| - Infection | - Fistula formation | - Mesh migration | - Irritation near by organ |
| - Recurrence | - Erosion | - Mesh deformation | - Obstruction |
| - Inflammation | - Seroma | - Allergic reaction | - Lymphocele |

PRECAUTIONS FOR USE

BIOMESH® CA.B.S.'Air® Composite prosthesis are delivered sterile (ethylene oxide sterilization).

Before any use, inspect the integrity of packaging and device (of which blister / peelable pouches). Do not use in the event of deterioration of the device and/or the packaging.

Do not use if the device is out of date.

COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned.

This device must be implanted only by a qualified surgeon and trained with the use of the product (knowledge of anatomy and visceral surgery).

IMPORTANT : DO NOT REUSE - DO NOT RESTERILIZE

THIS BALLOON MUST IMPERATIVELY BE WITHDRAWN

As specified on the product's labeling, the BIOMESH® C.A.B.S.'Air® Composite mesh is for single use only. It can not be re-used and/or re-sterilized (potential risks would be and are not limited to: loss of the product's sterility, risk of infection, loss of the product's efficiency, recurrence)

STORAGE

Store in a dry place away from sunlight and at room temperature in its original packaging.

IMPORTANT


For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

EXPLANTATION AND ELIMINATION OF DEVICES

Explantation and handling should be performed in accordance with the recommendations of ISO 12891-1:2015 "Implants for surgery - Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling". Any explanted device must be returned for analysis, following the current protocol. This protocol is available on demand from COUSIN BIOTECH. It is important to note that any implant that is not to be cleaned and disinfected before expedition must be contained in a sealed package. Explanted medical devices must be eliminated in accordance with the standards applicable in the country for the disposal of infectious hazardous waste. There are no specific recommendations for the elimination of a non-implanted device. The second paragraph is in case of an information request or complaint.

INFORMATION REQUEST AND COMPLAINT

Following its quality policy, COUSIN BIOTECH is committed to make every effort to produce and supply a high-quality medical device. However, if a health professional (client, user, prescriber, etc) has a complaint or cause for dissatisfaction with a product in terms of quality, safety or performance, he/she must inform COUSIN BIOTECH of this as soon as possible. In the case of failure of an implant or if it contributed to causing serious adverse effects for the patient, the health centre must follow the legal procedures in its country and inform COUSIN BIOTECH of this as quickly as possible. For any correspondence, please quote the reference and batch number, the reference coordinates and give a comprehensive description of the incident or complaint. Brochures, documentation and the surgical technique are available on request from COUSIN BIOTECH and its distributors.

SURGICAL PROCEDURE (except for CABST1R05T which comes without balloon)					
1-Installation: - The prosthesis equipped with its balloon crosses completely the hernial defect, -its positioning becomes then entirely retro parietal (intraperitoneal)		2-Deployment: - The expansion and the spreading of the prosthesis are obtained by inflating the balloon, - Parietal fixation is carried out by two sutures thanks to the pre-positioned threads		3- Prosthesis in place: -The prosthesis is fixed, -The balloon must imperatively be withdrawn	

Description of the prosthesis:

This prosthesis is in a shape of a flat plug with 2 fixation threads crimped with needles (4 needles) for the product references CABST1R05T and CABS AIRC7, and with 4 fixation threads crimped with needles (8 needles) for the CABS AIRC09 product reference.

The aim of its innovating expansion balloon concept is to make dissection easier and to secure the unfolding of the prosthesis. This balloon is located inside the plug; it is filled-up with 30cm³ of air. It acts as a pneumatic jack ensuring complete crease-free unfolding of the prosthesis.

Its double-sided design combines a smooth ePTFE side which reduces viscera adhesions with a knitted monofilament polypropylene side which allows better colonisation of the peritoneum.

Sizes of the prosthesis:

The medical device is a sterile kit.

The product range is composed of:

CABS AIR 5cm diameter: for small hernias and ventral hernias on trocar site

CABS AIR 7cm diameter: for umbilical and linea alba hernias, as well as small ventral hernias < 2cm

CABS AIR 9cm diameter: for umbilical and linea alba hernias, as well as small ventral hernias that measure between 2 and 5cm

Installation of the patient :

Patient lying on his back
Both legs kept in the body's axis
Each arm placed on an arm wrestler in a 90° angle from the body's axis

1-Incision :

The cutaneous incision (2 to 4cm according to the size of the defect) is centred on the hernia.

2-Dissection and resection of the peritoneal sac :

The hernial sac is identified and its content is re-introduced back inside the abdominal cavity. The sac is then left opened in order to allow the Intra-peritoneal prosthesis to be put in place.

Next step consists in liberating epiploico-peritoneal adhesions around the edges of the hernia so that the prosthesis can fully expand.

3-Positioning and fixating the prosthesis :

BIOMESH® CA.B.S.'Air® Composite is introduced through the hernial orifice after being rolled up into a tube shape or folded in an umbrella-like way.

The balloon is filled up in order to allow the expansion of the prosthesis and the placement of the prosthesis is checked so that it remains centred above the hernial orifice.

The prosthesis is fixed with the two pre-mounted threads that will transfix the peritoneum and the muscular wall.

Because of its Intraperitoneal placement, two additional points should be used to ensure secure fixation of the prosthesis and to avoid any epiploic tissue being caught in between the prosthesis and the wall.

Be careful not to suture the balloon to the implant when fixing the prosthesis.

The balloon is deflated and removed.

4-Closure of the incision :

Musculo-aponeurotic closure in front of the prosthesis with resorbable suture thread. Sub-cutaneous and then cutaneous closure, normally without drainage.

IMPLANT CARD

Be sure to scan your implant card immediately upon receipt so that you can keep track of it in case of loss.

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