SOFT LIFT® Kit

STERILE UROGENITAL IMPLANTS: SUBURETHRAL SUPPORT TAPE

en	Instructions for use	Page	2
	Surgical Procedure		37
fr	Notice d'instructions	Page	4
	Technique Opératoire		37
de	Gebrauchsanweisung	Seite	7
	Operationstechnik		37
it	Istruzioni per l'uso	Pagina	10
	Tècnica Operatoria		37
es	Instrucciones de uso	Pagina	12
	Técnica Operatoria		37
pt	Instruções de uso	Pagina	15
	Técnica Operatoria		37
ελ	Οδηγίες χρήσης	σελίδα	18
	Χειρουργική Διαδικασία		37
hu	Használati útmutató	Oldal	21
	Sebészeti eljárás		37
lt	Naudojimo instrukcija	Puslapis	23
	Chirurgine procedūra		37
nl	Gebruiksaanwijzing	Page	25
	Chirurgische procedure		37
pl	Instrukcja obsługi	Strona	28
	Technika operacyjna		37
sk	Návod na použitie	Strana	31
	Chirurgický postup		37
cs	Návod k použití	Strana	33
	Chirurgický postup		37
ro	Instrucțiuni de utilizare	Pagina	35
	Procedura chirurgicala		37

COUSIN BIOTECH

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59117 Wervicq-Sud -

FRANCE

Tél. : +33 (0) 3 20 14 41 20

Fax : +33 (0) 3 20 14 40 13

www.cousin-surgery.com

Made in France

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NOT220 / 231002

Version du 02/10/2023

This release is the last update of the instructions for use and replaces the previous edition

SOFT LIFT®Kit Suburethral support tape

DESCRIPTION

The SOFT LIFT® kit comprises a tape and 3 ancillary devices (1 pair for the transobturator technique and 1 for the retropubic and suprapubic approach) and a winged guide.

They are single use and delivered sterile.

The SOFT LIFT® prosthesis consists of a tape of knitted monofilament polypropylene, packaged in a plastic sheath.

It is prolonged by self-tightening loops for placing.

The ancillary devices have different shapes depending on the chosen approach: 1 curved needle for the retropubic and suprapubic approach, a pair of helical needles for the in-out and out-in transobturator approaches.

The end of the ancillary needles has a slot for fixing self-tightening loops.

A single-use, sterile guide is supplied with the implant to secure the passage of the ancillary device during the in-out transobturator approach.

When placing the self-tightening loop in the slot, always make sure it is as tight as possible and that the thread is in the slot designed for this purpose so that the prosthesis does not come apart from the ancillary when passing in the patient's tissues. (See figure 1 below).


Fig 1

MATERIALS

- **Prosthesis:** polypropylene (knitted monofilament)
 - **Thread:** Polyethylene terephthalate + biocompatible dye + dimethyl siloxane
 - **Ancillary device:** Stainless steel needle, polyoxymethylene handle.
 - **Guide:** Liquid crystal polymer + biocompatible dye.
- Not human nor animal origin – Non absorbable.**

INDICATIONS

The SOFT LIFT® prosthesis is designed for use as a suburethral support sling for the surgical treatment by vaginal approach in women.

PERFORMANCE

These biocompatible, not biodegradable prostheses cause reactive fibrosis that produces additional urethral support.

These prostheses are easy to insert due to their specially designed shape, excellent flexibility and optimum, rapid colonization and integration.

CONTRAINDICATIONS

Do not implant in the following cases:

- Allergy to any one of the components
- For women receiving anticoagulant therapy.
- Urinary tract infection
- Pregnant women
- Skeletally-immature children

Any patient who could possibly be pregnant should be warned that pregnancy may abolish the benefit of the procedure.

SIDE EFFECTS

As with all medical implants, this product may lead to side effects, such as:

Perforation or damage to blood vessels, bladder, urethra, rectum, intraperitoneal organs or nerves.

As with all implants, local irritation and/or a reaction to the foreign body may occur.

At tissue level, reactions to the implant may include vaginal extrusion, erosion through the urethra or surrounding tissues, displacement of the device, fistula or inflammation. If any such reaction occurs, it may be necessary to completely remove the sling.

As with all foreign bodies, the polypropylene sling may aggravate a pre-existing infection

Over-correction may lead to temporary or permanent obstruction of the lower urinary tract and to retention of urine.

The known risks of surgery for the treatment of incontinence include hematoma, acute or chronic pain, scarification of the vagina, retraction or release of the mesh, infection, erosion, de novo urgency or voiding difficulties, dyspareunia, migration of the device, the total failure of the intervention leading to relapse. One or several repeated surgeries could also be necessary to treat those complications, as some complications are not always completely corrected.

USAGE PRECAUTIONS

The SOFT LIFT®Kit device should be implanted by qualified surgeons only (knowledge of gynaecological and urological anatomy and surgery) who are extremely familiar with the device.

The SOFT LIFT®Kit implant, guide and ancillary devices are delivered sterile (ethylene oxide sterilization). Check that the packaging is intact before using (do not use if the labels and/or protective packaging are damaged). Do not use if the device is damaged.

The plastic sheath that covers the sling must be removed and must not be implanted.

Do not reuse or resterilize the implant, the guide and ancillary devices.

IMPORTANT: DO NOT REUSE – DO NOT RESTERILIZE

As indicated on the product labelling, the implant, guide and ancillary devices are designed for single use. They should in no event be reused and/or resterilized (potential risks include but are not limited to: loss of product sterility, risk of infection, loss of effectiveness, relapse).

STORAGE

To be stored in a dry place away from sunlight, at room temperature, in its original packaging.

RETRIEVAL AND DISPOSAL OF DEVICES

Devices should be retrieved and handled according to the recommendations of ISO 12891-1:2015 “Implants for surgery – Retrieval and analysis of surgical implants” Part 1: “Retrieval and Handling”.

Any explanted device must be sent back for analysis, following the current protocol. This protocol is available on request from COUSIN BIOTECH. It is important to note that any implant that must not be cleaned or disinfected prior to despatch must be contained in sealed packaging.

The removed medical device must be disposed of in accordance with the country’s standards for the disposal of infectious waste.

There are no specific recommendations for the disposal of a non-implanted device.

INFORMATION REQUESTS AND COMPLAINTS

Following its quality policy, COUSIN BIOTECH is committed to making every effort to produce and supply a high quality medical device. However, if a health professional (client, user, prescriber, etc.) has a complaint or cause for dissatisfaction with a product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In case of failure of an implant or if it has contributed to causing a serious adverse reaction in the patient, the health centre must follow the legal procedures in that country, and inform COUSIN BIOTECH immediately.

For any correspondence, please specify the reference, batch number, contact person’s details and a comprehensive description of the incident or complaint.

Brochures, documentation and surgical technique are available on request from COUSIN BIOTECH and its distributors.

IMPORTANT

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

IMPLANT CARD

Be sure to scan your implant card immediately upon receipt so that you can keep track of it in case of loss.

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