# NAJA® & ResCUBE

## NAJA® & ResCUBE ancillaries

en	Instructions for use of single use NAJA® & ResCUBE ancillary	Page	2
fr	Notice d'instructions de l'ancillaire NAJA® & ResCUBE à usage unique	Page	4
de	Anleitung zur Verwendung des NAJA® & ResCUBE Hilfsinstrumentes zur einmaligen Verwendung	Seite	6
it	Istruzioni per l'uso di uno strumentario NAJA® & ResCUBE monouso	Pagina	8
es	Instrucciones de uso del accesorio monouso NAJA® & ResCUBE	Página	10
pt	Instruções para utilização de acessório NAJA® & ResCUBE descartável	Página	12
ελ	Οδηγίες χρήσης για μία μόνο χρήση του βοηθητικού εξοπλισμού NAJA® & ResCUBE	Σελίδα	14
lv	Palīģierīču lietošanas instrukcija NAJA® & RESCUBE Vienreizējās spriegotājs	Puslapis	16
pl	Instrukcja użycia Przyrząd NAJA® & RESCUBE do użytku jednorazowego	Strona	18
ru	Инструкция по применению одноразового проводника NAJA® & ResCUBE	Стр.	20

**COUSIN BIOTECH**  
8 Rue Abbé Bonpain  
59117 Wervicq-Sud - FRANCE  
Tél. : +33 (0) 3 20 14 41 20  
Fax : +33 (0) 3 20 14 40 13  
[www.cousin-surgery.com](http://www.cousin-surgery.com)

Date de marquage CE: Mai 2017

**NOT247/230626**  
Version du 26/06/2023

Caution: Federal law (USA) restricts this device to sale, distribution and use by or on the order of a physician.

This release is the last update of the instructions of use and replace the previous edition

**Single use NAJA® & ResCUBE Ancillary**  
NAJA® & ResCUBE ancillaries

**DESCRIPTION**

This ancillary is intended for the placement of the NAJA® & ResCUBE ligament systems for spinal stabilisation. This ancillary is a Class Is medical device, intended for single use and is sterile supplied.

**MATERIALS**

The different elements of the ancillary are made from:

- IXEF®
- Stainless steel (for metallic components)

All the materials used are non-absorbable and are neither of human nor animal origin.

**PRECAUTIONS FOR USE**

The ancillary must stay stored in its original packaging, which ensures safety.

Carry out a visual examination of each instrument before each use in order to detect any possible deterioration. The instruments must not be used except with elements corresponding to the original COUSIN BIOTECH system.

Use of specific ancillary COUSIN BIOTECH system is required for implantation of the implant.

In addition, the instruments must not be used except by a surgeon with appropriate training in the placement of NAJA® & ResCUBE implants.

**Sterile device (Sterilized using Gamma irradiation)**

Store the ancillary in its original protective packaging, and do not remove them from the packaging until immediately before use.

Prior to use, check the product expiration date and verify the integrity of the sterile packaging. Do not use, if the package is damaged.

**Single-use device**

According to the labelling of this product, the ancillary is a single use device.

Products intended for single-use must not be re-used. Re-use or reprocessing (e.g. cleaning and resterilization) may compromise the structural integrity of the device and/or lead to device failure which may result in patient injury, illness or death. Furthermore, reuse or reprocessing of single-use devices may create a risk of contamination e.g. due to the transmission of infectious material from one patient to another. This could result in injury or death of the patient or user.

Contaminated ancillary must not be reprocessed. Any ancillary that has been contaminated by blood, tissue, and/or bodily fluids/matter should never be used again and should be handled according to hospital protocol. Even though it may appear undamaged, the ancillary may have small defects and internal stress patterns that may cause material fatigue.

**Storage precautions of the ancillary**

The device must be handled and stored with great care, in a dry place, protected from sunlight and at room temperature in its original packaging.

**INFORMATION REQUEST AND CLAIMS**

Following its quality policy, COUSIN BIOTECH is committed to make every effort to produce and supply a high quality medical device. However, if a health professional (client, user, prescriber ...) had a claim or cause of dissatisfaction with a product in terms of quality, safety or performance, he must inform COUSIN BIOTECH as soon as possible.

In case of failure of an implant or if it contributed to cause serious adverse effects for the patient, the health center must follow the legal procedures in his country, and inform COUSIN BIOTECH in the shortest time .

For any correspondence, thank you to specify the reference, batch number, the coordinates of a reference and a comprehensive description of the incident or claim.

Brochures, documentation and surgical technique are available on request from COUSIN BIOTECH and its distributors.

If further information is needed or required, please contact your COUSIN BIOTECH representative or distributor, or COUSIN BIOTECH at [spine@cousin-biotech.com](mailto:spine@cousin-biotech.com).

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