BIOLIG[®]

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NOT067 / 240531
Version du 31/05/2024

Instructions for use

BIOLIG® REINFORCEMENT ARTICULAR LIGAMENTS STERILE ANCHOR DEVICES STERILE PRODUCT FOR SINGLE USE

DESCRIPTION

Reinforcing ligaments and tendons for orthopedic surgery

MATERIALS

Polyethylene terephthalate (No resorbable)

Stainless steel

Neither human nor animal origin.

INDICATIONS

The ligaments are intended for the surgical reinforcement of instabilities in different joints. They can be located intra- or extra-articularly. The reinforcing ligament or tendinous should always be used to complement a suture or natural tissue graft.

CONTRAINDICATIONS

Do not use:

- In a septic environment or in an infected site,
- allergy to one of components
- growing child.

POSSIBLE ADVERSE SIDE EFFECTS

- Routine surgical complications
- superficial or deep infection
- anesthesia complications
- local hematoma, blood loss
- Postoperative loss of motion, frozen joint
- Inflammatory / allergic reactions / arthro fibrosis (around 1%) / Bone resorption (osteolysis)
- Ongoing instability episodes
- Neurovascular injury, peripheral neuropraxia
- Wound dehiscence, serous leak
- Post-surgical pain, persistent pain
- Ligament injury
- Permanent elongation of the ligament / Rerupture of the ligament

PRECAUTIONS FOR USE

Inspect the packaging to be sure it is intact (do not use if the labels and/or the protective packaging are damaged).

Do not use if the device is damaged.

DO NOT REUSE

DO NOT RESTERILIZE

COUSIN BIOTECH does not make any guarantees or recommendations as regards the use of a particular brand of fixation devices. The properties of these devices are the subject of modifications by the manufacturer and over which COUSIN BIOTECH has no control. Suture clips or sharp mounting screws must not be used.

This device should be implanted **only** by a qualified surgeon who has been trained in the use of this product.

Postoperative precautions

The device has not been tested for heating or migration in the MR environment.

The device has been evaluated for safety and compatibility in the MR (Magnetic Resonance) environment. According to medical literature and in view of the materials used in the implant, the device can be considered as MR safe.

PRECAUTIONS FOR STORAGE

To be stored in a dry place, protected from sunlight, at room temperature in the original packaging.

EXPLANATION AND DISPOSAL OF DEVICES

Explanation and handling should be carried out in accordance with ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling". Any explanted device must be returned for analysis as per the current protocol. This protocol is available from COUSIN BIOTECH on request. It is important to note that any implant that has not been cleaned and disinfected prior to dispatch must be placed in a sealed package. Explanted medical devices

must be disposed of in accordance with local regulations governing the disposal of infectious hazardous waste. The disposal of non-implanted devices is not subject to specific recommendations. The second paragraph covers information requests or complaints.

INFORMATION REQUESTS AND COMPLAINTS

In accordance with its quality policy, COUSIN BIOTECH is committed to producing and supplying high-quality medical devices. However, if health professionals (client, user, prescriber, etc.) wish to lodge a complaint or report dissatisfaction with a given product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In the event of implant failure or if the implant is implicated in serious adverse effects experienced by the patient, the health centre must follow the legal procedures in the country in question, and inform COUSIN BIOTECH without delay. The reference, batch number, reference contact details and a comprehensive description of the incident or complaint must be stated in all correspondence. Brochures, documentation and surgical techniques are available on request from COUSIN BIOTECH and its distributors.

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