# SACROMESH SOFT PROLAPS'<sup>®</sup>

## IMPLANT FOR THE TREATMENT OF PROLAPSE BY SACROCOLPOPEXY

en	<a href="#">Instructions for use</a>	Page	2
fr	<a href="#">Notice d'instructions</a>	Page	4
de	<a href="#">Gebrauchsanweisung</a>	Seite	7
it	<a href="#">Istruzioni per l'uso</a>	Pagina	10
es	<a href="#">Instrucciones de uso</a>	Página	13
pt	<a href="#">Nota de instruções</a>	Página	16
ελ	<a href="#">Οδηγίες χρήσης</a>	Σελίδα	19
lt	<a href="#">Naudojimo instrukcija</a>	Puslapis	22
fi	<a href="#">Käyttöohjeet</a>	Sivu	24
hu	<a href="#">Felhasználási útmutató</a>	Oldal	26
nl	<a href="#">Gebruiksaanwijzing</a>	Pagina	29
pl	<a href="#">Instrukcja obsługi</a>	Strona	32
ro	<a href="#">Instrucțiuni de utilizare</a>	Pagină	35
sv	<a href="#">Bruksanvisning</a>	Sida	38

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**Made in France**

**NOT123 / 240531**  
Version du 31/05/2024

This release is the last update of the instructions for use and replaces the previous edition

**SACROMESH SOFT PROLAPS®  
UROGENITAL IMPLANT**  
**Implant for the treatment of prolapse by sacrocolpopexy**

**DESCRIPTION**

Sacromesh Soft Prolaps® implants consist of knitted polypropylene monofilament. They are intended for sacrocolpopexy of the pelvic prolapse in women, by abdominal approach.

**MATERIALS**

100% polypropylene: knitted monofilament thread, weight 38 g/m2 .  
Marked thread : Polyethylene terephthalate tinted with a biocompatible dye (See on the package).  
Origin neither human, nor animal \_ Non resorbable.

**INDICATIONS**

Sacromesh Soft Prolaps® prosthesis are designed for use as reinforcement and support of tissues when treating genital prolapse surgically by sacrocolpopexy. They can be used in either conventional or laparoscopic surgery.

**PERFORMANCES**

The biocompatible and non-resorbable prosthesis have the advantage of being easy to install, show good resistance to sutures and traction, with excellent flexibility, optimal and rapid integration and colonization.

**CONTRAINDICATIONS**

Do not use the implant in the following circumstances:  
Allergy to any of the components  
In a place with progressive or latent infection or signs of tissue necrosis  
For women receiving anticoagulant therapy  
Urinary or genital tract infection or non-aseptic conditions  
Pregnancy  
Growing children  
Any patient who could possibly be pregnant should be warned that pregnancy may abolish the benefit of the procedure.

**SIDE EFFECTS**

As with all medical implants, this product may lead to side effects, such as:  
Perforation or damage to blood vessels, bladder, urethra, rectum, intraperitoneal organs or nerves.  
Obstruction/Occlusion  
Adhesion formation  
Allergic reaction  
Sexual dysfunction  
Adverse obstetric outcomes (preterm premature rupture of membranes)  
As with all implants, local irritation and/or a reaction to the foreign body may occur.  
At tissue level, reactions to the implant may include vaginal extrusion, erosion through the urethra or surrounding tissues, displacement of the device, fistula or inflammation. If any such reaction occurs, it may be necessary to completely remove the implant.  
As with all foreign bodies, the polypropylene mesh may aggravate a pre-existing infection  
Over-correction may lead to temporary or permanent obstruction of the lower urinary tract and to retention of urine.  
The known risks of surgery for the treatment of prolapse include hematoma, acute or chronic pain, scarification of the vagina, retraction or release of the mesh, infection, erosion, de novo urgency or voiding difficulties, dyspareunia, migration of the device, the total failure of the intervention leading to relapse. One or several repeated surgeries could also be necessary to treat those complications, as some complications are not always completely corrected.

**SURGICAL PROCEDURE**

Surgery under general anaesthetic.

Abdominal incision or installation of trocars.

Exposure of the lumbosacral angle.

Dissection of the prerectal fascia after incision of the anterior wall of Douglas' cul-de-sac. The dissection continues laterally down to the posterior elevator muscles, using a vaginal valve if necessary.

The mesh is fixed as distally as possible on the pelvic floor.

Dissection of the vagina's anterior wall.

The anterior mesh is fixed by a running stitch, avoiding transfixation. Depending on the mesh used by the surgeon and of his preference, the mesh is passed laterally to the uterine isthmus. The uterus is held in place, two meshes are passed on either side of it (for precut mesh or mesh cut longitudinally by the surgeon).

Fastening the anterior and posterior mesh to the Promontory (presacral ligament) without excessive tension .

Current recommendation is to do this stich on the Promontory with a non resorbable suture and avoid the use of staplers or tackers.

Peritonization must be as complete as possible to avoid contact between the prosthesis and the intraperitoneal constituents.

### **PRECAUTIONS FOR USE**

Sacromesh Soft Prolaps<sup>®</sup> devices must **always** be implanted by a qualified surgeon (familiar with the relevant anatomy and experienced in gynaecological and urological surgery) who has a good knowledge of the device, its intended use, the instruments and the surgical technique.

Sacromesh Soft Prolaps<sup>®</sup> implants are delivered sterile (sterilised using ethylene oxide), inspect the packaging to be sure it is intact (do not use if the labels and/or protective packaging are damaged).

Do not use if the device is damaged or out of date.

Do not reuse and do not resterilize the prosthesis.

### **STORAGE OF THE PROSTHESIS**

To be stored in a dry place away from sunlight, at room temperature, in its original packaging.

### **IMPORTANT**

**For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.**

### **EXPLANTATION AND DISPOSAL OF DEVICES**

Explantation and handling should be carried out in accordance with ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling". Any explanted device must be returned for analysis as per the current protocol. This protocol is available from COUSIN BIOTECH on request. It is important to note that any implant that has not been cleaned and disinfected prior to dispatch must be placed in a sealed package. Explanted medical devices must be disposed of in accordance with local regulations governing the disposal of infectious hazardous waste. The disposal of non-implanted devices is not subject to specific recommendations. The second paragraph covers information requests or complaints.

### **INFORMATION REQUESTS AND COMPLAINTS**

In accordance with its quality policy, COUSIN BIOTECH is committed to producing and supplying high-quality medical devices. However, if health professionals (client, user, prescriber, etc.) wish to lodge a complaint or report dissatisfaction with a given product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In the event of implant failure or if the implant is implicated in serious adverse effects experienced by the patient, the health centre must follow the legal procedures in the country in question, and inform COUSIN BIOTECH without delay. The reference, batch number, reference contact details and a comprehensive description of the incident or complaint must be stated in all correspondence. Brochures, documentation and surgical techniques are available on request from COUSIN BIOTECH and its distributors.

### **IMPLANT CARD**

**Be sure to scan your implant card immediately upon receipt so that you can keep track of it in case of loss.**

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