# 4D Mesh®

SEMI RESORBABLE PARIETAL REINFORCEMENT IMPLANT

		Page	
en	Instructions for use	2	
fr	Notice d'instructions	5	
de	Gebrauchsanweisung	9	
it	Istruzioni per l'uso	13	
es	Instrucciones de uso	16	
pt	Instruções de utilização	19	
ελ	Οδηγίες χρήσης	22	
cs	Návod k použití	26	
da	Brugsvejledning	29	
fi	Käyttöohjeet	32	
hu	Használati Útmutató	36	
nl	Gebruiksinstructies	39	
ro	Instrucțiuni de utilizare	43	
sk	Návod na použitie	46	
sv	Bruksanvisning	49	
tr	Kullanım Talimatları	52	
lv	Lietošanas instrukcijas	55	
et	Kasutusjuhend	58	
lt	Naudojimo instrukcijos	61	
en	Patient information leaflet (Australian requirements)	64	

**COUSIN BIOTECH**

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**NOT233\_240531**

Version du 31/05/2024

Caution : Federal law (USA) restricts this device to sale,  
distribution and use by or on the order of a physician.

This release is the last update of the instructions of use and replace the previous edition

## SEMI-RESORBABLE PARIETAL REINFORCEMENT IMPLANT STERILE PRODUCT FOR SINGLE USE

### **DESCRIPTION**

4DMESH® meshes are semi-resorbable parietal reinforcement implants.

### **MATERIALS**

Polypropylene (non resorbable) – Poly L Lactic Acid (resorbable). Blue thread (depending on references): monofilament polypropylene, non absorbable, colorant [phthalocyaninato(2-)] copper.

**No human nor animal origin. Semi-resorbable.**

### **INDICATIONS**

4DMESH® meshes are designed for the repair and reinforcement of inguinal and femoral hernias.

### **PERFORMANCES**

The device 4DMESH® is designed to be extra-peritoneally implanted. These products are biocompatible and semi-resorbable. Their main benefits are: can be cut to fit patient's anatomy, high resistance to suture and abdominal pressure, macroporous, transparent, lightweight after PLLA resorption.

### **CONTRAINDICATIONS**

Do not use in the following cases:

- Allergy to any of the components
- Infected site
- Pregnancy
- Growing children
- Anticoagulant therapy

### **UNDESIRABLE SIDE EFFECTS**

Like any implantable medical device, this implant is susceptible to generate possible undesirable side effects that may lead to reoperation:

- |                   |                      |                     |                           |
|-------------------|----------------------|---------------------|---------------------------|
| - Discomfort/Pain | - Adhesion formation | - Hematoma          | - Foreign Body reaction   |
| - Infection       | - Fitsula formation  | - Mesh migration    | - Irritation nearby organ |
| - Recurrence      | - Erosion            | - Mesh deformation  | - Lymphocele              |
| - Inflammation    | - Seroma             | - Allergic reaction |                           |
- A small additional risk of male infertility has been identified after bilateral groin hernia repair using meshes (open or laparoscopic)

### **PRECAUTIONS FOR USE**

4DMESH® must always be separated from abdominal cavity by peritoneum

4DMESH® prosthesis are delivered sterile (ethylene oxide sterilisation).

Before any use, inspect the integrity of packaging and device (of which blister / peelable pouches). Do not use in the event of deterioration of the device and/or the packaging. Do not use if the device is out of date. COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned.

This device must be implanted **only** by a qualified surgeon and trained with the use of the product (knowledge of anatomy and visceral surgery).

### **IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE**

As specified on the product's labelling, the 4DMESH® mesh is for single use only. It can not be re-used and/or re-sterilized (potential risks would be and are not limited to: loss of the product's sterility, risk of infection, loss of the product's efficiency, recurrence)

### **EXPLANTATION AND ELIMINATION OF DEVICES**

Devices should be retrieved and handled according to the recommendations of ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling".

Any explanted device must be sent back for analysis, following the current protocol. This protocol is available on request from COUSIN BIOTECH. It is important to note that any implant that must not be cleaned or disinfected prior to despatch must be contained in sealed packaging.

The removed medical device must be disposed of in accordance with the country's standards for the disposal of infectious waste.

There are no specific recommendations for the disposal of a non-implanted device.

### **STORAGE OF THE PROSTHESIS**

To be stored in a dry place away from light and at room temperature in its original packaging.

### **INFORMATION REQUESTS AND COMPLAINTS**

Following its quality policy, COUSIN BIOTECH is committed to making every effort to produce and supply a high quality medical device. However, if a health professional (client, user, prescriber, etc.) has a complaint or cause for dissatisfaction with a product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In case of failure of an implant or if it has contributed to causing a serious adverse reaction in the patient, the health centre must follow the legal procedures in that country, and inform COUSIN BIOTECH immediately.

For any correspondence, please specify the reference, batch number, contact person's details and a comprehensive description of the incident or complaint.

Brochures, documentation and surgical technique are available on request from COUSIN BIOTECH and its distributors.

### **IMPORTANT**

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

### **SURGICAL PROCEDURE**

#### **DESCRIPTION OF THE SURGICAL PROCEDURES:**

##### **Open Inguinal hernia repair : Lichtenstein**

- The patient is placed in the supine position
- Local or general anaesthesia
- Iliac incision of approximately 5 cm
- Dissection is performed accordingly to the Lichtenstein technique
- The hernia sac is identified and reduced
- The transversalis fascia can be tightened with a non-absorbable suture. With direct hernias it helps to avoid any immediate post-operative recurrence
- 4DMESH® mesh can be cut to fit the patient's anatomy
- The prosthesis is progressively put in place
  - begin by placing 4DMESH® on the pubic spine
  - the upper part of the prosthesis is placed on the conjoint tendon
  - carry placing the prosthesis towards the deep inguinal ring so that the first arm is placed behind
  - then the second arm is placed so that it encircles the spermatic cord
  - Fix the 4DMESH® mesh with non-resorbable sutures
    - apex is sutured to the pubic tubercle
    - the lower border of the mesh is sutured to the free edge of the inguinal ligament
    - continuous suture extends up just medial to the anterior superior iliac spine
    - suture the two tails of the mesh together around the spermatic cord
    - the infero-medial corner of the mesh is attached overlapping the pubic tubercle.
    - the mesh is anchored to the conjoined tendon by interrupted sutures
  - the cord then rests on the prosthesis
  - suction drain can be placed beneath the external oblique aponeurosis, especially in large inguinal hernias, where an extensive dissection was performed
  - the aponeurosis of external oblique is then closed with absorbable sutures
  - before closure of the surgical incision, its edges are infiltrated with a long-acting local anaesthetic, such as Naropain®
  - During peri-operative care of the patient, prophylactic antibiosis is usually given for 48 – 72 hours postoperatively

## Laparoscopic procedure (TAPP or TEP)

### 1/ Example of a TAPP repair

- General anaesthesia
- A pneumoperitoneum is created in the usual fashion
- intra-abdominal findings are reported [intra-abdominal pathology and inguinal hernia defects and sacs].
- The two additional trocars are inserted under direct vision.
- Dissection is initiated and carried out following classical TAPP technique
- The indirect inguinal hernia sac should be dissected carefully from the Spermatic Cord. Direct hernia sacs are easily dissected.
- The 4DMESH® mesh is inserted through the Trocar into the intra-abdominal cavity and deployed over the inguinal region. The blue thread on the oval and preshaped mesh is the anatomic landmark for the pubic bone. For the meshes with flaps, the cross shows the upper flap. The 4DMESH® mesh is attached to Cooper's Ligament, around and lateral to the Inferior Epigastric Vessels using tacks, staples or other mesh fixating device, according to the Surgeon's usual technique. The operator should be meticulous to avoid the triangle of doom, triangle of pain, iliac vessels and to place the fixations lateral to the inguinal ring. Mesh fixation points should be at least 1cm from the edge of the mesh with 1 cm spacing between fixation points.
- The operator should check the 4DMESH® mesh is well anchored to the surrounding structures.
- The peritoneum is closed meticulously and no defect between the peritoneum and the abdominal wall should be left open. In addition, it should cover the entire Mesh.
- The trocars are removed under direct vision. The fascia of the sub-umbilical trocar site is closed as needed.

### 2/ Example of a TEP repair

- General anaesthesia
- Create the pro-peritoneal space.
- The 10 mm trocar is inserted.
- Create the pneumo-pro-peritoneum.
- Two 5 mm ports are inserted under direct vision.
- If a peritoneal defect is created at this time it should be immediately closed.
- The anatomical landmarks are identical to when performing a TAPP repair.
- The indirect inguinal hernia sac should be dissected carefully from the Spermatic Cord. Direct hernia sacs are easily dissected. This is done by gentle traction with atraumatic graspers.
- The 4DMESH® mesh is inserted into the pre-peritoneal space and deployed over the inguinal region. The blue thread on the oval and preshaped mesh is the anatomic landmark for the pubic bone. For the meshes with flaps, the cross shows the upper flap.
- The 4DMESH® mesh is first attached or secured to Cooper's Ligament, and the superior aspect of the pubic ramus. The 4DMESH mesh is then fixed on the posterior aspect of the Linea Alba. The anchoring is continued around and lateral to the Inferior Epigastric Vessels. A few fixations are used to fix the 4DMESH mesh lateral to the internal ring. Mesh fixation points should be at least 1cm from the edge of the mesh with 1 cm spacing between fixation points.
- The 4DMESH® mesh should gently and generously cover the Iliac Vessels without major gaps.
- The operator should check the Mesh is well anchored to the surrounding structures. The Mesh should not migrate and remain in place.
- The pre-peritoneal space will be checked for any peritoneal defects. If any, they should be closed. The 5 mm trocars are removed under direct vision. The pre-peritoneal space will collapse.
- If there are any questions about a missed peritoneal defect, a completion laparoscopy should be performed.
- Trocars are removed and the fascial defect closed with the appropriate suture. The skin edges are approximated in the usual manner.

### Size of the prosthesis:

The 4DMESH® product range is composed of meshes of various sizes and shapes:

- Pre-cut meshes with keyhole for LICHTENSTEIN inguinal hernia repair
- Oval, pre-shaped and rectangular meshes for TAPP/TEP inguinal hernia repair

**IMPLANT DE RENFORCEMENT PARIETAL SEMI-RESORBABLE  
PRODUIT STERILE A USAGE UNIQUE**

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## Patient information leaflet 4DMESH®

Item	Patient information leaflet
1	a. <b>The name of the device :</b> 4DMESH® b. <b>The model of the device :</b> 4DMESH0715, 4DMESH1015, 4DMESH1215, 4DMESH1216, 4DMESH1317, 4DMESH1515u, 4DMESH1717, 4DMESHF812, 4DMESHM613, 4DMESHMBEL, 4DMESHMBES, 4DMESHPRL, 4DMESHPRLR, 4DMESHPRSL, 4DMESHPRSR, 4DMESHPRSR, 4DMESHPRXL, 4DMESHPRXR, 4DMESHRABA, FBOSRF128, FBOSRF952
2	a. <b>The intended purpose of the device :</b> 4DMESH® meshes are designed for the repair and reinforcement of inguinal and femoral hernias. b. <b>The kind of patient on whom the device is intended to be used :</b> Adult with inguinal and femoral hernias.
3	<b>Special operating instructions for the use of the device :</b> The device must be implanted only by a qualified surgeon and trained with the use of the product (knowledge of anatomy and visceral surgery).
4	a. <b>The intended performance of the device :</b> The device 4DMESH® is designed to be extra-peritoneally implanted. These products are biocompatible and semi-resorbable. Their main benefits are: can be cut to fit patient's anatomy, high resistance to suture and abdominal pressure, macroporous, transparent, lightweight after PLLA resorption. b. <b>Undesirable side effects :</b> Like any implantable medical device, this implant could generate possible undesirable effects : <ul style="list-style-type: none"><li>- Discomfort/pain</li><li>- Infection</li><li>- Recurrence</li><li>- Inflammation</li><li>- Adhesion formation</li><li>- Fistula formation</li><li>- Erosion</li><li>- Seroma</li><li>- Hematoma</li><li>- Mesh migration</li><li>- Mesh deformation</li><li>- Allergic reaction</li><li>- Foreign body reaction</li><li>- Irritation near by organ</li><li>- Lymphocele</li><li>- A small additional risk of male infertility has been identified after bilateral groin hernia repair using meshes (open or laparoscopic)</li></ul>

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5	<b>When 4DMESH® Should Not be Used :</b> Allergy to any of the components, infected site, pregnancy, growing children, anticoagulant therapy.
6	<b>MRI Information :</b> 4DMESH® implants are extraperitoneal partially resorbable parietal reinforcement implants. This surgical mesh is composed of Poly-L-Lactic Acid (PLLA) and Polypropylene (PP).
7	<b>Precautions :</b> 4DMESH® must always be separated from abdominal cavity by peritoneum  4DMESH® prostheses are delivered sterile (ethylene oxide sterilisation). Before any use, inspect the integrity of packaging and device (of which blister / peelable pouches). Do not use in the event of deterioration of the device and/or the packaging. Do not use if the device is out of date. COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned. This device must be implanted <b>only</b> by a qualified surgeon and trained with the use of the product (knowledge of anatomy and visceral surgery).
8	<b>Materials included in the device:</b> Polypropylene (non resorbable) – Poly L Lactic Acid (resorbable). Blue thread (depending on references): monofilament polypropylene, non absorbable, colorant [phthalocyaninato(2-)] copper.  No human nor animal origin. Semi-resorbable.
9	<b>Reporting a serious incident :</b> If you encounter any serious incident with your 4DMESH® implant you should report it immediately to your doctor. You should also report the incident to Cousin Biotech via the link <a href="mailto:serviceclients@cousin-surgery.com">serviceclients@cousin-surgery.com</a> or <a href="mailto:materiovigilance@cousin-biotech.com">materiovigilance@cousin-biotech.com</a> and to the Therapeutic Goods Administration using the following website: <a href="http://www.tga.gov">http://www.tga.gov</a>
10	<b>Legal Manufacturer :</b> Cousin Biotech, Allée des roses, 59117 Wervicq-Sud, France  <b>Manufacturing Facility :</b> Cousin Biotech, Allée des roses, 59117 Wervicq-Sud, France

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[www.cousin-biotech.com/ifu](http://www.cousin-biotech.com/ifu)  
IFU can be downloaded through the QR code.  
A hard copy can be sent within 7 days  
on request by email to  
[ifurequest@cousin-biotech.com](mailto:ifurequest@cousin-biotech.com)  
or by using the order form on our website.

LOT	en	Batch number	!	en	Caution (see instructions for use)
	fr	Numéro de lot		fr	Attention (voir notice d'instructions)
	de	Chargennummer		de	Gebrauchsanweisung beachten
	it	Número di lotto		it	Vedere manuale istruzioni
	es	Número de lote		es	Veanse las instrucciones de uso
	pt	Número do lote		pt	Atenção, consultar nota de instruções
	ελ	Αριθμός παρτίδας		ελ	Δείτε τις οδηγίες χρήσης
	cs	Číslo šarže		cs	Viz pokyny k použití
	da	Serienummer:		da	Se brugsvejledning
	nl	Eránumero		nl	Katso käyttöohjeet
	ro	Sarzs szám		ro	Olvassa el a használati leírást
	sk	Batchnummer		sk	Zie gebruksinstructies
	sv	Numär lot		sv	A se vedea instrucțiunile de utilizare
	tr	Číslo šarže		tr	Pozrite pokyny na použitie
	lv	Batch nummer		lv	Varning, se bruksanvisningen
	lt	Yama numarası		lt	Kullanım talimatlarına bakın
	fi	Partii number		fi	Ettevaatust (vaata kasutusjuhendit) Tähelepanu
	hu	Partijas numurs		hu	Skatit lietošanas instrukciju
	et	Partijos numeris		et	Atsargumo priemonės (žr. naudojimo instrukciją)
REF	en	Reference on the brochure	■	en	Manufacturer
	fr	Référence du catalogue		fr	Fabricant
	de	Verweis auf der Broschüre		de	Hersteller
	it	Referimento del catalogo		it	Produttore
	es	Referencia al folleto		es	Fabricante
	pt	Referência do catálogo		pt	Fabricante
	ελ	Αριθμός καταλόγου		ελ	Κατασκευαστής
	cs	Odkaz v brožuře		cs	Výrobce
	da	Reference på brochuren		da	Producent
	nl	Viittaus esitteessä		nl	Valmistaja
ro	Katalógus hivatkozás	ro	Gyártó		
sk	Referentie in de brochure	sk	Fabrikant		
sv	Referintä pe broşürä	sv	Producător		
tr	Odkaz v brožúre	tr	Výrobca		
lv	referens på broschyren	lv	Tillverkare		
lt	Broşürdeki referans	lt	Üretici		
fi	Viide brošüüris	fi	Tootja		
hu	Atsauce uz bukletu	hu	Ražotājs		
et	Nuoroda brošūroje	et	Gamintojas		
	en	Do not reuse	■	en	Use before: year and month
	fr	Ne pas réutiliser		fr	Utiliser jusque : année et mois
	de	Nicht wiederverwenden		de	Verbrauchen bis : Jahr und Monat
	it	Non riutilizzare		it	Utilizzare entro e non oltre: anno e mese
	es	No volver a utilizar		es	Úsese antes de: año y mes
	pt	Não reutilizar		pt	Utilizar até : ano e mês
	ελ	Να μην επαναχρησιμοποιείται		ελ	Ημερομηνία λήξης: έτος και μήνας
	cs	Nepoužívat znovu		cs	Spotřebovat do: rok a měsíc
	da	Må ikke genbruges		da	Bruges inden: år og måned
	nl	Älä käytä uudelleen		nl	Käytettävä ennen: vuosi ja kuukausi
ro	Ne használja újra	ro	Felhasználható: év és hónap		
sk	Niet hergebruiken	sk	Te gebruijen vóór: jaar en maand		
sv	A nu se reutiliza	sv	A se utiliza înainte de: anul și luna		
tr	Nepoužívat znova	tr	Spotrebovať do: rok a mesiac		
lv	återanvänd inte	lv	Använd före: år och månad		
lt	Tekrai kullanmayın	lt	Son kulanım: yıl ve ay		
fi	Mitte taaskasutada	fi	Kasutada enne: aasta ja kuu		
hu	Nesterilizēt atkārtoti	hu	Izlietot līdz (gads un mēnesis)		
et	Nenaudoti pakartotinai	et	Sunaudoti iki: metai ir mēnuo		

	<p>en fr de it es pt ελ cs da nl ro sk sv tr lv lt fi hu et</p> <p>Keep in a dry place A stocker dans un endroit sec Trocken lagern Conservare in un luogo asciutto Manténgase en un lugar seco Manter em lugar seco Φυλάξτε το προϊόν μακριά από την υγρασία Uchovajte na suchém místě Opbevares på et tørt sted Säilytä kuivassa Depozitați în loc uscat Op een droge plaats bewaren A se păstra într-un loc uscat Uchovajte na suchom mieste Förvara på ett torrt ställe Kuru bir yerde saklayın Hoida kuivas Uzglabāt sausā vietā Laikykite sausoje vietoje</p>		<p>en fr de it es pt ελ cs da nl ro sk sv tr lv lt fi hu et</p> <p>Keep away from sunlight A stocker à l'abri de la lumière du soleil Vor Sonnenlicht geschützt lagern Conservare al riparo dalla luce solare Mantener alejado de la luz solar Manter ao abrigo da luz solar Προφυλάξτε το προϊόν από το ηλιακό φως Uchovajte mimo dosah světla Opbevares mørkt Säilytä valolta suojaattuna Fénytől védendő Vrij van licht houden A se feri de lumină Uchovajte mimo dosahu svetla Håll dig borta från solljus Işıksız ortamda saklayın Hoida eemal päikesevalgusest Sargāt no gaismas Laikyti toliau nuo Saulės spindulių</p>
	<p>en fr de it es pt ελ cs da nl ro sk sv tr lv lt fi hu et</p> <p>Do not reuse Ne pas réutiliser Nicht wiederverwenden Non riutilizzare No volver a utilizar Não reutilizar Να μην επαναχρησιμοποιείται Nepoužívat znova Má ikke genbruges Älä käytä uudelleen Ne használja újra Niet hergebruiken A nu se reutiliza Nepoužívať znova återanvänd inte Tekrar kullanmayın Mitte taaskasutada Nesterilizēt atkārtoti Nenaudoti pakartotinai</p>		<p>en fr de it es pt ελ cs da nl ro sk sv tr lv lt fi hu et</p> <p>Do not use if the packaging is damaged Ne pas utiliser si l'emballage est endommagé Nicht verwenden, wenn Packung beschädigt ist Non utilizzare se la confezione è danneggiata No usar si el embalaje está dañado Não usar se a embalagem estiver danificada Να μην χρησιμοποιείται, εάν η συσκευασία έχει υποστεί κάποια βλάβη Nepoužívejte je-li obal poškozený Ne használja, ha sérült a csomagolás Må ikke bruges, hvis emballagen er beskadiget Älä käytä, mikäli pakkaus on vaurioitunut Niet gebruiken indien de verpakking beschadigd is A nu se utiliza dacă ambalajul este deteriorat Nepoužívajte ak je obal poškodený Använd inte om förpackningen är skadad Paket hasar görmüşse kullanmayın Mitte kasutada, kui pakend on kahjustatud Nelietot, ja iepakojums bojāts Nenaudoti, jei pakuotė pažeista</p>
	<p>en fr de it es pt ελ cs da nl ro sk sv tr lv lt fi hu et</p> <p>Double sterile barrier system Double barrière stérile Doppel-Sterilbarriere Doppia barriera sterile Doble barrera estéril Sistema duplo de barreira estéril Διπλό αποστειρωμένο φράγμα Dvojitá sterilní bariéra Dobbel steril skærm Dubbele steriele barrière Dublă barieră sterilă Dvojité sterilná bariéra Dubbel sterilbarriär Çifte steril bariyer Dubulta sterila barjera Dvigubas sterilus barjeras Valmistuspäivämäärä Gyártási dátum Kahekordne sterilne barjääär</p>		<p>en fr de it es pt ελ cs da nl ro sk sv tr lv lt fi hu et</p> <p>Medical device Dispositif medical Medizinprodukt Dispositivo medico Producto sanitario Dispositivo médico Ιατρική συσκευή Zdravotnický prostředek Medicinsk udstyr Medisch hulpmiddel Dispozitiv medical Zdravotnícka pomôcka Medicinteknisk produkt Tibbi dízenek Medicīniska ierīce Medicinos priemonē Lääkintälaitte Orvostechnikai eszköz Meditsiiniline seade</p>

	en fr de it es pt ελ cs da nl ro sk sv tr lv lt fi hu et	Manufacturing date Date de fabrication Herstellungsdatum Data di fabbricazione Fecha de fabricación Data de fabrico Ημερομηνία κατασκευής Datum výroby Fremstillingsdato Fabricagedatum Data fabricatie Dátum výroby Tillverkningsdatum Üretim tarihi Ražošanas datums Paganinimo data Valmistuspäivämäärä Gyártási dátum Tootmiskuupäev		en fr de it es pt ελ cs da nl ro sk sv tr lv lt fi hu et	MR safe Compatible IRM MRT-kompatibel Compatibile con MRI Compatible con IRM Compatível IRM Συμβατό με MRI Kompatibilní s MRI IRM kompatibel MRI-compatibel Compatibil RMN Kompatibilné s MRI MRT-kompatibel MR uyumu MR saderīgs Suderinamas su MRT Yhteensoviva MRI MRI-kompatibilis Kokkusobiv MRT
	en fr de it es pt ελ cs da nl ro sk sv tr lv lt fi hu et	Unique Device Identifier Numéro d'identifiant unique Einmalige Identifikationsnummer Identificazione unica del dispositivo Número de identificación única Número de identificação único Μοναδικός αριθμός αναγνώρισης Jedinečné identifikační číslo Enligt brugeridentifikationsnummer Uniek identificatienummer Număr de identificare unic Jedinečné identifikačné číslo Unikt identifieringsnummer Bireysel kimlik numarası Unikāls identifikācijas numurs Unikalus identifikavimo numeris Ainutkertainen tunnistenumero Egyedi azonosító szám Kordumatu tunnuskoodi number		en fr de it es pt ελ cs da nl ro sk sv tr lv lt fi hu et	Nombre d'implant dans l'emballage = Qty: Number of implants in the packaging = Qty : Anzahl der Implantate in der Packung = Anz.: Numero di impianti nella confezione = Qtà : Número de implantes en el embalaje = Cant.: Número de implantes na embalagem = Qtd.: Αριθμός εμφυτευμάτων στην συσκευασία = Ποσότητα: Počet implantátů v balení = Množ.: Antal implantater i pakken = antal: Aantal implantaten in de verpakking = aantal Numărul de implanturi în ambalaj = Cantitate Počet implantátů v balení = Množ.: Antal implantat i förpackningen = Kvantitet: Ambalajda bulunan implant sayısı=Miktar: Implantu skaits iepakojumā = daudz.: Implantu skaičius pakuotēje = Kiekis: Implantien lukumäärä pakauksessa = Määrä: A csomagolásban lévő implantátumok száma = Menny.: Implantaatide kogus pakendis = Arv:

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Cousin Biotech SAS au capital de 320 000 € - NAF 1396Z- N° TVA FR 34 398 460 261 - Siren 398 460 261 RCS Lille