QUICK ANCHOR BIOLIG®

NON RESORBABLE
REINFORCING JOINT LIGAMENT WITH AN ANCHORING SYSTEM

en	Instructions for use	Page	2
fr	Notice d'instructions	Page	4
de	Gebrauchsanweisung	Seite	6
it	Istruzioni per l'uso	Pagina	8
es	Instrucciones per l'uso	Pagina	10
pt	Nota de instruções	Pagina	12

COUSIN BIOTECH

Allée des Roses
59117 Wervicq-Sud -
FRANCE

Tél. : +33 (0) 3 20 14 41 20

Fax : +33 (0) 3 20 14 40 13

www.cousin-surgery.com

Made in France

CE
1639

NOT077/240712
Version du 12/07/2024

QUICK ANCHOR *BIOLIG*®
REINFORCING JOINT LIGAMENT WITH AN ANCHORING SYSTEM
STERILE PRODUCT FOR ONE USE

DESCRIPTION

The purpose of ligaments is to surgically reinforce instabilities in the various articulations. They may be fixed in an intra or extra-articular basis. Ligament reinforcement should always be used to complement suture or grafting of the natural ligament..

MATERIALS

Anchors: TA6V ELI ISO 5832-3

Ligament: Polyethylene terephthalate

Needles: Stainless steel (5832-1) (contact<24 hours – Limited contact)

Handle: Polyoxymethylene and Stainless steel (5832-1) (not implanted)

Nor human neither animal origin

INDICATIONS

Reinforcing ligaments for orthopedic surgery

POSSIBLE ADVERSE SIDE EFFECTS

- Inflammatory reactions
- Permanent elongation of the ligament
- Rupture of the ligament

WARNINGS

Inspect the packaging to be sure it is intact; Do not use if the labels and/or protective packaging are damaged

Do not use if the device is damaged.

DO NOT REUSE

DO NOT RESTERILIZE

Avoid using the device in a septic environment.

Do not implant in an infected site or in the case of allergy to the materials.

COUSIN BIOTECH advises the use of the provided ancillary. This one is of single use and is not re-sterilization.

COUSIN BIOTECH does not make any guarantees or recommendations as regards the use of a particular brand of fixation devices. The properties of these fastenings can be the subject to modifications from the manufacturer, which cannot be controlled by COUSIN BIOTECH. It is necessary to avoid fasten or the screws on sharp board.

It is not advised to implant the device in a child while he is still growing.

This device should be implanted **exclusively** by a qualified surgeon who has been trained in the use of this product.

Attention with the sterile handle of one use.

STORAGE

To be stored in a dry place away from sunlight, at room temperature, in its original packaging

EXPLANTATION AND DISPOSAL OF DEVICES

Explantation and handling should be carried out in accordance with ISO 12891-1:2015 “Implants for surgery – Retrieval and analysis of surgical implants” Part 1: “Retrieval and Handling”. Any explanted device must be returned for analysis as per the current protocol. This protocol is available from COUSIN BIOTECH on request. It is important to note that any implant that has not been cleaned and disinfected prior to dispatch must be placed in a sealed package. Explanted medical devices must be disposed of in accordance with local regulations governing the disposal of infectious hazardous waste. The disposal of non-implanted devices is not subject to specific recommendations. The second paragraph covers information requests or complaints.

INFORMATION REQUESTS AND COMPLAINTS

In accordance with its quality policy, COUSIN BIOTECH is committed to producing and supplying high-quality medical devices. However, if health professionals (client, user, prescriber, etc.) wish to lodge a complaint or report dissatisfaction with a given product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In the event of implant failure or if the implant is implicated in serious adverse effects experienced by the patient, the health centre must follow the legal procedures in the country in question, and inform COUSIN BIOTECH without delay. The reference, batch number, reference contact details and a comprehensive description of the incident or complaint must be stated in all correspondence. Brochures, documentation and surgical techniques are available on request from COUSIN BIOTECH and its distributors.

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