# ADHESIX® BIORING®

## ADJUSTABLE GASTRIC BANDING WITH ADHESIVE PORT

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Version du 11/07/2024

## ADHESIX® BIORING® DEVICE

### Adjustable gastric banding with adhesive port

#### **DESCRIPTION**

The **ADHESIX® BIORING®** adjustable gastric band is composed of:

- A radiopaque body comprising external reinforcement and an inflatable inner part. The extremities of the body form the closure system of the band after positioning it around the stomach.
- A radiopaque injection port, which is usually placed in the left subcostal or presternal region. This site allows the percutaneous adjustment of the diameter of the adjustable gastric band.
- A radiopaque catheter, which links the body of the adjustable gastric band to the injection port. The system is non-absorbable.
- A stiffener to secure the connection between the injection port and the catheter via a click-shut system.

#### **COMPOSITION**

The body comprises dimethylsiloxane. The external reinforcement is made from polyester which is fully impregnated with dimethylsiloxane, barium sulphate and a biocompatible dye.

The catheter is made from dimethylsiloxane loaded with barium sulphate, and printed with biocompatible ink.

The injection port is composed of polyoxymethylene with barium sulphate. The septum is made from dimethylsiloxane. A textile component is fixed to the port by means of a polyethylene terephthalate thread. The woven fabric/mesh comprises polypropylene coated with an adhesive resorbable hydrogel (polyvinylpyrrolidone – polyethylene glycol).

The stiffener comprises dimethylsiloxane and barium sulphate.

**No materials of human or animal origin were used. Semi-resorbable.**

#### **INDICATIONS**

The **ADHESIX® BIORING®** adjustable gastric band system is indicated for the reduction of excess body weight in persons presenting morbid obesity and weight loss resistance, i.e. persons whose Body Mass Index [BMI, formula = weight in kilograms/(height in metres)<sup>2</sup>] is greater than or equal to 40 kg/m<sup>2</sup>, or 35 to 40 kg/m<sup>2</sup> when combined with at least one comorbidity.

The surgeon can prescribe the adjustable gastric band at his/her own discretion.

#### **MRI SAFETY:**

The implant has not been tested for safety and compatibility with postoperative Magnetic Resonance Imaging (MRI). However, the implant mostly comprises dimethyl siloxane and non-metallic components deemed non-ferromagnetic in these quantities according to the literature. Moreover, the implant is presented as MRI safe in the mrisafety.com database with 1.5 and 3 Tesla. Consequently, the Adhesix® Bioring® device is considered MRI safe.

#### **CONTRAINDICATIONS**

- Inflammatory disease or congenital or acquired anomalies of the gastrointestinal tract,
- Severe organic disease,
- Potential bleeding disorder of the upper gastrointestinal tract or patients who regularly or systematically take aspirin or non-steroidal anti-inflammatory drugs,
- Diagnosis or pre-existing symptoms or family history of connective tissue autoimmune disease,
- Pregnancy,
- Patients with a drug or alcohol addiction,
- Patients who refuse to accept the dietary restrictions imposed by the procedure.
- Serious psychiatric disorders.
- Do not implant the device at the site of infection or in a septic environment.
- Do not implant the device if allergic to the components.
- Do not implant the device in growing children.

#### **POSSIBLE COMPLICATIONS**

Complications that may arise from the use of an adjustable gastric band include risks associated with the anaesthesia, surgical technique and the patient's level of intolerance to any foreign body implanted in the body.

#### **Peri-operative complications (due to the surgical technique):**

- Injury to the spleen
- Perforation of the stomach or oesophagus
- Other types of bleeding

#### **In the immediate postoperative period:**

Nausea, vomiting (especially in the first few days post-procedure and when the patient eats more than recommended).	Thrombo-embolic complications
Bloating due to gas	Infection-related complications
Dysphagia, dehydration	Postoperative bleeding

#### **Delayed postoperative complications:**

Band dilation	Band intolerance	Pyrosis
Band erosion	Pain at the port implantation site	Gastritis of the upper pouch
Oesophageal dilation	Band infection	Band migration
Pouch dilation	Band slippage	Displacement of the injection port
Band alteration	Gastroesophageal reflux	Infection of the injection port
Constipation		
Constrictive obstruction (due for example to oedema, food, inadequate calibration, band slippage or to a twisted or dilated upper pouch),		
Deflation of the band (due to leakage in the band, injection port or connection tube),		
Rapid weight loss may lead to symptoms of malnutrition, anaemia and related complications,		

#### **EXPLANTATION**

Explantations and handling should be carried out in accordance with ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling". Any explanted device must be returned for analysis as per the current protocol. This protocol is available from COUSIN BIOTECH on request. It is important to note that any implant that has not been cleaned and

disinfected before dispatch must be placed in a sealed package. Explanted medical devices must be disposed of in accordance with local regulations governing the disposal of infectious hazardous waste. The disposal of non-implanted devices is not subject to specific recommendations. The second paragraph covers information requests or complaints.

### **INFORMATION REQUEST AND COMPLAINTS**

In accordance with its quality policy, COUSIN BIOTECH is committed to producing and supplying high-quality medical devices. However, if health professionals (client, user, prescriber, etc.) wish to lodge a complaint or report dissatisfaction with a given product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In the event of implant failure or if the implant is implicated in serious adverse effects experienced by the patient, the health centre must follow the legal procedures in the country in question, and inform COUSIN BIOTECH without delay. The reference, batch number, reference contact details and a comprehensive description of the incident or complaint must be stated in all correspondence. Brochures, documentation and surgical techniques are available on request from COUSIN BIOTECH and its distributors.

### **WARNING**

The laparoscopic implantation technique is a state-of-the-art procedure. Surgeons who intend to perform this procedure should possess the necessary equipment and have extensive experience in laparoscopy.

### **PRECAUTIONS**

- Store in a dry place, away from light and at room temperature in the original packaging.
- Check the integrity of the packaging (do not use if the labels and/or protective packaging are damaged).
- Do not use if the device is damaged or has passed its expiry date.

**DO NOT REUSE**

**DO NOT RE-STERILIZE**

### **SURGICAL TECHNIQUE**

#### **Material provided**

- |  |   |
|--|---|
| <ul style="list-style-type: none"> <li>- One adjustable, radiopaque gastric band, with closure system</li> <li>- One syringe</li> <li>- One Huber needle</li> <li>- A gastric tube is available separately for small stomach calibration.</li> </ul> | <ul style="list-style-type: none"> <li>- One adhesive radiopaque injection port</li> <li>- One radiopaque catheter flushing device</li> <li>- One handling tag</li> <li>- One radiopaque stiffener</li> </ul> |
|--|---|

#### **Preparation of the material: All of the procedures should be performed in a sterile area.**

- Flush the band with physiological saline using the syringe and catheter rinsing device provided to completely remove the air from the device.
- Inspect the inflatable part of the band to detect any irregular inflation or leaks. When handling the band, make sure that it does not come into contact with sharp or pointed instruments.
- Clamp the catheter to prevent air from re-entering the device.
- Place the handling tag at the free end of the catheter.
- Flush the implantable injection port with physiological saline using the syringe and Huber needle.

Furthermore, to maximize the internal diameter of the band, it is advisable to fully aspirate the contents of the adjustable gastric band using the syringe and catheter-rinsing device, especially during procedures involving the "pars flaccida" technique, and in patients with a significant quantity of perigastric adipose tissue.


#### **Surgical technique:**

Prophylactic antibiotics: Antibiotic prophylaxis for skin and gut flora should be administered prior to beginning the procedure.

Various surgical techniques may be used, including: the perigastric technique, "pars flaccida" technique or a combination of both.

Adjustable gastric bands are usually placed laparoscopically, under general anaesthesia:

- The anaesthetist introduces the cuffed gastric tube into the stomach, for small stomach calibration.
- Dissect the lesser curvature or the right crus.
- 2nd window next to the His angle.
- Retro-gastric tunnelling.
- Introduce the adjustable gastric band into the abdomen. When passed through the trocar (12 mm minimum), the clip must be positioned on the reinforced part of the buckle.
- Inflate the gastric tube cuff with 15 cm<sup>3</sup> (ml) of physiological saline.
- Correctly position the band.
- Suture the band, using 3 to 6 sutures, to the gastric serosa, upstream and downstream (anterior gastro gastric suture enclosing the band in a "tunnel").
- Connect the tube to the injection port.
- General port placement recommendations:
  - Ensure that the catheter does not bend sharply.
  - When connecting the catheter to the port, avoid any twisting/spring effect on the catheter: check that the port is totally flat on a surface.
  - Ensure that the stiffener is properly clipped to the port (click-shut system).
  - In the abdominal wall, the path of the catheter must be oblique as opposed to perpendicular to the port. This prevents any lifting effect on the port, which could lead to port rotation over time:

**Incorrect catheter pathway:** the catheter may generate a lifting force on the port and lead to port rotation.

**Correct catheter pathway:** the catheter follows an oblique path through the muscle and generates no lifting force on the port which remains flat on the muscle.

The injection port is self-adhesive. Sutures are not required. Simply position it on the muscle:

- insert into the aponeurosis of the rectus abdominis muscle under the left costal edge,

- insert subcutaneously under the left costal edge,
- insert into the presternal region, on the xiphoid apophysis;
- insert in any other position at the surgeon's discretion.

The bands can also be implanted via laparotomy.

The band can be inflated immediately using physiological saline (0.9% NaCl ). ONLY use a **Huber needle** to inflate or deflate the band. The volume as well as the pre-inflation interval are at the discretion of the practitioner/surgeon (see the table below for maximum volume).

Band inflation: maximum volume of physiological saline according to band size

<b>Band sizes (diameter: mm)</b>	<b>Maximum volume of physiological saline (cm<sup>3</sup>)</b>
Ø 21/24	10
Ø 27/29	10
Ø 36/38	15

A typical band inflation algorithm is available on request.

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[Inhoudsopgave](#)

## Patient information leaflet ADHESIX® BIORING®

Item	Patient information leaflet
1	<p>a. <b>The name of the device</b> : ADHESIX® BIORING®</p> <p>b. <b>The model of the device</b> : BCBRI23ADU, BCBRI20ADU, BCBRIXLADU</p>
2	<p>a) The intended purpose of the device : ADHESIX® BIORING® is intended to be used as an adjustable gastric band for the treatment of obesity. It is indicated for the reduction of excess body weight in persons presenting morbid obesity and weight loss resistance.</p> <p>b) The kind of patient on whom the device is intended to be used : Persons presenting morbid obesity and weight loss resistance, i.e. persons whose Body Mass Index [BMI, formula = weight in kilograms/(height in metres)<sup>2</sup>] is greater than or equal to 40 kg/m<sup>2</sup>, or 35 to 40 kg/m<sup>2</sup> when combined with at least one comorbidity. The surgeon can prescribe the adjustable gastric band at his/her own discretion.</p>
3	<p><b>Special operating instructions for the use of the device</b> : No special operating instructions for the use of the device</p>
4	<p>a. <b>The intended performance of the device</b> :</p> <ul style="list-style-type: none"> <li>- The port is an X-Ray-opaque injection site and allows per-cutaneous adjustment of the gastric banding diameter.</li> <li>- The connection system clips on to secure the connection between the port and the catheter. The injection port is self-adhesive. Sutures are not required. Simply position it on the muscle as follows:</li> <li>- insert into the aponeurosis of the rectus abdominis muscle under the left costal edge,</li> <li>- insert subcutaneously under the left costal edge,</li> <li>- insert into the presternal region, on the xiphoid apophysis,</li> <li>- insert in any other position at the surgeon's discretion.</li> </ul> <p>Plug: In the event of an infection at the junction between the port and the catheter, the port is replaced with the plug until a new port is fitted. The plug is connected to the free end of the catheter of the ADHESIX® BIORING® adjustable gastric banding.</p> <p>Connector: When the catheter on the ADHESIX® BIORING® adjustable gastric banding is damaged (cut, leakage), the kit connector and catheter can be used to replace the damaged system.</p> <p>Always refer to the instructions for the ADHESIX® BIORING® adjustable gastric banding marketed by COUSIN BIOTECH.</p> <p>The product must be implanted only by a qualified surgeon experienced in using the ADHESIX® BIORING® System.</p>

This leaflet was revised in November 2022

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Item	Patient information leaflet
4	<p>b. <b>Undesirable side effects :</b>            Complications that may arise from the use of an adjustable gastric band include risks associated with the anaesthesia, the surgical technique, and the degree of intolerance of the patient to any foreign body implanted in the body.  <b>Peri-operative complications (due to the surgical technique) :</b> Injury to the spleen, Perforation of the stomach or oesophagus, Other types of bleeding  <b>In the immediate postoperative period :</b> Nausea, vomiting (especially in the first few days post-procedure and when the patient eats more than recommended), thrombo-embolic complication, bloating due to gas, dysphagia, dehydration, infection-related complications, postoperative bleeding  <b>Delayed postoperative complications :</b>            Band dilation, Band intolerance, Pyrosis, Band erosion, Pain at the port implantation site, Gastritis of the upper pouch, Oesophageal dilation, Band infection, Band migration, Pouch dilation Band slippage, Displacement of the injection port, Band alteration, Gastroesophageal reflux, Infection of the injection port, Constipation, Constrictive obstruction (due for example to oedema, food, inadequate calibration, band slippage or to a twisted or dilated upper pouch), Deflation of the band (due to leakage in the band, injection port or connection tube), Rapid weight loss may lead to symptoms of malnutrition, anaemia and related complications.</p>
5	<p><b>When the ADHESIX® BIORING® Should Not be Used :</b>            Do not use the implant in the following circumstances :</p> <ul style="list-style-type: none"> <li>- - Inflammatory disease or congenital or acquired anomalies of the gastrointestinal tract,</li> <li>- - Severe organic disease,</li> <li>- - Potential bleeding disorder of the upper gastrointestinal tract or patients who regularly or systematically take aspirin or non-steroidal anti-inflammatory drugs,</li> <li>- - Diagnosis or pre-existing symptoms or family history of connective tissue autoimmune disease,</li> <li>- - Pregnancy,</li> <li>- - Patients with a drug or alcohol addiction,</li> <li>- - Patients who refuse to accept the dietary restrictions imposed by the procedure.</li> <li>- - Serious psychiatric disorders.</li> <li>- - Do not implant the device at the site of infection or in a septic environment.</li> <li>- - Do not implant the device if allergic to the components.</li> <li>-- Do not implant the device in growing children.</li> </ul>

This leaflet was revised in November 2022

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Item	Patient information leaflet								
6	<p><b>MRI Information :</b> The implant has not been tested for safety and compatibility with postoperative Magnetic Resonance Imaging (MRI). However, the implant mostly comprises dimethyl siloxane and non-metallic components deemed non-ferromagnetic in these quantities according to the literature. Moreover, the implant is presented as MRI safe in the mrisafety.com database with 1.5 and 3 Tesla. Consequently, the Adhesix® Bioring® device is considered MRI safe..</p>								
7	<p><b>Precautions :</b> The band can be inflated immediately using physiological saline (0.9% NaCl). ONLY use a Huber needle to inflate or deflate the band. The volume and pre-inflation interval are at the discretion of the practitioner/surgeon (see the table below for maximum volume). Band inflation: maximum volume of physiological saline according to band size:</p> <table border="1" data-bbox="368 981 1302 1155"> <thead> <tr> <th>Band sizes (diameter : mm)</th> <th>Maximum volume of physiological saline (cm3)</th> </tr> </thead> <tbody> <tr> <td>Ø 21/24</td> <td>10</td> </tr> <tr> <td>Ø 27/29</td> <td>10</td> </tr> <tr> <td>Ø 36/38</td> <td>15</td> </tr> </tbody> </table> <ul style="list-style-type: none"> <li>- A typical band-inflation algorithm is available on request.</li> <li>- Store in a dry place, away from light and at room temperature in the original packaging.</li> <li>- Check the integrity of the packaging (do not use if the labels and/or protective packaging are damaged)</li> <li>- Do not use if the device is damaged or has passed its expiry date.</li> <li>- Do not re-use</li> <li>- Do no resterilize</li> </ul>	Band sizes (diameter : mm)	Maximum volume of physiological saline (cm3)	Ø 21/24	10	Ø 27/29	10	Ø 36/38	15
Band sizes (diameter : mm)	Maximum volume of physiological saline (cm3)								
Ø 21/24	10								
Ø 27/29	10								
Ø 36/38	15								
8	<p><b>Materials included in the device :</b> The body comprises dimethylsiloxane. The external reinforcement is made from polyester which is fully impregnated with dimethylsiloxane, barium sulphate and a biocompatible dye. The catheter is made from dimethylsiloxane loaded with barium sulphate, and printed with biocompatible ink. The injection port is composed of polyoxymethylene with barium sulphate. The septum is made from dimethylsiloxane. A textile component is fixed to the port by means of a polyethylene terephthalate thread. The woven fabric/mesh comprises polypropylene coated with an adhesive resorbable hydrogel (polyvinylpyrrolidone – polyethylene glycol). The stiffener comprises dimethylsiloxane and barium sulphate. No materials of human or animal origin were used. Semi-resorbable.</p>								

This leaflet was revised in November 2022

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Item	Patient information leaflet
9	<p><b>Reporting a serious incident :</b>            If you encounter any serious incident with your implant ADHESIX® BIORING® you should report it immediately to your doctor. You should also report the incident to Cousin Biotech via the link <a href="mailto:serviceclients@cousin-surgery.com">serviceclients@cousin-surgery.com</a> or <a href="mailto:materiovigilance@cousin-biotech.com">materiovigilance@cousin-biotech.com</a> and to the Therapeutic Goods Administration using the following website: <a href="http://www.tga.gov">http://www.tga.gov</a></p>
10	<p><b>Legal Manufacturer :</b> Cousin Biotech, Allée des roses , 59117 Wervicq-Sud, France</p> <p><b>Manufacturing Facility :</b> Cousin Biotech, Allée des roses, 59117 Wervicq-Sud, France</p>

This leaflet was revised in November 2022

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Cousin Biotech SAS au capital de 320 000 € - NAF 1396Z- N° TVA FR 34 398 460 261 - Siren 398 460 261 RCS Lille