# 4D Ventral®

SEMI RESORBABLE PARIETAL REINFORCEMENT IMPLANT  
EXTRAPERITONEAL

---

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**CE**  
0297

**NOT302\_240215**  
Version du 15/02/2024

**R<sub>x</sub>**  
only Caution : Federal law (USA) restricts this device to sale,  
distribution and use by or on the order of a physician.

This release is the last update of the instructions of use and replace the previous edition

Product under regulation 2017/745 (EU).

**SEMI RESORBABLE PARIETAL REINFORCEMENT IMPLANT  
EXTRAPERITONEAL  
SINGLE USE STERILE PRODUCT**

**MANUFACTURER IDENTIFICATION**

COUSIN BIOTECH  
ALLEE DES ROSES  
59117 WERVICQ SUD  
FRANCE

**IDENTIFICATION AND DESCRIPTION OF THE DEVICE**

4D Ventral® meshes are extraperitoneal semi-resorbable parietal reinforcement implants.

**INDICATIONS**

Repair of ventral hernia that require the addition of an extraperitoneal reinforcing or bridging material to obtain the desired surgical result.

**CONTRAINDICATIONS**

Do not use in the following cases:

- Allergy to any of the components
- Infected site
- Pregnancy
- Growing children

**UNDESIRABLE SIDE EFFECTS**

Like any implantable medical device, this implant is susceptible to generate possible undesirable side effects that may lead to reoperation:

- Discomfort/Pain
- Recurrence
- Hematoma
- Foreign body reaction
- Infection
- Adhesion
- Mesh migration
- Irritation nearby organ
- Inflammation
- Erosion/Extrusion
- Mesh deformation
- Obstruction
- Fistula
- Seroma/lymphocele
- Allergic reaction

**TARGET POPULATION**

Adult male or female suffering from a symptomatic ventral hernia necessitating a surgical repair with a method that includes an extraperitoneal implantation.

**IMPLANTED MATERIALS**

Polypropylene (PP) (non-resorbable) – Poly L Lactic Acid (PLLA) (resorbable)

	weight PP (g)	weight PLLA (g)
<b>4DVENT05RO</b>	0,12	0,18
<b>4DVENT07RO</b>	0,24	0,36
<b>4DVENT09RO</b>	0,39	0,59
<b>4DVENT12RO</b>	0,70	1,05
<b>4DVENT0715</b>	0,65	0,97
<b>4DVENT1015</b>	0,92	1,39
<b>4DVENT1515</b>	1,39	2,08
<b>4DVENT1520</b>	1,85	2,78
<b>4DVENT1530</b>	2,78	4,18
<b>4DVENT2020</b>	2,47	3,71
<b>4DVENT2025</b>	3,09	4,64
<b>4DVENT2535</b>	5,42	8,13
<b>4DVENT3030</b>	5,57	8,36
<b>4DVENT3040</b>	7,43	11,15

Origin neither human nor animal – Semi resorbable.

**PERFORMANCES**

The device 4D Ventral® contributes to reduce symptoms from hernia and contributes to avoid complications of hernia.

## **LINK TO SUMMARY OF SAFETY AND PERFORMANCE CHARACTERISTICS:**

[https://www.cousin-biotech.com/d-6318526d8b749\\_under-construction.pdf](https://www.cousin-biotech.com/d-6318526d8b749_under-construction.pdf)

## **MRI SAFETY**

The implant has not been tested for safety and compatibility in the case of postoperative investigation using Magnetic Resonance Imaging (MRI). However, the implant is mainly made with Polypropylene and Poly-L-Lactic Acid: non metallic components. Consequently, 4D Ventral® is considered MR safe.

## **PRECAUTIONS FOR USE**

4D Ventral® must always be separated from abdominal cavity by peritoneum

4D Ventral® implants are delivered sterile (ethylene oxide sterilisation).

When opening the carton box, the surgeon has to make sure that the lot number, the model and the size of the mesh are equal to the data on the label of the inner pouch. Otherwise, the use of the prosthesis is strictly prohibited.

Before any use, inspect the integrity of packaging and device (of which peelable pouches). Do not use in the event of deterioration of the device and/or the packaging. Do not use if the device is out of date. COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned. Mesh fixation points should be at least 1cm from the edge of the mesh with 1cm spacing between fixation points.

Unpacking and manipulation of the implant must be carried out following asepsis standards in order to guarantee the sterility.

Prevent the mesh from any contact with objects that may alter their surface.

This device must be implanted only by a qualified surgeon trained on the use of the product (knowledge of anatomy and visceral surgery).

## **IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE**

As specified on the product's labelling, the 4D Ventral® mesh is for single use only. It can not be re-used and/or re-sterilized (potential risks would be and are not limited to: loss of the product's sterility, risk of infection, loss of the product's efficiency, recurrence).

## **INSTRUCTIONS TO BE GIVEN TO THE PATIENTS BY THE SURGEON**

The surgeon should inform the patient of the potential physical and psychological restrictions and consequences of implanting the device. The patient must be informed of the surgical risks and possible side effects. The surgeon should invite the patient to return for further consultation if the patient presents symptoms that appear abnormal. An implant card for the patient is provided by the healthcare professional at the clinic. This card provides information to identify the device and traceability elements as well as the name, address and website of the manufacturer.

**The surgeon shall invite the patient to scan the implant card immediately upon receipt to keep track of it in case of loss.** The summary of the safety and performance characteristics of the device and the instructions for use can also be found on the COUSIN BIOTECH website.

## **STORAGE PRECAUTIONS**

To be stored in a dry place away from sunlight and at room temperature in its original packaging.


## **SURGICAL TECHNIQUE**

### **Extra-Peritoneal Hernia repair**

Example: Retrorectus mesh repair of a ventral hernia (Rives-Stoppa technique).

This technique utilizes the hernia sac to separate the mesh from the intra-abdominal contents.

- Superior to the umbilicus, dissection is performed above the posterior rectus fascia and under the rectus muscle.
- Below the umbilicus, dissection occurs in the preperitoneal space due to the lack of a posterior rectus sheath.
- The 4D Ventral® implant is placed in the newly formed pre-peritoneal space. The 4DVentral mesh has an oriented elasticity:

- Fixation of the mesh if needed, to the muscle layer above
- The peritoneum must be meticulously checked and no defect in the peritoneum should be left open to prevent any contact between the mesh and intra-abdominal contents.
- Redon's drains under aspiration are usually placed. The anterior sheath is closed. If tension is present, relaxing incisions may be used
- Skin closure.
- Using the Redon's drains, fill the pre-peritoneal space with an antiseptic solution and clamp the drains for 1 hour.
- A belt can be used to wrap the patient's abdomen for a few days to reinforce the abdominal wall.

## **EXPLANTATION AND DISPOSAL OF DEVICES**

Risks associated with the explantation: the surgeon who decides to remove the device should consider factors such as the risk of another surgery for the patient, and the difficulty of the explantation procedure. The removal of the implant is up to the surgeon's decision and must be subject to adequate postoperative follow-up.

Devices should be retrieved and handled according to the recommendations of ISO 12891-1:2015 “Implants for surgery – Retrieval and analysis of surgical implants” Part 1: “Retrieval and Handling”. Any explanted device must be sent back for analysis, following the current protocol. This protocol is available on request from COUSIN BIOTECH. It is important to note that any implant that must not be cleaned or disinfected prior to despatch must be contained in sealed packaging. The removed medical device must be disposed of in accordance with the country’s standards for the disposal of infectious waste.

There are no specific recommendations for the disposal of a non-implanted device.

### **IMPORTANT**

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

### **INFORMATION REQUEST AND COMPLAINTS**

Following its quality policy, COUSIN BIOTECH is committed to making every effort to produce and supply a high-quality medical device. However, if a health professional (client, user, prescriber, etc.) has a complaint or cause for dissatisfaction with a product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In case of failure of an implant or if it has contributed to causing a serious adverse reaction in the patient, the health centre must follow the legal procedures in that country and inform COUSIN BIOTECH immediately. For any correspondence, please specify the reference, batch number, contact person’s details and a comprehensive description of the incident or complaint.

Brochures, documentation and surgical technique are available on request from COUSIN BIOTECH and its distributors. If further information is needed or required, please contact your COUSIN BIOTECH representative or distributor.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

**IMPLANT DE RENFORCEMENT PARIETAL SEMI RESORBABLE  
EXTRAPERITONEAL  
PRODUIT STERILE A USAGE UNIQUE**

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