INTRASPINE®

DEVICE FOR THE INTER SPINOUS SPACE WITH LAMINAR SUPPORT

en	Instructions for use	Page	2
fr	Notice d'instructions	Page	6
de	Gebrauchsanweisung	Seite	11
it	Istruzioni per l'uso	Pagina	16
es	Instrucciones de uso	Página	21
pt	Instruções de uso	Página	26
ελ	Οδηγίες χρήσης	σελίδα	31
pl	Instrukcja obsługi	Strona	36
cs	Návod k použití	Strana	41
sl	Navodila za uporabo	Stran	46
ru	Инструкция по применению	Страница	51
kr	사용 지침	페이지	57
en	Patient information leaflet (Australian requirements)	Page	60

COUSIN BIOTECH

Allée des Roses
59117 Wervicq-Sud -
FRANCE

Tél. : +33 (0) 3 20 14 41 20

Fax : +33 (0) 3 20 14 40 13

www.cousin-surgery.com

Made in France

Ancillaires

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Date de marquage CE: mars 2007

NOT130 / 250416

Version du 16/04/2025

This release is the last update of the instructions of use and replace the previous edition

IntraSPINE®

DEVICE FOR THE INTER-SPINOUS SPACE WITH LAMINAR SUPPORT SINGLE USE STERILE PRODUCT

DESCRIPTION

The IntraSPINE® prosthesis is an interlaminar medical device that comprises a dimethyl siloxane (silicon) wedge (available in several sizes) adapted to the intervertebral space and a ligament in polyethylene terephthalate.

The wedge is covered with woven polyethylene terephthalate with an area covered by sheeting on the front face in the medullary zone.

Tension is applied to the ligament by a titanium loop (TA6V).

MATERIALS

IntraSPINE® prosthesis: Dimethyl siloxane, Polyethylene terephthalate, Titanium.

IntraSPINE® ancillaries and trial prostheses: Stainless steel, Acetal copolymer, Dimethyl siloxane.

Non-human and non-animal origin – Non-resorbable.

INDICATIONS

The IntraSPINE® implant is recommended for use in the lumbar spine including in the lumbosacral area (from L1 to S1) in the following cases:

- Arthropathic facet-syndrome
- Foraminal stenosis
- Degenerative discopathy
- Intervertebral ligament insufficiency

PERFORMANCE

The IntraSPINE® prosthesis is an elastic support device with a cushioning effect. Its interlaminar implantation enables it to be close to the centre of rotation, thereby producing the cushioning effect and providing tension on the interspinous ligamentary structures due to the elasticity of the supra-spinous wedge.

The principal function of the ligaments surrounding the upper and lower vertebrae in relation to the area receiving the prosthesis is to assist the posterior ligament during flexion.

IntraSPINE® has been validated for a long term implantation in the human body.

CONTRAINDICATIONS

1) Do not use the prosthesis in the following cases:

- Allergy to one of the components
- Growing child
- Infected site
- Pregnant woman

2) Factors likely to compromise successful implantation:

- Severe osteoporosis
- Major deformations of the spinal column
- Local bone tumors
- Systemic or metabolic disorders
- Infectious diseases
- Obesity
- Drug addictions
- Intense physical activity, for example competitive sports or hard labor
- Surgery at more than 2 levels with IntraSPINE Prosthesis

UNDESIRABLE SECONDARY EFFECTS

1) **Complications likely related to the device:**

- Infection of the prosthesis
- Device breakage (rupture of the ligament, fixed wing breakage, knitting tearing)
- Uncrimping of the needle
- Device loosening permanent elongation of the ligament
- Device migration displacement of the wedge
- Spondylolisthesis to the dynamically stabilized segment
- Persistent or worsening pain attributed to the dynamically stabilized level and required revision surgery with interbody fusion at that level
- Spinous process fractures intraoperative post-operative or during device removal (removal due to infection) and can be associated with dislocation of the implant
- Late reoperation due to absence of recovery
- Disability
- Discomfort
- Wrong level implantation, inappropriate or inaccurate positioning (posterior or wrong level)
- Seroma/swelling around the device
- Inflammatory reaction, allergy, foreign body reaction
- In certain cases the medical instrument had to be removed and pedicle screw fusion had to be performed.

2) **Complications unlikely to be related to the device, but related to surgery or to the patient's conditions:**

- Pseudomeningocele, fistula, dural tear, persistent CSF leakage, meningitis
- Loss of neurological function,
- Cauda equina syndrome, neuropathy, neurological deficits, arachnoiditis, and/or muscle loss, compression around nerves and/or pain
- Urinary retention or loss of bladder control or other types of urological system compromise
- Fracture, Bone loss or decrease in bone density
- Herniated nucleus pulposus, disc disruption or degeneration at, above, or below the level of surgery: Adjacent segment degeneration
- Non-union or pseudarthrosis, delayed union. Mal-union,
- Cessation of any potential growth of the operated portion of the spine
- Loss of or increase in spinal mobility or function
- Hemorrhage, hematoma, edema, embolism, stroke, excessive bleeding, phlebitis, wound necrosis, wound dehiscence, wound infection, damage to blood vessels, or other types of cardiovascular system compromise
- Development of respiratory problems

MRI SAFETY:

The implants are composed of non-ferromagnetic materials and present a geometry not likely to generate induced current. Moreover, as they are fixed to tissues, they are unlikely to be mobilised. A priori they can be considered compatible with an MRI scan. Their safety, in particular in terms

of heating and migration of the implant, has been evaluated through bibliographic data by comparison with data available on devices with similar composition, shape and use. This evaluation concluded a conditional compatibility for MRI scan up to 1.5 Tesla.

As a precautionary measure, it is recommended MRI scans be avoided within 48 hours of the implant placement, and to inform the person in charge of the scan of the recent implant placement, if such examination is essential.

It should be noted that devices presenting a high contrast with the biological environment can generate "artifacts" that must be taken into account for the perfect execution and interpretation of imaging exams.

For this purpose, the patient who has this implant must be advised to warn the health professionals concerned (radiologists and radiology operators) in as far as possible about the presence of this implant before these exams.

The Safety of the MR environment has not been tested.

POSTOPERATIVE PRECAUTIONS

The patient must be informed of the postoperative cares and of precautions.

The use of a brace could be necessary but is not mandatory. This decision is the surgeon one.

PRECAUTIONS FOR USE

Before usage, verify the integrity of the prosthesis, trial prostheses, ancillaries and packaging.

Do not use in the event of deterioration of the prosthesis and/or the labels, and/or the trial prostheses and/or the ancillaries and/or the packaging.

Do not use if the prosthesis has exceeded its use-by date.

Unpacking and manipulation of the implant must be carried out following asepsis standards in order to guarantee the sterility.

When opening the carton box, the surgeon has to make sure that the lot number, the model and the size of the prosthesis are equal to the data on the label of the inner pouch. Otherwise the use of the prosthesis is strictly prohibited.

Prevent the prosthesis, the trial prostheses and ancillaries from coming in to contact with objects that may alter their surface.

Operation conditions: Temperature from +32°C up to 42°C. Implants should be resistant to biologic fluids and the discharge of body tissues, which they contact during use. Therefore, the IntraSPINE® prosthesis has been validated for a long term implantation and, currently, COUSIN BIOTECH has 10 years of feedback concerning these products.

The IntraSPINE® prosthesis, trial prostheses and ancillaries must only be used by a qualified surgeon who is trained in spinal surgery and to the surgical technique of IntraSPINE®. The information contained in this notice is necessary but not sufficient to control the surgical technique.

The IntraSPINE® is made with soft materials; any contact with harder ones may alter its mechanical behaviour and lifetime. It is mandatory to avoid contact with other materials.

The prescription of the device is decided on by the surgeon, the only person qualified to do so.

Topping off is not contra-indicated but special attention is required. IntraSPINE® is made with soft materials; any contact with harder ones may alter its mechanical behaviour and lifetime. It is mandatory to avoid contact with other materials.

Use of the system of COUSIN BIOTECH ancillaries and specific trial prostheses is required for implantation of the IntraSPINE® prosthesis.

Implants must be handled and stored with great care, in a dry place, protected from sunlight and at room temperature.

Transportation: Transportation of the devices should be done in the closed vehicles in accordance with the rules for those transport vehicles.

There are no special transportation requirements for the device.

IMPORTANT : DO NOT REUSE AND DO NOT RESTERILISE THE INTRASPINE® PROSTHESIS


As indicated on the product labelling, the prosthesis is designed for single use. It should in no event be reused and/or resterilized (potential risks include but are not limited to: loss of product sterility, risk of infection, loss of effectiveness, relapse).

WARNING ABOUT ANCILLARIES AND TRIAL PROSTHESES

The IntraSPINE® ancillaries are Class I medical devices, intended for temporary use and are re-usable.

The IntraSPINE® trial prostheses are Class IIa medical devices, intended for temporary use and are re-usable.

Unlike the IntraSPINE® prosthesis, ancillaries and trial prostheses (included with the instrument kit/set) are furnished unsterile.

KITACINTRA / RCBININTRU

Before first use

Ancillaries and trial prostheses must be:

- Cleaned in the appropriate way
- Sterilised in an autoclave

DECONTAMINATION AND CLEANING OF ANCILLARIES

Decontamination and cleaning are performed under the responsibility of the healthcare centers.

The following prescribed methods and materials must be used to reduce the risk of transmission of UTA, (French health ministry - DGS/R13/2011/449 dated 01/12/2011).

This step has to take place before **the first utilization** and **immediately after the utilization** to avoid adhesion of particles or dry secretion to the instruments

The detachable instruments have to be dismantled

Products advised for the cleaning

- an enzymatic detergent neutral
- Or products adapted for the cleaning (neutral or soft alkaline) of surgical instruments in compliance with the regulations in effect.

- **Warning:** Don't use corrosive factor or caustic cleaning product
2 process are possible

A- Chart Automated cleaning in a disinfecter washer with manual Pre-cleaning:

	Step	Step Description	Step Instruction	Accessories	Duration
Cleaning Steps	1	Contamination Removal	Rinse product room temperature running tap water removing any visible organic material with assistance of a soft bristle brush	- Tap water - Soft bristle brush (Do not utilize metal cleaning brushes)	Until all visible soils removed
	2	Pre-soak	Prepare a detergent solution at ambient temperature (15-25°C / 59-77°F). Afterwards lumina, threads, joints and gaps of the instruments have to be flushed with a water jet pistol for a minimum time of 10 seconds for each position.	- Detergent - Water jet pistol	Until product is visually clean
	3	Ultrasonic cleaning	Clean the instrument in the ultrasonic tray including the solution. Be careful every unclench connector or screws and bolts by the vibrations will be reassured.	Ultrasonic tray	Minimum duration 10 minutes, this duration can vary according product
	4	Soak	Soak the instrument in tap water during 30 seconds. After, ultrasonic treatment the lumina, threads, joints and gaps of the instrument have to be flushed again with a water jet pistol for a minimum time of 10 seconds for each position.	Reverse osmosis water	Minimum 30 seconds of soak
	5	Automated Washer	Place entire disassembled device into the automated washer	- Automated Washer - Disassembled Device	Minimum total cycle time: 39 minutes

Instructions of cleaning in washer / disinfecter			
Cycles parameters	Time	Minimum temperature	Type of detergent/water
Pre-cleaning	2 minutes	<45°C / <113°F	Tap water
Cleaning	5 minutes	55°C / 131°F	Detergent
Draining	2 minutes	<45°C / <113°F	Tap water
Thermic draining	5 minutes	90°C / 194°F	Reverse osmosis water
Drying	25 minutes	> 70°C / > 158°F	Not applicable

Finishing Steps	6	Final Rinse	Rinse carefully with distilled water	Distilled water	Minimum 1 minute
	7	Final Drying	Dry device utilizing medical quality filtered air	Medical quality filtered air	Until product is visually dry
	8	Tidying	The instruments have to be placed in the right position inside the trays as described on the position map fixed on the tray for the sterilization.		

B- Chart Just manual cleaning :

	Step	Step Description	Step Instruction	Accessories	Duration
Cleaning Steps	1	Contamination Removal	Rinse product room temperature running tap water removing any visible organic material with assistance of a soft bristle brush	- Tap water - Soft bristle brush (Do not utilize metal cleaning brushes)	Until all visible soils removed
	2	Drying	Dry the device utilizing a dry non-shedding wipe. Medical quality filtered air may be utilized if available	- Non-shedding wipe - Medical quality filtered compressed air	Until product is visually dry
Disinfection Steps	3	Disinfection Application	Prepare a low-foaming neutral enzymatic detergent solution, using tap water (15-25°C / 59-77°F). Soak the instrument in the open position (if possible) during minimum 1 minute. During the immersion, actuate mobile parts of the instrument minimum 3 times (if possible), so that the detergent can access all parts of the instrument.	Detergent	Minimum duration 1 minute, this duration can vary according product
	4	Drying	After the soaking time, take out the instrument and wipe it with a disposable cloth. Then, place the instrument in a new bath of enzymatic detergent solution using warm tap water.	- Disposable cloth - Detergent	Until product is visually dry

Disinfection Steps	5	Manual Disinfection	Brush all surface of the instrument using a soft brush for 2 minutes. Adjust brushing time if needed. Actuate mobile parts of the instrument 3 times (when applicable) and use a syringe to access all difficult areas. Use a volume of 60 mL of the detergent solution.	- Soft bristle brush - Detergent - Syringe	Manual cleaning time duration is complete when the device's surface, joints, and crevices have been manually cleaned
	6	Rinse and finition	Soak the instrument in reverse osmosis water for 1 minute. Use a syringe and 60mL of reverse osmosis water to access all difficult areas. Repeat soaking 2 additional times using fresh water	- Reserve osmosis water - Syringe	Minimum 1 minute
	7	Final Rinse	Rinse carefully with distilled water	- Distilled water	Minimum 1 minute
	8	Final Drying	Dry devise utilizing medical quality filtered air	- Medical quality filtered air	Until product is visually dry
	9	Tidying	The instruments have to be placed in the right position inside the trays as described on the position map fixed on the tray for the sterilization.		

Note: In case of **suspection of Creutzfeld-Jakob Disease (CJD)**, incinerate the ancillaries and the fixation pliers

In case of no respect of the instructions mentioned above, the healthcare center will have to apply a validated equivalent cleaning process for which it will be fully responsible.

STERILIZATION OF ANCILLAIRES

Sterilization is performed under the responsibility of the healthcare centers.

The following prescribed methods and materials must be used to reduce the risk of transmission of UTA, (French health ministry - DGS/R13/2011/449 dated 01/12/2011).

1- The sterilization in autoclave is to make in a **specific container** (fenced basket, Tray, Plastic tray placed in a peelable pouches autoclavable).

NB: It isn't recommended to realize the sterilization with peelables pouches **autoclavable** without a container □ Ancillaries can pierce the films components the pouches.

2-A sterilization **BY AUTOCLAVE** is advised following this protocol:

Steam sterilization parameters	Values	Values
Cycle Type	Pre-vacuum	Pre-vacuum
Set Point Temperature	134°C / 273°F	132°C / 269.6°F
Cycle Time	Minimum 3 minutes	Minimum 4 minutes
Dry Time	Minimum 20 minutes	Minimum 20 minutes
Cool Time (inside and outside the autoclave chambers)	Minimum 40 minutes	Minimum 40 minutes

3 -Remark: After each cycle of cleaning/sterilization, insure of all the good working of the ancillaries: verify the integrality and the correct operation, of the locking systems (clipping, prehension) without excessive play.

4- If the ancillaries utilization is different, after sterilization stock them in a tray on a peelable pouches.

In case of no respect of the instructions mentioned above, the healthcare center will have to apply a validated equivalent sterilization process for which it will be fully responsible.

STORAGE PRECAUTIONS

Ancillaries and trial prostheses must be stored in their dedicated container or equivalent packaging to prevent any deterioration.

Keep in a dry place away from sunlight, at room temperature, in its original packaging.

EXPLANTATION AND DISPOSAL OF DEVICES

Explantation and handling should be carried out in accordance with ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling". Any explanted device must be returned for analysis as per the current protocol. This protocol is available from COUSIN BIOTECH on request. It is important to note that any implant that has not been cleaned and disinfected prior to dispatch must be placed in a sealed package. Explanted medical devices must be disposed of in accordance with local regulations governing the disposal of infectious hazardous waste. The disposal of non-implanted devices is not subject to specific recommendations. The second paragraph covers information requests or complaints.

INFORMATION REQUESTS AND COMPLAINTS

In accordance with its quality policy, COUSIN BIOTECH is committed to producing and supplying high-quality medical devices. However, if health professionals (client, user, prescriber, etc.) wish to lodge a complaint or report dissatisfaction with a given product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In the event of implant failure or if the implant is implicated in serious adverse effects experienced by the patient, the health centre must follow the legal procedures in the country in question, and inform COUSIN BIOTECH without delay. The reference, batch number, reference contact details and a comprehensive description of the incident or complaint must be stated in all correspondence. Brochures, documentation and surgical techniques are available on request from COUSIN BIOTECH and its distributors.

IMPLANT CARD

Be sure to scan your implant card immediately upon receipt so that you can keep trackt of it in case of loss.

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Patient information leaflet IntraSPINE®

Item	Patient information leaflet
1	a. The name of the device : IntraSPINE® b. The model of the device : INTRACAL08, INTRACAL10, INTRACAL12, INTRACAL14, INTRACAL16
2	a. The intended purpose of the device : The IntraSPINE® prosthesis is recommended in the lumbar spine including in the lumbosacral area (from L1 to S1 levels) in the following cases: - Arthropathic facet-syndrome - Foraminal stenosis - Degenerative discopathy - Intervertebral ligament insufficiency b. The kind of patient on whom the device is intended to be used : The target population is skeletally mature patients.
3	Special operating instructions for the use of the device : The IntraSPINE® prosthesis, trial prostheses and ancillaries must only be used by a qualified surgeon who is trained in spinal surgery and to the surgical technique of IntraSPINE®. The information contained in this notice is necessary but not sufficient to control the surgical technique.
4	a. The intended performance of the device : The IntraSPINE® prosthesis is an elastic support device with a cushioning effect. Its interlaminar implantation enables it to be close to the centre of rotation, thereby producing the cushioning effect and providing tension on the interspinous ligamentary structures due to the elasticity of the supra-spinous wedge. The principal function of the ligaments surrounding the upper and lower vertebrae in relation to the area receiving the prosthesis is to assist the posterior ligament during flexion. IntraSPINE® has been validated for a long term implantation in the human body. b. Undesirable side effects : • Complications likely related to the device - Infection of the prosthesis - Device breakage (rupture of the ligament, fixed wing breakage, knitting tearing) - Uncrimping of the needle - Device loosening permanent elongation of the ligament - Device migration displacement of the wedge - Spondylolisthesis to the dynamically stabilized segment - Persistent or worsening pain attributed to the dynamically stabilized level and required revision surgery with interbody fusion at that level - Spinous process fractures intraoperative post-operative or during device removal (removal due to infection) and can be associated with dislocation of the implant - Late reoperation due to absence of recovery

This leaflet was revised in March 2023

Page

63/72

Headquarters & Factory: **COUSIN BIOTECH** Allée des Roses - 59117 Wervicq-Sud – France **Tel** +33 (0) 3 20 14 41 20 –
Fax +33 (0) 3 20 14 40 13
www.cousin-surgery.com - **e-mail** f.pelletier@cousin-surgery.com

	- Disability - Discomfort - Wrong level implantation, inappropriate or inaccurate positioning (posterior or wrong level) - Seroma/swelling around the device - Inflammatory reaction, allergy, foreign body reaction - In certain cases the medical instrument had to be removed and pedicle screw fusion had to be performed. • Complications unlikely to be related to the device, but related to surgery or to the patient's conditions - Pseudomeningocele, fistula, dural tear, persistent CSF leakage, meningitis - Loss of neurological function, - Cauda equina syndrome, neuropathy, neurological deficits, arachnoiditis, and/or muscle loss, compression around nerves and/or pain - Urinary retention or loss of bladder control or other types of urological system compromise - Fracture, Bone loss or decrease in bone density - Herniated nucleus pulposus, disc disruption or degeneration at, above, or below the level of surgery: Adjacent segment degeneration - Non-union or pseudarthrosis, delayed union. Mal-union, - Cessation of any potential growth of the operated portion of the spine - Loss of or increase in spinal mobility or function - Hemorrhage, hematoma, edema, embolism, stroke, excessive bleeding, phlebitis, wound necrosis, wound dehiscence, wound infection, damage to blood vessels, or other types of cardiovascular system compromise - Development of respiratory problems
5	When IntraSPINE® Should Not be Used : 1) Do not use the prosthesis in the following cases: Allergy to one of the components, Growing child, Infected site, Pregnant woman 2) Factors likely to compromise successful implantation: Severe osteoporosis, Major deformations of the spinal column, Local bone tumors, Systemic or metabolic disorders, Infectious diseases, Obesity, Drug addictions, Intense physical activity, for example competitive sports or hard labor, Surgery at more than 2 levels with IntraSPINE® Prosthesis
6	MRI Information : The implants are composed of non-ferromagnetic materials and present a geometry not likely to generate induced current. Moreover, as they are fixed to tissues, they are unlikely to be mobilised. A priori they can be considered compatible with an MRI scan. Their safety, in particular in terms of heating and migration of the implant, has been evaluated through bibliographic data by comparison with data available on devices with similar composition, shape and use. This evaluation concluded a conditional compatibility for MRI scan up to 1.5 Tesla. As a precautionary measure, it is recommended MRI scans be avoided within 48 hours of the implant placement, and to inform the person in charge of the scan of the recent implant placement, if such examination is essential. It should be noted that devices presenting a high contrast with the biological environment can generate "artifacts" that must be taken into account for the perfect execution and interpretation of imaging exams.

This leaflet was revised in March 2023

Page

64/72

Headquarters & Factory: **COUSIN BIOTECH** Allée des Roses - 59117 Wervicq-Sud – France **Tel** +33 (0) 3 20 14 41 20 –

Fax +33 (0) 3 20 14 40 13

www.cousin-surgery.com - **e-mail** f.pelletier@cousin-surgery.com

6	For this purpose, the patient who has this implant must be advised to warn the health professionals concerned (radiologists and radiology operators) in as far as possible about the presence of this implant before these exams. The Safety of the MR environment has not been tested.
7	Precautions : Before usage, verify the integrity of the prosthesis, trial prostheses, ancillaries and packaging. Do not use in the event of deterioration of the prosthesis and/or the labels, and/or the trial prostheses and/or the ancillaries and/or the packaging. Do not use if the prosthesis has exceeded its use-by date. Unpacking and manipulation of the implant must be carried out following asepsis standards in order to guarantee the sterility. When opening the carton box, the surgeon has to make sure that the lot number, the model and the size of the prosthesis are equal to the data on the label of the inner pouch. Otherwise the use of the prosthesis is strictly prohibited. Prevent the prosthesis, the trial prostheses and ancillaries from coming in to contact with objects that may alter their surface. Operation conditions: Temperature from +32°C up to 42°C. Implants should be resistant to biologic fluids and the discharge of body tissues, which they contact during use. Therefore, the IntraSPINE® prosthesis has been validated for a long term implantation and, currently, COUSIN BIOTECH has 10 years of feedback concerning these products. The IntraSPINE® prosthesis, trial prostheses and ancillaries must only be used by a qualified surgeon who is trained in spinal surgery and to the surgical technique of IntraSPINE®. The information contained in this notice is necessary but not sufficient to control the surgical technique. The IntraSPINE® is made with soft materials; any contact with harder ones may alter its mechanical behaviour and lifetime. It is mandatory to avoid contact with other materials. The prescription of the device is decided on by the surgeon, the only person qualified to do so. Topping off is not contra-indicated but special attention is required. IntraSPINE® is made with soft materials; any contact with harder ones may alter its mechanical behaviour and lifetime. It is mandatory to avoid contact with other materials. Use of the system of COUSIN BIOTECH ancillaries and specific trial prostheses is required for implantation of the IntraSPINE® prosthesis. Implants must be handled and stored with great care, in a dry place, protected from sunlight and at room temperature. Transportation: Transportation of the devices should be done in the closed vehicles in accordance with the rules for those transport vehicles. There are no special transportation requirements for the device.
8	Materials included in the device: Dimethyl siloxane, Polyethylene terephthalate, Titanium.

This leaflet was revised in March 2023

Page

65/72

Headquarters & Factory: **COUSIN BIOTECH** Allée des Roses - 59117 Wervicq-Sud – France **Tel/** +33 (0) 3 20 14 41 20 –

Fax +33 (0) 3 20 14 40 13

www.cousin-surgery.com - **e-mail** f.pelletier@cousin-surgery.com

9	Reporting a serious incident : If you encounter any serious incident with your IntraSPINE® implant you should report it immediately to your doctor. You should also report the incident to COUSIN BIOTECH via the link serviceclients@cousin-surgery.com or materiovigilance@cousin-biotech.com and to the Therapeutic Goods Administration using the following website: http://www.tga.gov
10	Legal Manufacturer & Manufacturing Facility : COUSIN BIOTECH, Allée des Roses, 59117 Wervicq-Sud, France

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