## SEMI RESORBABLE PARIETAL REINFORCEMENT IMPLANT EXTRAPERITONEAL

|    |                                                                  |          |    |
|----|------------------------------------------------------------------|----------|----|
| en | <u>Instructions for use</u>                                      | Page     | 2  |
| fr | <u>Notice d'instructions</u>                                     | Page     | 4  |
| de | <u>Gebrauchsanweisung</u>                                        | Seite    | 6  |
| it | <u>Istruzioni per l'uso</u>                                      | Pagina   | 8  |
| es | <u>Instrucciones de uso</u>                                      | Pagina   | 10 |
| pt | <u>Instruções de utilização</u>                                  | Página   | 12 |
| ελ | <u>Οδηγίες χρήσης</u>                                            | Σελίδα   | 14 |
| bg | <u>Инструкции за употреба</u>                                    | страница | 16 |
| cs | <u>Návod k použití</u>                                           | Strana   | 18 |
| da | <u>Brugsvejledning</u>                                           | side     | 20 |
| fi | <u>Käyttöohjeet</u>                                              | Sivu     | 22 |
| nl | <u>Gebruiksaanwijzing</u>                                        | Pagina   | 24 |
| pl | <u>Instrukcja obsługi</u>                                        | Strona   | 26 |
| ro | <u>Instrucțiuni de utilizare</u>                                 | Pagină   | 28 |
| sk | <u>Návod na použitie</u>                                         | Strana   | 30 |
| lt | <u>Naudojimo instrukcijos</u>                                    | Puslapis | 32 |
| en | <u>Patient information leaflet<br/>(Australian requirements)</u> | Page     | 34 |

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Caution: Federal law (USA) restricts this device to sale,  
distribution and use by or on the order of a physician.

This release is the last update of the instructions of use and replace the previous edition

SEMI RESORBABLE PARIETAL REINFORCEMENT IMPLANT  
**EXTRAPERITONEAL**  
SINGLE USE STERILE PRODUCT

**DESCRIPTION**

4D Ventral® meshes are extraperitoneal partially resorbable parietal reinforcement implants.

**IMPLANTED MATERIALS**

Polypropylene (non-resorbable knitted monofilament) – Poly L Lactic Acid (resorbable knit)  
Origin neither human nor animal – Semi resorbable.

**INDICATIONS**

Repair of ventral hernia or other fascial defects that require the addition of an extraperitoneal reinforcing or bridging material to obtain the desired surgical result.

**PERFORMANCE**

The 4D Ventral® implant is designed for extraperitoneal implantation only.  
Meshes are biocompatible and partially resorbable. They have the advantage of being cut-able, of having a very high resistance to suture, of having a great flexibility.

**CONTRAINDICATIONS**

Do not implant intraperitoneally.

Do not use in the following cases:

- |                                    |                         |             |
|------------------------------------|-------------------------|-------------|
| - Allergy to any of the components | - Infected site         | - Pregnancy |
| - Growing children                 | - Anticoagulant therapy |             |

**UNDESIRABLE SIDE EFFECTS**

Like any implantable medical device, this implant is susceptible to generate possible undesirable side effects that may lead to reoperation:

- |                     |                      |                     |                           |
|---------------------|----------------------|---------------------|---------------------------|
| - Discomfort/Pain   | - Recurrence         | - Hematoma          | - Foreign body reaction   |
| - Infection         | - Adhesion formation | - Mesh migration    | - Irritation nearby organ |
| - Inflammation      | - Erosion            | - Mesh deformation  | - Lymphocele              |
| - Fistula formation | - Seroma             | - Allergic reaction |                           |

**USAGE PRECAUTIONS**

4D Ventral® implants are delivered sterile (ethylene oxide sterilisation).

Before any use, inspect the integrity of packaging and device (of which pelables sachets). Do not use in the event of deterioration of the device and/or the packaging.

Do not use if the device is out of date.

COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned. Mesh fixation points should be at least 1cm from the edge of the mesh with 1cm spacing between fixation points.

This device must be implanted only by a qualified surgeon trained on the use of the product (knowledge of anatomy and visceral surgery).

**IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE**

As specified on the product's labeling, the 4D Ventral® implants are for single use only. They can not be re-used and/or re-sterilised (potential risks would be and are not limited to: loss of the product's sterility, risk of infection, loss of the product's efficiency, recurrence)

**EXPLANTATION AND ELIMINATION OF DEVICES**

Devices should be retrieved and handled according to the recommendations of ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling".

Any explanted device must be sent back for analysis, following the current protocol. This protocol is available on request from COUSIN BIOTECH. It is important to note that any implant that must not be cleaned or disinfected prior to despatch must be contained in sealed packaging.

The removed medical device must be disposed of in accordance with the country's standards for the disposal of infectious waste.

There are no specific recommendations for the disposal of a non-implanted device.

**STORAGE**

To be stored in a dry place away from sunlight and at room temperature in its original packaging.

### **INFORMATION REQUESTS AND COMPLAINTS**

Following its quality policy, COUSIN BIOTECH is committed to making every effort to produce and supply a high quality medical device. However, if a health professional (client, user, prescriber, etc.) has a complaint or cause for dissatisfaction with a product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In case of failure of an implant or if it has contributed to causing a serious adverse reaction in the patient, the health centre must follow the legal procedures in that country, and inform COUSIN BIOTECH immediately.

For any correspondence, please specify the reference, batch number, contact person's details and a comprehensive description of the incident or complaint.

Brochures, documentation and surgical technique are available on request from COUSIN BIOTECH and its distributors.

### **IMPORTANT**

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

### **IMPLANT CARD**

**Be sure to scan your implant card immediately upon receipt so that you can keep track of it in case of loss.**



### **SURGICAL PROCEDURE**

#### **Extra-Peritoneal Hernia repair**

Example: Retrorectus mesh repair of a ventral hernia (Rives-Stoppa technique).

This technique utilizes the hernia sac to separate the mesh from the intra-abdominal contents.

- Superior to the umbilicus, dissection is performed above the posterior rectus fascia and under the rectus muscle.
- Below the umbilicus, dissection occurs in the preperitoneal space due to the lack of a posterior rectus sheath.
- The 4D Ventral® implant is placed in the newly formed pre-peritoneal space. The 4DVentral mesh has an oriented elasticity:

Direction of the  
elasticity

- Fixation of the mesh if needed, to the muscle layer above
- The peritoneum must be meticulously checked and no defect in the peritoneum should be left open to prevent any contact between the mesh and intra-abdominal contents.
- Redon's drains under aspiration are usually placed. The anterior sheath is closed. If tension is present, relaxing incisions may be used
- Skin closure.
- Using the Redon's drains, fill the pre-peritoneal space with an antiseptic solution and clamp the drains for 1 hour.
- A belt can be used to wrap the patient's abdomen for a few days to reinforce the abdominal wall.

**IMPLANT DE RENFORCEMENT PARIETAL SEMI RESORBABLE  
EXTRAPERITONEAL  
PRODUIT STERILE A USAGE UNIQUE**

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| Item | Patient information leaflet                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1    | <p>a. <b>The name of the device :</b> 4DVENTRAL®</p> <p>b. <b>The model of the device :</b> 4DVENT05RO, 4DVENT07RO, 4DVENT09RO, 4DVENT12RO, 4DVENT0715, 4DVENT1015, 4DVENT1515, 4DVENT1520, 4DVENT1530, 4DVENT2020, 4DVENT2025, 4DVENT2535, 4DVENT3030, 4DVENT3050, 4DVENT3040 and 4DVENT4545</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| 2    | <p>a. <b>The intended purpose of the device :</b><br/>Repair and extraperitoneal parietal reinforcement of ventral hernias.</p> <p>b. <b>The kind of patient on whom the device is intended to be used :</b><br/>Adult male or female suffering from a symptomatic ventral hernia necessitating a surgical mesh repair with an extraperitoneal implantation method.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| 3    | <p><b>Special operating instructions for the use of the device :</b><br/>This device must be implanted only by a qualified surgeon trained on the use of the product (knowledge of anatomy and visceral surgery). The 4DVENTRAL® implant is designed for extraperitoneal implantation only.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| 4    | <p>a. <b>The intended performance of the device :</b><br/>4DVENTRAL® are designed for parietal reinforcement and hernia repair by open or laparoscopic approach. Prostheses are used as hernia defect coverings/patches in conjunction with surgical repair to help stop visceral contents from passing through the hernia orifice. 4DVENTRAL® intend to cause some degree of reactional fibrosis. Semi-resorbable meshes are proposed as a potential solution which combines the advantages of an absorbable mesh to reduce inflammation, and of a non-absorbable mesh to provide long term strength.</p> <p>b. <b>Undesirable side effects :</b><br/>Like any implantable medical device, this implant is susceptible to generate possible undesirable side effects that may lead to reoperation: Discomfort/Pain, Recurrence, Adhesions, Erosions, Irritation of nearby organ, Infection, Inflammation, Fistula formation, Seroma, Mesh deformation, Hematoma, Mesh migration, Allergic reaction, Foreign body reaction. Seek medical advice if you experience any of these symptoms.</p> |
| 5    | <p><b>When 4DVENTRAL® Should Not be Used :</b><br/>The following are contraindications for 4DVENTRAL®: Allergy to any of the components, Infected site, Pregnancy, Growing children, Anticoagulant therapy</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |

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| Item | Patient information leaflets                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
|------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 6    | <b>MRI Information :</b><br>4DVENTRAL® are non-conducting and non-magnetic devices. They present no risk in Magnetic Resonance Imaging (MRI) environment.                                                                                                                                                                                                                                                                                                                                                                                             |
| 7    | <b>Precautions :</b><br>After surgery with 4DVENTRAL®, please thoroughly follow instructions given by your surgeon. The 4DVENTRAL® mesh is designed for permanent implantation. If you feel persistent pain or any undesirable effect, you should contact a healthcare professional related to your operation with the 4DVENTRAL device.                                                                                                                                                                                                              |
| 8    | <b>Materials included in the device :</b><br>4DVENTRAL® implants are extraperitoneal partially resorbable parietal reinforcement implants. This surgical mesh is composed of Poly-L-Lactic Acid (PLLA) and Polypropylene (PP)                                                                                                                                                                                                                                                                                                                         |
| 9    | <b>Reporting a serious incident :</b><br>If you encounter any serious incident with your 4DVENTRAL® implant you should report it immediately to your doctor. You should also report the incident to Cousin Biotech via the link <a href="mailto:serviceclients@cousin-surgery.com">serviceclients@cousin-surgery.com</a> or <a href="mailto:materiovigilance@cousin-biotech.com">materiovigilance@cousin-biotech.com</a> and to the Therapeutic Goods Administration using the following website: <a href="http://www.tga.gov">http://www.tga.gov</a> |
| 10   | <b>Legal Manufacturer :</b><br>Cousin Biotech, Allée des roses, 59117 Wervicq-Sud, France<br><b>Manufacturing Facility :</b><br>Cousin Biotech, Allée des roses, 59117 Wervicq-Sud, France                                                                                                                                                                                                                                                                                                                                                            |

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Cousin Biotech SAS au capital de 320 000 € - NAF 1396Z- N° TVA FR 34 398 460 261 - Siren 398 460 261 RCS Lille