Sterile Spinal Dynamic Posterior Stabilization Device

en	Instructions for use	Page	2	MDR CE marking date: 02/10/2025
fr	Notices d'instructions	Page	6	Date de marquage CE MDR : 02/10/2025
de	Gebrauchsanweisung	Seite	13	<i>Datum der MDR-CE-Kennzeichnung: 02/10/2025</i>
it	Istruzioni per l'uso	Pagina	19	<i>Data di marcatura CE MDR: 02/10/2025</i>
es	Instrucciones de uso	Pagina	25	<i>Fecha del marcado CE MDR: 02/10/2025</i>
pt	Nota de instruções	Pagina	31	<i>MDR Data da marcação CE: 02/10/2025</i>
nl	Gebruiksaanwijsties	Page	37	<i>Datum MDR CE-markering: 02/10/2025</i>
pl	Instrukcja obsługi	Strona	43	<i>Data oznakowania CE zgodnie z MDR: 02/10/2025</i>
ro	Instrucțiuni de utilizare	Pagină	49	<i>Data marcării CE MDR: 02/10/2025</i>
ελ	Οδηγίες χρήσης	Σελίδα	54	<i>Ημερομηνία σήμανσης MDR CE: 02/10/2025</i>
cs	Návod k použití	Strana	60	<i>Datum udělení označení MDR CE: 02/10/2025</i>
en	Patient information leaflet (Australian requirements)	Page	65	

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Made in France

NOT315_251010
Version du 10/10/2025

Caution : Federal law (USA) restricts this device to sale,
distribution and use by or on the order of a physician.

This release is the last update of the instructions of use and replace the previous edition
Product under regulation 2017/745 (EU).

Sterile Spinal Dynamic Posterior Stabilization Device

SINGLE USE STERILE PRODUCT

MANUFACTURER IDENTIFICATION

COUSIN BIOTECH
ALLEE DES ROSES
59117 WERVICQ SUD
FRANCE

IDENTIFICATION AND DESCRIPTION OF THE DEVICE

The sterile spinal dynamic posterior stabilization device BDyn is intended to restore the stabilization of the spine segment from T10 to S1 by preserving the anatomical lordosis and the deadening of the intervertebral joint. It is composed of the BDyn spinal shock absorber fixed on the vertebrae thanks to the dedicated screws.

The intended purpose of BDyn is to induce the load sharing, relieve the facets loads, prevent the adjacent syndrome, allow the combination with the fusion and maintain the mobility and ensure the absorption of mechanical loads of the spine functional unit from T10 to S1.

INDICATIONS

The sterile spinal dynamic posterior stabilization device BDyn is intended for posterior stabilization from thoracic vertebrae T10 to sacrum S1 with or without bone graft for the following indications:

- Degenerative intervertebral disc disease and/or articular facets confirmed by further examinations
- Spinal canal stenosis
- Degenerative spondylolisthesis grade 1
- Segmental instability

CONTRAINDICATIONS

- Active infectious process or significant risk of infection (immunocompromise)
- Signs of local inflammation
- Fever or leukocytosis
- Morbid obesity
- Pregnancy
- Mental illness
- Grossly distorted anatomy caused by congenital abnormalities
- Any other medical or surgical condition which would preclude the potential benefit of spinal implant surgery, such as the presence of congenital abnormalities, elevation of the sedimentation rate unexplained by other diseases, elevation of the white blood count
- Suspected or documented metal allergy or intolerance
- Any case where the implant components selected for use would be too large or too small to achieve a successful result
- Any patient having inadequate tissue coverage over the operative site or inadequate bone stock or quality
- Any patient in which implant utilization would interfere with anatomical structures or expected physiological performance
- Any patient unwilling to follow postoperative instructions
- Any case not describe in the indications
- Traumas (i.e. fracture or dislocation)
- Abnormal curvatures (i.e. scoliosis and/or hyper lordosis)
- Tumors.
- Spondylolisthesis grade 2 and more
- Pseudarthrosis and/or failed previous fusion
- Severe bone resorption, osteomalacia, severe osteoporosis

UNDESIRABLE SIDE EFFECTS

All of the possible adverse events associated with spinal surgery and without instrumentation are possible:

- Infection
- Pseudomeningocele, fistula, breach dura, persistent CSF leakage, meningitis
- Loss of neurological function, sensorial and/or motor, including complete or incomplete paralysis, dysesthesias, hyperesthesia, anesthesia, paresthesia, appearance of radiculopathy, and/or the development or continuation of pain, numbness, neuroma, spasms, sensory loss, tingling sensation, and/or visual deficits
- Cauda equina syndrome, neuropathy, transient or permanent neurological deficits, paraplegia, paraparesis, reflex deficits, irritation, arachnoiditis, and/or muscle loss
- Urinary retention or loss of bladder control or other types of urological system compromise
- Scar formation possibly causing by a neurological compromise or compression around nerves and/or pain
- Fracture, microfracture, resorption, damage or penetration of any spinal bone (including the sacrum, pedicles, and/or vertebral body) and/or bone graft or bone graft harvest site at, above and/or below the level of surgery
- Herniated nucleus pulposus, disc disruption or degeneration at, above, or below the level of surgery, canal adjacent stenosis
- Non-union or pseudarthrosis, delayed union. Mal union
- Cessation of any potential growth of the operated portion of the spine
- Loss of or increase in spinal mobility or function
- Inability to perform the activities of daily living
- Bone loss or decrease in bone density
- Graft donor site complications including pain, fracture, or wound healing problems
- Ileus, gastritis, bowel obstruction or loss of bowel control or other types of gastrointestinal system compromise
- Hemorrhage, hematoma, occlusion, seroma, oedema, hypertension, embolism, stroke, excessive bleeding, phlebitis, wound necrosis, wound dehiscence, damage to blood vessels, or other types of cardiovascular system compromise
- Reproduction system compromise, sterility, sexual dysfunction
- Development of respiratory problems, e.g. pulmonary embolism, atelectasis, bronchitis, pneumonia, etc
- Change in mental status
- Death

All of the possible adverse events associated with spinal surgery with instrumentation are possible. A listing of potential adverse events linked to the medical device includes, not limited to:

- Early or late loosening of any or all of the components
- Disassembly, bending and/or breakage of any or all of the components (screw breakage)
- Foreign body (allergic) reaction to implants, debris, corrosion products (from crevice, fretting, and/or general corrosion), including metallosis, tumor formation and/or autoimmune disease
- Pressure on the skin from components parts with inadequate tissue coverage over the implant possibly causing skin penetration, irritation, fibrosis, neurosis, and/or pain
- Tissue or nerve damage caused by improper positioning and placement of implants or instruments
- Post-operative change in spinal curvature, loss of correction, height, and/or reduction

TARGET POPULATION

The Sterile Spinal Dynamic Posterior Stabilization device BDyn is only intended to skeletally mature adults.

IMPLANTED MATERIALS

The BDyn device components are manufactured from medical grade titanium TA6V (ISO 5832-3; ASTM F 136), silicone, long-term implantable polyurethane. Titanium alloy Ti6Al4V ELI (ISO 5832-3; ASTM F 136), Polydimethyl siloxane (PDMS) MED 4770 unrestricted, Polycarbonate urethane (PCU) Bionate® II 80A. Origin neither human nor animal – Non resorbable.

Reference Material	Titanium alloy Ti6Al4V ELI	Polydimethyl siloxane	Polycarbonate urethane (PCU) Bionate® II 80A
RCBDYSD50U	6.74 g	0.22 g	0.19 g
RCBDYMD50U	11.94 g	0.22 g	0.19 g
RCBDYSD55U	7.33 g	0.22 g	0.19 g
RCBDYMD55U	13.62 g	0.22 g	0.19 g
RCBDYLD55U	22.38 g	0.22 g	0,19 g

Quantity mentioned in the table above are with a confidence interval of 2%.

ACCESSORIES

BDyn instruments: stainless steel, Propylux (Polypropylene) or POM (PolyOxyMethylen). Non-human and non-animal origin – Non-resorbable.

ACTION AND PERFORMANCES

The BDyn spinal shock absorber of the sterile spinal dynamic posterior stabilization device BDyn consists of a metallic cylindrical hollow part containing elastomer components made of silicone elastomers (polydimethyl siloxane (PDMS)) and long-term implantable polycarbonate urethane (PCU) which are bending out under the effect of a metallic piston rod connected with the vertebra of the treated segment by the pedicular screws tested and approved by the company COUSIN BIOTECH. The combination of rigid and flexible components helps to maintain the mobility and ensure the absorption of mechanical loads in flexion/extension, compression, axial rotation, and lateral bending. This ensures to reduce intradiscal pressure and relieve facet. Several configurations are available and allow the combination of dynamic and fusion solution.

SCREWS SYSTEM APPROVED

REFERENCES	Description	Screws system approved
RCBDYSD50U	Dynamic shock absorber on small mobile rod (Ø5.0mm) allowing the preservation of mobility on one level	G2S system from Neuro France Implants
RCBDYSD55U	Dynamic shock absorber on small mobile rod (Ø5.5mm) allowing the preservation of mobility on one level	G2S system from Neuro France Implants 3K from Fradis Diplomat from Signus
RCBDYMD50U	Dynamic shock absorber on medium-sized rod (Ø5.0mm) allowing preservation of mobility on one level and fusion on one or two lower levels	G2S system from Neuro France Implants
RCBDYMD55U	Dynamic shock absorber on medium-sized rod (Ø5.5mm) allowing preservation of mobility on one level and fusion on one or two lower levels	G2S system from Neuro France Implants 3K from Fradis Diplomat from Signus
RCBDYLD55U	Dynamic shock absorber on long-sized rod (Ø5.5mm) allowing preservation of mobility one level and fusion up to 4 lower levels	G2S system from Neuro France Implants 3K from Fradis Diplomat from Signus

LINK TO SUMMARY OF SAFETY AND PERFORMANCE CHARACTERISTICS :

https://www.cousin-biotech.com/d-68ece07c929b8_tf06sscpbdyn.pdf

MRI SAFETY

The BDyn implant has been tested for safety and compatibility in the case of postoperative investigation using Magnetic Resonance Imaging (MRI). Non-clinical testing has demonstrated that the BDyn orthopaedic implants are MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions. Failure to follow these conditions may result in injury.

Maximum Static Magnetic Field Strength	1.5T / 3.0T
Type of Nuclei	Hydrogen
Static Magnetic Field (B0) Orientation	Horizontal
Magnet Type	Cylindrical-bore
Maximum Spatial Field Gradient	87.1 T/m (8710 gauss/cm)

Maximum B0*ldB0/drl product	101.5 T2/m
RF Excitation	1.5T: Circularly Polarized (CP) 3T: Multichannel-2 (MC-2)
RF Transmit / Receive Coil Type	Integrated Whole Body Transmit and Receive Coil
RF Power	Normal Operating Mode
Maximum Whole-Body SAR	Whole body SAR ≤ 2 W/kg Note: Under the scan conditions defined above, the BDyn implants are expected to produce a maximum temperature rise of less than 6.1 °C at 1.5T and less than 3.4 °C at 3T after 15 minutes of continuous scanning.
Patient Characteristics	1.5T: Patients with uncompromised thermoregulation and under controlled conditions (a medical doctor or a dedicated trained person can respond instantly to heat induced physiological stress) 3T: Patients with uncompromised thermoregulation and under uncontrolled conditions, or patients with compromised thermoregulation (all persons with impaired systemic or reduced local thermoregulation) under controlled conditions.
MR Image Artifact	MR image quality may be compromised if the imaging area of interest is in the same area of the implant. Some manipulation of scan parameters may be required to compensate for the artifact. In non-clinical testing, the image artifact caused by the device extends approximately 36 mm from the implant when imaged with a spin echo pulse sequence and 22 mm with a gradient echo, both at 1.5T.

Consequently, the BDyn is considered MR conditional.

PRECAUTIONS FOR USE

The BDyn devices are delivered sterile, sterilized by ethylene oxide.

Before any use, inspect the integrity of the packaging and device (including peelable pouches).

Do not use in the event of deterioration of the labels and/or the device and/or the packaging.

Do not use if the device is out of date.

The installation of BDyn devices can only be made with pedicular screws tested and approved by the company COUSIN BIOTECH.

You should never use implants made of stainless steel and titanium alloy in the same construct.

Many factors must be considered during spinal surgery, which means that the results obtained are highly variable. This device, like all spinal implants, can only withstand the loads of the body if it is in addition to a bone structure. In the case of alteration of the bone support, loosening, disassembly and/or rupture of the device may occur. Preoperative and operating procedures, including precise knowledge of suitable surgical techniques, and proper selection of the good reference of the device adapted to the patient and its narrow setting up are important considerations in the successful use of the device by the surgeon.

The BDyn device must be implanted only by a qualified surgeon, having knowledge in the use of the product and who has the knowledge of the anatomy, spinal surgery, pedicle screws fixation technique and specific BDyn device surgical technique.

IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE

In accordance with the labelling of this product, BDyn implants are for single use; they must not in any case be reused or re-sterilised (potential risks include, but are not limited to, loss of sterility of the products, risk of infection, loss of product efficacy, recurrence, etc.)

INSTRUCTIONS TO BE GIVEN TO THE PATIENTS BY THE SURGEON

The surgeon should inform the patient of the potential physical and psychological restrictions and consequences of implanting the device. The patient must be informed of the surgical risks and possible side effects. The surgeon should invite the patient to return for further consultation if the patient presents symptoms that appear abnormal. An implant card for the patient is provided by the healthcare professional at the clinic. This implant card provides information, for the patient, to identify the device and traceability elements as well as the name, address and website of the manufacturer. The surgeon shall invite the patient to scan the implant card immediately upon receipt to keep track of it in case of loss. The surgeon shall also inform the patient that the summary of the safety and performance characteristics of the device and the instructions for use can be found on the COUSIN BIOTECH website.

INSTRUMENTS SYSTEM USED WITH BDYN DEVICES

INSTRUMENTS REFERENCE	USE	LEGAL MANUFACTURER
BDYNCONT00	Container	COUSIN BIOTECH
RCBANDH10U	Torque handle (10 N)	COUSIN BIOTECH
RCBANSOUNU	Reversible probe	COUSIN BIOTECH
RCBANRH50U	Rod holder	COUSIN BIOTECH
RCBANLPERU	Lumbar perforator	COUSIN BIOTECH
RCBANTPERU	Thoracic perforator	COUSIN BIOTECH
RCBANLRBEU	Bender	COUSIN BIOTECH
RCBANDISTU	Distraction Clamp	COUSIN BIOTECH
RCBANCOMPU	Contraction Clamp	COUSIN BIOTECH
RCBANTD50U	Trial implant	COUSIN BIOTECH
RCBANBDHOU	BDyn holder	COUSIN BIOTECH
RCBANRODP	Rod pusher	COUSIN BIOTECH
014112C	Square Tip	NEURO FRANCE Implants
013615	Cannulated Tap 5.0	NEURO FRANCE Implants
011300	Fork	NEURO FRANCE Implants
911537	Kirschner wire	NEURO FRANCE Implants
013001N	Straight ratchet cannulated handle	NEURO FRANCE Implants
010201C	Cannulated screwdriver	NEURO FRANCE Implants
010494	Nut holder	NEURO FRANCE Implants

010505C	Nut driver	NEURO FRANCE Implants
011308	Persuader	NEURO FRANCE Implants
010001N	Ratchet cannulated T-handle	NEURO FRANCE Implants
013616	Cannulated Tap 6.0	NEURO FRANCE Implants
010610C	Counter-torque	NEURO FRANCE Implants
016103	One-way probe	NEURO FRANCE Implants

WARNING ABOUT INSTRUMENTS

The BDyn instruments are intended for temporary use and are re-usable.

Unlike the BDyn implants, instruments are furnished unsterile. For more details on instruments, please refer to the instructions for use specific to information on the care, decontamination, cleaning, disinfection, maintenance and sterilisation of class I reusable surgical instruments manufactured and/or distributed by COUSIN BIOTECH.

DECONTAMINATION, CLEANING AND STERILIZATION OF INSTRUMENTS

Please refer to Instructions for Use NOT348 related to information on the care, decontamination, cleaning, disinfection, maintenance and sterilisation of class I reusable surgical instruments manufactured and/or distributed by COUSIN BIOTECH.

STORAGE PRECAUTIONS

Store in a dry place away from sunlight and at room temperature in its original packaging.

SURGICAL TECHNIQUE

The placement of the spinal dynamic posterior stabilization device BDyn begins by the preparation of the pedicles of the vertebrae and by the pedicle aim. Once the vertebra has been prepared, the fixation of the pedicle screws tested and approved by COUSIN BIOTECH can be placed laterally to the facet joint.

Before placing the BDyn rods, it is necessary to optimize the location of polyaxial screw in order to facilitate the positioning of BDyn device. The trial prosthesis RCBANTD50U enables to better estimate the bulk of the BDyn dynamic rod and optimise the implantation of screw. It must be perfectly flattened against the bottom of polyaxial screws heads without affecting articular facets.

We can then proceed to the installation of BDyn rods. The mobile piston rod of the BDyn spinal shock absorber must be fixed in the head of the cranial screw concentrically to the cylindrical body. A 2 mm space between the top of the cylindrical body of the BDyn spinal shock absorber and the lower part of the screw head need to be ensured in order to preserve the sufficient shock absorption amplitude necessary for the good working of the BDyn device. The laser marking "up" has to be positioned face to the surgeon, respecting the «up» direction and the arrow must be aligned with the two screw head opening. The surgery ends by the final tightening by using the dedicated torque limiting of the plugs on the head of the polyaxial pedicle screws by using systematically the suitable antirotator.

It is essential to position the implant according to the indications. The lower rigid rod under the shock absorbing cylinder of the implants must not be counter-bent in the opposite direction. These operations could effectuate the concentration of internal stresses which could become the place of an eventual failure of the implant.

Depending to the surgeon's decision, it is possible to cut the medium-sized or the long-sized fixed rod for a single level of fusion or a three level of fusion respectively.

EXPLANTATION AND DISPOSAL OF DEVICES

Risks associated with the explantation: the surgeon who decides to remove the device should consider factors such as the risk of another surgery for the patient, and the difficulty of the explantation procedure. The removal of the implant is up to the surgeon's decision and must be subject to adequate postoperative follow-up. Devices should be retrieved and handled according to the recommendations of ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling". Any explanted device must be sent back for analysis, following the current protocol. This protocol is available on request from COUSIN BIOTECH. It is important to note that any implant that must not be cleaned or disinfected prior to despatch must be contained in sealed packaging. The removed medical device must be disposed of in accordance with the country's standards for the disposal of infectious waste.

There are no specific recommendations for the disposal of a non-implanted device.

Explantation explanations are available in the surgical technique of BDyn.

IMPORTANT

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

INFORMATION REQUEST AND COMPLAINTS

Following its quality policy, COUSIN BIOTECH is committed to making every effort to produce and supply a high-quality medical device. However, if a health professional (client, user, prescriber, etc.) has a complaint or cause for dissatisfaction with a product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In case of failure of an implant or if it has contributed to causing a serious adverse reaction in the patient, the health centre must follow the legal procedures in that country and inform COUSIN BIOTECH immediately. For any correspondence, please specify the reference, batch number, contact person's details and a comprehensive description of the incident or complaint.

Brochures, documentation and surgical technique are available on request from COUSIN BIOTECH and its distributors. If further information is needed or required, please contact your COUSIN BIOTECH representative or distributor.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

IMPLANT CARD INFORMATION LEAFLET

How to complete the implant card ?

The implant card, for the patient, should be completed as follows:

Front side of the implant card

The user has to complete the implant card before transmitting the implant card to the patient. Please complete it as follows:

	Carte Implant - Implant Card A conserver / Keep it	Carte Implant - Implant Card
Patient name	 _____	 COUSIN SURGERY <i>We care for Surgery</i>
Date of surgical procedure	 . . / . . /	
Name of practitioner	 _____	

Back side of the implant card

No fields need to be completed by the user. The following information on the device is printed directly by the manufacturer on the back of the implant card. A QR code allowing the patient to consult the instruction for use, as well as the website address, are available on the implant card.

Information to identify the device	Traceability information
• Indication that the product is a medical device • Description of the device family • Trade name of the device • Reference of the device • Device is MR safe • Denomination name of the device in sales countries	• UDI number • Batch number • Name and address of the legal manufacturer

Explanation of the used symbols in the Implant Card

All symbols are commonly used and explained at the end of the instruction for use.

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CONTRAINDICATIONS

- Active infectious process or significant risk of infection (immunocompromise)
- Signs of local inflammation
- Fever or leukocytosis
- Morbid obesity
- Pregnancy
- Mental illness
- Grossly distorted anatomy caused by congenital abnormalities
- Any other medical or surgical condition which would preclude the potential benefit of spinal implant surgery, such as the presence of congenital abnormalities, elevation of the sedimentation rate unexplained by other diseases, elevation of the white blood count
- Suspected or documented metal allergy or intolerance
- Any case where the implant components selected for use would be too large or too small to achieve a successful result
- Any patient having inadequate tissue coverage over the operative site or inadequate bone stock or quality
- Any patient in which implant utilization would interfere with anatomical structures or expected physiological performance
- Any patient unwilling to follow postoperative instructions
- Any case not describe in the indications
- Traumas (i.e. fracture or dislocation)
- Abnormal curvatures (i.e. scoliosis and/or hyper lordosis)
- Tumors
- Spondylolisthesis grade 2 and more
- Pseudarthrosis and/or failed previous fusion
- Severe bone resorption, osteomalacia, severe osteoporosis

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- Early or late loosening of any or all of the components
- Disassembly, bending and/or breakage of any or all of the components (screw breakage)
- Foreign body (allergic) reaction to implants, debris, corrosion products (from crevice, fretting, and/or general corrosion), including metallosis, tumor formation and/or autoimmune disease

- Pressure on the skin from components parts with inadequate tissue coverage over the implant possibly causing skin penetration, irritation, fibrosis, neurosis, and/or pain
- Tissue or nerve damage caused by improper positioning and placement of implants or instruments
- Post-operative change in spinal curvature, loss of correction, height, and/or reduction

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ACCESSORIES

BDyn instruments: stainless steel, Propylux (Polypropylene) or POM (PolyOxyMethylen).
Non-human and non-animal origin – Non-resorbable.

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https://www.cousin-biotech.com/d-66d5a5fc142af_tf06sscpissue-3draft.pdf

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Maximum Static Magnetic Field Strength	1.5T / 3.0T
Type of Nuclei	Hydrogen
Static Magnetic Field (B0) Orientation	Horizontal
Magnet Type	Cylindrical-bore
Maximum Spatial Field Gradient	87.1 T/m (8710 gauss/cm)
Maximum B0*IdB0/drl product	101.5 T2/m
RF Excitation	1.5T: Circularly Polarized (CP)

	3T: Multichannel-2 (MC-2)
RF Transmit / Receive Coil Type	Integrated Whole Body Transmit and Receive Coil
RF Power	Normal Operating Mode
Maximum Whole-Body SAR	Whole body SAR ≤ 2 W/kg Note: Under the scan conditions defined above, the BDyn implants are expected to produce a maximum temperature rise of less than 6.1 °C at 1.5T and less than 3.4 °C at 3T after 15 minutes of continuous scanning.
Patient Characteristics	1.5T: Patients with uncompromised thermoregulation and under controlled conditions (a medical doctor or a dedicated trained person can respond instantly to heat induced physiological stress) 3T: Patients with uncompromised thermoregulation and under uncontrolled conditions, or patients with compromised thermoregulation (all persons with impaired systemic or reduced local thermoregulation) under controlled conditions.
MR Image Artifact	MR image quality may be compromised if the imaging area of interest is in the same area of the implant. Some manipulation of scan parameters may be required to compensate for the artifact. In non-clinical testing, the image artifact caused by the device extends approximately 36 mm from the implant when imaged with a spin echo pulse sequence and 22 mm with a gradient echo, both at 1.5T.

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INSTRUMENTS REFERENCE	USE	LEGAL MANUFACTURER
BDYNCONT00	Container	COUSIN BIOTECH
RCBANDH10U	Torque handle (10 N)	COUSIN BIOTECH
RCBANSOUNU	Reversible probe	COUSIN BIOTECH
RCBANRH50U	Rod holder	COUSIN BIOTECH
RCBANLPERU	Lumbar perforator	COUSIN BIOTECH
RCBANTPERU	Thoracic perforator	COUSIN BIOTECH
RCBANLRBEU	Bender	COUSIN BIOTECH
RCBANDISTU	Distraction Clamp	COUSIN BIOTECH
RCBANCOMPU	Contraction Clamp	COUSIN BIOTECH
RCBANTD50U	Trial implant	COUSIN BIOTECH
RCBANBDHOU	BDyn holder	COUSIN BIOTECH
RCBANRODPU	Rod pusher	COUSIN BIOTECH

014112C	Square Tip	NEURO FRANCE Implants
013615	Cannulated Tap 5.0	NEURO FRANCE Implants
011300	Fork	NEURO FRANCE Implants
911537	Kirschner wire	NEURO FRANCE Implants
013001N	Straight ratchet cannulated handle	NEURO FRANCE Implants
010201C	Cannulated screwdriver	NEURO FRANCE Implants
010494	Nut holder	NEURO FRANCE Implants
010505C	Nut driver	NEURO FRANCE Implants
011308	Persuader	NEURO FRANCE Implants
010001N	Ratchet cannulated T-handle	NEURO FRANCE Implants
013616	Cannulated Tap 6.0	NEURO FRANCE Implants
010610C	Counter-torque	NEURO FRANCE Implants
016103	One-way probe	NEURO FRANCE Implants

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Patient information leaflet BDyn

Item	Patient information leaflet
1	a. The name of the device : BDyn b. The model of the device : RCBDYMD55U, RCBDYSD55U
2	a. The intended purpose of the device : To restore the stabilization of the intervertebral thoracic, lumbar and sacral segments by preserving the anatomical lordosis and the deadening of the intervertebral joint. b. The kind of patient on whom the device is intended to be used : Adult human and female or male more than 40 kg.
3	Special operating instructions for the use of the device : The BDyn device must be implanted only by a qualified surgeon, having knowledge in the use of the product and who has the knowledge of the anatomy, spinal surgery, pedicle screws fixation technique and specific BDyn device surgical technique.
4	a. The intended performance of the device : The BDyn spinal shock absorber of the sterile spinal dynamic posterior stabilization device BDyn consists of a metallic cylindrical hollow part containing elastomer components made of silicone elastomers (polydimethyl siloxane (PDMS)) and long-term implantable polycarbonate urethane (PCU) which are bending out under the effect of a metallic piston rod connected with the vertebra of the treated segment by the pedicular screws tested and/or approved by the company COUSIN BIOTECH. The elastomer components assure the absorption of the mechanical loads applied on the intervertebral joint in compression, pulling, flexion-traction and lateral flexion. Some configurations enable fixation across several spinal segments. b. Undesirable side effects : Like any implantable medical device, this implant is susceptible to generate possible undesirable side effects that may lead to reoperation: Infection, Pseudomeningocele, fistula, breach dura, persistent CSF leakage, meningitis, loss of neurological function, sensorial and/or motor, including complete or incomplete paralysis, dysesthesias, hyperesthesia, anesthesia, paresthesia, appearance of radiculopathy, and/or the development or continuation of pain, numbness, neuroma, spasms, sensory loss, tingling sensation, and/or visual deficits, cauda equina syndrome, neuropathy, transient or permanent neurological deficits, paraplegia, paraparesis, reflex deficits, irritation, arachnoiditis, and/or muscle loss, urinary retention or loss of bladder control or other types of urological system compromise, scar formation possibly causing by a neurological compromise or compression around nerves and/or pain, fracture, microfracture, resorption, damage or penetration of any spinal bone (including the sacrum, pedicles, and/or vertebral body) and/or bone graft or bone graft harvest site at, above and/or below the level of surgery, herniated nucleus pulposus, disc disruption or degeneration at, above, or below the level of surgery, canal adjacent stenosis, Non-union or pseudarthrosis, delayed union, mal union, Cessation of any potential growth of the operated portion of the spine, Loss of or increase in spinal mobility or function, Inability to perform the activities of daily living, Bone loss or decrease in bone density, Graft donor site complications including pain, fracture, or wound healing problems, Ileus, gastritis, bowel obstruction or loss of bowel control or other types of gastrointestinal system compromise, hemorrhage, hematoma, occlusion, seroma, edema, hypertension, embolism, stroke, excessive bleeding, phlebitis, wound necrosis, wound dehiscence, damage to blood vessels, or other types of cardiovascular system compromise, Reproduction system compromise, sterility, sexual dysfunction, Development of respiratory problems, e.g. pulmonary embolism, atelectasis, bronchitis, pneumonia, etc, Change in mental status, Death
5	When BDyn Should Not be Used : Active infectious process or significant risk of infection (immunocompromise), signs of local inflammation, Fever or leukocytosis, Morbid obesity, Pregnancy, Mental illness, Grossly distorted anatomy caused by congenital abnormalities, Any other medical or surgical condition which would preclude the potential benefit of spinal implant surgery, such as the presence of congenital abnormalities, elevation of the sedimentation rate unexplained by other diseases, elevation of the white blood count Suspected or documented metal allergy or intolerance, Any case where the implant components selected for use would be too large or too small to achieve a successful result, Any patient having inadequate tissue coverage over the operative site or inadequate bone stock or quality, Any patient in which implant utilization would interfere with anatomical structures or expected physiological performance, Any patient unwilling to follow postoperative instructions, Any case not describe in the indications, Traumas (i.e. fracture or dislocation), Abnormal curvatures (i.e. scoliosis and/or hyper lordosis), Tumors, Spondylolisthesis grade 2 and more, Pseudarthrosis and/or failed previous fusion, Severe bone resorption, osteomalacia, severe osteoporosis

This leaflet was revised in March 2023

Page 56/60

Headquarters & Factory: **COUSIN BIOTECH** Allée des Roses - 59117 Wervicq-Sud – France **Tel** +33 (0) 3 20 14 41 20 –

Fax +33 (0) 3 20 14 40 13

www.cousin-surgery.com - **e-mail** f.pelletier@cousin-surgery.com

Item	Patient information leaflet
6	MRI Information : The implants are composed of non-ferromagnetic materials and present a geometry non susceptible to generate induced currents. Moreover, as they are fixed to bone /or/ tissues, they are unlikely to be mobilised. A priori they can be considered compatible with an MRI scan. Their safety, in particular in terms of heating and migration of implant has been evaluated through bibliographic data by comparison with data available on devices with similar composition, shape and use. This evaluation concluded to a safety use for MRI scan between 1.5 and 3 Tesla. As a precautionary measure, it is recommended to avoid MRI scans within the 48h of the implant placement, and to inform the person in charge of the scan of the recent implant placement, if such examination is essential. It has to be noted that the devices which present a high contrast with the biological environment can generate « artifacts » that has to be taken into account for the perfect execution and interpretation of imaging exams. For this purpose, it has to be recommended to the patient who has this implant to warn as far as possible the concerned health professionals (radiologists and radiology operators) about the presence of this implant before these exams.
7	Precautions : The BDyn devices are delivered sterile, sterilized by ethylene oxide. Before any use, inspect the integrity of the packaging and device (including peelable pouches). Do not use in the event of deterioration of the labels and/or the device and/or the packaging. Do not use if the device is out of date. The sterile spinal dynamic posterior stabilization device BDyn is only intended to adult human female or male more than 40 Kg. COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned. The installation of BDyn device can only be made with pedicular screws tested and/or approved by the company COUSIN BIOTECH. The eventual substitution of the BDyn spinal shock absorber has to be carried out only with the dedicated replacement rigid rods. You should never use implants made of stainless steel and titanium alloy in the same construct. A successful result is not always achieved in every surgical case. This fact is especially true in spinal surgery where many extenuating circumstances may compromise the results. This device system is not intended to be the sole means of spinal support. No spinal implant can withstand body loads without the support of bone. In this event, bending, loosening, disassembly and/or breakage of the device(s) will eventually occur. Preoperative and operating procedures, including precise knowledge of suitable surgical techniques, and proper selection of the good reference of the device adapted to the patient and its narrow setting up are important considerations in the successful use of the device by the surgeon. The BDyn device must be implanted only by a qualified surgeon, having knowledge in the use of the product and who has the knowledge of the anatomy, spinal surgery, pedicle screws fixation technique and specific BDyn device surgical technique Postoperative precautions
8	Materials included in the device : The BDyn Device components are manufactured from medical grade titanium TA6V (ISO 5832-3; ASTM F 136), silicone, long-term implantable polyurethane. BDyn implant: Titanium alloy Ti6Al4V ELI (ISO 5832-3; ASTM F 136), Polydimethyl siloxane (PDMS) MED 4770 unrestricted Polycarbonate urethane (PCU) Bionate® II 80A
9	Reporting a serious incident : If you encounter any serious incident with your BDyn implant you should report it immediately to your doctor. You should also report the incident to COUSIN BIOTECH via the link serviceclients@cousin-surgery.com or materiovigilance@cousin-biotech.com and to the Therapeutic Goods Administration using the following website: http://www.tga.gov
10	Legal Manufacturer & Manufacturing Facility : COUSIN BIOTECH, Allée des Roses, 59117 Wervicq-Sud, France

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