SEMI RESORBABLE PARIETAL REINFORCEMENT IMPLANT

en	Instructions for use	Page	2	MDR CE marking date: 30/09/2025
fr	Notices d'instructions	Page	7	Date de marquage CE MDR 20/09/2025
de	Gebrauchsanweisung	Seite		MDR CE-Kennzeichnungsdatum: 30.09.2025
it	Istruzioni per l'uso	Pagina		Data di marcatura CE MDR: 30/09/2025
es	Instrucciones de uso	Página		Fecha de marcado CE MDR: 30/09/2025
pt	Instruções de utilização	Página		Data de marcação CE MDR: 30/09/2025
ελ	Οδηγίες χρήσης	Σελίδα		Ημερομηνία σήμανσης CE MDR: 30/09/2025
cs	Návod k použití	Stránka		Datum označení CE podle MDR: 30/09/2025
da	Brugsvejledning	Side		MDR CE-mærkningsdato: 30/09/2025
fi	Käyttöohjeet	Sivu		MDR CE-merkinnän päivämäärä: 30/09/2025
hu	Használati Útmutató	Oldal		MDR CE jelölés dátuma: 2025/09/30.
nl	Gebruiksaanwijzing	Pagina		MDR CE-markeringdatum: 30/09/2025
ro	Instrucțiunile de utilizare	Pagină		Data marcatului CE MDR: 30/09/2025
sk	Návod na použitie	Strana		Dátum označenia CE podľa MDR: 30/09/2025
sv	Bruksanvisning	Sida		MDR CE-märkningsdatum: 30/09/2025
tr	Kullanım Talimatları	Sayfa		MDR CE işaretleme tarihi: 30/09/2025
lv	Lietošanas instrukcijas	Lappuse		MDR CE marķējuma datums: 30/09/2025
et	Kasutusjuhend	Lehekülg		MDR CE-märgistuse kuupäev: 30/09/2025
lt	Naudojimo instrukcijos	Puslapis		MDR CE žymėjimo data: 2025/09/30
en	Patient information leaflet (Australian requirements)	Page		

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NOT300_251017

Version du 17/10/2025

Caution : Federal law (USA) restricts this device to sale, distribution and use by or on the order of a physician.

This release is the last update of the instructions of use and replace the previous edition
Product under regulation 2017/745 (EU).

**SEMI-RESORBABLE PARIETAL REINFORCEMENT IMPLANT
STERILE PRODUCT FOR SINGLE USE**

MANUFACTURER IDENTIFICATION

COUSIN BIOTECH
ALLEE DES ROSES
59117 WERVICQ SUD
FRANCE

IDENTIFICATION AND DESCRIPTION OF THE DEVICE

4DMESH® meshes are semi-resorbable parietal reinforcement implants.

INTENDED USE

Intended to be used as parietal reinforcement implant.

INDICATIONS

4DMESH® meshes are designed for the repair and reinforcement of inguinal and femoral hernias.

CONTRAINDICATIONS

Do not use in the following cases:

- Allergy to any of the components
- Infected site
- Pregnancy
- Growing children

UNDESIRABLE SIDE EFFECTS

Like any implantable medical device, this implant is susceptible to generate possible undesirable side effects that may lead to reoperation:

- | | | | |
|-------------------|---------------------|---------------------|---------------------------|
| - Discomfort/Pain | - Adhesion | - Hematoma | - Foreign Body reaction |
| - Infection | - Fistula | - Mesh migration | - Irritation nearby organ |
| - Recurrence | - Erosion/extrusion | - Mesh deformation | - Obstruction |
| - Inflammation | - Seroma/Lymphocele | - Allergic reaction | |
- A small additional risk of male infertility has been identified after bilateral groin hernia repair using meshes (open or laparoscopic surgery)

TARGET POPULATION

Adult male or female suffering from groin hernias with an indication of surgical repair with extraperitoneal implantation. These include:

- Inguinal hernias
- Femoral hernias

IMPLANTED MATERIALS

Polypropylene (non resorbable) (PP) – Poly L Lactic Acid (resorbable) (PLLA). Blue thread (depending on references): monofilament polypropylene (non resorbable), dye: [phthalocyaninato(2-)] copper.

	Weight PLLA (g)	Weight PP (g)	Weight of the blue dye [Phthalocyaninato(2-)] copper (mg)
4DMESHPRSR	1,02	0,25	< 0,01
4DMESHPRSL	1,02	0,25	< 0,01
4DMESHPRLR	1,30	0,33	< 0,01
4DMESHPRLL	1,30	0,33	< 0,01
4DMESHPRXR	1,44	0,36	< 0,01
4DMESHPRXL	1,44	0,36	< 0,01
4DMESHRABA	1,34	0,33	< 0,01

4DMESH1215	1,60	0,40	< 0,01
4DMESH1317	1,39	0,35	< 0,01
4DMESH1216	1,29	0,32	< 0,01
4DMESH1717	2,60	0,65	NA (no landmark)
FBIOSRF128	0,67	0,17	NA (no landmark)
4DMESHM613	0,62	0,16	NA (no landmark)
4DMESHF812	0,85	0,21	NA (no landmark)
4DMESH1515	2,03	0,51	NA (no landmark)
4DMESH0715	0,95	0,24	NA (no landmark)
4DMESH0510	0,45	0,11	NA (no landmark)
4DMESH1015	1,16	0,29	< 0,01
4DMESHMBES	0,40	0,10	NA (no landmark)
4DMESHMBEL	0,62	0,16	NA (no landmark)
FBIOSRF952	0,38	0,09	NA (no landmark)

No human nor animal origin. Semi-resorbable.

PERFORMANCES

The device 4DMESH® contributes to reduce symptoms from hernia and contributes to avoid complications of hernia.

LINK TO SUMMARY OF SAFETY AND PERFORMANCE CHARACTERISTICS:

https://www.cousin-biotech.com/d-68f1fd7631392_tf07csscpissue-4.pdf

MRI SAFETY

The implant has not been tested for safety and compatibility in the case of postoperative investigation using Magnetic Resonance Imaging (MRI). However, the implant is mainly made with Polypropylene and Poly-L-Lactic Acid: non metallic components. Consequently, 4DMESH® is considered MR safe.

PRECAUTIONS FOR USE

4DMESH® must always be separated from abdominal cavity by peritoneum

4DMESH® prosthesis are delivered sterile (ethylene oxide sterilisation).

When opening the carton box, the surgeon has to make sure that the lot number, the model and the size of the mesh are equal to the data on the label of the inner pouch. Otherwise, the use of the prosthesis is strictly prohibited.

Before any use, inspect the integrity of packaging and device (of which blister / peelable pouches). Do not use in the event of deterioration of the device and/or the packaging. Do not use if the device is out of date. COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned.

Unpacking and manipulation of the implant must be carried out following asepsis standards in order to guarantee the sterility.

Prevent the mesh from any contact with objects that may alter their surface.

Meshes for laparoscopic surgery have a blue landmark that should be positioned towards the medial side of the patient (see the surgical technique below).

This device must be implanted by a qualified surgeon (familiar with the relevant anatomy and experienced in visceral surgery).

IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE

As specified on the product's labelling, the 4DMESH® mesh is for single use only. It can not be re-used and/or re-sterilized (potential risks would be and are not limited to: loss of the product's sterility, risk of infection, loss of the product's efficiency, recurrence)

INSTRUCTIONS TO BE GIVEN TO THE PATIENTS BY THE SURGEON

The surgeon should inform the patient of the potential physical and psychological restrictions and consequences of implanting the device. The patient must be informed of the surgical risks and possible side effects. The surgeon should invite the patient to return for further consultation if the patient presents symptoms that appear abnormal. An implant card for the patient is provided by the healthcare professional at the clinic. This card provides information to identify the device and traceability elements as well as the name, address and website of the manufacturer. Please refer to the section "Implant card information leaflet" for instructions on how to complete it. **The surgeon shall invite the patient to scan the implant card immediately upon receipt to keep track of it in case of loss.** The summary of the safety and performance characteristics of the device and the instructions for use can also be found on the COUSIN BIOTECH website.

STORAGE PRECAUTIONS

To be stored in a dry place away from light and at room temperature in its original packaging.

SURGICAL TECHNIQUE

DESCRIPTION OF THE SURGICAL PROCEDURES:

Open Inguinal hernia repair : Lichtenstein

- The patient is placed in the supine position
- Local or general anaesthesia
- Iliac incision of approximately 5 cm
- Dissection is performed accordingly to the Lichtenstein technique
- The hernia sac is identified and reduced
- The transversalis fascia can be tightened with a non-absorbable suture. With direct hernias it helps to avoid any immediate post-operative recurrence
- 4DMESH® mesh can be cut to fit the patient's anatomy
- The prosthesis is progressively put in place
 - begin by placing 4DMESH® on the pubic spine
 - the upper part of the prosthesis is placed on the conjoint tendon
 - carry placing the prosthesis towards the deep inguinal ring so that the first arm is placed behind
 - then the second arm is placed so that it encircles the spermatic cord
 - Fix the 4DMESH® mesh with non-resorbable sutures
 - apex is sutured to the pubic tubercle
 - the lower border of the mesh is sutured to the free edge of the inguinal ligament
 - continuous suture extends up just medial to the anterior superior iliac spine
 - suture the two tails of the mesh together around the spermatic cord
 - the infero-medial corner of the mesh is attached overlapping the pubic tubercle.
 - the mesh is anchored to the conjoint tendon by interrupted sutures
 - the cord then rests on the prosthesis
 - suction drain can be placed beneath the external oblique aponeurosis, especially in large inguinal hernias, where an extensive dissection was performed
 - the aponeurosis of external oblique is then closed with absorbable sutures
 - before closure of the surgical incision, its edges are infiltrated with a long-acting local anaesthetic, such as Naropein®
 - During peri-operative care of the patient, prophylactic antibiotics is usually given for 48 – 72 hours postoperatively

Laparoscopic procedure (TAPP or TEP)

1/ Example of a TAPP repair

- General anaesthesia
- A pneumoperitoneum is created in the usual fashion
- intra-abdominal findings are reported [intra-abdominal pathology and inguinal hernia defects and sacs].
- The two additional trocars are inserted under direct vision.
- Dissection is initiated and carried out following classical TAPP technique
- The indirect inguinal hernia sac should be dissected carefully from the Spermatic Cord. Direct hernia sacs are easily dissected.
- The 4DMESH® mesh is inserted through the Trocar into the intra-abdominal cavity and deployed over the inguinal region. The blue thread on the oval and preshaped mesh is the anatomic landmark for the pubic bone. For the meshes with flaps, the cross shows the upper flap. The 4DMESH® mesh is attached to Cooper's Ligament, around and lateral to the Inferior Epigastric Vessels using tacks, staples or other mesh fixing device, according to the Surgeon's usual technique. The operator should be meticulous to avoid the triangle of doom, triangle of pain, iliac vessels and to place the fixations lateral to the inguinal ring. Mesh fixation points should be at least 1cm from the edge of the mesh with 1 cm spacing between fixation points.
- The operator should check the 4DMESH® mesh is well anchored to the surrounding structures.
- The peritoneum is closed meticulously and no defect between the peritoneum and the abdominal wall should be left open. In addition, it should cover the entire Mesh.
- The trocars are removed under direct vision. The fascia of the sub-umbilical trocar site is closed as needed.

2/ Example of a TEP repair

- General anaesthesia
- Create the pro-peritoneal space.
- The 10 mm trocar is inserted.
- Create the pneumo-pro-peritoneum.
- Two 5 mm ports are inserted under direct vision.
- If a peritoneal defect is created at this time it should be immediately closed.
- The anatomical landmarks are identical to when performing a TAPP repair.
- The indirect inguinal hernia sac should be dissected carefully from the Spermatic Cord. Direct hernia sacs are easily dissected. This is done by gentle traction with atraumatic graspers.
- The 4DMESH® mesh is inserted into the pre-peritoneal space and deployed over the inguinal region. The blue thread on the oval and preshaped mesh is the anatomic landmark for the pubic bone. For the meshes with flaps, the cross shows the upper flap.
- The 4DMESH® mesh is first attached or secured to Cooper's Ligament, and the superior aspect of the pubic ramus. The 4DMESH mesh is then fixed on the posterior aspect of the Linea Alba. The anchoring is continued around and lateral to the Inferior Epigastric Vessels. A few fixations are used to fix the 4DMESH mesh lateral to the internal ring. Mesh fixation points should be at least 1cm from the edge of the mesh with 1 cm spacing between fixation points.
- The 4DMESH® mesh should gently and generously cover the Iliac Vessels without major gaps.

- The operator should check the Mesh is well anchored to the surrounding structures. The Mesh should not migrate and remain in place.
- The pre-peritoneal space will be checked for any peritoneal defects. If any, they should be closed. The 5 mm trocars are removed under direct vision. The pre-peritoneal space will collapse.
- If there are any questions about a missed peritoneal defect, a completion laparoscopy should be performed.
- Trocars are removed and the fascial defect closed with the appropriate suture. The skin edges are approximated in the usual manner.

Size of the prosthesis:

The 4DMESH® product range is composed of meshes of various sizes and shapes:

- Pre-cut meshes with keyhole for LICHTENSTEIN inguinal hernia repair
- Oval, pre-shaped and rectangular meshes for TAPP/TEP inguinal hernia repair

EXPLANTATION AND DISPOSAL OF DEVICES

Risks associated with the explantation: the surgeon who decides to remove the device should consider factors such as the risk of another surgery for the patient, and the difficulty of the explantation procedure. The removal of the implant is up to the surgeon's decision and must be subject to adequate postoperative follow-up.

Devices should be retrieved and handled according to the recommendations of ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling". Any explanted device must be sent back for analysis, following the current protocol. This protocol is available on request from COUSIN BIOTECH. It is important to note that any implant that must not be cleaned or disinfected prior to despatch must be contained in sealed packaging. The removed medical device must be disposed of in accordance with the country's standards for the disposal of infectious waste.

There are no specific recommendations for the disposal of a non-implanted device.

IMPORTANT

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

INFORMATION REQUEST AND COMPLAINTS

Following its quality policy, COUSIN BIOTECH is committed to making every effort to produce and supply a high-quality medical device. However, if a health professional (client, user, prescriber, etc.) has a complaint or cause for dissatisfaction with a product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In case of failure of an implant or if it has contributed to causing a serious adverse reaction in the patient, the health centre must follow the legal procedures in that country and inform COUSIN BIOTECH immediately. For any correspondence, please specify the reference, batch number, contact person's details and a comprehensive description of the incident or complaint.

Brochures, documentation and surgical technique are available on request from COUSIN BIOTECH and its distributors. If further information is needed or required, please contact your COUSIN BIOTECH representative or distributor.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

IMPLANT CARD INFORMATION LEAFLET

How to complete the implant card ?

The implant card, for the patient, should be completed as follows:

Front side of the implant card

The user has to complete the implant card before transmitting the implant card to the patient. Please complete it as follows:

	Carte Implant - Implant Card A conserver / Keep it	Carte Implant - Implant Card
Patient name	 _____	 COUSIN SURGERY <i>We care for Surgery</i>
Date of surgical procedure	 . . / . . /	
Name of practitioner	 _____	

Back side of the implant card

No fields need to be completed by the user. The following information on the device is printed directly by the manufacturer on the back of the implant card. A QR code allowing the patient to consult the instruction for use, as well as the website address, are available on the implant card.

Information to identify the device	Traceability information
• Indication that the product is a medical device • Description of the device family • Trade name of the device • Reference of the device • Device is MR safe • Denomination name of the device in sales countries	• UDI number • Batch number • Name and address of the legal manufacturer

Explanation of the used symbols in the Implant Card

All symbols are commonly used and explained at the end of the instruction for use.

**IMPLANT DE RENFORCEMENT PARIETAL SEMI-RESORBABLE
PRODUIT STERILE A USAGE UNIQUE**

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Patient information leaflet 4DMESH®

Item	Patient information leaflet
1	a. The name of the device : 4DMESH® b. The model of the device : 4DMESH0715, 4DMESH1015, 4DMESH1215, 4DMESH1216, 4DMESH1317, 4DMESH1515u, 4DMESH1717, 4DMESHF812, 4DMESHM613, 4DMESHMBEL, 4DMESHMBES, 4DMESHPRLL, 4DMESHPRLR, 4DMESHPRSL, 4DMESHPRSR, 4DMESHPRSR, 4DMESHPRXL, 4DMESHPRXR, 4DMESHHRABA, FBIOSRF128, FBIOSRF952
2	a. The intended purpose of the device : 4DMESH® meshes are designed for the repair and reinforcement of inguinal and femoral hernias. b. The kind of patient on whom the device is intended to be used : Adult with inguinal and femoral hernias.
3	Special operating instructions for the use of the device : The device must be implanted only by a qualified surgeon and trained with the use of the product (knowledge of anatomy and visceral surgery).
4	a. The intended performance of the device : The device 4DMESH® is designed to be extra-peritoneally implanted. These products are biocompatible and semi-resorbable. Their main benefits are: can be cut to fit patient's anatomy, high resistance to suture and abdominal pressure, macroporous, transparent, lightweight after PLLA resorption. b. Undesirable side effects : Like any implantable medical device, this implant could generate possible undesirable effects : - Discomfort/pain - Infection - Recurrence - Inflammation - Adhesion formation - Fistula formation - Erosion - Seroma - Hematoma - Mesh migration - Mesh deformation - Allergic reaction - Foreign body reaction - Irritation near by organ - Lymphocele - A small additional risk of male infertility has been identified after bilateral groin hernia repair using meshes (open or laparoscopic)

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5	When 4DMESH® Should Not be Used : Allergy to any of the components, infected site, pregnancy, growing children, anticoagulant therapy.
6	MRI Information : 4DMESH® implants are extraperitoneal partially resorbable parietal reinforcement implants. This surgical mesh is composed of Poly-L-Lactic Acid (PLLA) and Polypropylene (PP).
7	Precautions : 4DMESH® must always be separated from abdominal cavity by peritoneum 4DMESH® prosthesis are delivered sterile (ethylene oxide sterilisation). Before any use, inspect the integrity of packaging and device (of which blister / peelable pouches). Do not use in the event of deterioration of the device and/or the packaging. Do not use if the device is out of date. COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned. This device must be implanted only by a qualified surgeon and trained with the use of the product (knowledge of anatomy and visceral surgery).
8	Materials included in the device: Polypropylene (non resorbable) – Poly L Lactic Acid (resorbable). Blue thread (depending on references): monofilament polypropylene, non absorbable, colorant [phthalocyaninato(2-)] copper. No human nor animal origin. Semi-resorbable.
9	Reporting a serious incident : If you encounter any serious incident with your 4DMESH® implant you should report it immediately to your doctor. You should also report the incident to Cousin Biotech via the link serviceclients@cousin-surgery.com or materiovigilance@cousin-biotech.com and to the Therapeutic Goods Administration using the following website: http://www.tga.gov
10	Legal Manufacturer : Cousin Biotech, Allée des roses, 59117 Wervicq-Sud, France Manufacturing Facility : Cousin Biotech, Allée des roses, 59117 Wervicq-Sud, France

This image shows a full page of white paper with horizontal blue ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

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Cousin Biotech SAS au capital de 320 000 € - NAF 1396Z- N° TVA FR 34 398 460 261 - Siren 398 460 261 RCS Lille