## STERILE MINI-INVASIVE DEVICES FOR INTERSPINOUS SPACE WITH LAMINAR SUPPORT

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Made in France

**NOT296**  
Version du 09/09/2025

This release is the last update of the instructions of use and replace the previous edition  
**Product under regulation 2017/745 (EU).**

**STERILE MINI-INVASIVE DEVICES FOR INTERSPINOUS SPACE WITH LAMINAR SUPPORT  
SINGLE USE STERILE PRODUCT**

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**IDENTIFICATION AND DESCRIPTION OF THE DEVICE**

The BatSPINE prosthesis is an interlaminar medical device that comprises a dimethyl siloxane (silicon) wedge (available in several sizes) with a nitinol frame adapted to the intervertebral space. The wedge is covered with woven polyethylene terephthalate with an area covered by sheeting on the front face in the medullary zone.

**INDICATIONS**

BatSPINE prosthesis is an interspinous with laminar support for discal and/or facet assistance indicated for use from L1-S1 in skeletally mature patients with at least moderate impairment in function, who experience relief in flexion from their symptoms of leg/buttocks/groin pain, with or without back pain, and who have undergone at least 6 months of nonoperative treatment and indicated for use in:

- Arthropathic facet-syndrome
- Foraminal stenosis
- Degenerative discopathy

**CONTRAINDICATIONS**

**1) Do not use the prosthesis in the following cases:**

- Allergy to one of the components
- Growing child
- Infected site
- Pregnant woman

**2) Factors likely to compromise successful implantation:**

- Severe osteoporosis
- Major deformations of the spinal column
- Local bone tumors
- Systemic or metabolic disorders
- Infectious diseases
- Obesity
- Drug addictions
- Surgery at more than 2 levels with BatSPINE prosthesis
- Spondylolisthesis from grade 2 to 4
- Intervertebral disc with Pfirrmann classification from grade 5
- S1 Spinous Process length inferior to 18mm when implanted in L5/S1

**UNDESIRABLE SECONDARY EFFECTS**

**1) Complications likely related to the device:**

- Infection of the prosthesis
- Device breakage (fixed wing breakage, knitting tearing)
- Device migration displacement of the wedge
- Wrong level implantation, inappropriate or inaccurate positioning (posterior or wrong level)
- Seroma/swelling around the device
- Inflammatory reaction, allergy, foreign body reaction
- Spondylolisthesis to the dynamically stabilized segment
- Persistent or worsening pain attributed to the dynamically stabilized level and required revision surgery with interbody fusion at that level
- Spinous process fractures intraoperative post-operative or during device removal (removal due to infection) and can be associated with dislocation of the implant
- Late reoperation due to absence of recovery
- Disability
- Discomfort
- In certain cases, the medical instrument had to be removed and pedicle screw fusion had to be performed.

**2) Complications unlikely to be related to the device, but related to surgery or to the patient's conditions:**

- Pseudomeningocele, fistula, dural tear, persistent CSF leakage, meningitis
- Loss of neurological function,
- Cauda equina syndrome, neuropathy, neurological deficits, arachnoiditis, and/or muscle loss, compression around nerves and/or pain
- Urinary retention or loss of bladder control or other types of urological system compromise
- Fracture, bone loss or decrease in bone density
- Herniated nucleus pulposus, disc disruption or degeneration at, above, or below the level of surgery: Adjacent segment degeneration
- Non-union or pseudarthrosis, delayed union, mal-union
- Cessation of any potential growth of the operated portion of the spine

- Loss of or increase in spinal mobility or function
- Hemorrhage, hematoma, edema, embolism, stroke, excessive bleeding, phlebitis, wound necrosis, wound dehiscence, wound infection, damage to blood vessels, or other types of cardiovascular system compromise
- Development of respiratory problems

#### **TARGET POPULATION**

Skeletally mature patients.

#### **IMPLANTED MATERIALS**

BatSPINE prosthesis: Dimethyl siloxane, Polyethylene terephthalate, Nitinol.

| Reference Material         | RCBBATS08U      | RCBBATS10U      | RCBBATS12U      | RCBBATS14U      | RCBBATS16U      |
|----------------------------|-----------------|-----------------|-----------------|-----------------|-----------------|
| Dimethyl siloxane          | 2.61 +/- 0.5 g  | 3.05 +/- 0.5 g  | 3.47 +/- 0.5 g  | 3.96 +/- 0.5 g  | 4.44 +/- 0.5 g  |
| Polyethylene terephthalate | 0.67 +/- 0.01 g | 0.84 +/- 0.01 g | 1.00 +/- 0.01 g | 1.17 +/- 0.01 g | 1.34 +/- 0.01 g |
| Nitinol                    | 0.37 +/- 0.01 g | 0.40 +/- 0.01 g | 0.43 +/- 0.01 g | 0.46 +/- 0.01 g | 0.52 +/- 0.01 g |

#### **ACCESSORIES**

BatSPINE instruments: stainless steel, polyphenylsulfone, aluminium, polyamide, PTFE.

Non-human and non-animal origin – Non-resorbable.

#### **ACTION AND PERFORMANCES**

The BatSPINE prosthesis is an elastic support device with a cushioning effect. Its interlaminar implantation enables it to be close to the centre of rotation, thereby producing the cushioning effect and providing tension on the interspinous ligamentary structures due to the elasticity of the supra-spinous wedge.

#### **LINK TO SUMMARY OF SAFETY AND PERFORMANCE CHARACTERISTICS :**

[https://www.cousin-biotech.com/d-6318526d8b749\\_under-construction.pdf](https://www.cousin-biotech.com/d-6318526d8b749_under-construction.pdf)

#### **MRI SAFETY**

Non-clinical testing has demonstrated that the BatSPINE range is MR Conditional in accordance with the ASTM F2503-20 standard definition. A patient with this device can be safely scanned in an MR system meeting the following conditions. Failure to follow these conditions may result in injury.

|                                               |                                                                                                                                                                                                                                                                                                                                                                                         |
|-----------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Maximum Static Magnetic Field Strength</b> | 1.5T / 3.0T                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Type of Nuclei</b>                         | Hydrogen                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Static Magnetic Field (B0) Orientation</b> | Horizontal                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Magnet Type</b>                            | Cylindrical-bore                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Maximum Spatial Field Gradient</b>         | 1.5 T: 67.63 T/m (6763 G/cm)<br>3 T: 33.81 T/m (3381 G/cm)                                                                                                                                                                                                                                                                                                                              |
| <b>RF Excitation</b>                          | 1.5T: Circularly Polarized (CP)<br>3T: Multichannel-2 (MC-2)                                                                                                                                                                                                                                                                                                                            |
| <b>RF Transmit Coil Type</b>                  | Whole body and extremities local transmit coils can be used                                                                                                                                                                                                                                                                                                                             |
| <b>Receive Coil Type</b>                      | Any receive coil can be used                                                                                                                                                                                                                                                                                                                                                            |
| <b>RF Operating mode</b>                      | Normal Operating Mode                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Maximum Whole-Body SAR</b>                 | Whole body SAR ≤ 2 W/kg                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Scan duration and wait time</b>            | Scanning for 1h maximum without a cooling period                                                                                                                                                                                                                                                                                                                                        |
| <b>MR Image Artifact</b>                      | The presence of the item within the MR image field of view may produce an image artifact. Some manipulation of scan parameters may be needed to compensate for the artifact. In non-clinical testing, the image artifact caused by the device extends approximately 8 mm from the device edge when imaged with a spin echo pulse sequence and 8 mm with a gradient echo, both at 1.5 T. |

#### **PRECAUTIONS FOR USE**

Before usage, verify the integrity of the prosthesis, instruments and packaging.

Do not use in the event of deterioration of the prosthesis and/or the instruments, and/or the labels and/or the packaging.

Do not use if the prosthesis has exceeded its use-by date. Unpacking and manipulation of the implant must be carried out following asepsis standards in order to guarantee sterility.

When opening the carton box, the surgeon or the operating room nurse has to make sure that the lot number, the model and the size of the prosthesis are equal to the data on the label of the inner pouch. Otherwise, the use of the prosthesis is strictly prohibited.

Prevent the prosthesis and instruments from coming in to contact with objects that may alter their surface.

The BatSPINE prosthesis is made with soft materials; any contact with harder ones may alter its mechanical behaviour and lifetime. It is mandatory to avoid contact with other materials.

Implants must be resistant to biological fluids and secretions from the tissues with which they come into contact during use.

Implants must be handled and stored with great care, in a dry place, protected from sunlight and at room temperature.

Do not maintain the folding more than 30 minutes. Use of COUSIN BIOTECH instrument is required for implantation of the BatSPINE prosthesis.

The BatSPINE prosthesis and instruments must only be used by a qualified neurosurgeon, orthopaedic surgeon or spine surgeon who is trained in spinal surgery, with mini invasive experience and trained to the surgical technique of BatSPINE. The information contained in these instructions for use is necessary but not sufficient to control the surgical technique.

The prescription of the device is decided on by the surgeon, the only person qualified to do so. Topping off is not contra-indicated but special attention is required.

Transportation: Transportation of the devices should be done in the closed vehicles in accordance with the rules for those transport vehicles. There are no special transportation requirements for the device.

### **SURGICAL TECHNIQUE**

1. Place the patient in the prone knee position or prone position.
2. Perform an incision of 2.5 cm long at 2.5 cm away from midline and insert a generic sleeve by interlaminar approach with the use of dilators.
3. With the curator, prepare interspinous and interlaminar space making sure to penetrate the contra-lateral side for device position. Control the space created by using the space-checker.
4. Then the interlaminar space is distracted gradually by inserting the distractor-sizers, size by size with a clockwise movement. Measure the interlaminar space with the distractor-sizers to determine the implant size.
5. Once the size of the device has been chosen, use the ramp-distractor between the laminae to perform distraction.
6. Place the prosthesis in the compressor to compress it allowing the insertion into the BatSPINE vehicle instrument. Connect the BatSPINE vehicle to the gun and the ramp to the BatSPINE vehicle. Then press the trigger to release the prosthesis in the intervertebral space. Remove the ramp-distractor with anti-clockwise movement and use the positioner to maintain the prosthesis in the interlaminar space.
7. Control the position of the device.
8. Remove the sleeve and suture the skin.

### **IMPORTANT: DO NOT REUSE AND DO NOT RESTERILISE THE BATSPINE PROSTHESIS**

As indicated on the product labelling, the prosthesis is designed for single use. It should in no event be reused and/or resterilized (potential risks include but are not limited to: loss of product sterility, risk of infection, loss of effectiveness, relapse).

### **WARNING ABOUT INSTRUMENTS**

The BatSPINE instruments are intended for temporary use and are re-usable.

Unlike the BatSPINE prosthesis, instruments are furnished unsterile. For more details on instruments, please refer to the instructions for use specific to information on the care, decontamination, cleaning, disinfection, maintenance and sterilisation of class I reusable surgical instruments manufactured and/or distributed by COUSIN BIOTECH.

The Non-sterile instruments for BatSPINE are divided into 3 instruments system comprising :

- Articulated arm for MIS (RCBANARM1U)
- Clamp for articulated arm MIS (RCBANARM2U)
- Generic sleeve 50 mm (RCBANSL05U)
- Generic sleeve 70 mm (RCBANSL07U)
- Generic sleeve 90 mm (RCBANSL09U)
- Generic sleeve handle (RCBANSLHAU)
- Compression system (RCBANBCOMU)
- Small and large inlays (RCBANSUINU / RCBANSLINU / RCBANLUINU / RCBANLLINU)
- Small and large vehicles (RCBANSIMVU / RCBANLIMVU)
- Gun (RCBANBGUNU)
- 5 dilatation tubes for MIS (RCBANDILAU) – from 2 to 18 mm
- Space checker (RCBANECU)
- Spinous curator (RCBANCURAU)
- 5 Sizers distractor 8, 10, 12, 14, 16 mm – (RCBANSD08U / RCBANSD10U / RCBANSD12U / RCBANSD14U / RCBANSD16U)
- Small and large obturators (RCBANSOBTU / RCBANLOBTU)
- 2 sizers handle (RCBANSIHAU)
- Small and large distracting ramp (RCBANSRAMU/ RCBANLRAMU)
- Positioner (RCBANPOSU)

### **DECONTAMINATION, CLEANING AND STERILIZATION OF INSTRUMENTS**

Please refer to Instructions for Use NOT348 related to information on the care, decontamination, cleaning, disinfection, maintenance and sterilisation of class I reusable surgical instruments manufactured and/or distributed by COUSIN BIOTECH.

### **POST-OPERATIVE PRECAUTIONS**

The patient should be informed of the post-operative care and precautions to be observed.

The use of a brace may be necessary but is not mandatory. This decision can only be made by the surgeon.

### **INSTRUCTIONS TO BE GIVEN TO THE PATIENTS BY THE SURGEON**

The surgeon should inform the patient of the potential physical and psychological restrictions and consequences of implanting the device. The patient must be informed of the surgical risks and possible side effects. The surgeon should invite the patient to return for further consultation if the patient presents symptoms that appear abnormal. An implant card for the patient is provided by the healthcare professional at the clinic. This implant card provides information, for the patient, to identify the device and traceability elements as well as the name, address and website of the manufacturer. The surgeon shall invite the patient to scan the implant card immediately upon receipt to keep track of it in case of loss. The surgeon shall also inform the patient that the summary of the safety and performance characteristics of the device and the instructions for use can be found on the COUSIN BIOTECH website.

### **STORAGE PRECAUTIONS**

Keep in a dry place away from sunlight, at room temperature, in its original packaging. Implants must be handled and stored with great care, in a dry place, protected from sunlight and at room temperature in their original packaging.

### **EXPLANTATION AND DISPOSAL OF DEVICES**

Explantation and handling should be carried out in accordance with ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling". Any explanted device must be returned for analysis as per the current protocol. This protocol is available from COUSIN BIOTECH on request. It is important to note that any implant that has not been cleaned and disinfected prior to dispatch must be placed in a sealed package. Explanted medical devices must be disposed of in accordance with local regulations governing the disposal of infectious hazardous waste. The disposal of non-implanted devices is not subject to specific recommendations. The second paragraph covers information requests or complaints.

Risks associated with the explantation: the surgeon who decides to remove the device should consider factors such as the risk of another surgery for the patient, and the difficulty of the explantation procedure. The removal of the implant is up to the surgeon's decision and must be subject to adequate postoperative follow-up.

#### INFORMATION REQUESTS AND COMPLAINTS

In accordance with its quality policy, COUSIN BIOTECH is committed to producing and supplying high-quality medical devices. However, if health professionals (client, user, prescriber, etc.) wish to lodge a complaint or report dissatisfaction with a given product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In the event of implant failure or if the implant is implicated in serious adverse effects experienced by the patient, the health centre must follow the legal procedures in the country in question, and inform COUSIN BIOTECH without delay. The reference, batch number, reference contact details and a comprehensive description of the incident or complaint must be stated in all correspondence. Brochures, documentation and surgical techniques are available on request from COUSIN BIOTECH and its distributors.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

#### How to complete the implant card ?

The implant card, for the patient, should be completed as follows:

##### Front side of the implant card

The user has to complete the implant card before transmitting the implant card to the patient. Please complete it as follows:

|                            | Carte Implant - Implant Card                                                                          | Carte Implant - Implant Card                                                                                                               |
|----------------------------|-------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------|
|                            | <b>A conserver / Keep it</b>                                                                          |                                                                                                                                            |
| Patient name               |  _____               | <br><b>COUSIN SURGERY</b><br><i>We care for Surgery</i> |
| Date of surgical procedure |  . . / . . / . . . . |                                                                                                                                            |
| Name of practitioner       |  _____             |                                                                                                                                            |

##### Back side of the implant card

No fields need to be completed by the user. The following information on the device is printed directly by the manufacturer on the back of the implant card. A QR code allowing the patient to consult the instruction for use, as well as the website address, are available on the implant card.

| Information to identify the device                                                                                                                                                                                                                                                                                | Traceability information                                                                                                                     |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------|
| <ul style="list-style-type: none"> <li>• Indication that the product is a medical device</li> <li>• Description of the device family</li> <li>• Trade name of the device</li> <li>• Reference of the device</li> <li>• Device is MR safe</li> <li>• Denomination name of the device in sales countries</li> </ul> | <ul style="list-style-type: none"> <li>• UDI number</li> <li>• Batch number</li> <li>• Name and address of the legal manufacturer</li> </ul> |

#### Explanation of the used symbols in the Implant Card

All symbols are commonly used and explained at the end of the instruction for use.

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