BIOMESH®

CA.B.S'Air® Composite

STERILE INTRAPERITONEAL NON-RESORBABLE PARIETAL REINFORCEMENT IMPLANTS

en	Instructions for use	Page	2	CE marking date: 09/05/2025
fr	Notice d'instructions	Page	6	Date de marquage CE :En 09/05/2025
de	Gebrauchsanweisung	Seite	10	Datum der CE-Kennzeichnung: 09/05/2025
it	Istruzioni per l'uso	Pagina	14	Data di marcatura CE: 09/05/2025
es	Instrucciones de uso	Pagine	19	Fecha de marcado CE: 09/05/2025
pt	Instruções de utilização	Página	23	Data da Marcação CE: 09/05/2025
ελ	Οδηγίες χρήσης	Σελίδα	27	Ημερομηνία σήμανσης CE: 09/05/2025
nl	Gebruiksaanwijzing	Pagina	31	Datum CE-markering: 09/05/2025
da	Brugervejledning	Side	35	CE-mærkningsdato: 09/05/2025
fi	Käyttöohje	Sivu	39	CE-merkinnän päivämäärä: 09/05/2025
lv	Lietošanas instrukcija	Lapa	43	CE marķējuma datums: 09/05/2025
ro	Instrucțiuni de utilizare	Pagina	47	Data marcatului CE: 09/05/2025
sv	Bruksanvisning	Sida	51	CE-märkningsdatum: 09/05/2025
tr	Kullanım Talimatları	Sayfa	55	CE işareti tarihi: 09/05/2025
en	Patient information leaflet (Australian requirements)	Page	59	CE marking date: 09/05/2025

R_x
only

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CE
0476

NOT313_251103
Version du 03/11/2025

Caution : Federal law (USA) restricts this device to sale,
distribution and use by or on the order of a physician

This release is the last update of the instructions of use and replace the previous edition
Product under regulation 2017/745 (EU).

STERILE INTRAPERITONEAL NON-RESORBABLE PARIETAL REINFORCEMENT IMPLANTS**BIOMESH® CA.B.S.'AIR® COMPOSITE****Expandable parietal prosthesis with inflatable balloon for visceral surgery****SINGLE USE STERILE PRODUCT****MANUFACTURER IDENTIFICATION**

COUSIN BIOTECH
 ALLEE DES ROSES
 59117 WERVICQ SUD
 FRANCE

IDENTIFICATION AND DESCRIPTION OF THE DEVICE

The BIOMESH® CA.B.S.'AIR® COMPOSITE prosthesis consists of:

- A parietal prosthesis made of an ePTFE and polypropylene plug, two or four polyethylene terephthalate with biocompatible colorant (D&C Green #6) and ePTFE fixation threads crimped with stainless steel needles (see on package).
- A balloon which allows the deployment of the parietal prosthesis (except for CABST1R05T).

INDICATIONS

Repair and parietal reinforcement of ventral, umbilical and incisional hernias requiring reinforcement with non absorbable support material.

CONTRAINDICATIONS

Do not use in the following cases:

- Allergy to any of the components
- Infected site
- Pregnancy
- Growing children
- Anticoagulant therapy

UNDESIRABLE SIDE EFFECTS

Like any implantable medical device, this implant is susceptible to generate possible undesirable effects:

- | | | |
|---------------------|---------------------------|----------------------|
| - Discomfort/Pain | - Recurrence | - Adhesion formation |
| - Erosions | - Irritation nearby organ | - Infection |
| - Inflammation | - Fistula formation | - Seroma |
| - Mesh deformation | - Hematoma | - Mesh migration |
| - Allergic reaction | - Foreign body reaction | - Obstruction |

TARGET POPULATION

Adult with ventral hernia including umbilical and incisional hernia necessitating a surgical repair with a method that includes an intraperitoneal implantation.

IMPLANTED MATERIALS

Polypropylene, Polyethylene terephthalate, ePTFE, Biocompatible colorant (D&C Green #6)

Reference	Weight of ePTFE (g)	Weight of Polypropylene (g)	Weight of PET (g)
CABST1R05T	0,01 g	0,57 g	0,19 g
CABSAIRC7	0,02 g	1,13 g	0,19 g
CABSAIRC09	0,04 g	2,03 g	0,19 g

Materials in contact with the patient during the surgery :

Stainless steel, Poly vinyl chloride.

Origin neither human nor animal – Non-absorbable.

PERFORMANCES

The **BIOMESH® CA.B.S.'AIR® COMPOSITE** prosthesis has a non-adhesive side, which is designed to avoid or reduce visceral adhesions during intra peritoneal use or extra peritoneal use in the presence of a peritoneal defect for ventral, umbilical and incisional hernias.

Products, biocompatible and non resorbable, cause reactional fibrosis taking over the reinforcement after six months of implantation. They have the advantage of having a very high resistance to suture, have a great flexibility and allow a optimum and rapid integration and colonisation.

If necessary, the surgeon has the possibility of adding two additional points of fixing, according to his appreciation.

LINK TO THE SAFETY AND CLINICAL PERFORMANCE SUMMARY :
https://www.cousin-biotech.com/d-69087e49880cc_tf03sscpissue-g.pdf

MRI SAFETY

The implant has not been tested for safety and compatibility in the case of postoperative investigation using Magnetic Resonance Imaging (MRI). However, the implant is mainly made with polypropylene, polyethylene terephthalate and non metallic component. Consequently, **BIOMESH® CA.B.S.'AIR® COMPOSITE** is considered MR safe.

PRECAUTIONS FOR USE

BIOMESH® CA.B.S.'AIR® COMPOSITE prosthesis are delivered sterile (ethylene oxide sterilization).

Before any use, inspect the integrity of packaging and device (of which blister / peelable pouches). Do not use in the event of deterioration of the device and/or the packaging.

Do not use if the device is out of date.

COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned.

This device must be implanted only by a qualified surgeon and trained with the use of the product (knowledge of anatomy and visceral surgery).

IMPORTANT : DO NOT REUSE - DO NOT RESTERILIZE

THIS BALLOON MUST IMPERATIVELY BE WITHDRAWN

As specified on the product's labeling, the **BIOMESH® CA.B.S.'AIR® COMPOSITE** mesh is for single use only. It cannot be re-used and/or re-sterilized (potential risks would be and are not limited to: loss of the product's sterility, risk of infection, loss of the product's efficiency, recurrence)

INSTRUCTIONS TO BE GIVEN TO THE PATIENTS BY THE SURGEON

The surgeon should inform the patient of the potential physical and psychological restrictions and consequences of implanting the device. The patient must be informed of the surgical risks and possible side effects. The surgeon should invite the patient to return for further consultation if the patient presents symptoms that appear abnormal. An implant card for the patient is provided by the healthcare professional at the clinic. This implant card provides information, for the patient, to identify the device and traceability elements as well as the name, address and website of the manufacturer. The surgeon shall invite the patient to scan the implant card immediately upon receipt to keep track of it in case of loss. The surgeon shall also inform the patient that the summary of the safety and performance characteristics of the device and the instructions for use can be found on the COUSIN BIOTECH website.

STORAGE

To be stored in a dry place away from sunlight and at room temperature in its original packaging.

SURGICAL TECHNIQUE

SURGICAL PROCEDURE (except for CABST1R05T which comes without balloon)					
1-Installation: - The prosthesis equipped with its balloon crosses completely the hernial defect, -its positioning becomes then entirely retro parietal (intraperitoneal)		2-Deployment: - The expansion and the spreading of the prosthesis are obtained by inflating the balloon, - Parietal fixation is carried out by two sutures thanks to the pre-positioned threads		3- Prosthesis in place: -The prosthesis is fixed, -The balloon must imperatively be withdrawn	

Description of the prosthesis:

This prosthesis is in a shape of a flat plug with 2 fixation threads crimped with needles (4 needles) for the product references CABST1R05T and CABS AIRC7, and with 4 fixation threads crimped with needles (8 needles) for the CABS AIRC09 product reference.

The aim of its innovating expansion balloon concept is to make dissection easier and to secure the unfolding of the prosthesis. This balloon is located inside the plug; it is filled-up with 30cm³ of air. It acts as a pneumatic Jack ensuring complete crease-free unfolding of the prosthesis.

Its double-sided design combines a smooth ePTFE side which reduces viscera adhesions with a knitted monofilament polypropylene side which allows better colonisation of the peritoneum.

Sizes of the prosthesis:

The medical device is a sterile kit.

The product range is composed of:

CABS AIR 5cm diameter: for small hernias and ventral hernias on trocar site

CABS AIR 7cm diameter: for umbilical and linea alba hernias, as well as small ventral hernias < 2cm

CABSAIR 9cm diameter: for umbilical and linea alba hernias, as well as small ventral hernias that measure between 2 and 5cm

Installation of the patient :

Patient lying on his back

Both legs kept in the body's axis

Each arm placed on an arm wrestler in a 90° angle from the body's axis

1-Incision :

The cutaneous incision (2 to 4cm according to the size of the defect) is centred on the hernia.

2-Dissection and resection of the peritoneal sac :

The hernial sac is identified and its content is re-introduced back inside the abdominal cavity. The sac is then left opened in order to allow the Intra-peritoneal prosthesis to be put in place.

Next step consists in liberating epiploico-peritoneal adhesions around the edges of the hernia so that the prosthesis can fully expand.

3-Positioning and fixating the prosthesis :

BIOMESH® CA.B.S.'AIR® COMPOSITE is introduced through the hernial orifice after being rolled up into a tube shape or folded in an umbrella-like way.

The balloon is filled up in order to allow the expansion of the prosthesis and the placement of the prosthesis is checked so that it remains centred above the hernial orifice.

The prosthesis is fixed with the two pre-mounted threads that will transfix the peritoneum and the muscular wall.

Because of its Intraperitoneal placement, two additional points should be used to ensure secure fixation of the prosthesis and to avoid any epiploic tissue being caught in between the prosthesis and the wall.

Be careful not to suture the balloon to the implant when fixing the prosthesis.

The balloon is deflated and removed.

4-Closure of the incision :

Musculo-aponeurotic closure in front of the prosthesis with resorbable suture thread. Sub-cutaneous and then cutaneous closure, normally without drainage.

EXPLANTATION AND DISPOSAL OF DEVICES

Devices should be retrieved and handled according to the recommendations of ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling". Any explanted device must be sent back for analysis, following the current protocol. This protocol is available on request from COUSIN BIOTECH. It is important to note that any implant that must not be cleaned or disinfected prior to despatch must be contained in sealed packaging. The removed medical device must be disposed of in accordance with the country's standards for the disposal of infectious waste.

There are no specific recommendations for the disposal of a non-implanted device.

Risks associated with the explantation: the surgeon who decides to remove the device should consider factors such as the risk of another surgery for the patient, and the difficulty of the explantation procedure. The removal of the implant is up to the surgeon's decision and must be subject to adequate postoperative follow-up.

IMPORTANT

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

INFORMATION REQUEST AND COMPLAINTS

Following its quality policy, COUSIN BIOTECH is committed to making every effort to produce and supply a high-quality medical device. However, if a health professional (client, user, prescriber, etc.) has a complaint or cause for dissatisfaction with a product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In case of failure of an implant or if it has contributed to causing a serious adverse reaction in the patient, the health centre must follow the legal procedures in that country and inform COUSIN BIOTECH immediately. For any correspondence, please specify the reference, batch number, contact person's details and a comprehensive description of the incident or complaint.

Brochures, documentation and surgical technique are available on request from COUSIN BIOTECH and its distributors. If further information is needed or required, please contact your COUSIN BIOTECH representative or distributor.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

IMPLANT CARD INFORMATION LEAFLET

How to complete the implant card ?

The implant card, for the patient, should be completed as follows:

Front side of the implant card

The user has to complete the implant card before transmitting the implant card to the patient. Please complete it as follows:

	Carte Implant - Implant Card	Carte Implant - Implant Card
	A conserver / Keep it	
Patient name	 _____	 COUSIN SURGERY <i>We care for Surgery</i>
Date of surgical procedure	 . . / . . /	
Name of practitioner	 _____	

Back side of the implant card

No fields need to be completed by the user. The following information on the device is printed directly by the manufacturer on the back of the implant card. A QR code allowing the patient to consult the instruction for use, as well as the website address, are available on the implant card.

Information to identify the device	Traceability information
• Indication that the product is a medical device • Description of the device family • Trade name of the device • Reference of the device • Device is MR safe • Denomination name of the device in sales countries	• UDI number • Batch number • Name and address of the legal manufacturer

Explanation of the used symbols in the Implant Card

All symbols are commonly used and explained at the end of the instruction for use.

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Patient information leaflet

BIOMESH® CA.B.S'.Air® Composite

Item	Patient information leaflet
1	a. The name of the device : BIOMESH® CA.B.S'.Air® Composite b. The model of the device : CABST1R05T / CABS AIRC7 / CABS AIRC09
2	a. The intended purpose of the device : BIOMESH® CA.B.S'.Air® Composite is intended for the repair and parietal reinforcement of ventral, umbilical and incisional hernias requiring reinforcement with non-absorbable support material. b. The kind of patient on whom the device is intended to be used : Adult with ventral hernia including umbilical and incisional hernia necessitating a surgical repair with a method that includes an intraperitoneal implantation.
3	Special operating instructions for the use of the device : <i>Installation of the patient :</i> Patient lying on his back Both legs kept in the body's axis Each arm placed on an arm wrester in a 90° angle from the body's axis <i>1-Incision :</i> The cutaneous incision (2 to 4cm according to the size of the defect) is centred on the hernia. <i>2-Dissection and resection of the peritoneal sac :</i> The hernial sac is identified and its content is re-introduced back inside the abdominal cavity. The sac is then left opened in order to allow the Intra-peritoneal prosthesis to be put in place. Next step consists in liberating epiploico-peritoneal adhesions around the edges of the hernia so that the prosthesis can fully expand. <i>3-Positioning and fixating the prosthesis :</i> BIOMESH® CA.B.S'.AIR® COMPOSITE is introduced through the hernial orifice after being rolled up into a tube shape or folded in an umbrella-like way. The balloon is filled up in order to allow the expansion of the prosthesis and the placement of the prosthesis is checked so that it remains centred above the hernial orifice. The prosthesis is fixed with the two pre-mounted threads that will transfix the peritoneum and the muscular wall. Because of its Intraperitoneal placement, two additional points should be used to ensure secure fixation of the prosthesis and to avoid any epiploic tissue being caught in between the prosthesis and the wall. Be careful not to suture the balloon to the implant when fixing the prosthesis. The balloon is deflated and removed. <i>4-Closure of the incision :</i> Musculo-aponeurotic closure in front of the prosthesis with resorbable suture thread. Sub-cutaneous and then cutaneous closure, normally without drainage.

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Item	Patient information leaflet															
4	a. The intended performance of the device : The BIOMESH® CA.B.S.'AIR® Composite prosthesis has a non-adhesive side, which is designed to avoid or reduce visceral adherences during intra peritoneal use or extra peritoneal use in the presence of a peritoneal defect for ventral, umbilical and incisional hernias. Products, biocompatible and non resorbable, cause reactional fibrosis taking over the reinforcement after six months of implantation. They have the advantage of having a very high resistance to suture, have a great flexibility and allow a optimum and rapid integration and colonisation. If necessary, the surgeon has the possibility of adding two additional points of fixing, according to his appreciation. b. Undesirable side effects : Like any implantable medical device, this implant is susceptible to generate possible undesirable effects: <table><tr><td>- Discomfort/Pain</td><td>- Recurrence</td><td>- Adhesion formation</td></tr><tr><td>- Erosions</td><td>- Irritation nearby organ</td><td>- Infection</td></tr><tr><td>- Inflammation</td><td>- Fistula formation</td><td>- Seroma</td></tr><tr><td>- Mesh deformation</td><td>- Hematoma</td><td>- Mesh migration</td></tr><tr><td>- Allergic reaction</td><td>- Foreign body reaction</td><td>- Obstruction</td></tr></table>	- Discomfort/Pain	- Recurrence	- Adhesion formation	- Erosions	- Irritation nearby organ	- Infection	- Inflammation	- Fistula formation	- Seroma	- Mesh deformation	- Hematoma	- Mesh migration	- Allergic reaction	- Foreign body reaction	- Obstruction
- Discomfort/Pain	- Recurrence	- Adhesion formation														
- Erosions	- Irritation nearby organ	- Infection														
- Inflammation	- Fistula formation	- Seroma														
- Mesh deformation	- Hematoma	- Mesh migration														
- Allergic reaction	- Foreign body reaction	- Obstruction														
5	The potential risks that could arise if the instructions and the precautions for use are not met are (for example): loss of product sterility risk of infection loss of effectiveness recurrence															
6	MRI Information: The implant has not been tested for safety and compatibility in the case of postoperative investigation using Magnetic Resonance Imaging (MRI). However, the implant is mainly made with polypropylene, polyethylene terephthalate and non-metallic component. Consequently, BIOMESH® CA.B.S.'AIR® Composite is considered MR safe.															
7	Precautions : The device is intended to be permanently implanted in the body. No maintenance is therefore possible. If the patient experiences any symptoms that appear abnormal, the patient is invited to return to the surgeon for further consultation. An implant card for the patient is provided by the healthcare professional at the clinic. This implant card provides information, for the patient, to identify the device and traceability elements as well as the name, address and website of the manufacturer. The removal of the implant is up to the surgeon's decision and must be subject to adequate postoperative follow-up. Some risks are associated with the explantation. The surgeon who decides to remove the device should consider factors such as the risk of another surgery for the patient and the difficulty of the explantation procedure.															
8	Materials included in the device : The implanted materials are Polypropylene, Polyethylene Terephthalate (PET) and ePTFE for anti-adherent side. The biological safety of the materials and any manufacturing residuals has been demonstrated through chemical characterization and biological tests.															
9	Reporting a serious incident : If you encounter any serious incident with your BIOMESH® CA.B.S'.Air® Composite implant you should report it immediately to your doctor. You should also report the incident to Cousin Biotech via the link serviceclients@cousin-surgery.com or materiovigilance@cousinbiotech.com and to the Therapeutic Goods Administration using the following website: https://www.tga.gov.au															
10	Legal Manufacturer : Cousin Biotech, Allée des Roses, 59117 Wervicq-Sud, France Manufacturing Facility : Cousin Biotech, Allée des Roses, 59117 Wervicq-Sud, France															

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Cousin Biotech SAS au capital de 320 000 € - NAF 1396Z- N° TVA FR 34 398 460 261 - Siren 398 460 261 RCS Lille