# ADJUSTABLE ACL CORTICAL ATTACHMENT

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**NOT331\_251210**

Version du 10/12/2025

Caution: Federal law (USA) restricts this device to sale, distribution and use by or on the order of a physician.

This release is the last update of the instructions of use and replace the previous edition

**Product under regulation 2017/745 (EU).**

**ADJUSTABLE ACL CORTICAL ATTACHMENT  
STERILE CORTICAL FIXATION DEVICES  
SINGLE USE STERILE PRODUCT**

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**IDENTIFICATION AND DESCRIPTION OF THE DEVICE**

Cortical Fixation System for tendons and ligaments.

The adjustable fixation system comprises an adjustable loop for graft suspension and a titanium button for the cortical support.

The cortical fixation devices are intended to be used to anchor tendons and ligaments during Anterior Cruciate Ligament (ACL) reconstruction.

**INDICATIONS**

The Cortical Fixation Systems are used for the instability of the knee due to a traumatic rupture or tear of the anterior cruciate ligament (ACL).

**CONTRAINDICATIONS**

- Known hypersensitivity to the implant material. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Insufficient quantity or quality of bone.
- Blood supply and previous infections which may tend to retard healing.
- Active infection.
- Conditions which tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period.
- Not to be used for the bone – patella tendon – bone technique
- Do not implant in children while they are growing.

**UNDESIRABLE SIDE EFFECTS**

Possible complications linked to the implant, that may lead to :

- Surgery failure, instability, and reintervention related to
  - o Breaking or elongation of the loop,
  - o Migration, breaking, upturn or collapse of the button
- Reaction to a foreign body
- Allergic reaction or intolerance to a component of the implant
- Local and transient inflammatory reactions
- Infection

Possible complications linked to all the ACL surgery using the implant :

- Pain in the incision area or surgical site
- Discomfort
- Infection
- Hematoma
- Wound dehiscence
- Dynamic valgus of the knee causing pain and instability
- Cyclops syndrome limiting extension of knee and leading to additional surgery
- Rupture of the ACL graft leading to surgery failure, instability, and reintervention
- Meniscus lesion
- Osteolysis, localised demineralisation of the bone

**TARGET POPULATION**

Adult and minor (with mature skeleton) undergoing an ACL reconstruction procedure.

**IMPLANTED MATERIALS**

- **Materials implanted:**
  - o Cortical button: Non-ferromagnetic titanium alloy (Ti TA6V ELI)
  - o Adjustable loop: 100% polyethylene (UHMWPE)
  - o Extended plate for cortical fixation System (optional, not included in the packaging): Non-ferromagnetic titanium alloy (Ti TA6V ELI)

Neither animal nor human origin. Non-absorbable

- **Materials in contact with the patient during the operation:**
  - o Pulling and flipping braids: Polyethylene terephthalate (PET) + biocompatible colorant D&C green#6
  - o Silicone marking

Reference	Weight of PE (g)	Weight of Titanium TA6V4 compliant with ISO 5832-3 (g)
OCBGEFRB1U	0,2	0,19

## **PERFORMANCES**

This system facilitates ACL repair/reconstruction by graft fixation in the bone. It has been approved for long-term implantation in the human body.

The main benefit of the device is a contribution to the benefit of ACL surgery and depends directly on the global success of the surgery: regaining knee stability after a rupture of the ACL ligament, which allows to resume work and/ or sport activities.

## **LINK TO SUMMARY OF SAFETY AND PERFORMANCE CHARACTERISTICS:**

[https://www.cousin-biotech.com/d-61951ae4e7fce\\_page-en-construction.pdf](https://www.cousin-biotech.com/d-61951ae4e7fce_page-en-construction.pdf)

## **MRI SAFETY**

The implants are composed of lightly ferromagnetic materials (TA6V4 ELI compliant with ISO 5832-3) and non-ferromagnetic materials (PE, PET, silicone) and present a geometry not likely to generate induced current. Moreover, as they are fixed to tissues, they are unlikely to be mobilized. Furthermore, as devices are implanted more than 30mm apart, they are unlikely to generate heat. Their safety, in particular in terms of heating and migration of implant, has been evaluated through bibliographic data and by comparison with data available on devices with similar composition, shape and use tested for MRI compatibility. This evaluation concluded a conditional compatibility for MRI of 1,5 Tesla to 3 Tesla.

As a precautionary measure, it is recommended MRI scans be avoided within 48 hours of the implant placement, and to inform the person in charge of the scan of the recent implant placement if such examination is essential.

It should be noted that devices presenting a high contrast with the biological environment can generate "artifacts" that must be taken into account for the perfect execution and interpretation of imaging exams.

For this purpose, patients implanted with this device must be advised to warn the health professionals concerned (radiologists and radiology operators) to the best of their ability about the presence of this implant before these exams.

## **PRECAUTIONS FOR USE**

Before use, check the integrity of the implant and the packaging.

Do not use if the implant, labels and/or protective packaging is damaged.

Do not use the product after the expiry date indicated on the label. This could compromise the performance of the device and the patient's safety.

Unpacking and handling of the implant must be carried out in accordance with asepsis standards, to ensure sterility.

When opening the cardboard box, the surgeon must check that the batch number and model match to those indicated on the label on the inside cover. If they do not, use of the implant is strictly forbidden.

This implant can be used with standard instruments available on the market and specified for orthopaedic surgery for ligament or tendon anchorage.

Surgical instruments are subject to wear during normal use. Instruments that have been used many times or been subject to very high stress are likely to break. Instruments must only be used for their intended purpose. COUSIN BIOTECH recommends that instruments should be regularly inspected for wear and bending.

<b>DO NOT REUSE</b>	<b>DO NOT RESTERILIZE</b>
---------------------	---------------------------

The implant is for single use, as specified on the label. It must not on any account be reused or re-sterilized (potential risks include, but are not limited to: loss of sterility of the product, risk of infection, loss of effectiveness of the product).

COUSIN BIOTECH does not provide any guarantee or recommendations regarding the use of a specific make of fixation device. The properties of these devices are subject to alteration by the manufacturer, over which COUSIN BIOTECH has no control.

The surgeon shall ensure that the bone tunnel has the correct length and diameter.

This device must **only** be implanted by a qualified surgeon trained in the use of the product. Avoid using staples or sharp-edged screws.

## **INSTRUCTIONS TO BE GIVEN TO THE PATIENTS BY THE SURGEON**

The surgeon should inform the patient of the potential physical and psychological restrictions and consequences of implanting the device. The patient must be informed of the surgical risks and possible side effects. The patient must be informed of the limitations inherent in the nature and function of the implant, in particular regarding the weight and constraints placed by the body on the device before the bone has sufficiently healed. The surgeon should invite the patient to return for further consultation if the patient presents symptoms that appear abnormal. An implant card for the patient is provided by the healthcare professional at the clinic. This implant card provides information, for the patient, to identify the device and traceability elements as well as the name, address and website of the manufacturer. The surgeon shall invite the patient to scan the implant card immediately upon receipt to keep track of it in case of loss. The surgeon shall also inform the patient

that the summary of the safety and performance characteristics of the device and the instructions for use can be found on the COUSIN BIOTECH website.

### **STORAGE PRECAUTIONS**

Store in a dry place away from sunlight, at room temperature and in the original packaging.

### **LINK TO SURGICAL TECHNIQUE**

The surgical technique is available in electronic format on the website:

[https://www.cousin-biotech.com/d-61951ae4e7fce\\_page-en-construction.pdf](https://www.cousin-biotech.com/d-61951ae4e7fce_page-en-construction.pdf)

### **SURGICAL TECHNIQUE**

#### **Preliminary information**

The system is adapted for the use of the right internal and semitendinosus tendons with the all-inside and inside-out techniques as well as retrodrilling.

In case of retrodrilling, COUSIN BIOTECH recommends a tunnel diameter of 5 mm to insert the implant. In case of a full tunnel, it is preferable to use the additional extended plate.

Use of the device with a tunnel diameter of less than 5 mm has not been validated. If the tunnel diameter is greater than 5 mm, the surgeon can decide whether or not to use the extended plate.

#### **Technique**

1. Suture the graft and put it in the loop.
2. Create a bone tunnel to insert the implant with a diameter of 5 mm (recommended diameter) and a graft cavity suited to its diameter.
3. Insert the pulling and flipping braids and adjustment braids through the bone tunnel.
4. Pull on one of the pulling and flipping braids until the button comes out of the bone tunnel through the cortical bone whilst maintaining the other pulling and flipping braid and adjustment braids.
5. Pull on the second pulling and flipping braid in order to tip over the button and allow cortical contact.
6. Validate fixation by pulling the graft in the opposite direction.
7. While still applying counter tension to the graft, slide the adjustment braids one after another to adjust the loop to the required length, starting with the white braid. Do not apply tension to the black marked braid.
8. To lock the loop, continue to apply tension to the white adjustment braid and pull firmly on the black marked adjustment braid.
9. Fix the graft in the opposite bone tunnel using the required method.
10. If necessary, slide the adjustment braids alternately once again to obtain the required tension. And repeat step 8.
11. Remove the pulling and flipping braid(s) from the button.
12. Section the adjustment braids at a minimum distance of 10 mm from the cortical button.

#### **Use of the extended plate (not included in the packaging)**

It may be necessary to use an extended plate to increase the contact surface of the button on the cortical bone (full tibial or femoral tunnel).

Repeat steps 1 to 5 above.

6. Remove one of the two pulling and flipping braids.
7. Put the button on the extended plate. Ensure that the remaining suture is on the split side to avoid hindering the interlocking.
8. While still applying counter tension to the graft, slide the adjustment braids one after another to adjust the loop to the required length, starting with the white braid. Do not apply tension to the black braid.
9. To lock the loop, continue to apply tension to the white adjustment braid and pull firmly on the black adjustment braid.
10. Fix the graft in the opposite bone tunnel using the required method.
11. If necessary, slide the adjustment braids alternately once again to obtain the required tension. Repeat step 9.
12. Remove the remaining pulling and flipping braid from the button.
13. Section the adjustment braids at a minimum distance of 10 mm from the cortical button.

### **EXPLANTATION AND DISPOSAL OF DEVICES**

Risks associated with the explantation: the surgeon who decides to remove the device should consider factors such as the risk of another surgery for the patient, and the difficulty of the explantation procedure. The removal of the implant is up to the surgeon's decision and must be subject to adequate postoperative follow-up.

Devices should be retrieved and handled according to the recommendations of ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling". Any explanted device must be sent back for analysis, following the current protocol. This protocol is available on request from COUSIN BIOTECH. It is important to note that any implant that must not be cleaned or disinfected prior to despatch must be contained in sealed packaging. The removed medical device must be disposed of in accordance with the country's standards for the disposal of infectious waste. There are no specific recommendations for the disposal of a non-implanted device.

Risks associated with the explantation: the surgeon who decides to remove the device should consider factors such as the risk of another surgery for the patient, and the difficulty of the explantation procedure. The removal of the implant is up to the surgeon's decision and must be subject to adequate postoperative follow-up.

## IMPORTANT

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

## INFORMATION REQUEST AND COMPLAINTS

Following its quality policy, COUSIN BIOTECH is committed to making every effort to produce and supply a high-quality medical device. However, if a health professional (client, user, prescriber, etc.) has a complaint or cause for dissatisfaction with a product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In case of failure of an implant or if it has contributed to causing a serious adverse reaction in the patient, the health centre must follow the legal procedures in that country and inform COUSIN BIOTECH immediately. For any correspondence, please specify the reference, batch number, contact person's details and a comprehensive description of the incident or complaint. Brochures, documentation and surgical technique are available on request from COUSIN BIOTECH and its distributors. If further information is needed or required, please contact your COUSIN BIOTECH representative or distributor.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## IMPLANT CARD INFORMATION LEAFLET

## How to complete the implant card ?

The implant card, for the patient, should be completed as follows:

### Front side of the implant card

The user has to complete the implant card before transmitting the implant card to the patient. Please complete it as follows:

IMPLANT CARD

KEEP IT

Patient name

?

Date of surgical procedure

31

.

.

/

.

.

/

.

.

.

.

Name of practitioner

+

IMPLANT CARD

### Back side of the implant card

No fields need to be completed by the user. The following information on the device is printed directly by the manufacturer on the back of the implant card. A QR code allowing the patient to consult the instruction for use, as well as the website address, are available on the implant card.

Information to identify the device	Traceability information
<ul style="list-style-type: none"> <li>• Indication that the product is a medical device</li> <li>• Description of the device family</li> <li>• Trade name of the device</li> <li>• Reference of the device</li> <li>• Device is MR safe</li> <li>• Denomination name of the device in sales countries</li> </ul>	<ul style="list-style-type: none"> <li>• UDI number</li> <li>• Batch number</li> <li>• Name and address of the legal manufacturer</li> </ul>

### Explanation of the used symbols in the Implant Card

All symbols are commonly used and explained at the end of the instruction for use.

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