SUMMARY OF SAFETY AND CLINICAL PERFORMANCE of Sterile Intraperitoneal Non Resorbable Parietal Reinforcement Implants

ACTION	WRITTEN BY	VERIFIED AND APPROVED BY
NAME	Regulatory Affairs Specialist	Regulatory Affairs Director
NAME	O.MOREAUX	F.PELLETIER
DATE	03/12/2025	03/12/2025
VISA		

en	<u>SUMMARY OF SAFETY AND CLINICAL PERFORMANCE of Sterile Intraperitoneal Non Resorbable Parietal Reinforcement Implants</u>
de	<u>KURZBERICHT ÜBER SICHERHEIT UND KLINISCHE LEISTUNG steriler intraperitonealer nicht resorbierbarer parietaler Verstärkungsimplantate</u>
it	<u>SINTESI RELATIVA ALLA SICUREZZA E ALLA PRESTAZIONE CLINICA degli impianti di rinforzo parietale intraperitoneali sterili non riassorbibili</u>
es	<u>RESUMEN DE SEGURIDAD Y RENDIMIENTO CLÍNICO de Implantes de refuerzo parietal intraperitoneales no reabsorbibles estériles</u>
nl	<u>SAMENVATTING VAN DE VEILIGHEID EN KLINISCHE PRESTATIES van steriele intraperitoneale, niet-resorbeerbare pariëtale verstevigingsimplantaten</u>

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Abbreviations List

CD	Crohn's disease
CER	Clinical Evaluation Report
EC	European Commission
EU	European Union
ePTFE	Extracted polytetrafluoroethylene
IFU	Instructions For Use
MDR	Medical Device Regulation
PMCF	Post-Market Clinical Follow-up
PMS	Post-Market Surveillance
PP	Polypropylene
PSUR	Periodic Safety Update Report
QoL	Quality of Life
SRN	Single Registration Number
SSCP	Summary of Safety and Clinical Performance
UDI	Unique Device Identification
VL	Ventral laparoscopy

List of tables

Table 1: INTRAMESH® T1/ PROMESH® SURG INTRA and BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB and their commercial references	17
Table 2: Classification of INTRAMESH® T1/ PROMESH® SURG INTRA and BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB	19
Table 3: Models and sizes of products per family	21
Table 4: Residual risks after mitigation	23
Table 5: Cumulative field safety corrective actions	24
Table 6: Summary of clinical data obtained for INTRAMESH® T1/ PROMESH® SURG INTRA and BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB.....	25

Table of Contents

SECTION A: SSCP FOR USERS/HEALTHCARE PROFESSIONALS	17
1. IDENTIFICATION OF THE DEVICE AND GENERAL INFORMATION	17
1.1. DEVICE TRADE NAME(S)	17
1.2. MANUFACTURER'S NAME AND ADDRESS	18
1.3. MANUFACTURER'S SRN	18
1.4. BASIC UDI-DI	18
1.5. MEDICAL DEVICE NOMENCLATURE DESCRIPTION / TEXT	18
1.6. CLASS OF DEVICE	18
1.7. YEAR WHEN THE FIRST CERTIFICATE WAS ISSUED COVERING THE DEVICE	19
1.8. AUTHORISED REPRESENTATIVE IF APPLICABLE, NAME AND THE SRN	20
1.9. NOTIFIED BODY'S NAME AND SRN	20
2. INTENDED PURPOSE OF THE DEVICE	20
2.1. INTENDED PURPOSE	20
2.2. INDICATIONS AND TARGET POPULATION(S)	20
2.3. CONTRAINDICATIONS AND/LIMITATIONS	20
3. DEVICE DESCRIPTION	21
3.1. DESCRIPTION OF THE DEVICE	21
3.2. A REFERENCE TO PREVIOUS GENERATIONS OR VARIANTS IF SUCH EXIST AND A DESCRIPTION OF THE DIFFERENCES	22
3.3. DESCRIPTION OF ANY ACCESSORIES WHICH ARE INTENDED TO BE USED IN COMBINATION WITH THE DEVICE	22
3.4. DESCRIPTION OF ANY OTHER DEVICES and products WHICH ARE INTENDED TO BE USED IN COMBINATION WITH THE DEVICE	23
4. RISKS AND WARNINGS	23
4.1. RESIDUAL RISKS AND UNDESIRABLE EFFECTS	23
4.2. WARNINGS AND PRECAUTIONS	24
4.3. OTHER ASPECTS OF SAFETY INCLUDING A SUMMARY OF ANY FIELD SAFETY CORRECTIVE ACTION IF APPLICABLE	24
5. SUMMARY OF CLINICAL EVALUATION AND INFORMATION ON POST-MARKET CLINICAL FOLLOW-UP	25
5.1. SUMMARY OF OTHER CLINICAL DATA related to equivalent device, IF APPLICABLE	25
5.2. SUMMARY OF CLINICAL PERFORMANCE AND SAFETY	25
5.2.1. SUMMARY OF CLINICAL DATA FROM CONDUCTED INVESTIGATIONS OF THE DEVICE BEFORE THE CE-MARKING, IF APPLICABLE	25
5.2.2. SUMMARY OF CLINICAL DATA FROM OTHER SOURCES, IF APPLICABLE	25
5.2.3. CONCLUSIONS ON CLINICAL PERFORMANCE AND SAFETY REQUIREMENTS	33

5.3.	ONGOING POST MARKET CLINICAL FOLLOW-UP	34
6.	POSSIBLE DIAGNOSTIC OR THERAPEUTIC ALTERNATIVES	34
7.	SUGGESTED PROFILE AND TRAINING FOR USERS	35
8.	REFERENCE TO ANY HARMONISED STANDARDS AND COMMON SPECIFICATIONS APPLIED	35
9.	REVISION HISTORY	36
SECTION B: SSCP FOR PATIENTS.....		38
1.	DEVICE IDENTIFICATION AND GENERAL INFORMATION.....	38
1.1.	DEVICE AND TRADE NAME.....	38
1.2.	MANUFACTURER: NAME AND ADDRESS	38
1.3.	BASIC UDI-DI	38
1.4.	YEAR WHEN THE DEVICE WAS FIRST CE-MARKED.....	38
2.	INTENDED USE OF THE DEVICE.....	38
2.1.	INTENDED PURPOSE	38
2.2.	INDICATIONS AND INTENDED PATIENT GROUPS	39
2.3.	CONTRAINDICATIONS.....	39
3.	DEVICE DESCRIPTION.....	39
3.1.	DEVICE DESCRIPTION AND MATERIAL/ SUBSTANCES IN CONTACT WITH PATIENT TISSUES	39
3.2.	INFORMATION ABOUT MEDICINAL SUBSTANCES IN THE DEVICE, IF ANY.....	40
3.3.	DESCRIPTION OF HOW THE DEVICE IS ACHIEVING ITS INTENDED MODE OF ACTION.....	40
3.4.	DESCRIPTION OF ACCESSORIES, IF ANY.....	40
4.	RISKS AND WARNINGS.....	40
4.1.	HOW POTENTIAL RISKS HAVE BEEN CONTROLLED OR MANAGED	40
4.2.	REMAINING RISKS AND UNDESIRABLE EFFECTS.....	40
4.3.	WARNINGS AND PRECAUTIONS	40
4.4.	SUMMARY OF ANY FIELD SAFETY CORRECTIVE ACTION (FSCA INCLUDING FSN), IF APPLICABLE	41
5.	SUMMARY OF CLINICAL EVALUATION AND POST-MARKET CLINICAL FOLLOW-UP	41
5.1.	CLINICAL BACKGROUND OF THE DEVICE	41
5.2.	the clinical evidence for the ce-marking	42
5.3.	SAFETY.....	42
6.	POSSIBLE DIAGNOSTIC OR THERAPEUTIC ALTERNATIVES	42
6.1.	GENERAL DESCRIPTION OF THERAPEUTIC ALTERNATIVES	43
7.	SUGGESTED TRAINING FOR USERS	43
8.	REFERENCE TO ANY HARMONISED STANDARDS AND COMMON SPECIFICATIONS APPLIED	43
ABSCHNITT A: SSCP FÜR ANWENDER / MEDIZINISCHE FACHKRÄFTE.....		51
9.	IDENTIFIZIERUNG DES PRODUKTS UND ALLGEMEINE INFORMATIONEN	51

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INTRODUCTION

This Summary of Safety and Clinical Performance (SSCP) is produced in alignment with the requirements of COUSIN BIOTECH's Clinical Evaluation Report (TF03_CER_ISSUE-2) and REGULATION (EU) 2017/745 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL, referencing the MDCG 2019 – 9 - Summary of safety and clinical performance: A guide for manufacturers and notified bodies.

This SSCP is intended to provide public access to an updated summary of the main aspects of the safety and clinical performance of the device. The SSCP is not intended to replace the Instructions For Use (IFU) as the main document to ensure the safe use of the device, nor is it intended to provide diagnostic or therapeutic suggestions to intended users or patients.

This SSCP is an important source of information for intended users – both healthcare professionals and if relevant for patients. It is one of several means intended to fulfil the objectives of the Medical Device Regulation (MDR) to enhance transparency and provide adequate access to information. The SSCP for Sterile Intraperitoneal Non Resorbable Reinforcement Implants composed of INTRAMESH® T1/ PROMESH® SURG INTRA and BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB includes both information intended for healthcare professionals (Section A) as well as for the patient (Section B).

The SSCP for INTRAMESH® T1/ PROMESH® SURG INTRA and BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB is based on the following documentation pertaining to the INTRAMESH® T1/ PROMESH® SURG INTRA and BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB devices:

- Clinical Evaluation Report (CER)
- Post-Marketing Surveillance (PMS) plan and regular reports
- Post-Marketing Clinical Follow-up (PMCF) plan and regular reports
- Design verification/validation reports
- If appropriate, the IFU

The SSCP shall be reviewed on an annual basis and updated if needed to ensure that any clinical and/or safety information in the SSCP remains correct and complete after implementation of the PMCF reports and of the Periodic Safety Update Report (PSUR). Any pertinent modification in the documentation pertaining to INTRAMESH® T1/ PROMESH® SURG INTRA and BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB shall be implemented in all sections of the SSCP if needed.

The SSCP of INTRAMESH® T1/ PROMESH® SURG INTRA and BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB must contain all the information detailed in this template. However, the manufacturer may add further information from the Technical Documentation (TD) of the device to enhance the comprehension of the mandatory information if it does not affect the readability of the SSCP and if it excludes any element of a promotional nature.

For the sake of objectivity and impartiality to the patients, information in the SSCP shall adequately summarize both favorable and unfavorable data.

The IFU shall provide all relevant information needed to directly access the SSCP in Eudamed. The IFU shall be updated regularly on this database. The following applies to the IFU:

- It shall state that the SSCP is available in the European database on medical devices (Eudamed), where it is linked to the Basic UDI-DI.
- It should provide the URL to the Eudamed public website: <https://ec.europa.eu/tools/eudamed>
- It should state the value of the Basic UDI-DI. Alternatively, another metadata can be stated provided it can be used to unambiguously search and find the intended SSCP in Eudamed.

The MDR requires that the SSCP shall be written in a way that is clear to the intended user and, if relevant, to the patient. The SSCP should therefore be translated into the languages accepted in the Member States where the device is sold/envisaged to be sold. By analogy, this is also a requirement for an IFU.

SECTION A: SSCP FOR USERS/HEALTHCARE PROFESSIONALS

The following information is intended for users/healthcare professionals. A summary intended for patients can be found in Section B.

1. IDENTIFICATION OF THE DEVICE AND GENERAL INFORMATION

INTRAMESH® T1/ PROMESH® SURG INTRA implants are intended to be used for the repair and parietal reinforcement of ventral, umbilical and incisional hernias requiring reinforcement with non-absorbable support material.

BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB implants are intended for repair and parietal reinforcement of ventral, umbilical and incisional hernias requiring reinforcement with non-absorbable support material.

1.1.DEVICE TRADE NAME(S)

Intraperitoneal non resorbable parietal reinforcement implants include two sub-families, the *INTRAMESH® T1/ PROMESH® SURG INTRA* and *BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB* implants.

The references *PROMESH® SURG* are strictly identical to those initially CE marked (*INTRAMESH® T1* and *BIOMESH® CA.B.S.'Air® Composite*). They are only different in terms of labels, brand name, external design of the cardboard box and they have their own IFU. The intended purpose and performances, materials, suppliers, manufacturing processes, packaging of the devices and sterilization method are identical. All information that refers to COUSIN BIOTECH product references is applicable to *PROMESH® SURG* devices.

The devices trade names are *INTRAMESH® T1/ PROMESH® SURG INTRA* and *BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB* and the associated commercial references, which differ in terms of model and size, are presented below in Table 1:

Table 1: *INTRAMESH® T1/ PROMESH® SURG INTRA* and *BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB* and their commercial references

FAMILY	REFERENCE
<u>INTRAMESH® T1 / PROMESH® SURG INTRA</u>	<u>FBIOT11015/ PSIN1015RS</u>
	<u>FBIOT11515/ PSIN1515SQ</u>
	<u>FBIOT11520/ PSIN1520RS</u>
	<u>FBIOT12020/ PSIN2020SQ</u>
	<u>FBIOT12025/ PSIN2025RS</u>
	<u>FBIOT13030/</u>

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Note: *PROMESH® SURG INTRA* is identical to *INTRAMESH® T1* and *PROMESH® SURG UMB* is identical to *BIOMESH® C.A.B.S.'Air® Composite*.

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1.3. MANUFACTURER'S SRN

The Single Registration Number (SRN) of the manufacturer is FR-MF-000001179.

1.4. BASIC UDI-DI

INTRAMESH® T1/PROMESH® SURG INTRA : 376018557INTRAMESHT1HQ
BIOMESH® CA.B.S.'Air Composite / PROMESH® SURG UMB : 376018557CABSAIRCOMPOKG

The present SSCP encompasses the two following sub-families (**INTRAMESH® T1** and **CA.B.S.'Air Composite** devices) for the reasons below:

- They all are Class III medical devices
- They all are non-resorbable
- Their mechanism of action is the same
- They have the same indications
- The design relies on the same technological platform
- The Benefit/Risk analysis is detailed in the CER by sub-family

1.5. MEDICAL DEVICE NOMENCLATURE DESCRIPTION / TEXT

Not applicable yet.

1.6. CLASS OF DEVICE

COUSIN BIOTECH PROPERTY - REPRODUCTION FORBIDDEN
CONFIDENTIAL

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The principal characteristics of *INTRAMESH® T1/ PROMESH® SURG INTRA* and *BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB* are listed below in Table 2:

Table 2: Classification of *INTRAMESH® T1/ PROMESH® SURG INTRA* and *BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB*

CATEGORY	CLASSIFICATION
Duration	Permanent/ Long-term
Invasive device	Yes
Body contact	Implant (Interaction with the abdominal wall)
Biological effect	No
Usage	Single use
Sterile	The devices are sterile (ethylene oxide sterilization)*
Active medical device	No
Software	No
Medicinal substance	No medicinal products are incorporated in the devices
Biological material	No material derived from animal tissue or human blood is incorporated in the devices**
Radioactive device	The devices are not radioactive

According to Rule 8 of Annex VIII of the EU MDR 2017/745, *INTRAMESH® T1/ PROMESH® SURG INTRA* and *BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB* devices are considered to be Class III medical devices.

* The implants are designed for long-term implantation. They are sterilized by ethylene oxide and designed for single-use.

** The implanted materials of *INTRAMESH® T1*, consisting of polypropylene (PP) and polytetrafluoroethylene (ePTFE) of the anti-adherent side and the dimethyl siloxane-biocompatible black ink, are considered non resorbable and do not contain products derived from human or animal origin.

The implanted materials of *BIOMESH® CA.B.S.'Air® Composite* are PP, polyethylene terephthalate, ePTFE and a biocompatible colorant (D and C Green #6). The materials that are in contact with the patient during the surgery are stainless steel and poly vinyl chloride and they are non-absorbable. The duration of contact of these materials during the surgery is less than 60 minutes. No material of human or animal origin is in contact with the patient.

1.7. YEAR WHEN THE FIRST CERTIFICATE WAS ISSUED COVERING THE DEVICE

The CE marking certificate number of *INTRAMESH® T1* and *BIOMESH® CA.B.S.'Air® Composite* and *PROMESH® SURG* was FR19/81843462 under the MDD 92/42/CEE. The *INTRAMESH® T1* devices are CE marked since 2002, the *BIOMESH® CA.B.S.'Air® Composite* devices are CE marked since 2004 and the *PROMESH® SURG* devices were first marked in March 2020.

The CE marking number is MDR 00019-A for the EU Quality Management System Certificate and MDR 00019-B for the EU Technical Documentation Assessment Certificate, obtained on 2025-05-09 in accordance with Regulation (EU) 2017/745.

1.8. AUTHORISED REPRESENTATIVE IF APPLICABLE, NAME AND THE SRN

The authorized representative of this company is COUSIN BIOTECH.

The Single Registration Number (SRN) of this authorized representative is N° FR-MF-000001179.

1.9. NOTIFIED BODY'S NAME AND SRN

The Notified Body (NB) for the *INTRAMESH® T1/ PROMESH® SURG INTRA* and *BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB* and the validation of this SSCP is the KIWA.

Its single identification number is CE 0476.

2. INTENDED PURPOSE OF THE DEVICE

2.1. INTENDED PURPOSE

INTRAMESH® T1/ PROMESH® SURG INTRA and *BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB* implants are intended to be used for the repair and parietal reinforcement of ventral, umbilical and incisional hernias requiring reinforcement with non-absorbable support material.

2.2. INDICATIONS AND TARGET POPULATION(S)

Indications for *INTRAMESH® T1/ PROMESH® SURG INTRA* implants:

This device is indicated for the repair and parietal reinforcement of ventral, umbilical and incisional hernias requiring reinforcement with non-absorbable support material.

Indications for *BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB* implants:

This balloon-deployed device is indicated for the repair and parietal reinforcement of ventral, umbilical and incisional hernias requiring reinforcement with non-absorbable support material.

Patient populations:

Adult with ventral** hernia including umbilical and incisional hernia necessitating a surgical repair with a method that includes an intraperitoneal implantation.

***Ventral hernia is caused by weakness of the anterior abdominal wall due to midline defects, previous incisions, or increased intra-abdominal pressure. Ventral hernias typically include umbilical hernia, incisional (post surgery) hernias, and hernias at the junction of the abdominal muscles.*

2.3. CONTRAINDICATIONS AND/LIMITATIONS

Contra indications:

- Allergy to any of the components
- Infection at the implantation site
- Pregnancy
- Growing children
- Patients undergoing anticoagulant therapy


3. DEVICE DESCRIPTION

3.1. DESCRIPTION OF THE DEVICE

Intraperitoneal non resorbable parietal reinforcement implants include two sub-families, the *INTRAMESH® T1/ PROMESH® SURG INTRA* and *BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB* implants.

The references *PROMESH® SURG* are strictly identical to those initially CE marked (*INTRAMESH® T1* and *BIOMESH® CA.B.S.'Air® Composite*). They are only different in terms of labels, brand name, external design of the cardboard box and they have their own IFU. The intended purpose and performances, materials, suppliers, manufacturing processes, packaging of the devices and sterilization method are identical. All information that refers to COUSIN BIOTECH product references is applicable to *PROMESH® SURG* devices. The main components and models with their corresponding sizes and presentations are presented below in Table 3.

Table 3: Models and sizes of products per family

FAMILY	REFERENCE	DESCRIPTION AND SIZE (mm)	PRODUCT PRESENTATION
INTRAMESH® T1 / <u>PROMESH® SURG INTRA</u>	FBIOT11015/ PSIN1015RS	Mesh 10 x 15	
	FBIOT11515/ PSIN1515SQ	Mesh 15 x 15	
	FBIOT11520/ PSIN1520RS	Mesh 15 x 20	
	FBIOT12020/ PSIN2020SQ	Mesh 20 x 20	
	FBIOT12025/ PSIN2025RS	Mesh 20 x 25	
	FBIOT13030/ PSIN3030SQ	Mesh 30 x 30	
	FBIOT13050/ PSIN3050RT	Mesh 30 x 50	
	FBIOT1R12X/ PSIN12DIAM	Mesh Ø 12	
BIOMESH® CA.B.S.'Air® Composite/ <u>PROMESH® SURG UMB</u>	CABST1R05T/ PSUMBIN05D	1 Mesh Ø 5 and 2 sutures threads	
	CABSAIRC7/ PSUMBIN07D	Kit with mesh Ø 7, inflatable balloon, syringe and 2 sutures	
	CABSAIRC09/ PSUMBIN09D	Kit with mesh Ø 9 inflatable balloon, syringe and 4 sutures	

INTRAMESH® T1 weighs 120 ± 10 g/m² and has two sides. Side A is made up of a smooth microperforated ePTFE material to minimize adhesions, while side B is composed of a macroporous PP matrix for fast tissue ingrowth. These implants can be used in laparoscopic surgery or open surgery. They have a non-adhesive

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side designed to avoid visceral adhesions during intra peritoneal use or extra peritoneal use in the presence of a peritoneal defect. The marking in black ink is present on the non-adhesive side. The meshes have a shape memory and are flexible, biocompatible and non-resorbable and are considered to be thin, biocompatible and non-resorbable.

Technique: During the surgery, the defect is measured, and the relevant mesh size is selected by the surgeon to ensure enough overlap. The mesh should overlap the defect(s) 3–5 cm on all sides. For larger defects, the overlap should be at least 5 cm. The mesh is rolled up ePTFE side inwards and introduced into the abdomen. The correct orientation of *INTRAMESH® T1* is essential. The white anti-adherent side (ePTFE) is designed to minimize bowel adhesions to the prosthesis. This side should be placed towards the bowels. The landmarks on the ePTFE side allow correct centering of the mesh over the defect. The PP side should face the abdominal wall. The PP mesh surface should never be placed against bowel or other visceral structures. *INTRAMESH® T1* should always be adequately fixated. Fixation devices or nonabsorbable monofilament sutures are recommended to secure the prosthesis. The mesh must lay flat.

BIOMESH® CA.B.S.'Air® Composite weighs $326.9 \pm 2 \text{ g/m}^2$ and is made up of ePTFE and PP. The device comprises a parietal prosthesis with two or four fixation threads crimped with stainless steel needles and a balloon permitting the deployment of the parietal prosthesis (except for CABST1R05T). The mesh can be deployed by the balloon system which is present inside the device and is inflated with air. The implant, which is a ready-to-use kit, can be quickly and safely fixed by the pre-attached sutures crimped with needles. At the end of the procedure, the balloon must be imperatively withdrawn.

The aim of this expansion balloon concept is to facilitate dissection and secure the unfolding of the prosthesis. This balloon is located inside the plug; it is filled-up with 30cm³ of air. It acts as a pneumatic jack ensuring complete crease-free unfolding of the prosthesis. Its double-sided design combines a smooth ePTFE side, designed to avoid visceral adhesions, with a knitted monofilament PP side, which allows better settlement of the peritoneum.

Technique: The cutaneous incision (2 to 4 cm according to the size of the defect) is centered on the hernia. The hernial sac is identified and its content is re-introduced back inside the abdominal cavity. The sac is then left open in order to allow the intra-peritoneal prosthesis to be settled in place. The epiploic-peritoneal adhesions around the hernia edges are then liberated so that the prosthesis can fully expand. Next, *BIOMESH® CA.B.S.'Air® Composite* is introduced through the hernial orifice after being rolled up into a tube shape or folded in an umbrella-like way. The balloon is filled up to permit the prosthesis expansion and the placement of the prosthesis is checked so that it stays centered above the hernial orifice. The prosthesis is fixed with the two pre-mounted threads that will transfix the peritoneum and the muscular wall. Because of its intraperitoneal placement, two additional points should be used to ensure secure fixation of the prosthesis and avoid any epiploic tissue being caught in between the prosthesis and the wall. The balloon is deflated and removed. At the end, musculo-aponeurotic closure in front of the prosthesis is performed with resorbable suture thread. Sub-cutaneous and then cutaneous closure is done, normally without drainage.

3.2.A REFERENCE TO PREVIOUS GENERATIONS OR VARIANTS IF SUCH EXIST AND A DESCRIPTION OF THE DIFFERENCES

Not applicable, as there are no previous versions of the devices.

3.3.DESCRPTION OF ANY ACCESSORIES WHICH ARE INTENDED TO BE USED IN COMBINATION WITH THE DEVICE

No specific accessories are intended to be used with *INTRAMESH® T1/ PROMESH® SURG INTRA* and *BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB* implants.

3.4. DESCRIPTION OF ANY OTHER DEVICES AND PRODUCTS WHICH ARE INTENDED TO BE USED IN COMBINATION WITH THE DEVICE

No other devices are intended to be used in combination with *INTRAMESH® T1/ PROMESH® SURG INTRA* and *BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB* implants in addition to generic surgical instruments.

4. RISKS AND WARNINGS

4.1. RESIDUAL RISKS AND UNDESIRABLE EFFECTS

In the risk analysis report, risks have been classified into different categories:

- Risks related to device properties and characteristics (Raw material, energy, packaging type, sterilization, sterility maintaining)
- Risks related to biological, chemical and non-biocompatibility properties (Material characterization, leachable, *in vivo* condition and microbiological safety)
- Risks related to formation (Labels, IFU, cautions)
- Risks related to device use (Implantation)
- Risks related to usability (Transport, storage condition, user/device interface, human factors)
- Risks related to implanted device functioning (Failure, maintenance, aging, expected early or late adverse events)

Table 4: Residual risks after mitigation

Risks	Before reduction	After reduction
High risk	5	0
Medium risk	57	16
Weak risk	1	47

Based on the output of the risk analysis, 63 risks have been identified after mitigation from each risk category and each phase. High risks have been reduced significantly and the residual risks are considered acceptable with respect to the effect. After risk assessment and control, no high risks were identified. The medium risks have been reduced as much as possible. After reduction, the remaining sixteen medium risks have been evaluated individually and are deemed acceptable in terms of patient safety. Overall, 47 low residual risks have been evaluated and considered as acceptable by COUSIN BIOTECH.

The 63 risks identified manifest a favorable risk/benefit analysis.

The rate of adverse events related to *INTRAMESH® T1/ PROMESH® SURG INTRA* and *BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB* implants has been demonstrated as acceptable in view of the benefits if the devices are used according to the instructions for use. Paying attention to the sterility of the implants can majorly avoid complications.

4.2. WARNINGS AND PRECAUTIONS

Warnings about undesirable side effects:

Like any implantable medical device, this implant is susceptible to generate undesirable side effects:

- Discomfort/ pain
- Infection
- Recurrence
- Inflammation
- Adhesion formation
- Fistula formation
- Erosion
- Seroma
- Hematoma
- Mesh migration
- Mesh deformation
- Allergic reaction
- Foreign body reaction
- Irritation of nearby organs
- Obstruction

Precautions required by the manufacturer:

- Prostheses are delivered sterile (ethylene oxide sterilization).
- The implants should not be reused and/ or resterilized as potential risks related to such acts include loss of product sterility, risk of infection, loss of effectiveness, recurrence and others.
- Before use, inspect the integrity of the packaging and the device (especially the blister / peelable pouches). Do not use in the event of deterioration of the device and/or the packaging.
- Do not use if the device is out of date.
- Cousin BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned.
- This device must be implanted only by a qualified surgeon and trained with the use of the product (knowledge of anatomy and visceral surgery).

4.3. OTHER ASPECTS OF SAFETY INCLUDING A SUMMARY OF ANY FIELD SAFETY CORRECTIVE ACTION IF APPLICABLE

Table 5: Cumulative field safety corrective actions

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The bowel adhesions are risks that can occur after meshes implant and are mentioned in the IFU. They cannot be avoided sometimes when a foreign object such as prosthesis is inserted intraperitoneally. Products with defects are not recommended to be used as indicated in the precautions. These risks are mentioned in the risk analysis and considered as acceptable after risk mitigation.

These field safety corrective actions have been integrated in the PMS activities and reports. No new field action has been implemented since 2008.

5. SUMMARY OF CLINICAL EVALUATION AND INFORMATION ON POST-MARKET CLINICAL FOLLOW-UP

5.1. SUMMARY OF OTHER CLINICAL DATA RELATED TO EQUIVALENT DEVICE, IF APPLICABLE

Not applicable.

5.2. SUMMARY OF CLINICAL PERFORMANCE AND SAFETY

5.2.1. SUMMARY OF CLINICAL DATA FROM CONDUCTED INVESTIGATIONS OF THE DEVICE BEFORE THE CE-MARKING, IF APPLICABLE

Not applicable.

5.2.2. SUMMARY OF CLINICAL DATA FROM OTHER SOURCES, IF APPLICABLE

In order to assess the performance and safety of the *INTRAMESH® T1/ PROMESH® SURG INTRA* and *BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB* implants, analysis of clinical data from the literature, clinical data from the manufacturer, PMS and other clinical data was performed. In total, five articles, three PMCF studies, one investigator initiated study, one users survey report, one preclinical study as PMS activities were included.

5.2.2.1. CLINICAL DATA

Table 6: Summary of clinical data obtained for *INTRAMESH® T1/ PROMESH® SURG INTRA* and *BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB*

DATA SOURCE	PRODUCT	DESCRIPTION (Authors, population...)	MAIN RESULTS-ADDRESS THE PERFORMANCE AND/OR THE SAFETY
Clinical data from the manufacturer			
PMCF study: Cohort follow-up Bressolette, 2004 The study was conducted from	<i>INTRAMESH® T1</i>	A prospective observational cohort « Suivi de cohorte de patient ayant eu une cure d'événtration	The medical and surgical conditions of <i>INTRAMESH® T1</i> use were highlighted. The procedure was performed under general anesthesia and the prosthesis was placed intraperitoneally (80% of the patients were operated on by laparoscopy through the trans-peritoneal route, 12% by laparotomy and 8% had a conversion). A clinical follow-up was performed at 1, 2 and 6 months.

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25/05/2002 to 19/10/2004.		traitée par INTRAMESH® T1 » N=51 patients Gender : 57% males and 43% females; Age: mean: 56 years	Performance: The device was used for parietal repair and reinforcement. The success of the surgery was evaluated at 97% as one case of recurrence was declared at 18 months. Safety: In total, 37% patients had postoperative complications at 1 month, 21% at 2 months and 29% at 6 months. It has been reported that about 50% of patients did not face any complications. In total, 30% of patients suffered from moderate pain. This complication was declared immediately after the procedure. Overall, it was indicated that 12% declared seromas and 3% had a hematoma. Only 3% of patients had recurrence at 18 months after the operation. The safety conclusions were: - No major perioperative and postoperative complications were witnessed. - No unexpected adverse events or deficiencies were reported. - Hernia recurrence rate was 3% (at 18 months). - Short and midterm (18 months) safety. The study demonstrated accepted efficacy and safety of INTRAMESH® T1.
PMCF study: Cohort follow-up Friis Andersen, 2006 The study was conducted from 01/03/2006 to 31/08/2006	INTRAMESH® T1	A prospective observational cohort « Suivi de cohorte de patients ayant subis une cure d'événtration ou d'hernie traitée par l'implant de réfection de paroi INTRAMESH T1 » N= 41 patients Gender: 54% males and 46% females; Age: median: 50 years, range: 28-87 years	In addition to studying the medical and surgical conditions of the implant use, all events related to the operation and the device were demonstrated in this cohort. The procedure was performed by ventral laparoscopy (VL) during an average period of 70 minutes. Patients were hospitalized during two days and followed up at 1, 3, 12 and 24 months. Performance: Parietal repair was attained during the process. The success of the surgery was considered to be 100% in the short-term at 3 months with absence of recurrence. A suspected case of recurrence was reported at 5 months. Safety: A short term safety was documented. It has been indicated that no unexpected adverse events or deficiencies have been noted. The minor complications reported were three infections, four seromas and one perforation of the intestines. Major complications were: an ileus, an exploratory laparotomy without result and a hematoma.

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PMCF study : Hernia Club, 2021 The study was conducted from 08/01/2020 to June 2021.	BIOMESH® CA.B.S.'Air® Composite	A prospective, longitudinal and multicentric register “Efficacy and safety of the CA.B.S.'air® Composite for umbilical, incisional hernia repair. Interim report of data collected to date in Club Hernie registry for CA.B.S.'air® Composite” N= 152 patients Gender; 85 males (56%) and 67 females (44%); Age (mean ± SD) : (58.5 ± 13.1 years)	The study demonstrated acceptable results in terms of safety and efficacy but only in the short-term. This study was designed to assess the performance of the device and ensure its maintenance over time in terms of the patient's quality of life and recurrence. Furthermore, the intraoperative, postoperative complications and reinterventions were observed to evaluate the safety. In total, 152 patients with clinically diagnosed ventral hernia who underwent surgery using the Sterile Intraperitoneal Non-Resorbable Parietal Reinforcement Implants BIOMESH® CA.B.S.'Air® Composite . Mean follow-up was realized at 36.1 ± 14.9 months (minimum 13 months, maximum 69 months) Follow-up has been realized for 100% of patients, with at least one long term visit. Regarding surgery background of the included patients, 71.1% have primary ventral hernias whereas 28.9% have secondary ventral hernia (incisional hernia). 98% had open surgery with most of the intraperitoneal approach (80%). Among them, 3 patients (60%) had COO conversion (open to open conversion) and 2 patients (40%) had CLO conversion (conversion from laparoscopic to open surgery). Mesh was fixed with separate non-absorbable stitches in 82.9% of cases. <u>Performance and safety results:</u> Analysis of pain intensity demonstrated that patients' general pain was improved at one month after surgery (85.2%) and during follow-up after surgery. The rate of painful patients goes from 73.7% preoperatively to 14.8% at D30 Seven patients (4.6%) had intraoperative complications. None of patients experienced post operative complications. These complications are not major and expected. No recurrence was observed to date. The results of the medical devices studied are consistent with published results for ventral hernia cures, with a very low recurrence rate and a surgical outcome, evaluated by the patients themselves (PROM) objectively excellent (61%) or good (20%). This present study highlights good safety and performance after BIOMESH® CA.B.S.'Air®
--	---	--	---

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			Composite implantation for the treatment of ventral hernia. Given the low rate of per operative or postoperative complications, the absence of major complications BIOMESH® CA.B.S.'Air® Composite – Club Hernie - Interim Report and recurrence hernia to date, the benefit-risk ratio seems favorable. Note: This is the second interim report (17th March 2023) The first report was signed on 14th October 2021 This study is ongoing and shall be terminated in 2026 after the 5-year follow-up of the last patient included in the study.
Investigator initiated study: LeBlanc, 2014 The study was conducted from May 2012 to September 2013 and a follow-up until 3 months after surgery.	BIOMESH® CA.B.S.'Air® Composite	A prospective, observational, non-comparative, non-randomized study “Study report on performance of Bio-mesh® CA.B.S Air Composite Mesh” N= 22 patients Gender : 95% males and 5% females	The performance of the device and its effectiveness in reducing complications when used on optimal implantation techniques were evaluated. Patients were then followed-up for Quality of life (QoL) at 1 week, 1 and 3 months after surgery. Performance: Parietal repair was attained during the process. The success of the surgery was evaluated at 100% in the short-term. The QoL of all patients was improved. In total, 22 (100%) patients recovered prior to hospital discharge. Other benefits that were addressed were patients satisfaction and the ease of use of the implant for surgeons. In total, 17 surgeons rated the implant use as excellent. Furthermore, 60% of patients had medications prescribed for pain at one week after surgery, 28% at 1 month and 11% at 3 months. Safety : Very limited data on safety was presented in this study. It was indicated that only few and minor postoperative complications were documented. A dense and vascular adhesion at the hernia site was observed. No major perioperative complications from the device were reported. Note: This study was not completed due to administrative reasons and not safety reasons. The results obtained were limited to umbilical hernia in men.

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Clinical data from the literature

A. Hussain, H. Mahmood, J. Nicholls, S. El-Hasani. Int J Surg. 2008 Feb;6(1):15-9. The study was conducted from February 2006 to July 2007.	INTRAMESH® T1	An observational prospective study “Laparoscopic ventral hernia repair. Our experience of 61 consecutive series: Prospective study” N=61 patients Gender: 32 males (52.45%) and 29 females (47.54%); Age: mean: 53.42 range: 39-80 years	ePTFE mesh was inserted in 58 patients. Overall, the patient classification was 42% incisional hernias following laparotomy (17 midline incisions, 4 grid iron incisions and 5 gynecological operations), 39.34% umbilical and paraumbilical hernias, 6.55% spigelian hernias, 6.55% epigastric hernias and 4.9% port site incisional hernias. The mean follow-up period was 8 months. Performance: Parietal repair was accomplished and the success of the surgery was calculated at 98% in the short-term as one recurrence case (1.6%) was reported at 8 months. Safety: Post-operative: Ten complications were reported, resulting in an overall complication rate of 16.39%: - Two infections (3.2%) took place, one infected mesh and one infected seroma (It is indicated that seromas are commonly reported after hernia operations) - Five patients (8%) presented severe post-operative pain, which was related to the technique. Follow-up at 8 months: The postoperative pain continued in 8% of the population. Lung consolidation and pleura effusion was a complication of laparoscopy that had been noted. The general morbidity rate was 16% in the short-term. No major complications from the device were reported. The laparoscopic ventral hernia repair was effective.
Wen W, Majerus B, Van De Moortel M, Lobue S, Fobe D, Philippart P, Berwouts L, Coteur J, Gabriels K, Van der Speeten K. Acta Chir Belg. 2017 Apr 25:1-8.	INTRAMESH® T1	A multicentre prospective registry “Laparoscopic ventral hernia repair using a composite mesh with polypropylene and expanded polytetrafluoroethylene: a prospective, multicentre registry” N= 90 patients	For intraperitoneal use in primary and incisional ventral hernia repair, the safety and efficacy of INTRAMESH® T1 was evaluated. Overall, 61.1% of patients presented primary ventral hernias and 38.9% had incisional ventral hernias. The follow-up was performed at 1 and 12 months . Performance: The efficacy of the device was evaluated in terms of absence of recurrence. The surgery success was evaluated at 100% due to absence of recurrence (0%) at 12 months. Therefore, parietal repair was attained.

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The study was conducted between January 2013 and September 2014.		Gender: 60 males (66.6%) and 30 females (33.3%); Age: mean: 56.4 years, range: 25-80 years	Safety: The use of <i>INTRAMESH® T1</i> showed short and mid-term safety assessed at > 1 year. No recurrence occurred during the follow-up (median of 408 days). In total, 10 complications were reported. Among the complications that were reported at 12 months: - 4.4% were considered major (1 bowel obstruction and 1 re-laparoscopy) - 21% were considered to be minor (seroma, local numbness, pain and swelling). <i>INTRAMESH® T1</i> demonstrated acceptable efficiency with favourable short and midterm outcomes and morbidity.
Bensaadi H, Paolino L, Valenti A, Polliand C, Barrat C, Champault G. Am Surg. 2014 Jan;80(1):57-65 The study was conducted between January 2007 and August 2011.	<i>BIOMESH® CA.B.S.'Air® Composite</i> compared to Ventralex Hernia Patch	A prospective, comparative, randomized and controlled study "Intraperitoneal Tension-free Repair of a Small Midline Ventral Abdominal Wall Hernia: Randomized Study with a Mean Follow-up of 3 Years" N= 42 patients (Information on gender and age: Not available)	The aim of the study was to determine the long-term recurrence and complication rates after small abdominal wall hernia repair with two different bilayer prostheses. Patients were divided into two groups, one group, in which <i>BIOMESH® CA.B.S.'Air® Composite</i> was used and the other, in which Ventralex Hernia Patch was used. Performance: When using <i>BIOMESH® CA.B.S.'Air® Composite</i>, the success of the surgery was reported at 100% at a mean follow-up of 3 years. This success percentage was based on absence of recurrence in the group where <i>BIOMESH® CA.B.S.'Air® Composite</i> was employed, while 11.7% patients had recurrence in the second arm group. Parietal repair was attained. <i>BIOMESH® CA.B.S.'Air® Composite</i> demonstrated better performance than the comparator. It was easier to work with as the associated balloon in the device facilitated superior deployment of the prosthesis allowing for good fixation with four sutures. It had been indicated that 94% of patients in the <i>BIOMESH® CA.B.S.'Air® Composite</i> group were satisfied compared to 73% in the comparator group. The quality of life in both groups, before and after (3-months) surgery was assessed using the SF-12 questionnaire. Safety: No reintervention due to late complications was required in the group where <i>BIOMESH® CA.B.S.'Air® Composite</i> was employed. Absence of complications was equally reported in this group. Less pain and less discomfort were observed in the Cabs'Air® group compared to the other group.

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			<i>BIOMESH® CA.B.S.'Air® Composite</i> use showed better performance and safety compared to the comparator.
GETAID, Beyer-Berjot L, Moszkowicz D, Bridoux V, Schneider L, Theuil L, et al. World J Surg. 2020 Jul;44(7):2394–400 The study was conducted between January 2000 and February 2018.	<i>BIOMESH® CA.B.S.'Air® Composite</i>	A multicenter retrospective comparative study “Mesh Repair in Crohn’s Disease: A Case-Matched Multicenter Study in 234 Patients” N= Mixed data (112 Crohn 120 Control) (Information on gender and age: Not available)	All consecutive Crohn’s Disease (CD) patients that underwent operations and mesh repair for ventral hernia (primary or incisional) between January 2000 and December 2017 were retrospectively included. The main outcomes of the study focused on the impact of CD on the septic morbidity after mesh repair. Performance: Parietal repair was observed. The success of the surgery/ procedure was valued at 86% for patients with CD and 91% for patients in the control group. This success percentage was based on recurrence rate which was 14% versus 8.3% for CD and control groups respectively. Safety: The process represented medium-term safety. CD patients manifested a higher risk of postoperative septic morbidity (18.4% vs. 5%, $p = 0.001$), entero-prosthetic fistula (7% vs. 0, $p < 0.01$) and mesh withdrawals (5.3% vs. 0, $p = 0.011$). CD appears to be a risk factor for septic morbidity after mesh repair

Other datasets of clinical interest

Preclinical study: Schreinemacher 2013	<i>INTRAMESH® T1</i>	“Coated meshes for hernia repair provide comparable intraperitoneal adhesion prevention.” N= NA 36 rats	The objective of this study was to examine the efficacy of several coated meshes for adhesion reduction (adhesions, incorporation and tissue reaction). Five commercially available meshes with a layered coating were placed intraperitoneally in rats and followed up for 90 days: polypropylene and polyester meshes, both coated with absorbable collagen (Parietene Composite and Parietex Composite, respectively), and three polypropylene meshes respectively coated with absorbable omega-3 fatty acids (C-Qur Edge), absorbable cellulose (Sepramesh IP), and nonabsorbable expanded polytetrafluoroethylene (Intramesh T1). Uncoated polypropylene and collagen meshes (Parietene and Permacol, respectively) served as the control condition Performance: Limited data on integration. Safety:
---	----------------------	---	--

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Ching SS, Sarela AI, Dexter SPL, Hayden JD, McMahon MJ. Surg Endosc. 2008 Oct;22(10):2244–50. The study was conducted between December 2002 and August 2007.	INTRAMESH® T1 Used for 21% of patients	A retrospective comparative study “Comparison of early outcomes for laparoscopic ventral hernia repair between nonobese and morbidly obese patient populations.” N= 168 patients (Intramesh T1 estimated around 37 Patients) Gender: 87 males (52.4%) and 79 females (47.5%); Age: median: 55 years, range: 24 – 92 years	Reduction of adhesions: The efficacy in terms of adhesion reduction was similar in all coated meshes The study included patients in which several meshes were used. The PP plus ePTFE, composite mesh (INTRAMESH® T1) was used for 21% of patients, the Proceed mesh (PP + polydioxanone polymer film + oxidized regenerated cellulose) was used for 13% of patients, the Parietex composite mesh (multifilament polyethylene terephthalate with a collagen-based hydrogel component) was used for 55% of patients and polyethylene terephthalate with dimethyl siloxane was used for 7% of patients. Note: This study has limitations as the article provides pooled information without distinction in product types. Performance: The success of the surgery was evaluated at 88% in the long-term as it was associated with 12% of recurrence at 19 months. The incidence of recurrence was significantly correlated with the size of the fascial defect and the size of the mesh. Safety: In total, 9% of patients had postoperative infections, 10% chronic pain, 9% seroma, 5% hematoma, 4% enterotomy, 2% chronic wound sinus and 2% enterocutaneous fistula. Infective complications occurred in 7% of patients. No significant difference in the incidence of perioperative complications or recurrence after LVHR was observed between the morbidly obese patients and the non– morbidly obese patients.
Users survey: Cousin Biotech Survey 2021	INTRAMESH® T1	« PMCF Survey Report STERILE INTRAPERITONEAL NON-RESORBABLE PARIETAL REINFORCEMENT IMPLANT / INTRAMESH T1 » N= NA	Performance: Contribution to general performance. Safety: Contribution to general safety. In the meantime, this survey has no final results on the safety and performance of the devices.

5.2.2.2. PMS DATA

In total, COUSIN BIOTECH has received 93 complaints between January 2006 and March 2020. Among the 93 customer complaints, 76 were considered as non-serious incidents and were not reported to competent

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authorities between January 2006 and December 2021. The number of materiovigilance events is very low compared to the number of devices sold since 2006 (0.011%).

These incidents were then classified into different typologies. Balloon defect, device breakage/ tearing, adhesion/ obstruction and device breakage/ needle were the four typologies that could be analysed.

Concerning materiovigilance events, 17 cases were reported between January 2006 and December 2021. Among the cases, 7 were related to balloon defect, 5 to adhesion/obstruction, 2 to device breakage, 1 to assembly, 1 to internal packaging and 1 to migration of prosthesis in the small intestine. It was indicated that the number of materiovigilance events was very low compared to number of sold devices since 2006 (0.014%).

Furthermore, 10 events, related to *BIOMESH® CA.B.S.'Air® Composite*, with no consequences on patients have been notified through the materiovigilance system. Overall, 7 were linked to balloon expansion, 2 problems related to threads use and 1 de-insertion of the biface at the interface. Among the other incidents were 5 cases of adhesion causing pain or intestinal occlusion/ obstruction.

Adhesion cases were evaluated and considered not to be related to the device. They are well-known and expected undesirable side effects. They cannot be avoided when a foreign object such as prosthesis is inserted intraperitoneally.

The complaint ratio for the devices was calculated:

BIOMESH® CA.B.S.'Air® Composite: $15 / 94\,372 = 1.5 / 10\,000$

INTRAMESH® T1: $2 / 51\,701 = 0.38 / 10\,000$

Based on the low frequency of complaints, the safety profile of the implants is confirmed.

5.2.3. CONCLUSIONS ON CLINICAL PERFORMANCE AND SAFETY REQUIREMENTS

Based on all the clinical and PMS data demonstrated previously, the following conclusions can be made:

INTRAMESH® T1/ PROMESH® SURG INTRA implants were studied in 280 patients with primary and incisional ventral, umbilical, paraumbilical, spigelian, epigastric, hypochondriac, subcostal, iliac, through the linea alba and port site incisional hernias.

BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB implants were studied in 216 patients with ventral including primary, incisional, parietal, umbilical, through the linea alba, epigastric and small incisional hernias.

Taking into consideration the safety data collected from the literature, PMCF studies, PMS and other datasets, one study (Wen, 2017) exposed the occurrence of two major events (1 bowel obstruction and 1 re-laparoscopy) related to the use of *INTRAMESH® T1/ PROMESH® SURG INTRA*. No new severe adverse events related to *BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB* implants have been identified. Generally, the frequency of hernia recurrence ranged between 1 and 10%. Taking into account the total number of patients treated, the safety profile of the devices was positive. No new complications, other than those identified in the technical documentation of the manufacturer, were reported.

The adverse events that were mostly reported were expected and mentioned in the IFU. Some of them were predictable as they are related to the surgery and not directly to the device. Considering the lack of new complications presented on *INTRAMESH® T1/ PROMESH® SURG INTRA* and *BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB* implants, it may be concluded that *INTRAMESH® T1/ PROMESH® SURG INTRA* and *BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB* implants have a good safety profile in all tested indications for a follow-up mean of 12 months.

Implant sterility, device biocompatibility, limitation of bowel adhesions and reduction of visceral adherences during intraperitoneal use in presence of a peritoneal defect are the safety claims made by COUSIN BIOTECH.

Globally, the findings are considered to be consistent with the safety claimed by COUSIN BIOTECH.

Furthermore, the results obtained demonstrate the performance of *INTRAMESH® T1/ PROMESH® SURG INTRA* and *BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB* implants.

Positive outcomes were observed in terms of parietal repair and reinforcement in addition to hernia treatment, manifested through surgery success in all indicated types of hernia repair. Based on percentage of recurrence, the surgery success was evaluated between 86 and 100%.

Performance data for *INTRAMESH® T1/ PROMESH® SURG INTRA* and *BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB* implants are coherent with those issued from the technical documentation of the products.

The available clinical data supports the use of *INTRAMESH® T1/ PROMESH® SURG INTRA* in inguinal, umbilical, femoral and ventral hernias repair and parietal reinforcement. The use of *BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB* implants in all forms of hernia repair, which requires reinforcement with non-absorbable support material, is similarly exposed through the clinical data.

The performance claims made by COUSIN BIOTECH are parietal repair, parietal reinforcement and absence of resorption of meshes.

The findings are generally considered to be consistent with the performances claimed by COUSIN SURGERY. The associated benefit-risk ratio seems favourable if the products are used following the recommendations of use.

5.3. ONGOING POST MARKET CLINICAL FOLLOW-UP

One PMCF study (Hernia Club, 2021) is still ongoing and will be terminated by 2026 after the 5-year follow-up of the last patient included in the study.

The available clinical data shows that there is no need for additional PMCF studies for the *INTRAMESH® T1/ PROMESH® SURG INTRA* and *BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB* implants. This consideration is based on the historic use of sterile intraperitoneal non-resorbable parietal reinforcement implants, the conformance to all applicable harmonized standards, the safety profile and the acceptable benefit/ risk ratio demonstrated in the clinical evaluation report on the devices.

All actions taken are described in the manufacturers PMCF plan and are evaluated in an annual PMCF report, the conclusions of which are included in the clinical evaluation report and the present SSCP.

6. POSSIBLE DIAGNOSTIC OR THERAPEUTIC ALTERNATIVES

Knowing that the surgery is indicated only when conservative measures have failed and watchful waiting is not appropriate, the available alternatives in terms of intervention are all surgical and associated with a mesh choice that possesses:

- Comparable surgical risks and complications rates
- Comparable short term success rates, none of them having demonstrated any superiority

- Comparable recurrence and reintervention rates, none of them having demonstrated any superiority.

7. SUGGESTED PROFILE AND TRAINING FOR USERS

The devices should be implanted only by qualified surgeons who have knowledge of abdominal and visceral anatomy and surgery and are familiar with the device.

8. REFERENCE TO ANY HARMONISED STANDARDS AND COMMON SPECIFICATIONS APPLIED

Hereafter, the list of the standards claimed for the Sterile Intraperitoneal Non Resorbable Parietal Reinforcement Implants, showing the standards harmonised according to Regulation 2017/745 and the state of the art of standards.

Harmonized standards following Regulation (EU) 2017/745 are identified with “*” in the present document. No common specifications have been identified for the moment.

Harmonised Standard Title	Date
*EN ISO 13485:2016/A11:2021 - Medical Devices - Quality Management Systems - Requirements for Regulatory purposes	March 2016/A11:September 2021
NF EN ISO 9001:2015 - Quality Management Systems - Requirements	October 2015
NF EN ISO/IEC 17025:2017 - General requirements for the competence to carry out tests and/or calibrations, including sampling. It covers testing and calibration performed using standard methods, nonstandard methods, and laboratory-developed methods	December 2017
NF EN ISO 14630:2013 - Non-active surgical implants - General requirements	February 2013
*EN ISO 15223-2021 - Medical devices – Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements	September 2021
*EN ISO 14971: 2019 - Medical devices – Application of risk management to medical devices	December 2019
NF EN 62366-1:2015 Amd 1:2020 - Medical Devices - Application of usability engineering to medical devices	August 2020
EN 556-1:2001/AC:2006 - Sterilization of medical devices. Requirements for medical devices to be designated "STERILE" – Part 1: Requirements for terminally sterilized medical devices	November 2006
*EN ISO 11135-1:2014 (+A1:2019) - Sterilization of health care products – Ethylene oxide – Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices	July 2014/ A1: November 2019
*EN ISO 11737-1:2018/A1 2021 - Sterilization of health care products — Microbiological methods — Part 1: Determination of a population of microorganisms on products	January 2018/ A1:June 2021
*EN ISO 11737- 2:2020 - Sterilization of medical devices – Microbiological methods - Part 2: Tests of sterility performed in the definition, validation and maintenance of a sterilization process	May 2020
*EN ISO 11607-1:2020 - Packaging for terminally sterilized medical devices – Part 1: Requirements for materials, sterile barrier systems and packaging systems	January 2020
*EN ISO 11607-2:2020 - Packaging for terminally sterilized medical devices – Part 2: Validation requirements for forming, sealing and assembly processes	January 2020
EN 14937:2009 - Sterilization of health care products - General requirements for characterization of a sterilizing agent and the development, validation and routine control of a sterilization process for medical devices (ISO 14937:2009)	October 2009
NF EN ISO 10993-1:2020 Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management process	December 2020
NF EN ISO 10993- 3:2014 - Biological evaluation of medical devices – Part 3: Tests for genotoxicity, carcinogenicity and reproductive toxicity	December 2014

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NF EN ISO 10993- 5:2010 - Biological evaluation of medical devices – Part 5: Tests for in vitro cytotoxicity	July 2010
NF EN ISO 10993-6:2017 - Biological evaluation of medical devices – Part 6: Tests for local effects after implantation	April 2017
ISO 10993-7/A1:2019 - Biological evaluation of medical devices - Part 7: Ethylene oxide sterilization residuals	December 2019
ISO 10993-10:2021 - Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization	November 2021
NF EN ISO 10993-11 :2018 - Biological evaluation of medical devices – Part 11: Tests for systemic toxicity (ISO 10993-11:2017)	May 2018
*EN ISO 10993-12:2021- Biological evaluation of medical devices – Part 12: Sample preparation and reference materials	January 2021
NF EN ISO 10993-17:2009 - Biological evaluation of medical devices – Part 17: Establishment of allowable limits for leachable substances	December 2009
NF EN ISO 10993-18:2020 - Biological evaluation of medical devices. Part 18: Chemical Characterization of materials	May 2020
NF EN ISO 14155:2020 - Clinical investigation of medical devices for human subjects - Good clinical practice	August 2020
NF EN ISO 14644-1:2016 - Cleanrooms and associated controlled environments—Part 1: Classification of air cleanliness by particle	February 2016
NF EN ISO 14644-2:2016 - Cleanrooms and associated controlled environments— Part 2: Monitoring to provide evidence of cleanroom performance related to air cleanliness by particle concentration	February 2016
NF EN ISO 14644-3:2019 - Cleanrooms and associated controlled environments—Part 3: Tests methods	August 2019
NF EN 17141:2020 - Cleanrooms and associated controlled environments — Biocontamination control	August 2020
NF 868-5:2018 - Packaging for terminally sterilized medical devices – Part 5: Sealable pouches and reels of porous materials and plastic film construction - Requirements and test methods	December 2018
ISO/TR 15499:2016 - Biological evaluation of medical devices—Guidance on the conduct of biological evaluation within a risk management process	December 2016
ISO 12891-1:2015 - Retrieval and analysis of surgical implants - Part 1 : Retrieval and handling	July 2015
NF ISO 2859-1:2000/A1:2011 - Sampling procedures for inspection by attributes – Part 1: Sampling plans indexed by acceptable quality level (AQL) for lot-by-lot inspection	April 2000/A1: September 2011
ISO TR 20416:2020 - Medical devices — Post-market surveillance for manufacturers	July 2020
ISO 20417:2021 - Medical devices — Information to be supplied by the manufacturer	April 2021
ISO 24971:2021 - Medical devices — Guidance on the application of ISO 14971	April 2021
ISO 16142:2016 -Medical devices – Recognized essential principles of safety and performance of medical devices – Part 1: general essential principles and additional specific essential principles for all non-IVD medical devices and guidance on the selection of standards	March 2016

9. REVISION HISTORY

SSCP revision number	Date issued	Change description	Revision validated by the Notified Body
A	01/03/2022	Creation of the document for MDR submission	☐ Yes Validation language : ☐ No (only applicable for class IIa or some IIb implantable devices (MDR, Article 52 (4) 2 nd

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			paragraph) for which the SSCP is not yet validated by the NB
B	26/07/2022	Addition of harmonised standards	☐ Yes Validation language : ☐ No
C	22/12/2022	Deletion of NF EN 1041+A1:2013. Modification of the BASIC UDI : 375018557MESHINTRARJ	☐ Yes Validation language : ☐ No
D	30/01/2023	According to KIWA clinical review : - Update of the BASIC-UDI. 375018557INTRAMESHT1FT for INTRAMESH® T1/PROMESH® SURG INTRA and 375018557CABSAIRCOMPOHD for CA.B.S'Air® Composite and PROMESH® SURG UMB According to KIWA technical file review : - Update of the list of standards (deletion of NF EN 1041) and addition of revision dates.	☐ Yes Validation language : ☐ No
E	14/04/2023	Addition of a rationale for using different BASIC-UDI	☐ Yes Validation language : ☐ No
F	01/12/2023	Update according to the last interim report of Club Hernie Registry	☐ Yes Validation language : ☐ No
G	27/02/2024	Update according to the new headquarters address	☐ Yes Validation language : ☐ No
H	03/12/2025	Update following the correction of a numerical error identified in the company prefix of the Basic UDI-DI. Update following the certification of CE marking in compliance with Regulation (EU) 2017/745 on medical devices (MDR).	☐ Yes Validation language : ☐ No

SECTION B: SSCP FOR PATIENTS

The following information is intended for patients.

1. DEVICE IDENTIFICATION AND GENERAL INFORMATION

1.1. DEVICE AND TRADE NAME

INTRAMESH® T1/ PROMESH® SURG INTRA and *BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB* implants.

1.2. MANUFACTURER: NAME AND ADDRESS

Manufacturer's name : COUSIN BIOTECH
Address : Allée des Roses
59117 Wervicq-Sud- France

1.3. BASIC UDI-DI

INTRAMESH® T1/PROMESH® SURG INTRA : 376018557INTRAMESHT1HQ
BIOMESH® CA.B.S.'Air Composite / PROMESH® SURG UMB : 376018557CABS AIRCOMPOKG

The present Technical File encompasses the two following sub-families (**INTRAMESH® T1** and **CA.B.S.'Air Composite** devices) for the reasons below:

- They all are Class III medical devices
- They all are non-resorbable
- Their mechanism of action is the same
- They have the same indications
- The design relies on the same technological platform
- The Benefit/Risk analysis is detailed in the CER by sub-family

1.4. YEAR WHEN THE DEVICE WAS FIRST CE-MARKED

The *INTRAMESH® T1* was first CE-marked in 2002, the *BIOMESH® CA.B.S.'Air® Composite* was first marked in 2004 and the two types of *PROMESH® SURG* devices were marked recently in 2020.

2. INTENDED USE OF THE DEVICE

2.1. INTENDED PURPOSE

The implants are generally used in hernia repair as they support the damaged tissue during the healing process.

INTRAMESH® T1/ PROMESH® SURG INTRA implants are used for the repair and reinforcement of the wall of abdominal cavity (parietal structure) of ventral, umbilical and incisional hernias.

BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB implants are used for the repair and reinforcement of the wall of abdominal cavity (parietal structure) of ventral, umbilical and incisional hernias.

2.2. INDICATIONS AND INTENDED PATIENT GROUPS

INTRAMESH® T1/ PROMESH® SURG INTRA and *BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB implants* are used to repair and reinforce the wall of abdominal cavity (parietal structure) of ventral, umbilical and incisional hernias that require reinforcement with non-absorbable support material.

The implants are intended to be used in adults who suffer from a symptomatic ventral hernia that requires a surgical repair with a method that includes an implantation in the abdominal cavity (intraperitoneal implantation).

2.3. CONTRAINDICATIONS

INTRAMESH® T1/ PROMESH® SURG INTRA and *BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB implants* should not be employed in:

- Patients that have allergies to any of the components of the device
- Patients with infections at the implantation sites
- Pregnant patients
- Growing children
- Patients that have anticoagulant treatments.

3. DEVICE DESCRIPTION

3.1. DEVICE DESCRIPTION AND MATERIAL/ SUBSTANCES IN CONTACT WITH PATIENT TISSUES

The *INTRAMESH® T1/ PROMESH® SURG INTRA* and *BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB devices* belong to the group of intraperitoneal non-resorbable parietal reinforcement implants. Meshes are designed for parietal reinforcement and hernia repair by open or laparoscopic approach (insertion of short, narrow tubes (trocars) via small (less than one centimeter) incisions).

INTRAMESH® T1 comes in various shapes and sizes (in cm).

Side A is made up of a smooth microperforated ePTFE material to reduce adhesions (scar tissue that forms between abdominal tissues and organs causing tissues and organs to stick together), while side B is composed of a macroporous PP matrix for fast tissue ingrowth. The implanted materials are considered non-resorbable. They are not derived from human or animal origin. They have a shape memory and are very flexible.

BIOMESH® CA.B.S.'Air® Composite is mainly made up of ePTFE and PP. The device is made up of a parietal prosthesis with two or four fixation threads compressed with stainless steel needles and a balloon permitting the deployment of the parietal prosthesis. The implanted materials are polypropylene, polyethylene terephthalate, ePTFE and a biocompatible colorant (D and C Green #6). The materials that are in contact with the patient during the surgery are stainless steel and poly vinyl chloride and they are non-resorbable. No material of human or animal origin is in contact with the patient.

The implants are single-use devices that are intended to be used on one individual during a single procedure. They are designed for long-term implantation and sterilized by ethylene oxide.

3.2. INFORMATION ABOUT MEDICINAL SUBSTANCES IN THE DEVICE, IF ANY

The devices do not contain any medicinal substances.

3.3. DESCRIPTION OF HOW THE DEVICE IS ACHIEVING ITS INTENDED MODE OF ACTION

Protheses are employed as coverings that prevent the internal or visceral contents from passing through the hernia opening in the abdominal wall.

The implants are always fixed to the peritoneum, which is the membrane that covers the intra-abdominal organs, as they need to remain centered above the hernia opening. The fixation method is variable based on the surgeon's choice.

Both *INTRAMESH® T1/ PROMESH® SURG INTRA* and *BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB* prostheses have a non-adhesive side (non-tending to stick), which is designed to avoid visceral or internal adherences during intraperitoneal use, or in the presence of a problem in the peritoneum due to the implant in addition to precautions taken by the surgeon to avoid internal adhesions during the surgery.

3.4. DESCRIPTION OF ACCESSORIES, IF ANY

No specific accessories are intended to be used with *INTRAMESH® T1/ PROMESH® SURG INTRA* and *BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB* implants.

4. RISKS AND WARNINGS

Please contact your healthcare professional if you believe that you are experiencing side effects related to the device or its use or if you are concerned about risks. This document does not replace a consultation with your healthcare professional, if needed.

4.1. HOW POTENTIAL RISKS HAVE BEEN CONTROLLED OR MANAGED

The manufacturer has performed a risk assessment of the devices in which all the risks are identified. Subsequent actions are taken by the manufacturer to reduce any potential risk without impacting the benefit/risk ratio for the patient and his/her safety. The main techniques for risk reduction are presented in the documentation provided to the user (instructions for use, labels,...).

4.2. REMAINING RISKS AND UNDESIRABLE EFFECTS

The risk management report shows that the residual risks associated with *INTRAMESH® T1/ PROMESH® SURG INTRA* and *BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB* implants are acceptable for the patient's safety once all the measures to reduce the risk have been applied.

The rate of adverse events related to the use of *INTRAMESH® T1/ PROMESH® SURG INTRA* and *BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB* implants was rare and acceptable when devices are used according the instructions for use.

4.3. WARNINGS AND PRECAUTIONS

The warnings associated with the undesirable side effects generated by the implantable devices are:

- Discomfort/pain
- Infection
- Recurrence (Reappearance of hernia)
- Inflammation

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- Adhesion formation (Formation of scar tissue between abdominal tissues and organs causing tissues and organs to stick together)
- Fistula formation (Formation of abnormal connection between two body parts, due to surgery for example)
- Erosion (Removal or displacement of the meshes)
- Seroma (Collection of body fluids in the place where tissues were removed)
- Hematoma
- Mesh migration
- Mesh deformation
- Allergic reaction
- Foreign body reaction
- Irritation of nearby organs
- Obstruction (Abdominal blocking)

The precautions taken by the manufacturer are:

- The prostheses are sterilized by ethylene oxide and delivered sterilized.
- It is highly recommended not to reuse or resterilize the implants.
- It is forbidden to inspect the integrity of the packaging and the device before use.
- It is recommended not to use the device in case of device and/ or packaging damage.
- It is recommended not to use the device if the expiration date is over.
- The device must be implanted and handled only by a qualified surgeon who is trained to use the device and has knowledge in anatomy and visceral surgery. However, the type of implant and implant fixation means that should be used is left to the surgeon.

4.4. SUMMARY OF ANY FIELD SAFETY CORRECTIVE ACTION (FSCA INCLUDING FSN), IF APPLICABLE

Three field safety corrective actions have been reported since the device has been on the market. They concern bowel adhesions (Formation of scar tissue inside the abdomen) and product textile defects. Such adhesions can occur after meshes implant and cannot be avoided sometimes when a foreign object is inserted in the abdominal cavity.

Products with defects are not recommended to be used as indicated in the precautions.

These risks are mentioned in the risk analysis and considered acceptable after risk reduction. These actions were integrated in the post-market surveillance documentation of the manufacturer. No new field action has been implemented since 2015.

5. SUMMARY OF CLINICAL EVALUATION AND POST-MARKET CLINICAL FOLLOW-UP

5.1. CLINICAL BACKGROUND OF THE DEVICE

Abdominal wall hernias are a challenging health problem. In order to repair the hernia, conventional and traditional open suture techniques were previously applied. Such techniques were associated with high hernia recurrence rates. The introduction of meshes significantly reduced the hernia recurrence rate. Such implants aim to restore the organ layers.

They are designed with one surface, mostly made of polypropylene or polyester, that is in contact with the abdominal wall to promote tissue in-growth or healing. The other surface acts as a long-term anti-adhesive barrier for the intestine to avoid abdominal adhesences.

The *INTRAMESH® T1/ PROMESH® SURG INTRA* is thin, has a great flexibility when being manipulated and possess shape memory characteristics.

The *BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB* with its balloon and pre-attached sutures permits quick and safe fixation of the implant. It simplifies the running of surgical procedures and ensures deep tissue reconstruction.

5.2. THE CLINICAL EVIDENCE FOR THE CE-MARKING

Based on information retrieved from clinical studies and post-marketing surveillance, *INTRAMESH® T1/ PROMESH® SURG INTRA* and *BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB* implants were used in a total of 591 patients with all types of ventral hernia, primary of resulting from a surgical incision. Positive outcomes in terms of surgery success (86 to 100%, depending on the initial risk factors of the patient), hernia repair and low frequency of hernia recurrence (1 to 10%) were witnessed in all indications and types of hernia mentioned. Two major complications were documented when *INTRAMESH® T1* was used. These complications were abdominal blocking and re-laparoscopy (repetition of surgical intervention). Patients were generally satisfied when questioned during their postoperative follow-up visits.

This outcome highlights positive clinical evidence associated with absence of new serious adverse events or complications when employing *INTRAMESH® T1/ PROMESH® SURG INTRA* and *BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB* implants.

5.3. SAFETY

The adverse events found in the clinical and post-market surveillance data have been addressed in the instructions for use and are largely related to the surgical intervention rather than the device:

- Infection
- Recurrence
- Adhesion
- Seroma
- Discomfort
- Regular pain
- Hematoma
- Obstruction

These complications have been included in the manufacturer's risk management process and were mentioned in the technical documentation. They are considered acceptable under normal conditions of use when weighed against benefits to the patient.

Taking into account the total number of patients treated, the safety profile of the devices was positive. No new complications, other than those identified in the technical documentation of the manufacturer, were reported.

The benefit/risk ratio is favorable and acceptable for the use of *INTRAMESH® T1/ PROMESH® SURG INTRA* and *BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB* implants in hernia repair and parietal reinforcement.

6. POSSIBLE DIAGNOSTIC OR THERAPEUTIC ALTERNATIVES

When considering alternative treatments, it is recommended to contact your healthcare professional who can evaluate your individual situation.

6.1. GENERAL DESCRIPTION OF THERAPEUTIC ALTERNATIVES

In case of failure of conservative measures, the therapeutic alternatives are surgical procedures associated with a mesh choice that has comparable surgical risks, complication rates, short term success rates, recurrence rates and reintervention rates than that on *INTRAMESH® T1/ PROMESH® SURG INTRA* and *BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB* implants.

7. SUGGESTED TRAINING FOR USERS

INTRAMESH® T1/ PROMESH® SURG INTRA and *BIOMESH® CA.B.S.'Air® Composite/ PROMESH® SURG UMB* implants are devices intended to be used by a qualified surgeon with knowledge in anatomy and visceral surgery. No specific training is required from the patients.

8. REFERENCE TO ANY HARMONISED STANDARDS AND COMMON SPECIFICATIONS APPLIED

Hereafter, the list of the standards claimed for the Sterile Intraperitoneal Non Resorbable Parietal Reinforcement Implants, showing the standards harmonised according to Regulation 2017/745 and the state of the art of standards.

Harmonized standards following Regulation (EU) 2017/745 are identified with “*” in the present document. No common specifications have been identified for the moment.

Harmonised Standard Title	Date
*EN ISO 13485:2016/A11:2021 - Medical Devices - Quality Management Systems - Requirements for Regulatory purposes	March 2016/A11:September 2021
NF EN ISO 9001:2015 - Quality Management Systems - Requirements	October 2015
NF EN ISO/IEC 17025:2017 - General requirements for the competence to carry out tests and/or calibrations, including sampling. It covers testing and calibration performed using standard methods, nonstandard methods, and laboratory-developed methods	December 2017
NF EN ISO 14630:2013 - Non-active surgical implants - General requirements	February 2013
*EN ISO 15223-2021 - Medical devices – Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements	September 2021
*EN ISO 14971: 2019 - Medical devices – Application of risk management to medical devices	December 2019
NF EN 62366-1:2015 Amd 1:2020 - Medical Devices - Application of usability engineering to medical devices	August 2020
EN 556-1:2001/AC:2006 - Sterilization of medical devices. Requirements for medical devices to be designated "STERILE" – Part 1: Requirements for terminally sterilized medical devices	November 2006
*EN ISO 11135-1:2014 (+A1:2019) - Sterilization of health care products – Ethylene oxide – Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices	July 2014/ A1: November 2019
*EN ISO 11737-1:2018/A1 2021 - Sterilization of health care products — Microbiological methods — Part 1: Determination of a population of microorganisms on products	January 2018/ A1:June 2021
*EN ISO 11737- 2:2020 - Sterilization of medical devices – Microbiological methods - Part 2: Tests of sterility performed in the definition, validation and maintenance of a sterilization process	May 2020
*EN ISO 11607-1:2020 - Packaging for terminally sterilized medical devices – Part 1: Requirements for materials, sterile barrier systems and packaging systems	January 2020
*EN ISO 11607-2:2020 - Packaging for terminally sterilized medical devices – Part 2: Validation requirements for forming, sealing and assembly processes	January 2020

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EN 14937:2009 - Sterilization of health care products - General requirements for characterization of a sterilizing agent and the development, validation and routine control of a sterilization process for medical devices (ISO 14937:2009)	October 2009
NF EN ISO 10993-1:2020 Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management process	December 2020
NF EN ISO 10993- 3:2014 - Biological evaluation of medical devices – Part 3: Tests for genotoxicity, carcinogenicity and reproductive toxicity	December 2014
NF EN ISO 10993- 5:2010 - Biological evaluation of medical devices – Part 5: Tests for in vitro cytotoxicity	July 2010
NF EN ISO 10993-6:2017 - Biological evaluation of medical devices – Part 6: Tests for local effects after implantation	April 2017
ISO 10993-7/A1:2019 - Biological evaluation of medical devices - Part 7: Ethylene oxide sterilization residuals	December 2019
ISO 10993-10:2021 - Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization	November 2021
NF EN ISO 10993-11 :2018 - Biological evaluation of medical devices – Part 11: Tests for systemic toxicity (ISO 10993-11:2017)	May 2018
*EN ISO 10993-12:2021- Biological evaluation of medical devices – Part 12: Sample preparation and reference materials	January 2021
NF EN ISO 10993-17:2009 - Biological evaluation of medical devices – Part 17: Establishment of allowable limits for leachable substances	December 2009
NF EN ISO 10993-18:2020 - Biological evaluation of medical devices. Part 18: Chemical Characterization of materials	May 2020
NF EN ISO 14155:2020 - Clinical investigation of medical devices for human subjects - Good clinical practice	August 2020
NF EN ISO 14644-1:2016 - Cleanrooms and associated controlled environments—Part 1: Classification of air cleanliness by particle	February 2016
NF EN ISO 14644-2:2016 - Cleanrooms and associated controlled environments— Part 2: Monitoring to provide evidence of cleanroom performance related to air cleanliness by particle concentration	February 2016
NF EN ISO 14644-3:2019 - Cleanrooms and associated controlled environments—Part 3: Tests methods	August 2019
NF EN 17141:2020 - Cleanrooms and associated controlled environments — Biocontamination control	August 2020
NF 868-5:2018 - Packaging for terminally sterilized medical devices – Part 5: Sealable pouches and reels of porous materials and plastic film construction - Requirements and test methods	December 2018
ISO/TR 15499:2016 - Biological evaluation of medical devices—Guidance on the conduct of biological evaluation within a risk management process	December 2016
ISO 12891-1:2015 - Retrieval and analysis of surgical implants - Part 1 : Retrieval and handling	July 2015
NF ISO 2859-1:2000/A1:2011 - Sampling procedures for inspection by attributes – Part 1: Sampling plans indexed by acceptable quality level (AQL) for lot-by-lot inspection	April 2000/A1: September 2011
ISO TR 20416:2020 - Medical devices — Post-market surveillance for manufacturers	July 2020
ISO 20417:2021 - Medical devices — Information to be supplied by the manufacturer	April 2021
ISO 24971:2021 - Medical devices — Guidance on the application of ISO 14971	April 2021
ISO 16142:2016 -Medical devices – Recognized essential principles of safety and performance of medical devices – Part 1: general essential principles and additional specific essential principles for all non-IVD medical devices and guidance on the selection of standards	March 2016

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