INTRAMESH® T1

STERILE INTRAPERITONEAL NON-RESORBABLE PARIETAL REINFORCEMENT IMPLANTS

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Version du 11/02/2026

only

Caution : Federal law (USA) restricts this device to sale, distribution and use by or on the order of a physician.

This release is the last update of the instructions of use and replace the previous edition

Product under regulation 2017/745 (EU).

IMPLANT INTRAMESH® T1
STERILE INTRAPERITONEAL NON-RESORBABLE PARIETAL REINFORCEMENT IMPLANTS
 SINGLE USE STERILE PRODUCT

MANUFACTURER IDENTIFICATION

COUSIN BIOTECH
 ALLEE DES ROSES
 59117 WERVICQ SUD
 FRANCE

IDENTIFICATION AND DESCRIPTION OF THE DEVICE

INTRAMESH® T1 implants are parietal reinforcement implants.

INDICATIONS

Repair and parietal reinforcement of ventral, umbilical and incisional hernias requiring reinforcement with non absorbable support material.

CONTRAINDICATIONS

Do not use in the following cases:

- Allergy to any of the components
- Infected site
- Pregnancy
- Growing children
- Anticoagulant therapy

UNDESIRABLE SIDE EFFECTS

Like any implantable medical device, this implant is susceptible to generate possible undesirable effects:

- | | | | |
|-------------------|----------------------|---------------------|----------------------------|
| - Discomfort/Pain | - Adhesion formation | - Hematoma | - Foreign Body Reaction |
| - Infection | - Fistula formation | - Mesh migration | - Irritation near by organ |
| - Recurrence | - Erosion | - Mesh deformation | - Obstruction |
| - Inflammation | - Seroma | - Allergic reaction | |

TARGET POPULATION

Adult with ventral hernia including umbilical and incisional hernia necessitating a surgical repair with a method that includes an intraperitoneal implantation.

IMPLANTED MATERIALS

Polypropylene and ePTFE for anti-adherent side.

Biocompatible ink with dimethyl siloxane depending on references (See on the package).

References	Weight of ePTFE (g)	Weight of Polypropylene (g)	Weight of Biocompatible ink (µg)
FBIOT11015	0.11 g	1.63 g	< 12
FBIOT11515	0.17 g	2.45 g	< 12
FBIOT11520	0.23 g	3.26 g	< 12
FBIOT12020	0.30 g	4.35 g	< 12
FBIOT12025	0.38 g	5.44 g	< 12
FBIOT13030	0.68 g	9.79 g	< 12
FBIOT13050	1.13 g	16.31 g	Non applicable for this reference
FBIOT1R12X	0.08 g	1.23 g	< 12

Origin neither human nor animal – Non resorbable.

PERFORMANCES

These products can be used in laparoscopic surgery or open surgery. **INTRAMESH® T1** products have a no-adhesive side, which is designed to avoid or reduce visceral adherences during intra peritoneal use or extra peritoneal use in the presence of a peritoneal defect. Marking in black ink, composed of a cross and dots is on the non-adhesive side. The meshes, biocompatible and non resorbable have the advantage of having a shape memory and a great flexibility.

LINK TO SUMMARY OF SAFETY AND PERFORMANCE CHARACTERISTICS :

https://www.cousin-biotech.com/d-69452bf5e54ae_tf03sscpissue.pdf

MRI SAFETY

The implant has not been tested for safety and compatibility in the case of postoperative investigation using Magnetic Resonance Imaging (MRI). However, the implant is mainly made with polypropylene and non metallic component. Consequently, **INTRAMESH® T1** is considered MR safe

PRECAUTIONS FOR USE

INTRAMESH® T1 implants are delivered sterile. Before any use, inspect the integrity of packaging and device (of which blister / pelables pouches). Do not use in the event of deterioration of the device and/or the packaging. Do not use if the device is out of date.

COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned.

This device must be implanted only by a qualified surgeon and trained to the use of the product (knowledge of anatomy and visceral surgery).

As indicated on the product labelling, the implant, is designed for single use. They should in no event be reused and/or reesterilized (potential risks include but are not limited to: loss of product sterility, risk of infection, loss of effectiveness, recurrence).

IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE

Do not cut out the mesh and suture it at least 1cm of the edge. The trocar to be used must be chosen according to the size of the implant to be used. This is to allow passage without excessive force. Always roll the implant with the ePTFE face inside for the passage in the trocar. For use in keyhole surgery, we recommend the use of non-notched atraumatic forceps.

INSTRUCTIONS TO BE GIVEN TO THE PATIENTS BY THE SURGEON

The surgeon should inform the patient of the potential physical and psychological restrictions and consequences of implanting the device. The patient must be informed of the surgical risks and possible side effects. The surgeon should invite the patient to return for further consultation if the patient presents symptoms that appear abnormal. An implant card for the patient is provided by the healthcare professional at the clinic. This implant card provides information, for the patient, to identify the device and traceability elements as well as the name, address and website of the manufacturer. The surgeon shall invite the patient to scan the implant card immediately upon receipt to keep track of it in case of loss. The surgeon shall also inform the patient that the summary of the safety and performance characteristics of the device and the instructions for use can be found on the COUSIN BIOTECH website.

STORAGE PRECAUTIONS

Store in a dry place away from sunlight and at room temperature in its original packaging.

SURGICAL TECHNIQUE

Ports are placed according to surgeon's preference and to hernia's location. Use appropriate trocar size to allow the prosthesis to slide down the trocar with minimal force.

The content of the hernia sac is reduced. Once all the content of the hernia sac has been reduced, it is critical that the surgeon evaluates what has been removed to ensure that there is no bowel injury and that all adhesions that could be lysed have been lysed.

The defect is measured, and the relevant size of mesh is selected ensuring enough overlap. The mesh should overlap the defect(s) 3–5 cm on all sides. For larger defects, the overlap should be at least 5 cm.

The mesh is rolled up ePTFE side inwards and introduced into the abdomen through the trocar: it is very important that **INTRAMESH® T1** is oriented correctly. The white anti-adherent side (ePTFE) is designed to minimize bowel adhesions to the prosthesis. Place this side of the mesh towards bowels. It is still recommended to pull down omentum wherever possible beneath the prosthesis to reduce the risk of visceral adhesion. The landmarks on the ePTFE side allow correct centering of the mesh over the defect.

The polypropylene side should face the abdominal wall. The polypropylene mesh surface should never be placed against bowel or other visceral structures.

Mesh Fixation

INTRAMESH® T1 should always be fixated. Fixation devices or nonabsorbable monofilament sutures are recommended to secure the prosthesis. If other fixation devices are used, they must be indicated for use in hernia repair. Care should be taken to ensure that the prosthesis is adequately fixated to the abdominal wall. If necessary additional fasteners and/or sutures should be used.

Once the mesh is secured, remove all ports under direct vision to ensure there is no bleeding. Evaluate the mesh as the pneumoperitoneum is released. The mesh must lay flat.

EXPLANTATION AND DISPOSAL OF DEVICES

Devices should be retrieved and handled according to the recommendations of ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling". Any explanted device must be sent back for analysis, following the current protocol. This protocol is available on request from COUSIN BIOTECH. It is important to note that any implant that must not be cleaned or disinfected prior to despatch must be contained in sealed packaging. The removed medical device must be disposed of in accordance with the country's standards for the disposal of infectious waste.

There are no specific recommendations for the disposal of a non-implanted device.

Risks associated with the explantation: the surgeon who decides to remove the device should consider factors such as the risk of another surgery for the patient, and the difficulty of the explantation procedure. The removal of the implant is up to the surgeon's decision and must be subject to adequate postoperative follow-up.

IMPORTANT

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

INFORMATION REQUEST AND COMPLAINTS

Following its quality policy, COUSIN BIOTECH is committed to making every effort to produce and supply a high-quality medical device. However, if a health professional (client, user, prescriber, etc.) has a complaint or cause for dissatisfaction with a product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In case of failure of an implant or if it has contributed to causing a serious adverse reaction in the patient, the health centre must follow the legal procedures in that country and inform COUSIN BIOTECH immediately. For any correspondence, please specify the reference, batch number, contact person's details and a comprehensive description of the incident or complaint.

Brochures, documentation and surgical technique are available on request from COUSIN BIOTECH and its distributors. If further information is needed or required, please contact your COUSIN BIOTECH representative or distributor.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

IMPLANT CARD INFORMATION LEAFLET

How to complete the implant card ?

The implant card, for the patient, should be completed as follows:

Front side of the implant card

The user has to complete the implant card before transmitting the implant card to the patient. Please complete it as follows:

	<div style="background-color: #0056b3; color: white; padding: 2px; font-weight: bold;">Carte Implant - Implant Card</div> <div style="text-align: center; font-weight: bold; color: #0056b3;">A conserver / Keep it</div>		<div style="background-color: #0056b3; color: white; padding: 2px; font-weight: bold;">Carte Implant - Implant Card</div>
Patient name	_____ _____		
Date of surgical procedure	. . / . . /		<div style="font-weight: bold; color: #0056b3; font-size: 1.2em;">COUSIN SURGERY</div> <div style="color: #0056b3; font-style: italic; font-size: 0.8em;">We care for Surgery</div>
Name of practitioner	_____ _____		

Back side of the implant card

No fields need to be completed by the user. The following information on the device is printed directly by the manufacturer on the back of the implant card. A QR code allowing the patient to consult the instruction for use, as well as the website address, are available on the implant card.

Information to identify the device	Traceability information
<ul style="list-style-type: none"> Indication that the product is a medical device Description of the device family Trade name of the device Reference of the device Device is MR safe Denomination name of the device in sales countries 	<ul style="list-style-type: none"> UDI number Batch number Name and address of the legal manufacturer

Explanation of the used symbols in the Implant Card

All symbols are commonly used and explained at the end of the instruction for use.

IMPLANT INTRAMESH® T1
IMPLANTS DE RENFORCEMENT PARIETAL INTRAPERITONEAL STERILE NON RESORBABLE
PRODUIT STERILE A USAGE UNIQUE

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Patient information leaflet INTRAMESH® T1

Item	Patient information leaflet
1	<p>a. The name of the device : INTRAMESH® T1</p> <p>b. The model of the device : FBIOT11015 / FBIOT11515 / FBIOT11520 / FBIOT12020 / FBIOT12025 / FBIOT13030 / FBIOT13050 / FBIOT1R12X</p>
2	<p>a. The intended purpose of the device : INTRAMESH® T1 implants are intended for the repair and parietal reinforcement of ventral, umbilical and incisional hernias requiring reinforcement with non-absorbable support material.</p> <p>b. The kind of patient on whom the device is intended to be used : Adult with ventral hernia including umbilical and incisional hernia necessitating a surgical repair with a method that includes an intraperitoneal implantation.</p>
3	<p>Special operating instructions for the use of the device :</p> <p>Ports are placed according to surgeon's preference and to hernia's location. Use appropriate trocar size to allow the prosthesis to slide down the trocar with minimal force.</p> <p>The content of the hernia sac is reduced. Once all the content of the hernia sac has been reduced, it is critical that the surgeon evaluates what has been removed to ensure that there is no bowel injury and that all adhesions that could be lysed have been lysed.</p> <p>The defect is measured, and the relevant size of mesh is selected ensuring enough overlap. The mesh should overlap the defect(s) 3–5 cm on all sides. For larger defects, the overlap should be at least 5 cm.</p> <p>The polypropylene side should face the abdominal wall. The polypropylene mesh surface should never be placed against bowel or other visceral structures.</p> <p>Mesh Fixation</p> <p>INTRAMESH® T1 should always be fixated. Fixation devices or nonabsorbable monofilament sutures are recommended to secure the prosthesis. If other fixation devices are used, they must be indicated for use in hernia repair. Care should be taken to ensure that the prosthesis is adequately fixated to the abdominal wall. If necessary additional fasteners and/or sutures should be used.</p> <p>Once the mesh is secured, remove all ports under direct vision to ensure there is no bleeding. Evaluate the mesh as the pneumoperitoneum is released. The mesh must lay flat.</p>

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Item	Patient information leaflet
4	<p>a. The intended performance of the device : These products can be used in laparoscopic surgery or open surgery. INTRAMESH® T1 products have a non-adhesive side, which is designed to avoid or reduce visceral adherences during intra peritoneal use or extra peritoneal use in the presence of a peritoneal defect. Marking in black ink, composed of a cross and dots is on the non-adhesive side. The meshes, biocompatible and non resorbable have the advantage of having a shape memory and a great flexibility</p> <p>b. Undesirable side effects : Like any implantable medical device, this implant is susceptible to generate possible undesirable effects:</p> <ul style="list-style-type: none"> -Discomfort/Pain - Adhesion formation - Hematoma - Foreign Body Reaction - Infection - Fistula formation - Mesh migration - Irritation near by organ - Recurrence - Erosion - Mesh deformation - Obstruction - Inflammation - Seroma - Allergic reaction
5	<p>The potential risks that could arise if the instructions and the precautions for use are not met are (for example): loss of product sterility risk of infection loss of effectiveness recurrence</p>
6	<p>MRI Information:</p> <p>The implant has not been tested for safety and compatibility in the case of postoperative investigation using Magnetic Resonance Imaging (MRI). However, the implant is mainly made with polypropylene and non-metallic component. Consequently, INTRAMESH® T1 is considered MR safe.</p>
7	<p>Precautions :</p> <p>The device is intended to be permanently implanted in the body. No maintenance is therefore possible. If the patient experiences any symptoms that appear abnormal, the patient is invited to return to the surgeon for further consultation. An implant card for the patient is provided by the healthcare professional at the clinic. This implant card provides information, for the patient, to identify the device and traceability elements as well as the name, address and website of the manufacturer.</p> <p>The removal of the implant is up to the surgeon's decision and must be subject to adequate postoperative follow-up. Some risks are associated with the explantation. The surgeon who decides to remove the device should consider factors such as the risk of another surgery for the patient and the difficulty of the explantation procedure.</p>

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Item	Patient information leaflet
8	<p>Materials included in the device :</p> <p>The implanted materials are Polypropylene and ePTFE for anti-adherent side.</p> <p>The biological safety of the materials and any manufacturing residuals has been demonstrated through chemical characterization and biological tests.</p>
9	<p>Reporting a serious incident :</p> <p>If you encounter any serious incident with your INTRAMESH® T1 implant you should report it immediately to your doctor. You should also report the incident to Cousin Biotech via the link serviceclients@cousin-surgery.com or materiovigilance@cousinbiotech.com and to the Therapeutic Goods Administration using the following website: https://www.tga.gov.au</p>
10	<p>Legal Manufacturer : Cousin Biotech, Allée des Roses, 59117 Wervicq-Sud, France</p> <p>Manufacturing Facility : Cousin Biotech, Allée des Roses, 59117 Wervicq-Sud, France</p>

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Cousin Biotech SAS au capital de 320 000 € - NAF 1396Z- N° TVA FR 34 398 460 261 - Siren 398 460 261 RCS Lille