# 4D DOME<sup>®</sup>

## SEMI-RESORBABLE PARIETAL REINFORCEMENT IMPLANT

en	<a href="#">Instructions for use</a>	Page	2	MDR CE marking date: 30/09/2025
fr	<a href="#">Notices d'instructions</a>	Page	6	Date de marquage CE MDR : 30/09/2025
de	<a href="#">Gebrauchsanweisung</a>	Seite	10	Datum der MDR-CE-Kennzeichnung: 30/09/2025
it	<a href="#">Istruzioni per l'uso</a>	Pagina	14	Data di marcatura CE MDR: 30/09/2025
es	<a href="#">Instrucciones de uso</a>	Página	18	Fecha del marcado CE MDR: 30/09/2025
pt	<a href="#">Instruções de utilização</a>	Página	22	MDR Data da marcação CE: 30/09/2025
ελ	<a href="#">Οδηγίες χρήσης</a>	Σελίδα	26	Ημερομηνία σήμανσης MDR CE: 30/09/2025
nl	<a href="#">Gebruiksaanwijzing</a>	Pagina	30	Datum MDR CE-markering: 30/09/2025
ro	<a href="#">Instrucțiuni de utilizare</a>	Pagina	34	Data marcajului MDR CE: 30/09/2025
tr	<a href="#">Kullanım talimatları</a>	Sayfa	38	MDR CE işareti tarihi: 30/09/2025
cs	<a href="#">Návod k použití</a>	Strana	42	Datum označení MDR CE: 30. 9. 2025
en	<a href="#">Patient information leaflet (Australian requirements)</a>	Page	46	MDR CE marking date: 30/09/2025

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**CE**  
0297  
**NOT301\_260216**  
Version du 16/02/2026

Rx  
only

Caution : Federal law (USA) restricts this device to sale distribution and use by or on the order of a physician

This release is the last update of the instructions of use and replace the previous edition Product under regulation 2017/745 (EU).

# 4DDOME®

## SEMI-RESORBABLE PARIETAL REINFORCEMENT IMPLANT SINGLE USE STERILE PRODUCT

### MANUFACTURER IDENTIFICATION

COUSIN BIOTECH  
ALLEE DES ROSES  
59117 WERVICQ SUD  
FRANCE

### IDENTIFICATION AND DESCRIPTION OF THE DEVICE

The 4DDOME® device consists of a semi-resorbable dome and a semi-resorbable reinforcement mesh for anterior and posterior parietal reinforcement of the groin.

### INTENDED USE

Intended to be used as parietal reinforcement implant.

### INDICATIONS

Repair and surgical reinforcement of inguinal and femoral hernias.

### CONTRAINDICATIONS

Do not implant in the following cases:

- Allergy to any of the components
- Infected site
- Pregnancy
- Growing children

### UNDESIRABLE SIDE EFFECTS

Like any implantable medical device, this implant is susceptible to generate possible undesirable side effects that may lead to reoperation:

- Discomfort/Pain
  - Infection
  - Inflammation
  - Fistula
  - Recurrence
  - Adhesion
  - Erosion/Extrusion
  - Seroma/lymphocele
  - Hematoma
  - Mesh migration
  - Mesh deformation
  - Allergic reaction
  - Foreign body reaction
  - Irritation nearby organ
  - Obstruction
- A small additional risk of male infertility has been identified after bilateral groin hernia repair using meshes

### TARGET POPULATION

Adult male or female suffering from groin hernias with an indication of surgical repair with extraperitoneal implantation. These include:

- Inguinal hernias
- Femoral hernias

### IMPLANTED MATERIALS

Polypropylene (PP) (non resorbable) - Poly-L-Lactic Acid (PLLA) (resorbable)

		Weight PP (g)	Weight PLLA (g)
4DDOM24GSR	Dome	0,03	0,26
	Flat mesh	0,21	0,64
4DDOM30GSR	Dome	0,05	0,41
	Flat mesh	0,21	0,64
4DDOM38GSR	Dome	0,08	0,60
	Flat mesh	0,21	0,64
4DDOME24SR	Dome	0,03	0,26
	Flat mesh	0,12	0,37
4DDOME30SR	Dome	0,05	0,41
	Flat mesh	0,12	0,37
4DDOME38SR	Dome	0,08	0,60
	Flat mesh	0,12	0,37

Origin neither human nor animal origin. Semi-resorbable

### PERFORMANCES

The device 4DDOME® contributes to reduce symptoms from hernia and contributes to avoid complications of hernia.

## **LINK TO SUMMARY OF SAFETY AND PERFORMANCE CHARACTERISTICS:**

[https://www.cousin-biotech.com/d-68e9066a67fce\\_tf07sscp4d.pdf](https://www.cousin-biotech.com/d-68e9066a67fce_tf07sscp4d.pdf)

### **MRI SAFETY**

The implant has not been tested for safety and compatibility in the case of postoperative investigation using Magnetic Resonance Imaging (MRI). However, the implant is made with Polypropylene and Poly-L-Lactic Acid: non metallic components. Consequently, 4DDOME® is considered MR safe.

### **PRECAUTIONS FOR USE**

The dome should always be installed with an anterior reinforcement.

The dome should not be cut.

4DDOME® devices must always be separated from abdominal cavity by peritoneum

4DDOME® devices are delivered sterile (ethylene oxide sterilisation).

When opening the carton box, the surgeon has to make sure that the lot number, the model and the size of the mesh are equal to the data on the label of the inner pouch. Otherwise, the use of the prosthesis is strictly prohibited.

Before any use, inspect the integrity of packaging and device (of which blister / peel-away pouches). Do not use in the event of deterioration of the device and/or the packaging.

Do not use if the device is out of date. COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned.

Unpacking and manipulation of the implant must be carried out following asepsis standards in order to guarantee the sterility.

Prevent the mesh from any contact with objects that may alter their surface.

This device must be implanted by a qualified surgeon (familiar with the relevant anatomy and experienced in visceral surgery)

### **IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE**

As indicated on the product labelling, the 4DDOME® implant is designed for single use. It should in no event be reused and/or resterilized (potential risks include but are not limited to: loss of product sterility, risk of infection, loss of effectiveness, recurrence).

### **INSTRUCTIONS TO BE GIVEN TO THE PATIENTS BY THE SURGEON**

The surgeon should inform the patient of the potential physical and psychological restrictions and consequences of implanting the device. The patient must be informed of the surgical risks and possible side effects. The surgeon should invite the patient to return for further consultation if the patient presents symptoms that appear abnormal. An implant card for the patient is provided by the healthcare professional at the clinic. This card provides information to identify the device and traceability elements as well as the name, address and website of the manufacturer. Please refer to the section "Implant card information leaflet" for instructions on how to complete it. **The surgeon shall invite the patient to scan the implant card immediately upon receipt to keep track of it in case of loss.** The summary of the safety and performance characteristics of the device and the instructions for use can also be found on the COUSIN BIOTECH website.


### **STORAGE PRECAUTIONS**

To be stored in a dry place away from sunlight and at room temperature in its original packaging.

### **SURGICAL TECHNIQUE**

#### **4DDOME® Surgical Technique (1)**

##### **Indirect Hernia**

① In indirect inguinal hernia, the sac is dissected up to the internal inguinal ring

② The 4DDOME® is placed in the defect maintaining the sac deeply. The 4DDOME® is then fixed by 1 or 4 stitches with non-resorbable suture to the resistant edge of the internal ring

③ Positioning the split onlay mesh around the spermatic cord and fixation with non-resorbable suture as in the Lichtenstein technique.

## 4DDOME® Surgical Technique (2)

### Direct Hernia

① In direct hernia the sac is dissected up to the transversalis fascia.

② The 4DDOME® is placed in the opening whilst holding the sac deeply. The 4DDOME® is then fixed by 7 or 8 stitches with non-resorbable suture to the resistant edge of the transversalis fascia.

③ Positioning the split onlay mesh around the spermatic cord and fixation with non-resorbable suture as in Lichtenstein's technique.

### **EXPLANATION AND DISPOSAL OF DEVICES**

Risks associated with the explantation: the surgeon who decides to remove the device should consider factors such as the risk of another surgery for the patient, and the difficulty of the explantation procedure. The removal of the implant is up to the surgeon's decision and must be subject to adequate postoperative follow-up.

Devices should be retrieved and handled according to the recommendations of ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling". Any explanted device must be sent back for analysis, following the current protocol. This protocol is available on request from COUSIN BIOTECH. It is important to note that any implant that must not be cleaned or disinfected prior to despatch must be contained in sealed packaging. The removed medical device must be disposed of in accordance with the country's standards for the disposal of infectious waste.

There are no specific recommendations for the disposal of a non-implanted device.

### **IMPORTANT**

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

### **INFORMATION REQUEST AND COMPLAINTS**

Following its quality policy, COUSIN BIOTECH is committed to making every effort to produce and supply a high-quality medical device. However, if a health professional (client, user, prescriber, etc.) has a complaint or cause for dissatisfaction with a product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In case of failure of an implant or if it has contributed to causing a serious adverse reaction in the patient, the health centre must follow the legal procedures in that country and inform COUSIN BIOTECH immediately. For any correspondence, please specify the reference, batch number, contact person's details and a comprehensive description of the incident or complaint.

Brochures, documentation and surgical technique are available on request from COUSIN BIOTECH and its distributors. If further information is needed or required, please contact your COUSIN BIOTECH representative or distributor.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## IMPLANT CARD INFORMATION LEAFLET

### How to complete the implant card ?

The implant card, for the patient, should be completed as follows:

#### **Front side of the implant card**

The user has to complete the implant card before transmitting the implant card to the patient. Please complete it as follows:

	Carte Implant - Implant Card A conserver / Keep it	Carte Implant - Implant Card
Patient name	 _____ _____	 <b>COUSIN SURGERY</b> <i>We care for Surgery</i>
Date of surgical procedure	 . . / . . / . . . . .	
Name of practitioner	 _____ _____	

#### **Back side of the implant card**

No fields need to be completed by the user. The following information on the device is printed directly by the manufacturer on the back of the implant card. A QR code allowing the patient to consult the instruction for use, as well as the website address, are available on the implant card.

Information to identify the device	Traceability information
<ul style="list-style-type: none"><li>• Indication that the product is a medical device</li><li>• Description of the device family</li><li>• Trade name of the device</li><li>• Reference of the device</li><li>• Device is MR safe</li><li>• Denomination name of the device in sales countries</li></ul>	<ul style="list-style-type: none"><li>• UDI number</li><li>• Batch number</li><li>• Name and address of the legal manufacturer</li></ul>

### Explanation of the used symbols in the Implant Card

All symbols are commonly used and explained at the end of the instruction for use.

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## Patient information leaflet 4DDOME®

Item	Patient information leaflet
1	<p>a. <b>The name of the device</b> : 4DDOME®</p> <p>b. <b>The model of the device</b> : 4DDOM24GSR, 4DDOM30GSR, 4DDOM38GSR, 4DDOME24SR, 4DDOME30SR, 4DDOME38SR</p>
2	<p>a. <b>The intended purpose of the device</b> : Repair and surgical reinforcement of inguinal and femoral hernias.</p> <p>b. <b>The kind of patient on whom the device is intended to be used</b> : Adult male or female suffering from a symptomatic ventral hernia necessitating a surgical mesh repair with an extraperitoneal implantation method.</p>
3	<p><b>Special operating instructions for the use of the device</b> : The device must be implanted only by a qualified surgeon and trained on the use of the product (knowledge of anatomy and visceral surgery).</p>
4	<p>a. <b>The intended performance of the device</b> :</p> <p>These implants, biocompatible and semi-resorbable make a reactional fibrosis which replaces the reinforcement after an implantation of 6 months. They have the advantages to have a very resistance to the suture, a quick and optimum integration and colonization and to have a resistance to the abdominal pressure.</p> <p>The 4D DOME® device is particularly designed to be implanted extra-peritoneally.</p> <p>b. <b>Undesirable side effects</b> :</p> <p>Like any implantable medical device, this implant could generate possible undesirable effects :</p> <ul style="list-style-type: none"> <li>- Discomfort/pain</li> <li>- Infection</li> <li>- Recurrence</li> <li>- Inflammation</li> <li>- Adhesions</li> <li>- Fistula formation</li> <li>- Erosion</li> <li>- Seroma</li> <li>- Hematoma</li> <li>- Mesh migration</li> <li>- Mesh deformation</li> <li>- Allergic reaction</li> <li>- Foreign body reaction</li> <li>- Irritation near by organ</li> <li>- Lymphocele</li> <li>- A small additional risk of male infertility has been identified after bilateral groin hernia repair using meshes (open or laparoscopic)</li> </ul>

This leaflet was revised January 2022

Page 46/52

**Headquarters:** Allée des roses - 59117 Wervicq-Sud – France **Tel** +33 (0) 3 20 14 41 19 - **Fax** +33 (0) 3 20 14 40 13  
**Factory:** Allée des roses - 59117 Wervicq-Sud – France **Tel** +33 (0) 3 20 14 41 20 - **Fax** +33 (0) 3 20 14 40 45  
[www.cousin-biotech.com](http://www.cousin-biotech.com) - **e-mail** f.pelletier@cousin-biotech.com  
 Cousin Biotech SAS au capital de 320 000 € - NAF 1396Z- N° TVA FR 34 398 460 261 - Siren 398 460 261 RCS Lille


5	<p><b>When 4DDOME® Should Not be Used :</b> Do not implant in the following cases: Allergy to any of the components, infected site, pregnancy, growing children, anticoagulant therapy</p>
6	<p><b>MRI Information :</b> 4DDOME ® are non-conducting and non-magnetic devices. They present no risk in Magnetic Resonance Imaging (MRI) environment.</p>
7	<p><b>Precautions :</b> The dome should always be installed with an anterior reinforcement. The dome should not be cut. 4D DOME® devices are delivered sterile (ethylene oxide sterilisation). Before any use, inspect the integrity of packaging and device (of which blister / peel-away pouches). Do not use in the event of deterioration of the device and/or the packaging. Do not use if the device is out of date. COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned. This device must be implanted only by a qualified surgeon and trained on the use of the product (knowledge of anatomy and visceral surgery).</p>
8	<p><b>Materials included in the device:</b> This surgical mesh is composed of Polypropylene (non resorbable) Poly-L-Lactic Acid (resorbable) Origin neither human nor animal origin. Semi-resorbable</p>
9	<p><b>Reporting a serious incident :</b> If you encounter any serious incident with your 4DDOME® implant you should report it immediately to your doctor. You should also report the incident to Cousin Biotech via the link <a href="mailto:serviceclients@cousin-surgery.com">serviceclients@cousin-surgery.com</a> or <a href="mailto:materiovigilance@cousin-surgery.com">materiovigilance@cousin-surgery.com</a> and to the Therapeutic Goods Administration using the following website: <a href="http://www.tga.gov">http://www.tga.gov</a></p>
10	<p><b>Legal Manufacturer :</b> Cousin Biotech, Allée des roses, 59117 Wervicq-Sud, France <b>Manufacturing Facility :</b> Cousin Biotech, Allée des roses, 59117 Wervicq-Sud, France</p>

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