Instruments for BatSPINE / Instruments pour BatSPINE / Instrumente für BatSPINE / Strumenti per BatSPINE / Instrumental para BatSPINE / Instrumentos para BatSPINE / Όργανα για BatSPINE / Narzędzia do BatSPINE / Instrumente pentru BatSPINE / Nástroje pro BatSPINE / Instrumente für BatSPINE

Instruments for IntraSPINE / Instruments pour IntraSPINE / Instrumente für IntraSPINE / Strumenti per IntraSPINE / Instrumental para IntraSPINE / Instrumentos para IntraSPINE / Όργανα για IntraSPINE / Narzędzia do IntraSPINE / Instrumente pentru IntraSPINE / Nástroje pro IntraSPINE / Instrumente für IntraSPINE

Instruments for BDyn / Instruments pour BDyn / Instrumente für BDyn / Strumenti per BDyn / Instrumental para BDyn / Instrumentos para BDyn / Όργανα για BDyn / Narzędzia do BDyn / Instrumente pentru Bdyn / Nástroje pro Bdyn / Instrumente für BDyn

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NOT348 / 260427

Version du 27/04/2026

This release is the last update of the instructions of use and replace the previous edition.

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PURPOSE

This manual provides information on the care, decontamination, cleaning, disinfection, maintenance and sterilisation of Class Ir reusable surgical instruments manufactured and/or distributed by COUSIN BIOTECH.

This document has been created to assist healthcare personal in the safe handling, effective maintenance and reprocessing of reusable instruments manufactured and/or distributed by COUSIN BIOTECH.

GLOSSARY

Contaminated : state of an item that has actually or potentially been in contact with micro-organisms.

Decontamination : use of physical or chemical methods to eliminate, deactivate or destroy blood-borne pathogens on a surface or element until they can no longer transmit infectious particles and until the surface or elements are safe for handling or destruction.

Disinfection : process used to reduce the number of viable micro-organisms on a product to a level previously specified as suitable for handling or use.

Maintenance : visual and/or functional check

Packaging : metal case for storing instruments, designed to ensure correct sterilisation of the instruments.

Washer/Disinfector : machine designed to clean and disinfect medical instruments and other items used in medical, dental, pharmaceutical and veterinary practices.

Cleaning : removal of all contamination from an item, necessary before any other treatment.

- manual: cleaning without the use of a washer or automatic washer/disinfector.
- automated: cleaning using a washer or automatic washer/disinfector.

Sterile : free from viable micro-organisms.

Sterilisation : validated process used to eliminate all forms of viable micro-organisms from an instrument.

Treatment : A process including decontamination, cleaning, disinfection, conditioning and sterilisation, necessary to prepare a new or used surgical instrument for use.

Note : In a sterilisation process, the nature of microbiological death is described by an exponential function. Therefore, the presence of micro-organisms on an item can be expressed in terms of probability. While this probability can be reduced to a very small number, it can never be reduced to zero. This probability can only be assured for validated processes.

CONSIDERATIONS

This instruction manual relates to all COUSIN BIOTECH reusable surgical instruments and should be read carefully. These instruments are specifically dedicated to the placement of spinal and orthopaedic surgery implants (see implant table page 9).

Spinal and orthopaedic surgery may require heavy instruments with several components, articulated or rotating parts and removable handles. Instruments are usually supplied in containers to optimise space requirements and follow the operating procedure as closely as possible.

During surgery, instruments are contaminated by blood, bone and tissue debris. They may also be contaminated by body fluids containing pathogens. All healthcare personnel must be familiar with the universal precautions necessary to avoid injuries caused by handling sharp and/or pointed instruments during and after surgical procedures and during reprocessing.

These instruments must only be used by medical care staff such as surgeons, IBODEs, qualified instrument technicians, reconditioning staff (disinfection, washing and sterilisation). These staff are trained to handle them and to master their use, as well as the operating technique for the implant in question.

New and used instruments must be carefully treated in accordance with these instructions before being used.

Healthcare establishments are responsible for cleaning, disinfecting, packaging and sterilising all instruments before they are used again or transported.

The main instruments must be complete and in good condition in order to be used correctly.

PROCESSING INSTRUCTIONS

Surgical instruments are supplied non-sterile and must be cleaned, disinfected and sterilised before use. After use, these instruments must be decontaminated, cleaned, stored, packaged and sterilised. The following information describes the correct steps, based on the organisations' recommendations (see references), for reprocessing surgical instruments to ensure a safe, efficient and sustainable life cycle.

Only instruments manufactured and/or distributed by COUSIN BIOTECH should be included in the dedicated containers. These instructions do not apply to containers containing instruments not manufactured and/or distributed by COUSIN BIOTECH.

The instructions below describe the necessary processing steps that new and used instruments must undergo to be sterile.

WARNINGS AND PRECAUTIONS

Protection of reprocessing personal

- **Universal precautions must be observed** by hospital staff working with contaminated or potentially contaminated medical instruments. Sharp instruments must be handled with great care.
- **Personal protective equipment (PPE) must be worn** when handling or working with contaminated or potentially contaminated equipment, instruments, and materials. PPE includes a laboratory coat, mask, goggles or face shield, gloves and overshoes.

Instrument protection

- **Do not use wire brushes or scouring pads** during manual cleaning procedures. These materials may damage the surface and finish of the instruments. Use soft bristle nylon brushes and swabs.
- **Do not place heavy instruments on top of fragile ones.**

Other precautions

- Do not allow contaminated instruments to dry before reprocessing. All subsequent cleaning and sterilisation steps are facilitated by not allowing blood, body fluids, tissue and bone debris, saline or disinfectant to dry on the instruments used.
- It is important to take care when manually cleaning instruments used for drilling, milling or piercing to ensure that all bone residues are removed.

Note: it is important to select enzyme solutions designed to dissolve blood, tissues and body fluids. Some enzyme solutions are specifically designed to dissolve faeces or other organic contaminants and may not be suitable for use with spinal and orthopaedic surgery instruments.

- In the event of breakage of the instrument in situ, a fragment of the instrument can be located by imaging methods adapted to the material.

CLEANING AGENTS

- Neutral enzymatic detergents or cleaning products suitable (neutral pH or mild alkaline) for surgical instruments that comply with current regulations must be used during manual cleaning procedures.

Manual scrubbing must always be carried out with the instrument immersed in the cleaning solution to avoid aerosol formation and splashing, which can spread contaminants. Cleaning agents must be easily and completely rinsed from instrument surfaces to prevent the build-up of detergent residues.

- **Alkaline agents with a pH of 12 or less** may be used to clean stainless steel and certain polymer instruments in countries where this is required by regulation or where there is a risk of prion diseases such as Transmissible Spongiform Encephalopathy (TSE) and Creutzfeld-Jakob Disease (CJD). **It is extremely important to neutralise and rinse alkaline cleaning agents completely and thoroughly from instruments.**

Note: If Creutzfeld-Jakob disease (CJD) is suspected, incinerate the affected instrument.

In the event of non-compliance with the above instructions, the healthcare centre must apply an equivalent validated cleaning process for which it will be entirely responsible.

Note: sharp instruments must be carefully inspected after treatment with alkaline detergents to ensure that the sharp edges are ready for use.

UNAUTHORISED AGENTS

- Corrosive cleaning/disinfecting agents **must not** be used.
- Instruments must not come into contact with fluorinated or chlorinated products, or with greasy detergents. Instruments with parts made of synthetic materials (plastics) must never be brought into contact with strong acids (pH < 4), bases, organic or ammoniacal solvents, oxidising chemicals or any other product likely to alter the material.
- **Stripping agents containing morpholine** must not be used in steam sterilisers. These agents leave residues that can damage polymer instruments.

INSPECTION ON RECEIPT

- Upon receipt, instrument sets should be inspected for completeness. In order to check the presence of the instruments on arrival and on departure from the place of use, refer to the plate or laser marking in the bottom of the container illustrating the positioning of each instrument in the container.
- Spinal and orthopaedic surgery procedures follow a precise order for the use of instruments. All instruments and sizes are required, except in special cases. Contact your distributor if requested instruments have been omitted but are required for the operation.

SHELF LIFE

- Repeated treatment, in accordance with the instructions in this leaflet, does not lead to any deterioration in the expected performance of COUSIN BIOTECH reusable surgical instruments unless otherwise stated.
- COUSIN BIOTECH does not indicate the maximum number of acceptable reuses for instruments. The lifespan of reusable surgical instruments depends on a number of factors, including the type of handling undergone and the treatment between uses. The best way to determine when a device should no longer be reused is to carry out a visual inspection and check that it is working properly before each use.

- If polymer surfaces become 'chalky', show visible surface wear (e.g. cracking or delamination) or if polymer instruments show any visual (e.g. discolouration) or functional abnormality they should be replaced.

CLEANING

WARNING :

Automatic cleaning using a washer/disinfector alone may not be sufficiently effective for spinal surgery instruments with lumens, slots, cannulae, blind holes, contacting surfaces and other complex features. A manual pre-cleaning process is recommended. Manual pre-cleaning should be combined with automatic cleaning, but automatic cleaning alone is not sufficient and is not recommended.

- For cleaning, all instruments (including removable handles) should be separated from each other where possible. Disassembly, when necessary, is generally straightforward. Take care not to lose any parts. If a part is lost, inform COUSIN BIOTECH as soon as possible.
- Cleaning must be carried out in accordance with the facility's established practices. Instruments must be removed from metal containers for cleaning procedures. Do not clean instruments in metal containers. Containers and lids where applicable should be cleaned separately from instruments.
- Instruments containing a polymer stopper or protector should be removed.
- Instruments with a standard cleaning connector must be cleaned with a suitable system, available in the care centre.
- The use of hard water should be avoided. Softened tap water may be used for the initial rinse. The final rinse should be carried out with purified (distilled) water to remove mineral deposits from the instruments. One or more of the following processes may be used to purify the water: ultrafiltration (UF), reverse osmosis (RO), deionisation or equivalent.

STERILISATION METHOD

Autoclave sterilisation should be carried out in the specific container supplied by COUSIN BIOTECH, or other container (wire basket, tray, plastic tray placed in an autoclavable peel-off bag) used by the care centre.

- Steam/moist heat is the recommended sterilisation method for spinal and orthopaedic surgical instruments manufactured and/or distributed by COUSIN BIOTECH.
- Ethylene oxide (EO), plasma and dry heat sterilisation methods are not recommended for the sterilisation of COUSIN BIOTECH reusable instruments.
- Instruments with removable sleeves must be disassembled for sterilisation. Instruments containing a polymer stopper or protective element must be removed.
- Avoid the use of acidic anticorrosion agents.

POINT OF USE - PREPARATION FOR TREATMENT

- Remove excess liquid and organic tissue from instruments with a disposable lint-free pad. Place instruments in a tray of distilled water or in a basket covered with a damp cloth. Do not allow saline solution, blood, body fluids, bone and tissue fragments or other organic debris to dry on the instruments before cleaning.

Note: immersion in proteolytic enzyme solutions facilitates cleaning, particularly for complex instruments with hard-to-reach areas (e.g. tubular designs, etc.). These enzymatic solutions dissolve proteinaceous material and prevent blood and protein-based materials from drying on the instruments. The manufacturer's instructions for the preparation and use of these solutions must be explicitly followed.

- Instruments should be cleaned within 30 minutes of use to minimise the possibility of drying prior to cleaning.
- Used instruments should be transported closed or covered to avoid any risk of further contamination.

PREPARATION BEFORE CLEANING

- Multi-component instruments should be disassembled for proper cleaning. Take care not to lose any component. If a part is lost, notify your sales representative as soon as possible.
- Implant instructions for use and surgical techniques may provide an additional source to illustrate assembly/disassembly instructions for specific instruments.

PREPARATION OF CLEANING AGENTS

- All cleaning agents should be prepared to the manufacturer's recommended dilution and temperature. Softened tap water may be used to prepare cleaning agents. Use of the recommended temperatures is important for optimum performance of the cleaning agents.
- Cleaning agents in dry powder form must be completely dissolved before use to avoid oxidation or corrosion of the instruments.
- Fresh cleaning solutions should be prepared when existing solutions become visibly contaminated (bloody and/or cloudy).

CLEANING PRODUCTS	
Water	Use cold distilled or reverse osmosis filtered water
Detergent	Prepare a detergent according to the manufacturer's instructions.
Enzymatic cleaner	Prepare an enzymatic cleaner according to the manufacturer's instructions.
Accessories for manual cleaning	Nylon soft bristle brushes and/or brush head, syringes, gloves, absorbent disposable cloth.
Ultrasonic cleaner	Ultrasonic cleaners should be checked regularly to ensure that they are operating correctly.

Manual cleaning/disinfection instructions :

	Step	Step description	Step instructions	Accessories	Duration
Cleaning steps	1	Contamination removal	Rinse product room temperature running tap water removing any visible organic material with assistance of a soft bristle brush	<ul style="list-style-type: none"> - Tap water - Soft bristle brush (do not utilize metal cleaning) 	Until all visible soil is removed
	2	Drying	Dry the device utilizing a dry non-shedding wipe. Medical quality filtered air may be utilized if available	<ul style="list-style-type: none"> - Non-scheduling wipe - Medical quality filtered compressed air 	Until product is visually dry
Disinfection steps	3	Disinfection Application	Prepare a low-foaming neutral enzymatic detergent solution, using tap water (15-25°C / 59-77°F). Soak the instrument in the open position (if possible) during minimum 1 minute. During the immersion, actuate mobile parts of the instrument minimum 3 times (if possible), so that the detergent can access all parts of the instrument.	<ul style="list-style-type: none"> - Detergent 	Minimum duration 1 minute, this duration can vary according product
	4	Drying	After the soaking time, take out the instrument and wipe it with a disposable cloth. Then, place the instrument in a new bath of enzymatic detergent solution using warm tap water.	<ul style="list-style-type: none"> - Disposable cloth - Detergent 	Until product is visually dry
	5	Manual Desinfection	Brush all surface of the instrument using a soft brush for 2 minutes. Adjust brushing time if needed. Actuate mobile parts of the instrument 3 times (when applicable) and use a syringe to access all difficult areas. Use a volume of 60 mL of the detergent solution.	<ul style="list-style-type: none"> - Soft bristle brush - Detergent - Syringe 	Manual cleaning time duration is complete when the device's surface, joints, and crevices have been manually cleaned
	6	Rinse and finition	Soak the instrument in reverse osmosis water for 1 minute. Use a syringe and 60mL of reverse osmosis water to access all difficult areas. Repeat soaking 2 additional times using fresh water.	<ul style="list-style-type: none"> - Reserve osmosis - Syringe 	Minimum 1 minute
	7	Final rinse	Rincer soigneusement avec de l'eau distillée	<ul style="list-style-type: none"> - Distilled water 	Minimum 1 minute
	8	Final drying	Dry devise utilizing medical quality filtered air	<ul style="list-style-type: none"> - Medical quality filtered air 	Until product is visually dry
	9	Tidying	The instruments have to be placed in the right position inside the trays as described on the position map fixed on the tray for the sterilization.		

These manual cleaning and disinfection instructions can be followed by automatic cleaning:

	Step	Description	Step Instruction	Accessories	Duration
Cleaning Steps	1	Contamination Removal	Rinse product room temperature running tap water removing any visible organic material with assistance of a soft bristle brush	- Tap water - Soft bristle brush (Do not utilize metal cleaning brushes)	Until all visible soil is removed
	2	Pre-soak	Prepare a detergent solution at ambient temperature (15-25°C / 59-77°F). Afterwards lumina, threads, joints and gaps of the instruments have to be flushed with a water jet pistol for a minimum time of 10 seconds for each position.	- Detergent - Water jet pistol	Until product is visually clean
	3	Ultrasonic cleaning	Clean the instrument in the ultrasonic tray including the solution. Be careful every unclench connector or screws and bolts by the vibrations will be reassure.	- Ultrasonic tray	Minimum duration 10 minutes, this duration can vary according product
	4	Soak	Soak the instrument in tap water during 30 seconds. After, ultrasonic treatment the lumina, threads, joints and gaps of the instrument have to be flushed again with a water jet pistol for a minimum time of 10 seconds for each position.	- Reserve osmosis water	Minimum 30 seconds of soak
	5	Automated Washer	Place entire disassembled device into the automated washer	- Automated Washer - Disassembled Device	Minimum total cycle time: 39 minutes

Instructions of cleaning in washer / disinfectant			
Cycles parameters	Time	Minimum temperature	Type of detergent/water
Pre-cleaning	2 minutes	<45°C / <113°F	Tap water
Cleaning	5 minutes	55°C / 131°F	Detergent
Draining	2 minutes	<45°C / <113°F	Tap water
Thermic draining	5 minutes	90°C / 194°F	Reserve osmosis water
Drying	25 minutes	> 70°C / > 158°F	Not applicable

INSPECTION, MAINTENANCE, TESTING AND LUBRICATION

- Carefully inspect each instrument to ensure that all visible contamination has been removed. If contamination is present, repeat the cleaning/disinfection process.
- Visually check for completeness, damage and/or excessive wear.
- Check that moving parts are in good working order. (e.g. joints, hinges, sliding parts, etc.).
- Hinged, rotating or articulated instruments must be lubricated with a water-soluble product (e.g. medical paraffin oil or equivalent lubricant) intended for surgical instruments to be sterilised.

Note: Do not use mineral oil or silicone lubricants as they :

- **Coat micro-organisms ;**
 - **Prevent direct surface-vapour contact;**
 - **Are difficult to remove.**
- Check that elongated instruments are not deformed.
 - When instruments are to be assembled, check that they fit correctly: It is not recommended to assemble instruments or components from different suppliers for reasons of compatibility, mechanics and design. COUSIN BIOTECH declines all responsibility in the event of incompatibility when using instruments or components from different sources.

PACKAGING

Individually packed instrument sets :

- Cleaned and disinfected instruments should be packed in their unassembled state.
- When packaging individual instruments, use single-use, medical-grade bags suitable for steam sterilisation, of an appropriate size, and double-wrap the instruments.

Pre-configured instrument storage containers should be used under the following conditions :

- Cleaned and disinfected instruments must be packed in their disassembled state.
- When storing instruments in sterilisation baskets with lids, use medical grade paper suitable for steam sterilisation and wrap the baskets using the double-wrap method. Baskets with lids can also be placed in an approved sterilisation container, fitted with a sealing lid for sterilisation.
- Wrapping and packaging must comply with EN ISO 11607, and be suitable for steam sterilisation (temperature resistance to at least 138°C (280°F), with sufficient steam permeability).
- Areas designated for specific instruments must contain only instruments specifically intended for these areas.
- Only instruments manufactured and/or distributed by COUSIN BIOTECH should be included in COUSIN BIOTECH instrument containers.

STERILISATION INSTRUCTIONS

Sterilisation is carried out under the responsibility of healthcare centres.

Refer to the table below for the minimum recommended sterilisation parameters which have been validated by COUSIN BIOTECH to provide a Sterility Assurance Level (SAL) of 10⁻⁶.

AUTOCLAVE sterilisation is recommended according to the following protocol :

Steam sterilization parameters	Values	Values
Cycle Type	Pre-vacuum	Pre-vacuum
Set Point Temperature	134°C / 273°F	132°C / 269.6°F
Cycle Time	Minimum 3 minutes	Minimum 4 minutes
Dry Time	Minimum 20 minutes	Minimum 20 minutes
Cool Time (inside and outside the autoclave chambers)	Minimum 40 minutes	Minimum 40 minutes

Note: the steriliser manufacturer's instructions for operation and load configuration must be followed explicitly.

- The centre is responsible for internal procedures for the reassembly, inspection and packaging of instruments once they have been thoroughly cleaned in a manner which ensures sterilising steam penetration and appropriate drying. Arrangements for the protection of sharp or potentially dangerous edges of instruments should also be recommended by the health care organisation. In any case, due to the differences between sterilisation equipment and the microbial load of the equipment when used in the clinic, it is generally the responsibility of each user to check and validate that the sterilisation cycle used will provide the correct level of assurance of appropriate instrument sterilisation.
- Steam/moist heat sterilisation is the recommended method for spinal and orthopaedic surgical instruments.
- Instruments must be correctly prepared and packed in COUSIN BIOTECH containers designed to allow steam to penetrate and come into direct contact with all surfaces.
- Ethylene oxide or plasma sterilisation methods must not be used.

Note: COUSIN BIOTECH recommends the use of its dedicated containers for instrument sterilisation, as they have been designed to hold each instrument in the right place for easy handling in the operating room.

STORAGE INSTRUCTIONS

- Containers of sterile and/or double-bagged instruments should be stored in a designated restricted area which is well ventilated and protected from dust, moisture, insects and extremes of temperature/humidity.
- Sterile instrument packs should be carefully examined prior to opening to ensure that the integrity of the packaging has not been compromised.

Note: maintaining the integrity of the sterile packaging is generally event-related. If a sterile pack is torn, perforated, shows signs of tampering or has been exposed to moisture, the instrument set must be repackaged and sterilised.

RESPONSIBILITIES OF CENTRES FOR INSTRUMENTS ON DEPOSIT

- Spinal and orthopaedic surgery instruments generally have a long service life; however, improper handling or inadequate protection can quickly shorten it. Instruments that no longer perform well due to long use, mishandling or improper maintenance should be returned to COUSIN BIOTECH. Report any problems with the instruments to COUSIN BIOTECH.

- The instructions provided in this manual have been validated by COUSIN BIOTECH in the laboratory and are suitable for the preparation of spinal and orthopaedic surgical instruments for use. It is the responsibility of the centre to ensure that reprocessing is carried out using the appropriate equipment and materials and that reprocessing staff are properly trained to achieve the desired result. Equipment and processes must be validated and monitored regularly. Any deviation from these instructions by the processor must be properly evaluated in terms of effectiveness in order to avoid any potential adverse consequences.

REQUESTS FOR INFORMATION AND COMPLAINTS

In accordance with its quality policy, COUSIN BIOTECH undertakes to do its utmost to produce and supply a quality medical device. If, however, a healthcare professional (customer, user, prescriber, etc.) has a complaint or reason for dissatisfaction concerning a product, in terms of quality, safety or performance, they should inform COUSIN BIOTECH as soon as possible. If an implant malfunctions or has contributed to a serious adverse effect for the patient, the care centre should follow the legal procedures in force in its country and inform COUSIN BIOTECH as soon as possible.

For all correspondence, please specify the reference, batch number, contact details of a referee, as well as a full description of the incident or complaint.

Brochures, documentation and operating techniques are available on request from COUSIN BIOTECH and its distributors.

REFERENCES

AAMI TIR 30 - Technical Information report: a compendium of processes, materials, test methods, and acceptance criteria for cleaning reusable medical devices

ANSI /AAMI ST79 and A1 - Steam sterilization and sterility assurance in health care facilities

AAMI TIR 12 - Designing, testing, and labeling reusable medical devices for reprocessing in health care facilities: A guide for medical device manufacturers

ISO 17664 -1 Sterilization of medical devices - Informations to supply by the manufacturer for the process of re-sterilization of medical devices

NF EN ISO 17 665-1 Sterilization of health care products — Moist heat — Part 1: Requirements for the development, validation and routine control of a sterilization process for medical devices

NOMENCLATURE

Family of product – Instruments for IntraSPINE	Product Codes	Family of product – Instruments for BDyn	Product Codes
Non Sterile Instruments for Devices for Interspinous Space with Laminar Support	ANCINTRALF	Non Sterile Instrument for Spinal Posterior Stabilization Devices	RCBANDH10U
	ANCINTRARG		RCBANSOUNU
	FANINTRA08		RCBANRH50U
	FANINTRA10		RCBANLPERU
	FANINTRA12		RCBANTPERU
	FANINTRA14		RCBANLRBEU
	FANINTRA16		RCBANDISTU
	RCBINTF08U		RCBANCOMPU
	RCBINTF10U		RCBANTD50U
	RCBINTF12U		RCBANBDHOU
	RCBINTF14U		RCBANRODPU
	RCBINTF16U		BDYNCONT00
	IMPACINTRA		RCBKITDY50
	PCALEINTRA		014112C
	PRTFAINTRA		013615
	RCBINTPFAU		011300
	FBIOMC2003		911537
	PINCEINTRA		013001N
	RCBANTENDU		010201C
	ROGNINTRAL		010494
	INTRACONT5		010505C
	INTRACONT6		011308
	RCBIS2BOXU		010001N
	RCBININTRU		013616
RCBISNORGU	010610C		
RCBINBIOIU	016103		
RCBINBIO2U			
RCBINSYSTU			

Family of product – Instruments for BatSPINE	Product Codes	
Non Sterile Instruments for Devices for Mini Invasive Interspinous Devices	RCBANSL05U RCBANSL07U RCBANSL09U RCBANLRAMU RCBANSRAMU RCBANSIHAU RCBANSLHAU RCBANS08U RCBANS10U RCBANS12U RCBANS14U RCBANS16U RCBANSOBTU RCBANLOBTU RCBANSUINU RCBANSLINU	RCBANLUINU RCBANLLINU RCBANSIMVU RCBANLIMVU RCBANARM1U RCBANARM2U RCBANBCOMU RCBANDILAU RCBANCHECU RCBANCURAU RCBANBGUNU RCBANPOSIU RCBANEBAAU RCBANECOMU RCBANBCONU

SURGICAL TECHNIQUE

IntraSPINE : https://www.cousin-biotech.com/d-6655aed293402_tech-op-intra-fr.pdf

BDyn : https://www.cousin-biotech.com/d-6655ae349037d_tech-op-bdynfr.pdf

BatSPINE : https://www.cousin-biotech.com/d-6657180a4a3d3_draft-tech-op-batspine15052024.pdf

UNDESIRABLE SIDE EFFECTS*

Undesirable side effects probably related to the device (which may be related to the implant or device):

- Pain, discomfort, numbness
- Bone/spine fracture(s)
- Infection / collection
- Dural tear

*only applicable to Instruments for IntraSPINE and Instruments for BDyn.

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