SEMI RESORBABLE PARIETAL REINFORCEMENT IMPLANT EXTRAPERITONEAL

en	Instructions for use	Page	2	MDR CE marking date: 30.09.2025
fr	Notice d'instructions	Page	5	Date de marquage CE MDR 30.09.2025
de	Gebrauchsanweisung	Seite	8	MDR CE-Kennzeichnungsdatum: 30.09.2025
it	Istruzioni per l'uso	Pagina	11	Data di marcatura CE MDR: 30.09.2025
es	Instrucciones de uso	Página	14	Fecha de marcado CE MDR: 30.09.2025
pt	Instruções de utilização	Página	17	Data de marcação CE MDR: 30.09.2025
ελ	Οδηγίες χρήσης	Σελίδα	20	Ημερομηνία σήμανσης CE MDR: 30.09.2025
bg	Инструкции за употреба	страница	24	Дата на маркировка CE по MDR: 30.09.2025
cs	Návod k použití	Strana	27	Datum označení CE podle MDR: 30.09.2025
da	Brugsvejledning	side	30	MDR CE-mærkningsdato: 30.09.2025
fi	Käyttöohjeet	Sivu	33	MDR CE-merkinnän päivämäärä: 30.09.2025
nl	Gebruiksaanwijzing	Pagina	36	MDR CE-markeringdatum: 30-09-2025
pl	Instrukcja obsługi	Strona	39	Data oznakowania CE MDR: 30.09.2025
ro	Instrucțiuni de utilizare	Pagină	42	Data marcajului CE MDR: 30.09.2025
sk	Návod na použitie	Strana	45	Dátum označenia CE podľa MDR: 30.09.2025
lt	Naudojimo instrukcijos	Puslapis	48	MDR CE žymėjimo data: 2025-09-30
no	Bruksanvisninger	Side	51	MDR CE merkingsdato: 30.09.2025
tr	Kullanım talimatları	Sayfa	54	MDR CE işareti tarihi: 30.09.2025
hu	Használati utasítás	Oldal	57	Az MDR szerinti CE-jelölés dátuma: 2025. szeptember 30.
lt	Naudojimo instrukcija	Puslapis	60	MDR CE ženkinimo data: 2025 09 30
en	Patient information leaflet (Australian requirements)	Page	63	

Rx
only

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Caution : Federal law (USA) restricts this device to sale,
distribution and use by or on the order of a physician.

NOT302_260518
Version du 18/05/2026

This release is the last update of the instructions of use and replace the previous edition

Product under regulation 2017/745 (EU).

**SEMI RESORBABLE PARIETAL REINFORCEMENT IMPLANT
EXTRAPERITONEAL
SINGLE USE STERILE PRODUCT**

MANUFACTURER IDENTIFICATION

COUSIN BIOTECH
ALLEE DES ROSES
59117 WERVICQ SUD
FRANCE

IDENTIFICATION AND DESCRIPTION OF THE DEVICE

4D Ventral® meshes are extraperitoneal semi-resorbable parietal reinforcement implants.

INTENDED USE

Intended to be used as parietal reinforcement implant.

INDICATIONS

Repair of ventral hernia that require the addition of an extraperitoneal reinforcing or bridging material to obtain the desired surgical result.

CONTRAINDICATIONS

Do not use in the following cases:

- Allergy to any of the components
- Infected site
- Pregnancy
- Growing children

UNDESIRABLE SIDE EFFECTS

Like any implantable medical device, this implant is susceptible to generate possible undesirable side effects that may lead to reoperation:

- Discomfort/Pain
- Infection
- Inflammation
- Fistula
- Recurrence
- Adhesion
- Erosion/Extrusion
- Seroma/lymphocele
- Hematoma
- Mesh migration
- Mesh deformation
- Allergic reaction
- Foreign body reaction
- Irritation nearby organ
- Obstruction

TARGET POPULATION

Adult male or female suffering from a symptomatic ventral hernia necessitating a surgical repair with a method that includes an extraperitoneal implantation.

IMPLANTED MATERIALS

Polypropylene (PP) (non-resorbable) – Poly L Lactic Acid (PLLA) (resorbable)

	weight PP (g)	weight PLLA (g)
4DVENT05RO	0,12	0,18
4DVENT07RO	0,24	0,36
4DVENT09RO	0,39	0,59
4DVENT12RO	0,70	1,05
4DVENT0715	0,65	0,97
4DVENT1015	0,92	1,39
4DVENT1515	1,39	2,08
4DVENT1520	1,85	2,78
4DVENT1530	2,78	4,18
4DVENT2020	2,47	3,71
4DVENT2025	3,09	4,64
4DVENT2535	5,42	8,13
4DVENT3030	5,57	8,36
4DVENT3040	7,43	11,15

Origin neither human nor animal – Semi resorbable.

PERFORMANCES

The device 4D Ventral® contributes to reduce symptoms from hernia and contributes to avoid complications of hernia.

LINK TO SUMMARY OF SAFETY AND PERFORMANCE CHARACTERISTICS:

https://www.cousin-biotech.com/d-68f1fd7631392_tf07csscpissue-4.pdf

MRI SAFETY

The implant has not been tested for safety and compatibility in the case of postoperative investigation using Magnetic Resonance Imaging (MRI). However, the implant is mainly made with Polypropylene and Poly-L-Lactic Acid: non metallic components. Consequently, 4D Ventral® is considered MR safe.

PRECAUTIONS FOR USE

4D Ventral® must always be separated from abdominal cavity by peritoneum

4D Ventral® implants are delivered sterile (ethylene oxide sterilisation).

When opening the carton box, the surgeon has to make sure that the lot number, the model and the size of the mesh are equal to the data on the label of the inner pouch. Otherwise, the use of the prosthesis is strictly prohibited.

Before any use, inspect the integrity of packaging and device (of which peelable pouches). Do not use in the event of deterioration of the device and/or the packaging. Do not use if the device is out of date. COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned. Mesh fixation points should be at least 1cm from the edge of the mesh with 1cm spacing between fixation points.

Unpacking and manipulation of the implant must be carried out following asepsis standards in order to guarantee the sterility.

Prevent the mesh from any contact with objects that may alter their surface.

This device must be implanted by a qualified surgeon (familiar with the relevant anatomy and experienced in visceral surgery).

IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE

As specified on the product's labelling, the 4D Ventral® mesh is for single use only. It can not be re-used and/or re-sterilized (potential risks would be and are not limited to: loss of the product's sterility, risk of infection, loss of the product's efficiency, recurrence).

INSTRUCTIONS TO BE GIVEN TO THE PATIENTS BY THE SURGEON

The surgeon should inform the patient of the potential physical and psychological restrictions and consequences of implanting the device. The patient must be informed of the surgical risks and possible side effects. The surgeon should invite the patient to return for further consultation if the patient presents symptoms that appear abnormal. An implant card for the patient is provided by the healthcare professional at the clinic. This card provides information to identify the device and traceability elements as well as the name, address and website of the manufacturer. Please refer to the section "Implant card information leaflet" for instructions on how to complete it. **The surgeon shall invite the patient to scan the implant card immediately upon receipt to keep track of it in case of loss.** The summary of the safety and performance characteristics of the device and the instructions for use can also be found on the COUSIN BIOTECH website.

STORAGE PRECAUTIONS

To be stored in a dry place away from sunlight and at room temperature in its original packaging.



SURGICAL TECHNIQUE

Extra-Peritoneal Hernia repair

Example: Retrorectus mesh repair of a ventral hernia (Rives-Stoppa technique).

This technique utilizes the hernia sac to separate the mesh from the intra-abdominal contents.

- Superior to the umbilicus, dissection is performed above the posterior rectus fascia and under the rectus muscle.
- Below the umbilicus, dissection occurs in the preperitoneal space due to the lack of a posterior rectus sheath.
- The 4D Ventral® implant is placed in the newly formed pre-peritoneal space. The 4DVentral mesh has an oriented elasticity:

- Fixation of the mesh if needed, to the muscle layer above
- The peritoneum must be meticulously checked and no defect in the peritoneum should be left open to prevent any contact between the mesh and intra-abdominal contents.
- Redon's drains under aspiration are usually placed. The anterior sheath is closed. If tension is present, relaxing incisions may be used
- Skin closure.
- Using the Redon's drains, fill the pre-peritoneal space with an antiseptic solution and clamp the drains for 1 hour.
- A belt can be used to wrap the patient's abdomen for a few days to reinforce the abdominal wall.

EXPLANTATION AND DISPOSAL OF DEVICES

Risks associated with the explantation: the surgeon who decides to remove the device should consider factors such as the risk of another surgery for the patient, and the difficulty of the explantation procedure. The removal of the implant is up to the surgeon's decision and must be subject to adequate postoperative follow-up.

Devices should be retrieved and handled according to the recommendations of ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling". Any explanted device must be sent back for analysis, following the current protocol. This protocol is available on request from COUSIN BIOTECH. It is important to note that any implant that must not be cleaned or disinfected prior to despatch must be contained in sealed packaging. The removed medical device must be disposed of in accordance with the country's standards for the disposal of infectious waste.

There are no specific recommendations for the disposal of a non-implanted device.

IMPORTANT

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

INFORMATION REQUEST AND COMPLAINTS

Following its quality policy, COUSIN BIOTECH is committed to making every effort to produce and supply a high-quality medical device. However, if a health professional (client, user, prescriber, etc.) has a complaint or cause for dissatisfaction with a product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In case of failure of an implant or if it has contributed to causing a serious adverse reaction in the patient, the health centre must follow the legal procedures in that country and inform COUSIN BIOTECH immediately. For any correspondence, please specify the reference, batch number, contact person's details and a comprehensive description of the incident or complaint.

Brochures, documentation and surgical technique are available on request from COUSIN BIOTECH and its distributors. If further information is needed or required, please contact your COUSIN BIOTECH representative or distributor.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

IMPLANT CARD INFORMATION LEAFLET

How to complete the implant card ?

The implant card, for the patient, should be completed as follows:

Front side of the implant card

The user has to complete the implant card before transmitting the implant card to the patient. Please complete it as follows:

	Carte Implant - Implant Card A conserver / Keep it	Carte Implant - Implant Card
Patient name	 _____	 COUSIN SURGERY <i>We care for Surgery</i>
Date of surgical procedure	 . . / . . /	
Name of practitioner	 _____	

Back side of the implant card

No fields need to be completed by the user. The following information on the device is printed directly by the manufacturer on the back of the implant card. A QR code allowing the patient to consult the instruction for use, as well as the website address, are available on the implant card.

Information to identify the device	Traceability information
• Indication that the product is a medical device • Description of the device family • Trade name of the device • Reference of the device • Device is MR safe • Denomination name of the device in sales countries	• UDI number • Batch number • Name and address of the legal manufacturer

Explanation of the used symbols in the Implant Card

All symbols are commonly used and explained at the end of the instruction for use.

**IMPLANT DE RENFORCEMENT PARIETAL SEMI RESORBABLE
EXTRAPERITONEAL
PRODUIT STERILE A USAGE UNIQUE**

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Patient information leaflet 4DVENTRAL®

Item	Patient information leaflet
1	a. The name of the device : 4DVENTRAL® b. The model of the device : 4DVENT05RO, 4DVENT07RO, 4DVENT09RO, 4DVENT12RO, 4DVENT0715, 4DVENT1015, 4DVENT1515, 4DVENT1520, 4DVENT1530, 4DVENT2020, 4DVENT2025, 4DVENT2535, 4DVENT3030, 4DVENT3050 and 4DVENT4545
2	a. The intended purpose of the device : Repair and extraperitoneal parietal reinforcement of ventral hernias. b. The kind of patient on whom the device is intended to be used : Adult male or female suffering from a symptomatic ventral hernia necessitating a surgical mesh repair with an extraperitoneal implantation method.
3	Special operating instructions for the use of the device : This device must be implanted only by a qualified surgeon trained on the use of the product (knowledge of anatomy and visceral surgery). The 4DVENTRAL® implant is designed for extraperitoneal implantation only.
4	a. The intended performance of the device : 4DVENTRAL® are designed for parietal reinforcement and hernia repair by open or laparoscopic approach. Prostheses are used as hernia defect coverings/patches in conjunction with surgical repair to help stop visceral contents from passing through the hernia orifice. 4DVENTRAL® intend to cause some degree of reactional fibrosis. Semi-resorbable meshes are proposed as a potential solution which combines the advantages of an absorbable mesh to reduce inflammation, and of a non-absorbable mesh to provide long term strength. b. Undesirable side effects : Like any implantable medical device, this implant is susceptible to generate possible undesirable side effects that may lead to reoperation: Discomfort/Pain, Recurrence, Adhesions, Erosions, Irritation of nearby organ, Infection, Inflammation, Fistula formation, Seroma, Mesh deformation, Hematoma, Mesh migration, Allergic reaction, Foreign body reaction. Seek medical advice if you experience any of these symptoms.
5	When 4DVENTRAL® Should Not be Used : The following are contraindications for 4DVENTRAL®: Allergy to any of the components, Infected site, Pregnancy, Growing children, Anticoagulant therapy

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Item	Patient information leaflets
6	MRI Information : 4DVENTRAL® are non-conducting and non-magnetic devices. They present no risk in Magnetic Resonance Imaging (MRI) environment.
7	Precautions : After surgery with 4DVENTRAL®, please thoroughly follow instructions given by your surgeon. The 4DVENTRAL® mesh is designed for permanent implantation. If you feel persistent pain or any undesirable effect, you should contact a healthcare professional related to your operation with the 4DVENTRAL device.
8	Materials included in the device : 4DVENTRAL® implants are extraperitoneal partially resorbable parietal reinforcement implants. This surgical mesh is composed of Poly-L-Lactic Acid (PLLA) and Polypropylene (PP)
9	Reporting a serious incident : If you encounter any serious incident with your 4DVENTRAL® implant you should report it immediately to your doctor. You should also report the incident to Cousin Biotech via the link serviceclients@cousin-surgery.com or materiovigilance@cousin-biotech.com and to the Therapeutic Goods Administration using the following website: http://www.tga.gov
10	Legal Manufacturer : Cousin Biotech, Allée des roses, 59117 Wervicq-Sud, France Manufacturing Facility : Cousin Biotech, Allée des roses, 59117 Wervicq-Sud, France

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Cousin Biotech SAS au capital de 320 000 € - NAF 1396Z- N° TVA FR 34 398 460 261 - Siren 398 460 261 RCS Lille