C DISC

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fr	<u>Notices d'instructions</u>	Page	2	<i>Date de marquage CE MDR : En attente de marquage CE</i>
en	<u>Instructions for use</u>	Page	10	<i>Date of CE marking MDR: Waiting for CE marking</i>

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IDENTIFICATION AND DESCRIPTION OF THE DEVICE

The C DISC sterile cervical disc device is an implantable medical device consisting of three elements: a lower plate, an upper plate and a mobile core positioned between these two plates.

These elements are supplied pre-assembled on a grasping forceps to facilitate the handling and placement of the device during surgery.

The C DISC sterile cervical disc device is manufactured using controlled manufacturing processes. The mechanical and functional characteristics of the device can only be guaranteed if the instructions described in this manual are followed.

INDICATIONS

The indications for the C DISC sterile cervical disc device are symptomatic degenerative disc disease of the cervical spine from C3 to T1 in adults, which may be accompanied by radiculopathy, moderate myelopathy, functional/neurological deficit in the arm and/or neck, at one or more levels (maximum 3 levels), with at least one of the following pathologies:

- Disc herniation,
- Disc protrusion with root canal stenosis,
- Mild to moderate spondylarthrosis

The mobility of the level to be operated on must be assessed beforehand.

CONTRAINDICATIONS

The implantation of a C DISC sterile cervical disc device is contraindicated in the following cases:

Contraindications related to the pathology

- Severe cervi-osteoarthritis or ossification myelopathy (OPLL),
- Severe facet osteoarthritis,
- Predominant posterior compression spinal canal stenosis,
- Spondylolisthesis, retrolisthesis,
- Significant segmental instability,
- Posture of the cervical spine in segmental kyphosis,
- Fracture,
- Surgical history of laminectomy and/or fusion by posterior approach of the segment concerned.

Patient-related contraindications

- Osteoporosis or severe osteopenia,
- Metabolic Bone Diseases,
- Systemic Inflammatory Osteoarticular Diseases,
- Rheumatoid arthritis,
- Bone tumour in the affected area,
- Any acute, chronic or local infection,
- Severe alcoholism,
- Pregnancy,
- Morbid (class III) obesity,
- Lack of patient cooperation,
- Senility,

Implant-related contraindications

- Allergy to at least one of the implant materials,

ADVERSE REACTIONS

The list of side effects listed below is not exhaustive.

Potential or actual adverse effects of cervical disc prostheses are:

- Migration of one of the components or implant,
- Insertion of prosthetic plates into the vertebral bodies,
- Loss of mobility,
- Failure of the medical device (wear, breakage, disassembly, etc.),
- Inflammation, Tissue reaction to prosthesis materials
- Fracture of bone structures near the medical device
- Allergies to the materials of the implant

Potential or actual adverse effects related to the surgical procedure are:

- Dysphagia, hoarseness/dysfunction of the vocal cords,
- Visceral wound (esophagus, airway, etc.),
- Airway obstruction,
- Pneumonia, Atelectasis,
- Phlebitis, Thrombotic accidents
- Cerebrovascular accident (CVA),
- Myocardial infarction,
- Haemorrhage requiring a blood transfusion,
- Hematoma and/or poor healing,
- Infection, Sepsis,
- Death,
- Dysesthesia, numbness, Paresthesia,
- Damage to neurovascular structures,
- Soft tissue injury, spinal cord injury, nerve root injury,
- Dural tear, cerebrospinal-fluid leak, compressive anterior meningocele,
- Change in spinal alignment, sagging, spinal instability,
- Implant loosening and migration leading to instability
- Increase of operating time, Delayed surgery
- Ectopic ossification

Other potential or actual adverse reactions:

- Adjacent level degeneration/disease
- Heterotopic ossification of the cervical spine
- Persistence or postoperative recurrence of pain or functional/neurological deficit,
- Re-intervention due to one of the set of effects listed above.

TARGET POPULATION

Adult with a mature skeleton.

IMPLANTED MATERIALS

Main components:

- The inferior endplate is made of cobalt-chromium-molybdenum alloy (CoCrMo – NF ISO 5832-12 standard). Its external surface is coated with a non-alloy titanium powder sprayed by plasma torch (Ti – NF ISO 13179-1 and ASTM F1580 standards).
- The superior endplate is made of cobalt-chromium-molybdenum alloy (CoCrMo – NF ISO 5832-12 standard). Its external surface is coated with a non-alloy titanium powder sprayed by plasma torch (Ti – NF ISO 13179-1 and ASTM F1580 standards).
- The mobile nucleus is made of very high molecular weight polyethylene (UHMWPE) (NF ISO 5834-1, NF ISO 5834-2 and ASTM F648 standards).


Non-implantable component:

- The grasping forceps on which the three prosthetic elements are pre-assembled is made of titanium alloy TiAl6V4 (NF EN ISO 5832-3 and ASTM F136 standards).

The C DISC sterile cervical disc device is designed for cement-free implantation.

Reference	1900	19100	19010	19020	19110	19260
Mass (g)	51	52	52	52,8	53,3	53,5

Reference	19120	19270	19130	19280	19290
Mass (g)	55,6	56	57,5	58	60,7

ACTIONS AND PERFORMANCE

Cervical disc arthroplasty consists of replacing one or more degenerative intervertebral discs anteriorly. Its objective is to restore a disc height of the level considered, while offering a degree of functional mobility compatible with the physiological kinematics of the spinal segment. It also contributes to the maintenance of reduced or absent pain following satisfactory decompression.

The purpose of cervical disc arthroplasty is to:

- Restore the mobility of the operated floor, with a minimum mobility of 3° in flexion-extension for at least 80% of patients at 2 years of age;
- Decrease pain according to the VAS (Visual Analog Scale) score by at least 20 points at 2 years for neck pain and radicular pain (Arm Pain);
- And/or reduce pain according to the NDI (Neck Disability Index) score by at least 15% at 2 years.

LINK TO THE SECURITY FEATURES AND PERFORMANCE SUMMARY:

https://www.cousin-biotech.com/d-6a69a8401895b_tf17asscpissue-a.pdf

MRI SAFETY

The components of the C DISC prosthesis are non-ferromagnetic and have a geometry that does not generate induced currents. They are attached to the bone or implant, which reduces the risk of mobilization. Non-clinical trials and literature reviews have shown that the device does not exhibit hazardous heating or implant migration when examined under appropriate MRI conditions. The evaluation concluded that MRI scans up to 3 Tesla are conditionally compatible.

MRI examinations are authorized only under the following conditions:

- Static magnetic field of 3 Tesla or less;
- Maximum spatial gradient of the magnetic field of 720 Gauss.cm⁻¹ or less;
- Maximum Specific Absorption Rate (SAR) of 3 W.kg⁻¹ for 15 minutes of acquisition.

The MRI system should be used in accordance with standard clinical practice and should not be used outside of the conditions specified above.

- The patient should be positioned according to standard MRI procedures.
- Continuous monitoring of the patient is recommended during the examination. If the patient experiences an unusual sensation of warmth or discomfort, the examination should be discontinued immediately.

As a precaution, it is recommended:

- Avoid MRI scans within 48 hours of implant placement;
- Inform the staff in charge of the examination of the recent placement of the implant, if an MRI is essential.

It should be noted that any device with a high contrast with the biological medium can generate artifacts that can affect the quality of the imaging. It is therefore advisable to recommend that the patient warn the health professionals (radiologists and radiology technicians) of the presence of the implant before the examination.

PRECAUTIONS FOR USE

Pre-operative precautions

- The effective mobility of the floor considered for arthroplasty,
- The patient's activity or profession: if the patient has to lift weights or make significant muscular efforts, the resulting forces may compromise the biological fixation and/or the implantable device,
- Cognitive disorders, mental disorders or use of toxic substances (alcoholism, for example): these situations may lead the patient to not comply with certain precautions or limitations, which can compromise the integrity of the implant or lead to other complications,
- Certain degenerative bone diseases: the evolution of degenerative bone diseases at the time of implantation can reduce the lifespan of the device, especially in the case of significant osteoporosis or osteomalacia.

Intraoperative precautions

- Reusable instruments for the implantation of the C DISC prosthesis must be cleaned and sterilized before use, in accordance with the instruction leaflet accompanying the instruments.
- The C DISC prosthesis is delivered pre-assembled and sterile on a grasping forceps specially designed for this purpose. Instructions for the handling and use of this assembly are described in the surgical technique.
- The surgical environment and the handling of the instruments and device must maintain a sterile state. The implant should be handled with care, without shock or impact, and should not meet metal surfaces or any surface that could damage its surfaces. Implants with scratches or alterations should not be used,
- The titanium powder-based coating must not meet a cloth or any other material that may contaminate it,
- Medical devices must not undergo any treatment or modification,
- When preparing the disc space, introducing the jigs, and implanting the device, any stress, vibration, shock, or impaction applied to the cervical segment should be moderate. The accessibility of the disc space must be sufficient, thanks to the distraction instruments provided,
- Before releasing the implant from the grasping forceps, a tabletop X-ray should be taken to confirm the correct positioning of the implant.

Post-operative precautions

- The patient should be informed of the limitations imposed by the C DISC prosthesis and warn him of the authorized level of activity,
- Sports with high energy or exposing the head to shocks, such as combat sports or acrobatics, are not recommended, as they can cause wear, dislocation or migration of the implant,
- Post-operative care must be organized to maintain joint mobility and prevent any dislocation,
- Any patient with the C DISC prosthesis must be monitored by an orthopaedic surgeon or neurosurgeon to detect signs of wear, dislocation or migration before the appearance of functional symptoms,
- Regular radiographic examinations should be performed and compared to post-operative images to detect any displacement, component fracture, bone disassembly or remodeling,
- The patient should be encouraged to report any unusual changes in the operated joint immediately.

INTENDED USER

The C DISC STERILE CERVICAL DISC DEVICE is intended exclusively for use by a surgeon (orthopaedic surgeon or neurosurgeon).

WARNINGS

Pre-operative precautions

- A thorough knowledge of the surgical technique of cervical disc prostheses is essential to ensure the correct placement of the C DISC cervical disc prosthesis.
- The appropriate use of instruments, especially those for determining the size and height of the implant (width and height templates) as well as for implanting the prosthesis (grasping forceps), is

imperative. Damaged or deformed instruments can lead to poor implant positioning and compromise the success of the procedure.

- An adequate stock of implants of various sizes should be available during surgery, including larger and smaller sizes than expected. It is important to verify the integrity of all packaging and implants before use.
- It is necessary to discuss with the patient the aspects related to the surgery and the implant. This discussion should address the limitations of joint reconstruction, patient-specific restrictions, the possible consequences of these limitations, and the importance of following the physician's preoperative instructions.
- Although rare, any allergies or reactions to the materials that make up the implant should be considered and eliminated.
- Even if implants consist only of standardised inert materials for surgical implants, the potential interaction with other implantable devices or active drug substances must be considered.
- Implants should not be reused under any circumstances. If an implant is removed after implantation, it must be destroyed or returned to the manufacturer. Although a removed implant may appear intact, it may have small defects or signs of internal stresses compromising its integrity. Use only new implants and do not change their configuration before implantation.
- The grasping forceps should not be reused. After implantation, it must be destroyed or returned to the manufacturer. The C DISC prosthesis is delivered mounted on its grasping forceps.

Do not use prostheses whose assembly on the clamp appears to be defective.

Intraoperative precautions

- Discectomy, decompression and correct placement of the implant are decisive for clinical outcomes.
- Soft tissue retraction during the anterior cervical approach should be moderate and performed using an atraumatic retractor to avoid any damage to the vocal cords, or any visceral or vascular wound.
- The preparation of the disc space, in particular the vertebral plates by curettage or milling, should not be excessive. The subchondral bone, which is necessary for proper support of the prosthetic plates, must be preserved.
- The selection of the appropriate implant is necessary to respect the function of the operated floor. The pruning must be as covering as possible on the vertebral plates. The height of the implant should be chosen after estimation using the height templates, considering the height of the healthy adjacent discs.
- If the implant is not the optimal size or height, or if its primary stability is not ensured, migration, dislocation or fracture of the components may occur. It is particularly important to ensure good bone/prosthesis support conditions.
- COUSIN BIOTECH components and instruments should not be used with those of any other manufacturer. Dimensional compatibility cannot be guaranteed, so the use of components from other manufacturers may lead to unsatisfactory results.
- Before closing the surgical site, thoroughly clean the area to remove any bone fragments, ectopic tissue, or foreign particles. The presence of foreign objects at the metal/plastic interface can cause excessive wear or friction. Ectopic bone tissue can lead to dislocation or limited mobility. Check joint mobility for early contact or instability

PATIENT INFORMATION LEAFLET

See the patient information leaflet referenced by LPRCDRGB01.

INFORMATION COMMUNICATED BY THE SURGEON TO THE PATIENT

The surgeon should inform the patient of the potential physical and psychological restrictions and consequences of implanting the device. The patient must be informed of the surgical risks and possible side effects. The surgeon should invite the patient to return for further consultation if the patient presents symptoms that appear abnormal. An implant card for the patient is provided by the healthcare professional at the clinic. This card provides information to identify the device and traceability elements as well as the name, address and website of the manufacturer. Please refer to the section "Implant card information leaflet" for instructions on how to complete it. **The surgeon shall invite the patient to scan the implant card immediately upon receipt to keep track of it in case of loss.** The summary of the safety and performance characteristics of the device and the instructions for use can also be found on the COUSIN BIOTECH website.

IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE

PRECAUTIONS FOR STORING THE PROSTHESIS

- Prostheses must be kept in their original packaging.
- The various components are for single use. The prosthesis is packaged pre-assembled on its gripper and is delivered sterile.
- Dentures must not undergo any treatment, modification or intervention not authorized by the manufacturer.
- Dentures that are damaged or have been subjected to inappropriate treatment should not be implanted.
- Prostheses with damaged packaging should not be implanted.

SURGICAL TECHNIQUE

For the implantation of the device, please refer to the surgical technique TORCDRGB01_GB

EXPLANTATION AND DISPOSAL OF DEVICES

Explantation and disposal of devices

1. Decision to expel

The decision to remove the device is the responsibility of the surgeon, based on an assessment of the patient's clinical condition, imaging tests and the benefit-risk balance of a new procedure.

In particular, the surgeon should take into account the difficulty of explantation, the time since implantation, the degree of osteointegration, bone quality, the proximity of adjacent anatomical structures and the reconstruction strategy envisaged. Explantation should be performed by a surgeon experienced in cervical spine surgery and familiar with the device and its instrumentation.

2. Risks associated with explantation

Explantation may be complicated by osteointegration of the prosthetic trays and formation of bone around them.

Potential risks include:

- injury or fracture of the vertebral plates;
- loss of bone material;
- damage to neurological, vascular or soft tissue structures;
- bleeding, hematoma, infection, or cerebrospinal fluid leak;
- uncontrolled mobilization, deterioration or incomplete removal of the device;
- the need to use complementary instruments or to use arthrodesis or another reconstruction solution.

Excessive force or distraction should be avoided. After explantation, the surgeon should check that all components have been removed, assess the condition of the surrounding structures and arrange for appropriate post-operative follow-up.

3. Recovery and handling of explanted devices

The explanted devices should be retrieved and handled in accordance with ISO 12891 recommendations.

Any device intended for analysis should be returned to COUSIN BIOTECH in accordance with the current protocol, available upon request.

Unless otherwise instructed by COUSIN BIOTECH, the explanted device must not be cleaned, disinfected, sterilized, disassembled or modified prior to shipment. It must be stored with all of its components, identified, sealed and sealed, and transported in accordance with the requirements applicable to potentially contaminated materials.

4. Elimination

Where no return for analysis is required, the explanted device and any associated contaminated material should be disposed of in accordance with the procedures of the health facility and the applicable national and local requirements for infectious waste.

An explanted device should never be reused, reprocessed or reimplanted.

No specific recommendations apply to an unimplanted and uncontaminated device. Any device that has been in contact with the patient, body fluids, or the surgical site should be considered potentially contaminated and disposed of according to applicable procedures

IMPORTANT

For any further information relating to the use of this product, please contact your representative or distributor.

REQUESTS FOR INFORMATION AND COMPLAINTS

In accordance with its quality policy, COUSIN BIOTECH is committed to doing everything possible to produce and supply a quality medical device. If, however, a healthcare professional (customer, user, prescriber, etc.) has a complaint or a reason for dissatisfaction with a product, in terms of quality, safety or performance, he or she must inform COUSIN BIOTECH as soon as possible. In the event of an implant malfunction or if it has

contributed to causing a serious adverse reaction for the patient, the health center must follow the legal procedures in force in its country, and inform COUSIN BIOTECH as soon as possible.
For all correspondence, please specify the reference, the batch number, the contact details of a referent, as well as an exhaustive description of the incident or complaint.
Brochures, documentation and operating techniques are available on request from COUSIN BIOTECH and its distributors.

Any serious incident that occurs in connection with the device must be notified to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

IMPLANT CARD INFORMATION LEAFLET

How to complete the implant card?

The implant card, for the patient, should be completed as follows:

Front side of the implant card

The user has to complete the implant card before transmitting the implant card to the patient. Please complete it as follows:

Carte Implant - Implant Card		Carte Implant - Implant Card	
A conserver / Keep it			
Patient name	 _____	 COUSIN SURGERY <i>We care for Surgery</i>	
Date of surgical procedure	 . . / . . /		
Name of practitioner	 _____		

Back side of the implant card

No fields need to be completed by the user. The following information on the device is printed directly by the manufacturer on the back of the implant card. A QR code allowing the patient to consult the instruction for use, as well as the website address, are available on the implant card.

Information to identify the device	Information de traçabilité
• Indication that the product is a medical device • Description of the family device • Trade name of the device • Reference of the device • Device is MR Conditional • Denomination name of the device in sales countries	• UDI Number • Batch number • Name and address of the legal manufacturer

Explanation of the used symbols in the Implant Card

All symbols are commonly used and explained at the end of the instruction for use

SYMBOLS ON LABELS

Please refer to the end of the operating instruction for a description of all symbols on the labels.

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