STERILE NON-RESORBABLE ARTICULAR LIGAMENTS

en [Instructions for use](#)

Page 2

*MDR CE marking date: CE
pending*

fr [Notices
d'instructions](#)

Page 4

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CE
1639

NOT319_260826
Version du 26/08/2026

This release is the last update of the instructions of use and replace the previous edition

Product under regulation 2017/745 (EU).

STERILE NON-RESORBABLE ARTICULAR LIGAMENTS
STERILE PRODUCT FOR SINGLE USE

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IDENTIFICATION AND DESCRIPTION OF THE DEVICE

Sterile non resorbable articular ligaments are surgically reinforcing ligaments and tendons for orthopaedic surgery. The articular ligament device is composed of one braid made of polyethylene. There are different types of ligaments:

- Round braid set with stainless steel suture needles at each end
- Flat braid (tape) set with stainless steel suture needle at one end
- Round braid with a loop and set with stainless steel suture needle at the end
- Flat braid (tape) with a loop set with stainless steel suture needle at the end
- Round-flat braid (tape)

INDICATIONS

Sterile non resorbable articular ligaments are intended for the surgical reinforcement of instabilities in different joints. They can be located intra- or extra-articularly. The reinforcing ligament or tendinous should always be used to complement a suture or natural tissue graft.

CONTRAINDICATIONS

Sterile non resorbable articular ligament should not be used in following conditions:

- In a septic environment or in an infected site,
- Allergy to one of components
- Growing child
- Pregnant woman

UNDESIRABLE SIDE EFFECTS

- 1) Known undesirable effects concerning sterile non resorbable articular ligaments unlikely related to the device, but related to the surgery or patient condition include:
 - Routine surgical complications
 - Infection
 - Anaesthesia complications
 - Local hematoma
 - Blood loss
 - Inflammatory
 - Arthrofibrosis, postoperative loss of motion, frozen joint
 - Bone resorption (osteolysis)
 - Neurovascular injury, peripheral neuropraxia
 - Wound dehiscence
 - Loss of mobility after surgery, locked joint, ongoing instability episodes
- 2) Undesirable secondary effects related to the device include:
 - Reaction to a foreign body
 - Allergic reactions
 - Discomfort
 - Post-surgical pain, persistent pain
 - Ligament injury
 - Surgery failure, recurrence of instability and reintervention related to:
 - Breaking of the ligament
 - Permanent elongation of the ligament

TARGET POPULATION

Sterile non resorbable articular ligaments are intended for adult (or grown up) male or female suffering from joint instability and needing a surgical reinforcement.

IMPLANTED MATERIALS

- **Materials implanted:**

Polyethylene and Cobalt chromite blue green spinel (blue dye).

References	Weight of polyethylene (g)
OCBLIA050U	0.323
OCBLIA020U, OCBLIA020X	0.192
OCBLIB140U, OCBLIB140X	0.201
OCBLIB200U	0.390
OCBLIC000U	0.129
OCBLIC200U	0.106
OCBLIC020U	0.201
OCBLID200U	0.525
OCBLID140U	0.258
OCBLIE200U, OCBLIE200X	0.411
OCBLIE140U	0.202

Ligaments are set with suture needles made of stainless steel which are not implanted in the patient, they only have a transitory contact with the patient.

Neither human nor animal origin materials are contained in the device.

PERFORMANCES

Sterile non resorbable articular ligaments are designed for surgical reinforcement of instabilities in different joints and contributes to the success of the surgery and to the joint function recovery. These performances are demonstrated by mechanical and technical specifications, such as crimping resistance and breaking strength.

LINK TO SUMMARY OF SAFETY AND PERFORMANCE CHARACTERISTICS:

https://www.cousin-biotech.com/d-6318526d8b749_under-construction.pdf

MRI SAFETY

The implant has not been tested for safety and compatibility in the case of postoperative investigation using Magnetic Resonance Imaging (MRI). However, the implant is made with polyethylene. The needles are not implanted and only used for the implantation of the articular ligament. Consequently, the articular ligament is considered MR safe.

PRECAUTIONS FOR USE

Sterile non resorbable articular ligament is delivered sterile (ethylene oxide sterilization).

Before any use, inspect the integrity of packaging and device (of which peelable pouches). Do not use if the device is out of date.

Do not use in the event of deterioration of the device and/or the packaging.

Unpacking and manipulation of the implant must be carried out following asepsis standards in order to guarantee the sterility.

It is recommended that both the braid and the needle are pulled delicately, not just the needle.

When opening the carton box, the surgeon has to make sure that the lot number, the model and the size are equal to the data on the label of the inner pouch. Otherwise, the use of the prosthesis is strictly prohibited.

Prevent the sterile non resorbable articular ligament from any contact with objects that may alter their surface.

The articular ligaments must always be implanted by a qualified surgeon (familiar with the relevant anatomy and experienced in orthopaedic surgery) who has a good knowledge of the device, its intended use, the instruments, and the surgical technique.

The surgical technique for sterile non resorbable articular ligaments depends on the indication, the tissue involved, the patient's circumstances and other factors.

Do not tension the ligament by pulling on the needle, as this could uncrimping the needle.

COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned.

IMPORTANT : DO NOT REUSE – DO NOT RESTERILIZE

As specified on the product's labelling, sterile non resorbable articular ligament is for single use only. It cannot be re-used and/or re-sterilized (potential risks would be and are not limited to: loss of the product's sterility, risk of infection, loss of the product's efficiency, recurrence).

INSTRUCTIONS TO BE GIVEN TO THE PATIENTS BY THE SURGEON

The surgeon should inform the patient of the potential physical and psychological restrictions and consequences of implanting the device. The patient must be informed of the surgical risks and possible side effects. The patient must be informed of the limitations inherent in the nature and function of the implant, in particular regarding the weight and constraints placed by the body on the device before the bone has sufficiently healed. The surgeon should invite the patient to return for further consultation if the patient presents symptoms that appear abnormal. An implant card for the patient is provided by the healthcare professional at the clinic. This implant card provides information, for the patient, to identify the device and traceability elements as well as the name, address and website of the manufacturer. Please refer to the section "Implant card information leaflet" for instructions on how to complete it. **The surgeon shall invite the patient to scan the implant card immediately upon receipt to keep track of it in case of loss.** The surgeon shall also inform the patient that the summary of the safety and performance characteristics of the device and the instructions for use can be found on the COUSIN BIOTECH website.

STORAGE PRECAUTIONS

Store in a dry place away from sunlight, at room temperature and in the original packaging.

SURGICAL TECHNIQUE

The surgical technique for sterile non resorbable articular ligaments depends on the indication, the tissue involved, the patient's circumstances and other factors.

Do not tension the ligament by pulling on the needle, as this could uncrimping the needle.

EXPLANTATION AND DISPOSAL OF DEVICES

The surgeon who decides to remove the device should consider factors such as the risk of another surgery for the patient, and the difficulty of the explantation procedure. The removal of the implant is up to the surgeon's decision and must be subject to adequate postoperative follow-up.

Devices should be retrieved and handled according to the recommendations of ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling". Any explanted device must be sent back for analysis, following the current protocol. This protocol is available on request from COUSIN BIOTECH. It is important to note that any implant that is not cleaned or disinfected prior to despatch must be contained in sealed packaging. The removed medical device must be disposed of in accordance with the country's standards relative to the disposal of infectious waste.

There are no specific recommendations for the disposal of a non-implanted device.

IMPORTANT

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

INFORMATION REQUEST AND COMPLAINTS

Following its quality policy, COUSIN BIOTECH is committed to making every effort to produce and supply a high-quality medical device. However, if a health professional (client, user, prescriber, etc.) has a complaint or cause for dissatisfaction with a product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In case of failure of an implant or if it has contributed to causing a serious adverse reaction in the patient, the health centre must follow the legal procedures in that country and inform COUSIN BIOTECH immediately. For any correspondence, please specify the reference, batch number, contact person's details and a comprehensive description of the incident or complaint.

If further information is needed or required, please contact your COUSIN BIOTECH representative or distributor.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

IMPLANT CARD INFORMATION LEAFLET

How to complete the implant card ?

The implant card, for the patient, should be completed as follows:

Front side of the implant card

The user has to complete the implant card before transmitting the implant card to the patient. Please complete it as follows:

	Carte Implant - Implant Card A conserver / Keep it	Carte Implant - Implant Card
Patient name	 _____	 COUSIN SURGERY <i>We care for Surgery</i>
Date of surgical procedure	 . . / . . /	
Name of practitioner	 _____	

Back side of the implant card

No fields need to be completed by the user. The following information on the device is printed directly by the manufacturer on the back of the implant card. A QR code allowing the patient to consult the instruction for use, as well as the website address, are available on the implant card.

Information to identify the device	Traceability information
• Indication that the product is a medical device • Description of the device family • Trade name of the device • Reference of the device • Device is MR safe • Denomination name of the device in sales countries	• UDI number • Batch number • Name and address of the legal manufacturer

Explanation of the used symbols in the Implant Card

All symbols are commonly used and explained at the end of the instruction for use.

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